

Product Catalogue 2020

Manufacturing Medical Furniture & Equipment Since 1953











Welcome

Supplying Medical Furniture & Equipment Since 1953

Why Choose **Bristol Maid?**

With over 65 years of experience, Bristol Maid is a leading manufacturer and supplier of medical furniture & equipment. From design and manufacture to supply and installation, Bristol Maid has a reputation for providing high product quality, excellent customer service and value for money. We take pride in working in close partnership with our customers to deliver effective solutions based on innovative design and technology that will meet the demands of today while anticipating the solutions for tomorrow.

Our highly-trained staff are experienced in every aspect of customer support from simple product selection to bespoke technical solutions, from on-site audits to complex project management arrangements, according to our customers' needs.

Bristol Maid's commitment to quality is based on years of experience working with healthcare professionals in an exacting and regulated environment, and is supported by an ISO 9001:2015 guality management system that is registered with a UKAS accredited certification body.

If you require further assistance, please do not hesitate to contact our Customer Service team who are here to help; or alternatively can arrange for your Regional Sales Representative to visit and help find a solution that best suits your requirements.





















Family business

Hospital Metalcraft Ltd is a family company. We own our factory. production lines and delivery lorries. We continually invest into our staff and the development of new product ranges.

Healthcare focused

With an extensive portfolio that covers almost every discipline, we have solutions to cater for the demands of the healthcare environment.

The best quality

With over 65 years' experience, we are focussed on the design and manufacture of high quality products which will create value for our customers. With a worldwide customer base, the Bristol Maid brand is recognised by many thousands of satisfied customers.

Manufactured in the UK

Hospital Metalcraft is proud to be a British manufacturing success story. Strict quality criteria are applied to both the Bristol Maid product range designed and produced within our own manufacturing facility and complementary equipment from carefully selected partners.

Delivered using our own lorries

Wherever possible we use our own fleet of lorries to ensure that your equipment arrives safe and secure.

Approved supplier

We are an accredited NHS & public sector accredited supplier. This allows you to buy with confidence knowing your requirements for compliant and competitive tendering has already been met.

Support you can rely on

Our external sales representative, sales support engineers and our Internal Customer Service Team are here to help you at every stage, from advice with product selection, managing the delivery of your order through to full service and repair support.

How To Order

How To Order

Custom Design & Manufacture

Whether it is a trollev that is slightly lower, a work surface that is slightly larger, a cabinet that must fit into a specific space or a unique colour scheme on a bedside cabinet, we recognise that we may not always to be able to provide a product that meets your specific requirements.

If you require a unique product, please contact us so we can work together to develop a solution that will exactly match your needs. Our design capability allows us to provide everything from minor product specification changes through to a full design & development service.

Tenders & Contracts

Through many years' experience we can support you with the facilitation and project management of tenders and contracts

+44 (0)1258 484455 Sales@bristolmaid.com www.bristolmaid.com

Monday to Friday - 8:30am - 5:00pm O Blandford Heights, Blandford Forum Dorset, DT11 7TG, UK

- Site audits and equipment recommendations
- Scheduled deliveries
- Product placement
- Removal of packaging
- Installation
- Commissioning
- Training

NHS & Public Sector Accredited Supplier

Public sector organisations can take advantage of our nationally agreed discounts through the NHS Supply Chain and Shared Business Services. Buy with confidence knowing your requirement for compliant and competitive tendering has already been met.

NHS Supply Chain

Medical Furniture Equipment, Accessories, Sackholders & Bins - Ref: 2019/S 137-337185 Ward & Residential Furniture - Ref: 2017/S 207-427709

Shared Business Services

Design, Furniture & Appliances - Ref: SBS/19/MA/TYX/9363/23

How to Order

- manufacture & delivery)
- Your order should contain the following information...
- Item No.
- Description
- Quantity
- Delivery address
- Invoice address
- Full contact details, including name, email address, phone & fax numbers
- If applicable, please ensure our quotation reference is detailed on your order to ensure any agreed discounts are applied
- · Item codes & descriptions in this catalogue supersede all previous editions
- A full price list can be requested from our website
- · Payments can be made by cheque, bank transfer, credit or debit card (with the exception of American Express)
- a credit account may be applied for
- For our terms and conditions refer to www.bristolmaid.com

· Orders are required in writing, either by fax, email or post. Should you wish to order over the phone we would need to take full payment at the time of order by debit or credit card. Alternatively, we can send you a proforma invoice for payment (required to be paid in full prior to

Orders should be submitted as an 'official purchase order' on company headed paper with all information clearly displayed


Payment Options

BACS

- Payments to: HSBC Bank PLC
- Sort Code: 40-12-18
- Account Number: 11044052

Credit or Debit Card

Major credit cards accepted



• If you are a new customer, the first three orders will be processed on a proforma basis, paid in full prior to manufacture & delivery, after which

• For further information and assistance, please contact our customer services team on +44 (0)1258 484455





Monday to Friday - 8:30am - 5:00pm O Blandford Heights, Blandford Forum Dorset, DT11 7TG, UK

Bariatric Equipment

| Folding Wheelchairs | 8 |
|------------------------------|----|
| Portering Chairs | 9 |
| Mobile Shower/Commode Chairs | 11 |
| Bariatric Shower Chairs | 12 |
| Shower/Commode Chairs | 14 |
| Chairs | 15 |
| Rollators & Training Suits | 17 |

Bedside Cabinets & Overbed Tables

| Polymer Bedside Cabinets | 20 |
|---------------------------------------|----|
| Bedside Cabinets | 22 |
| Bedside Cabinet Options | 25 |
| Bedside Cabinet/Wardrobe Combinations | 26 |
| Economy Overbed Table | 30 |
| Polymer Overbed & Overchair Tables | 31 |
| Laminate Overbed & Overchair Tables | 32 |

Caretray Trolleys & Racks

| Features | 36 |
|--|----|
| Disassembly & Cleaning | 38 |
| Stainless Steel, Low Level Trolleys | 40 |
| Stainless Steel, Standard Level Trolleys | 42 |
| Stainless Steel, High Level Trolleys | 44 |
| Stainless Steel, Mobile Racks | 46 |
| Stainless Steel, Static Racks | 48 |
| Mild Steel, Low Level Trolleys | 50 |
| Mild Steel, Standard Level Trolleys | 52 |
| Mild Steel, Standard Level Trolleys | 54 |
| Mild Steel, Mobile Racks | 56 |
| Mild Steel, Static Racks | 58 |
| Caretray Trolley Options | 60 |

Chart Trolleys, Workstations Room Dividers & Screens

| Fixed Height Trolleys | 72 |
|------------------------------|----|
| Fixed Height Workstations | 74 |
| Variable Height Workstations | 75 |
| Chart Boards | 76 |
| Document Holders | 77 |
| Room Dividers | 78 |
| Curtain Screens | 80 |
| Solid Screens | 81 |

Concealment Trolleys

| Variable Height Trolleys | 84 |
|--------------------------|----|
| Fixed Height Trolleys | 85 |
| Body Trays & Covers | 85 |

Cylinder Storage & Transportation

| Cylinder Trolleys - F/G Size | 88 |
|--------------------------------------|----|
| Cylinder Trolleys - J Size | 89 |
| Bulk Transportation Trolleys | 89 |
| Cylinder Trolleys - E, CD/D, AF Size | 90 |
| Vall Brackets & Cylinder Holders | 91 |

Dispensers, Holders & Racks

| Apron Dispensers | 94 |
|----------------------|----|
| Pocket Mask Holders | 94 |
| Glove Box Holders | 95 |
| Alcohol Gel Holders | 96 |
| Syringe Dispenser | 96 |
| Sharps Box Holder | 97 |
| Urine Bottle Holder | 97 |
| Catheter Bag Holders | 97 |
| Catheter Racks | 98 |
| Drainage Racks | 99 |
| | |

Dressing / Instrument Trolleys

| Removable Shelves, Stainless Steel | 102 |
|------------------------------------|-----|
| Removable Shelves, Mild Steel | 103 |
| Fixed Shelves, Sides Up | 104 |
| Fixed Shelves, Sides Down | 105 |
| Drawer Units | 107 |
| Options | 108 |
| | |

Drug & Medicine Dispensing

| Dispensing Trolleys | 112 |
|---|-----|
| CLIQ [®] & RFID Solutions | 115 |
| Patient Administration Trolleys | 116 |
| Drawer Trolleys | 118 |
| Tray Trolleys | 119 |
| Blister Pack (MDS) Trolleys | 120 |
| Nomad [®] Cassette Trolleys | 122 |
| Biodose [®] & MultiMeds Trolleys | 124 |
| Medinoxx Trolleys | 125 |
| Laptop Holders | 126 |
| Security Clamp | 128 |
| Distribution Trolleys | 129 |
| | |

Drug & Medicine Storage

| Drug & Medicine Cabinets | 132 | |
|--|-----|--|
| Controlled Drug Cabinets | 137 | |
| CLIQ™ | 140 | |
| Drug & Medicine Cupboards | 142 | |
| Options | 144 | |
| Patients / Residents Own Medication Cabinets | 146 | |
| Tablet Counter & Storage Boxes | 150 | |
| Pharmacy & Vaccine Refrigerators | 151 | |

Examination / Treatment **Couches & Chairs**

| Examination & Treatment Couches | 154 |
|---|-----|
| Leg & Arm Rests | 162 |
| Couch Steps | 163 |
| Bariatric Examination & Treatment & Couches | 164 |
| Therapy Couches | 170 |
| Examination & Treatment Chairs | 174 |
| Bariatric Examination & Treatment Chairs | 176 |
| Phlebotomy Chairs | 178 |
| Options | 181 |

General Purpose Trolleys

| Stainless Steel Trolleys | 186 |
|--------------------------|-----|
| Polymer Trolleys | 188 |
| Mild Steel Trolleys | 189 |

HTM71 Storage Solutions

| ScanCell® Carts | 192 |
|--|-----|
| Distribution Trolleys | 197 |
| Open Transfer/Transport Trolleys | 199 |
| HTM71 Racking | 200 |
| Trays, Baskets, Dividers & Accessories | 202 |

Infusion Stands & Accessories

| Nobile Pump Stands | 208 |
|-------------------------------|-----|
| Nobile Infusion Stands | 210 |
| Clamping Infusion Stands | 215 |
| Vall Mounted Infusion Holders | 215 |
| | |

Linen Trolleys & Bags

| Soiled Linen Trolleys | 218 |
|-----------------------|-----|
| Linen Trolleys | 222 |
| Hamper Trolleys | 223 |
| Soiled Linen Bags | 223 |
| Clean Linen Trolleys | 224 |
| | |

Maternity Equipment

| Variable Height Baby Cribs | 230 |
|----------------------------|-----|
| Fixed Height Baby Cribs | 232 |
| Variable Height Baby Bath | 233 |

Medical Records

| Medical Records Trolleys | 23 |
|------------------------------------|----|
| Medical Record Workstations | 24 |
| Laptop/Motion C5 Holders | 24 |
| Carefile Filing Pockets | 24 |
| CLIQ [®] & RFID Solutions | 24 |
| Medical Record Filing Trolleys | 24 |
| Medical Records Transfer Trolleys | 24 |
| Medical Records Cabinet | 25 |
| | |

Phlebotomy

| Caretray Trolley - Phlebotomy | 254 |
|-------------------------------|-----|
| Phlebotomy Chairs | 256 |

Portering Chairs

| Portering Chair - Rear Steer | 262 |
|--|------------|
| Portering Chair - Front Steer | 263 |
| Portering Chair - Bariatric, Rear Steer Options | 265 266 |

Procedure Trolleys & Carts

| Plaster Trolleys | 270 | |
|---------------------|-----|--|
| Mobile Workstations | 271 | |
| Procedure Carts | 272 | |

Resuscitation Trolleys & Carts

| Resuscitation Trolleys | 280 |
|--|-----|
| Caretray Trolley - Resuscitation | 282 |
| Resuscitation Carts - Polymer | 284 |
| Airways / Difficult Intubation Trolley | 287 |

Sackholders & Waste Bins

| Removable Body Bins | 29 |
|------------------------------------|----|
| Free Standing Sackholders | 29 |
| Fixed Body Plastic Sackholders | 29 |
| Removable Body Plastic Sackholders | 29 |

Seating Solutions

| Office, Task & Operators Chairs | 300 |
|--|-----|
| 24 Hour Task & Operator Chairs | 302 |
| Bariatric Chairs | 304 |
| Saddle Stools | 306 |
| Low Back Patient Chairs | 308 |
| Medium & High Back Patient Armchairs | 309 |
| Reclining Armchairs | 310 |
| Bariatric Reclining Armchairs | 320 |
| Reclining Porters Armchairs | 322 |
| Rise & Recline Porters Armchairs | 324 |
| Seat, Back Cushion & Headrest Options | 326 |
| Overnight Chair | 328 |
| Soft Seating | 330 |
| Visitor & Waiting Room Seating | 332 |
| Waiting Room Seating | 334 |
| Tech Chairs & Stools | 340 |
| PU Tech Chairs & Stools | 342 |
| Static Safe Tech Chairs & Stools | 344 |
| Static Safe & Sterile Tech Chairs & Stools | 346 |
| Medical Grade Clean Room Tech Chairs | 348 |

Sepsis Trolleys

| Sepsis Trolleys | 352 |
|-----------------|-----|
| Options | 353 |

Shelving & Racking Solutions

| Polymer Racking - MetroMax I® | 360 |
|---|-----|
| Polymer Racking - MetroMax Q® | 362 |
| Zinc Chrome Racking - Olympic [®] | 364 |
| High Density Storage - Top-Track [®] | 366 |

Sterile Services Equipment

| Preparation Tables - Variable Height | 370 |
|--------------------------------------|-----|
| Preparation Tables - Fixed Height | 371 |
| Sterile Supplies Trolleys | 372 |
| Trolley Washer Compatible | 374 |

Storage Cupboards, Cabinets & Lockers

| Storage Cupboards | 378 |
|---------------------------------------|-----|
| Storage Cabinets | 380 |
| First Aid & Clinical Reagent Cabinets | 381 |
| Key Cabinets | 382 |
| Storage / Personal Effects Lockers | 384 |

Theatre Equipment

| Bowl Stands | 388 |
|-----------------------------------|-----|
| Bucket Stands | 391 |
| Instrument Trolleys | 392 |
| Mayo Tables & Operating Platforms | 394 |
| Suture Storage Trolleys | 395 |
| Sterile Container Storage | 396 |

Workstations & Speciality Storage

| Workstations | 400 |
|------------------------|-----|
| Storage Units | 401 |
| Catheter Storage Units | 402 |





- Folding Wheelchairs 8
 - Portering Chairs 9
- Mobile Shower/Commode Chairs 11
 - Bariatric Shower Chairs 12
 - Shower/Commode Chairs 14
 - Chairs 15
 - Rollators & Training Suits 17



Bariatric Folding Wheelchair - Minimaxx & Disc

Features

- Maximum capacity 325kg
- Powder coated chromium molybdenum steel frame
- Rigid folding mechanism
- Low seat height to assist access
- Viscoelastic foam cushion with two way stretch vapour permeable incontinence cover and non-slip base
- Padded push handles to assist the attendant/porter
- Puncture proof wheels & castors
- Parking brakes as standard/optional disc brake for improved control
- Side/width adjustable polyurethane moulded footrests, 100kg capacity

Dimensions (w x d x h)

- Seat Width 610/710mm
- Seat Height 430mm
- Seat Depth 460mm
- Backrest Height 510mm

Options

- 5X-0100-099-005 Butterfly Armrests, Pair
- 5X-0100-099-006 Angled Armrests, Pair
- 5X-0156-099-030 Elevating Calf Support, Right Hand
- 5X-0156-099-031 Elevating Calf Support, Left Hand
- 5X-0100-099-001 Universal Table
- 5X-0100-099-020 IV & O2 Holder, 100mm Dia

Spares

- 5X-0230-061-045 Viscoelastic Cushion with Incontinence Cover. 610mm
- 5X-0230-071-045 Viscoelastic Cushion with Incontinence Cover, 710mm



Universal Table





| Item | Description | Overall Dimensions (w x d x h) | Weight |
|-------------------|----------------------|--------------------------------|--------|
| 5X-0100-061-000/K | Minimaxx, 610mm | 850 x 1100 x 990mm | 27.0kg |
| 5X-0100-071-000/K | Minimaxx, 710mm | 950 x 1100 x 990mm | 27.5kg |
| 5X-0109-061-000/K | Minimaxx Disc, 610mm | 850 x 1100 x 990mm | 28.0kg |
| 5X-0109-071-000/K | Minimaxx Disc, 710mm | 950 x 1100 x 990mm | 28.5kg |

Bariatric Folding Wheelchair - Minimaxx Push Motor

Features

- Maximum capacity 325kg
- Powder coated chromium molybdenum steel frame
- Rigid folding mechanism
- Low seat height to assist access
- Viscoelastic foam cushion with two way stretch vapour permeable incontinence cover and non-slip base
- Padded push handles to assist the attendant/porter
- Puncture proof wheels & castors
- Parking brakes as standard/optional disc brake for improved control
- Side/width adjustable polyurethane moulded footrests, 100kg capacity • Push motor to assist the attendant/porter, reduces the risk of sore shoulders & stress on wrists
- Compatible on flat surfaces up to 5 degrees
- Operated using a control panel fitted to the handle, offering drive & reverse modes with speed control from 0.6-3.1mph
- One charge provides a range of 2.75miles
- Supplied with one removable battery and charger (battery has to be removed from the unit when charging)

Dimensions (w x d x h)

- Seat Width 610/710mm
- Seat Height 430mm
- Seat Depth 460mm
- Backrest Height 510mm

Options

- 5X-0100-099-005 Butterfly Armrests, Pair
- 5X-0100-099-006 Angled Armrests, Pair
- 5X-0156-099-030 Elevating Calf Support, Right Hand
- 5X-0156-099-031 Elevating Calf Support, Left Hand
- 5X-0100-099-001 Universal Table
- 5X-0100-099-020 IV & O2 Holder, 100mm Dia

Spares

- 5X-0230-061-045 Viscoelastic Cushion with Incontinence Cover, 610mm
- 5X-0230-071-045 Viscoelastic Cushion with Incontinence Cover, 710mm

| Item | Description | Overall Dimensions (w x d x h) | Weight |
|-------------------|-----------------------------|--------------------------------|--------|
| 5X-0110-061-000/K | Minimaxx, Push Motor, 610mm | 850 x 1100 x 990mm | 50.0kg |
| 5X-0110-071-000/K | Minimaxx, Push Motor, 710mm | 950 x 1100 x 990mm | 50.5kg |

5X-0100-061-000/K

IV & O² Holder

Elevating Calf Support





Push Motor



Handle Mounted Control Panel

Portering Chair - Bariatric, Rear Steer, Sliding Foot Rest



Features

- Maximum capacity 380kg
- Seamless & contoured upholstery using flame retardant vinyl
- Moulded drop down arm rests allows side access for patients and helpers
- Sliding rubber foot rest aids access as well as supporting feet during transit
- X-ray/notes pocket on rear of seat
- Two hooks for patients belongings
- Fittings to accept IV poles mounted to rear uprights
- Retractable diagonal seatbelt
- 200mm fixed wheels/200mm total stop braking castors (cushioned)
- We recommend that two members of staff are in attendance when used with a patient in excess of 190kg

Dimensions (w x d x h)

- Overall Dimensions 930 x 1000 x 1090mm
- Seat Width 740mm
- Seat Height 530mm
- Seat Depth 670mm
- Backrest Height 480mm

Standard Colours

- Paint Grey
- Upholstery Bristol Blue

Options

- G/500/CYL Cylinder Holder, CD/D
- G/CLS/500 Calf/Leg Supports (pair)
- G/DP Infusion Pole
- G/CARD/STATION Head of Station
- G/CARD/LOCK Card/Coin Lock, Factory Fitted
- 2LOKE011 Release Card

Bariatric Mobile Shower / Commode Chairs

Features

- Suitable for showering or toileting (supplied with a bucket & lid)
- Powder coated chromium molybdenum steel frame
- Available in two widths 610/710mm
- Rotating & removable straight or butterfly armrests to assist access
- · Open back to allow easy sling access & assist positioning
- Polyurethane moulded seat & back
- Seat supplied with a removable lid
- Seat panel removable for cleaning
- Long open hole in seat (220mm at the widest point)
- Deep bucket with lid is removable from the front or back of the chair
 - Side/width adjustable polyurethane moulded footrests, 100kg capacity
 - 125mm swivel castors, 2 off braking
 - Maximum capacity 325kg

Dimensions (w x d x h)

- Seat Width 610/710mm
- Seat Height 550mm
- Seat Depth 560mm
- Backrest Height 480mm

Standard Colours

- Paint White
- Polyurethane Blue


Options

- 5X-0150-999-022 Seat Panel 610mm
- 5X-0150-999-023 Seat Panel 660/710mm
- 5X-0140-999-001 Commode Bucket & Lid
- 5X-0156-099-030 Elevating Calf Support, Right Hand
- 5X-0156-099-031 Elevating Calf Support, Left Hand
- 5X-0156-099-010 Butterfly Armrests, Pair
- 5X-0150-099-010/K Straight Armrests, Pair

| ltem | Description | Overall Dimensions (w x d x h) | Weight |
|---------------------------|---|--------------------------------|--------|
| Standard Armrest | | | |
| 5X-0150-061-000 | Bariatric Mobile Shower / Commode Chair - Standard Armrests, 610mm | 670 x 1120 x 990mm | 27.0kg |
| 5X-0150-071-000 | Bariatric Mobile Shower / Commode Chair - Standard Armrests, 710mm | 770 x 1120 x 990mm | 27.5kg |
| Butterfly Armrests | | | |
| 5X-0151-061-000 | Bariatric Mobile Shower / Commode Chair - Butterfly Armrests, 610mm | 840 x 1120 x 990mm | 27.5kg |
| 5X-0151-071-000 | Bariatric Mobile Shower / Commode Chair - Butterfly Armrests, 710mm | 940 x 1120 x 990mm | 28.0kg |

Alternative colours are available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Rear Steer Portering Chairs - Bariatric, Sliding Foot Rest in Aqua Marine - G/500/RS/VAQ

Vinyl - Ambla New Richmond in Colour



| sales@bristolmaid.com + | +44 (0)1258 4 | 484455 | www.bristoln |
|-------------------------|---------------|--------|--------------|
|-------------------------|---------------|--------|--------------|

G/500/RS/VBB Rear Steer, Sliding Foot Rest, Bristol Blue

43.8kg

| ltem | Description |
|------|-------------|





Bariatric Shower Chair - Tilt

Features

- Suitable for showering or toileting (supplied with a bucket & lid)
- Powder coated chromium molybdenum steel frame
- 32-degree tilt positioning (battery operation up to 21 cycles from one charge)
- Supplied with one removable battery & charger (battery must be removed from the unit when charging)
- Rotating & removable straight armrests to assist access (optional butterfly armrests)
- Open back to allow easy sling access & assist positioning
- Polvurethane moulded seat, back & adjustable headrest
- Seat supplied with a removable lid
- Seat panel removable for cleaning
- Long open hole in seat (220mm at the widest point)
- Optional deep bucket is removable from the front or back of the chair
- Side/width adjustable polyurethane moulded footrests, 100kg capacity (optional elevating calf supports)
- 150mm swivel castors, 2 off braking
- Maximum capacity 325kg

Dimensions (w x d x h)

- Seat Width 660mm
- Seat Height 510mm
- Seat Depth 420-480-540mm
- Backrest Height 530mm

Standard Colours

- Paint White
- Polyurethane Blue

Options

- 5X-0150-999-023 Seat Panel 660/710mm
- 5X-0140-999-001 Commode Bucket & Lid
- 5X-0156-099-030 Elevating Calf Support, Right Hand
- 5X-0156-099-031 Elevating Calf Support, Left Hand
- 5X-0156-099-010 Butterfly Armrests, Pair
- 5X-0150-099-010/K Straight Armrests, Pair

Bariatric Shower Chair - Rise & Tilt

Features

- Suitable for showering or toileting (supplied with a bucket & lid) Powder coated chromium molybdenum steel frame
- Height adjustable from 450-1000mm whilst tilting the user. 20 degrees in the highest position
- Battery operation up to 21 cycles from one charge
- Supplied with one removable battery & charger (battery must be removed from the unit when charging)
- Rotating & removable straight armrests to assist access (optional butterfly armrests)
- Open back to allow easy sling access & assist positioning
- Polyurethane moulded seat, back & adjustable headrest
- Seat supplied with a removable lid
- Seat panel removable for cleaning
- Long open hole in seat (220mm at the widest point)
- Optional deep bucket is removable from the front or back of the chair
- Side/width adjustable polyurethane moulded footrests, 100kg capacity (optional elevating calf supports)
- 150mm swivel castors, all braking
- Maximum capacity 325kg

Dimensions (w x d x h)

- Seat Width 660mm
- Seat Height 450-1000mm
- Seat Depth 420-480-540mm
- Backrest Height 240-530mm

Standard Colours

- Paint White
- Polyurethane Blue

Options

59.0kg

- 5X-0150-999-023 Seat Panel 660/710mm
- 5X-0140-999-001 Commode Bucket & Lid
- 5X-0156-099-030 Elevating Calf Support, Right Hand
- 5X-0156-099-031 Elevating Calf Support, Left Hand • 5X-0156-099-010 - Butterfly Armrests, Pair
- 5X-0150-099-010/K Straight Armrests, Pair

| Item | Description | Overall Dimensions (w |
|-----------------|---|-----------------------|
| 5X-0157-066-000 | Bariatric Shower Chair Rise & Tilt - Standard Armrests, 660mm | 740 x 1400 x 1000mm |



Handset to Tilt Chair



46.0kg

sales@bristolmaid.com +44 (0)1258 484455 www.bristolmaid.com



Bariatric Shower / Commode Chair

Features

- Suitable for showering or toileting (supplied with a bucket & lid) Powder coated chromium molybdenum steel frame
- Available in two widths 610/710mm
- Seat height adjustable 420-550mm
- Long armrests to assist access
- Open back to allow easy sling access & assist positioning
- Polyurethane moulded seat & back
- Seat supplied with a removable lid
- Seat panel removable for cleaning
- Long open hole in seat (220mm at the widest point)
- Optional deep bucket is removable from the front or back of the chair
- Maximum capacity 325kg

Dimensions (w x d x h) Standard Colours

- Seat Width 610/710mm Seat Height - 420-550mm
- Paint White Polyurethane - Blue
- Seat Depth 400-520mm

Options

- 5X-0150-999-022 Seat Panel 610mm
- 5X-0150-999-023 Seat Panel 660/710mm
- 5X-0140-999-001 Commode Bucket & Lid



Features Deep CMHR foam

- Vinyl or fabric upholstery
- Vinvl Fire retardant & antimicrobial
- · Fabric Fire retardant
- Black heavy duty steel five star base
- Black heavy duty gas lift
- Black height adjustable soft touch armrests

Bariatric High Back Chairs

- Easy lift ratchet height adjustable back
- Independent seat & back mechanism with side tension control
- Maximum capacity 315kg

Replace suffix on the product code with any of the following to select your colour i.e. Bariatric High Back in Lime - 5660XL/VLI

Vinyl - Ambla New Richmond in Colour





| ltem | Description | Overall Dimensions (w x |
|--------|-------------------|-------------------------|
| 5660XL | High Back, Swivel | 750 x 750 x 1080-1190mr |

| Item | Description | Overall Dimensions (w x d x h) | Weight |
|-----------------|--|--------------------------------|--------|
| 5X-0142-061-000 | Bariatric Shower / Commode Chair - 610mm | 730 x 1010 x 790 - 940mm | 15.0kg |
| 5X-0142-071-000 | Bariatric Shower / Commode Chair - 710mm | 830 x 1010 x 790 - 940mm | 16.0kg |







Bariatric Chairs

Features

- Deep CMHR foam
- Vinyl or fabric upholstery
- Vinyl Fire retardant & antimicrobial
- Fabric Fire retardant
- Steel or beech four leg frames
- Steel Powder coated
- Beech Varnished
- Maximum capacity 315kg

Alternative colours are available at extra cost. Replace suffix on the product code with any of the following to select your colour i.e. Bariatric, Steel Frame in Lime - 5670/VLI

Vinyl - Ambla New Richmond in Colour



| ltem | Description | Overall Dimensions (w x d x h) | Seat Height | Weight |
|------|------------------------|--------------------------------|-------------|--------|
| 5670 | Bariatric, Steel Frame | 920 x 600 x 885mm | 500mm | 15.8kg |
| 5671 | Bariatric, Beech Frame | 844 x 580 x 890mm | 490mm | 16.0kg |

Bariatric Rollator - King

Features

- Maximum capacity 295kg
- Heavy duty powder coated steel frame
- Rigid folding mechanism for effcient storage
 Height adjustable 775-990mm
- Large padded seat with back arch acting as a positive stop (355 x 508mm)
- 200mm cushioned wheels with easy to use hand brakes

Dimensions (w x d x h)

- Overall Dimensions 710 x 820 x 775-990mm
- Seat 508 x 355 x 500mm

Standard Colours

- Paint Red
- Vinyl Black



| Item | Description | Weight |
|-----------------|-----------------|--------|
| 5X-0170-050-000 | Rollator - King | 15.7kg |

Bariatric Walker - Stand Tall

Features

- Maximum capacity 325kg
- Heavy duty powder coated steel frame
- Height adjustable 1000-1300mm
- Low level handgrips to assist the user to stand up
- Full length width adjustable underarm pads
- 100mm swivel castors/125mm braking wheel
- Handlebar operated brakes with locking function

Standard Colours

- Paint Blue / Silver
- Vinyl Black

Dimensions (w x d x h)

• Overall Dimensions - 760 x 720 x 1000-1300mm



Bariatric Training Suit

Features

- Helps educate care staff to the movement difficulties
 experienced by obese patients
- Foam fitting inner suit
- Unisex design
- Ordinary clothing



| Item | Description | Weight |
|-----------------|---------------------|--------|
| 5X-0180-076-000 | Walker - Stand Tall | 23.0kg |



Bedside Cabinets & Overbed Tables

- Polymer Bedside Cabinets 20
 - Bedside Cabinets 22
- Bedside Cabinet Options 25
- Bedside Cabinet/Wardrobe Combinations 26
 - Economy Overbed Table 30
 - Polymer Overbed & Overchair Tables 31
 - Laminate Overbed & Overchair Tables 32

19

eClean™



Polymer Bedside Cabinet - Easy Clean

Features

- BS 3962 Part 6:1980 Assessment of resistance to mechanical damage
- BS EN 12720:2009 Furniture assessment of surface resistance to cold liquids
- BS EN 12721:2009 Furniture assessment of surface resistance to wet heat
- BS EN 12722:2009 Furniture assessment of surface resistance to dry heat
- · BS 5852:2006 Methods of the test for assessment of the ignitability
- BS 4875-7:2006 Domestic & contract storage furniture
- Developed as part of Design Bugs Out & evaluated by the Department of Health
- Manufactured from recyclable polyethylene offering good surface protection against liquid stains & heat marks
- Open backless design ensures that all areas are accessible
- Smooth surfaces allow for easy & effective cleaning
- Distinct upper, middle & lower zones
- Upper Nurse zone with integral rim to contain any spillages. Ideally positioned for staff to use while standing
- Middle Patient zone provides storage area for patient use
- Lower Patient zone provides shelf for washbowl, footwear etc.
- Drawers can be inserted from either side enabling the cabinet to be positioned on either side of the bed
- Upper & middle drawers can be fitted with an RFID lock, providing secure storage for patient medicines/property*
- Large 100mm twin wheel castors with easily accessible brakes

*cabinets with locks are supplied with one patient card - a programming card will be required per ward/directorate as per hospital policy.

Dimensions (w x d x h)

Drawers (internal)

Standard Colours

- Overall Dimensions 480 x 665 x 1237mm
- Body Grey White Drawers - Spearmint
- Upper 149 x 222 x 217mm
- Middle 338 x 491 x 197mm
- Lower 338 x 491 x 220mm

Options

- EBC900 Patient Card (Blue) Cabinets with locks are supplied with one patient card
- EBC901 Staff Card (Green)
- EBC902 Master Card, Staff (Red)
- EBC910 Programming Card Staff Pair (Black)
- 2LOKE134 Key Fob

| ltem | Description | Weight |
|--------|---|--------|
| EBC015 | Polymer Bedside Cabinet - No Locks | 29.0kg |
| EBC020 | Polymer Bedside Cabinet - Upper Drawer with Lock | 29.1kg |
| EBC025 | Polymer Bedside Cabinet - Middle Drawer with Lock | 29.1kg |
| EBC030 | Polymer Bedside Cabinet - Upper & Middle Drawers with Locks | 29.2kg |



Lock Cards

Features

- Compatible with battery powered RFID controlled electronic locks
- Credit card size plastic key cards with an embedded contactless RFID chip
- Available as user/staff/master access cards
- A programming card is required for each hospital or ward as per hospital procedure

Overall Dimensions (w x d x h)

• Card - 86 x 54 x 2mm



| Item | Description |
|----------|--|
| EBC900 | User Card - Patient (Blue) |
| EBC901 | User Card - Staff (Green) |
| EBC902 | Master Card - Staff (Red) |
| EBC910 | Programming Card - Staff Pair, (Black) |
| 2LOKE134 | Key Fob |





NC

BA

CC

CA

BF

Bedside Cabinet - Single Upper Drawer

Features

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients (with exception of door hinges which are more robust & practical)
- Constructed from melamine faced chipboard & laminate
- Upper drawer suitable for the storage of patient's medication or personal effects (fitted with a CAM lock & supplied with 2 keys)
- Optional personal effects drawer provides secure storage for personal possessions (fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button & RFID locking solutions available
- for the upper & personal effects drawer (see page 25)
- Lower storage offered through...
- Cupboard with double doors (adjustable shelf model dependant)
- Drawer (recessed shelf model dependant)
- Choice of two colour schemes...
- Two Tone Body grey white. Top, drawer & door fronts available in a range of colours
- Single Colour Beech throughout
- 50mm twin-wheel polypropylene castors

Dimensions (w x d x h)

- Overall Dimensions: 485 x 485 x 905mm
- Internal Dimensions...
- Upper Drawer: 413 x 381 x 182mm
- Personal Drawer: 413 x 381 x 102mm
- Lower Drawer: 413 x 381 x 505/297mm
- Cupboard: 445 x 443 x 537/394mm

Standard Colours

- Two Tone • PL - PS Solo RFID Electronic Lock (upper drawer)
- Body Grey White • TL - TESA RFID Electronic Lock (upper drawer)
 - Top / Drawer / Door Fronts Spearmint
 - Top / Drawer / Door Fronts Beech

| Grey White | Beech | Description | Weight |
|------------|-----------|--|--------|
| BC1C/GW | BC1C/BE | Single Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 34.0kg |
| BC1CPD/GW | BC1CPD/BE | Single Upper Drawer, Personal Drawer & Lower Cupboard with Double Doors | 35.8kg |
| BC1D/GW | BC1D/BE | Single Upper Drawer, Lower Drawer with Recessed Shelf | 40.4kg |
| BC1DPD/GW | BC1DPD/BE | Single Upper Drawer, Personal Drawer & Lower Drawer | 42.1kg |

Bedside Cabinet - Double Upper Drawer

Features

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients (with exception of door hinges which are more robust & practical)
- Constructed from melamine faced chipboard & laminate
- Upper drawers suitable for the storage of patient's medication or personal effects (1 x fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button & RFID locking solutions available for either drawer (see page 25)
- Lower storage offered through...
- Cupboard with double doors (adjustable shelf)
- Drawer (recessed shelf)
- · Choice of two colour schemes... • Two Tone - Body grey white. Top, drawer & door fronts available in a range of colours
- Single Colour Beech throughout
- 50mm twin-wheel polypropylene castors

Dimensions (w x d x h)

- Overall Dimensions 485 x 485 x 905mm
- Internal Dimensions...
- Upper Drawer: 180 x 381 x 182mm
- Lower Drawer: 413 x 381 x 505mm
- Cupboard: 445 x 443 x 537mm

Standard Colours

Two Tone

- Body Grev White
- Top / Drawer / Door Fronts Spearmint
- Single Colour
- Body Beech
- Top / Drawer / Door Fronts Beech

Options

- PL PS Solo RFID Electronic Lock (upper drawer)
- TL TESA RFID Electronic Lock (upper drawer)
- CL Electronic Push Button Lock (upper drawer)
- Lock Cards See page 25
- RPT Removable Plastic Top
- TH Towel Holder

| Grey White | Beech | Description | Weight |
|------------|---------|--|--------|
| BC2C/GW | BC2C/BE | Double Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 34.2kg |
| BC2D/GW | BC2D/BE | Double Upper Drawer, Lower Drawer with Recessed Shelf | 40.5kg |



22

NO

+44 (0)1258 484455 sales@bristolmaid.com www.bristolmaid.com

BC1C/GW/RPT

 Lock Cards - See page 25 RPT - Removable Plastic Top TH - Towel Holder

personal effects drawer)

Options

BC1DPD/BE

man -----

the second



Alternative Colours Available, Contact Us For Details



Bedside Cabinet - Upper Section & Rear Side Doors



BC3C/GW

Alternative Colours Available, Contact Us For Details

BC3CPD/BE

CC BF NC BA

Features

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients (with exception of door hinges which are more robust & practical)
- Constructed from melamine faced chipboard & laminate
- Upper section suitable for the storage of patient's medication or personal effects (fitted with a CAM lock & supplied with 2 keys)
- Rear side doors allow access to contents
- Optional personal effects drawer provides secure storage for personal possessions (fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button & RFID locking solutions available for the upper section & personal effects drawer (see page 25)
- Lower storage offered through...
- Cupboard with double doors (adjustable shelf model dependant)
- Drawer (recessed shelf model dependant)
- · Choice of two colour schemes...
- Two Tone Body grey white. Top, drawer & door fronts available in a range of colours
- Single Colour Beech throughout
- 50mm twin-wheel polypropylene castors

Dimensions (w x d x h)

- Overall Dimensions: 485 x 485 x 905mm
- Internal Dimensions...
- Upper Section: 445 x 135 x 249mm
- Rear Side Section: 290 x 445 x 249mm
- Personal Drawer: 413 x 381 x 102mm
- Lower Drawer: 413 x 381 x 505/297mm
- Cupboard: 445 x 443 x 537/394mm

Standard Colours

- Two Tone
- Body Grev White
- Top / Drawer / Door Fronts Spearmint
- Single Colour
- Body Beech
- Top / Drawer / Door Fronts Beech

TH - Towel Holder

personal effects drawer)

Lock Cards - See page 25

RPT - Removable Plastic Top

• PL - PS Solo RFID Electronic Lock (upper drawer)

• TL - TESA RFID Electronic Lock (upper drawer)

CL - Electronic Push Button Lock (upper &

| Grey White | Beech | Description | Weight |
|------------|-----------|--|--------|
| BC3C/GW | BC3C/BE | Upper Section, Lower Cupboard with Double Doors & Adjustable Shelf | 31.8kg |
| BC3CPD/GW | BC3CPD/BE | Upper Section, Personal Drawer & Lower Cupboard with Double Doors | 33.5kg |
| BC3D/GW | BC3D/BE | Upper Section, Lower Drawer with Recessed Shelf | 38.1kg |
| BC3DPD/GW | BC3DPD/BE | Upper Section, Personal Drawer & Lower Drawer | 39.8kg |

Electronic Push Button & RFID Locking Solutions

Through either electronic push button or contactless RFID, we are able to offer a secure locking solution to suit your specific requirements. The table below highlights the key features of each system - for more information please contact our customer services team.



| | Electronic Push Button Lock | TESA RFID Electronic Lock | PS Solo RFID Electronic Lock |
|--|--|---|---|
| RFID Compatibility | - | Mifare - 13.56 mHz | EM4200, EM4450, Hitag - 125 kHz |
| Levels of User Administration | 4 - Master / Nurse / User / Technician | 3 - Master / Nurse / User | 3 - Master / Nurse / User |
| Operation Modes | Private Use Public Use ^a | Private Use Public Use⁵ | Private Use Public Use⁵ |
| Audit Trail | No | No | No |
| Retro Fit Into Existing Cylinders | Yes | Yes - Dependent on cabinet type | No |
| Programming | 10 Button Keypad & # Button | Programming Card | Programming Card |
| penings On One Battery | 15,000 | 15,000 | 20,000 |
| Low Battery Warning | Yes | Yes | Yes |
| Batteries | 2 x AAA (1.5V) | 3 x AA (1.5v) | 1 x CR2 (3V) |
| Low Battery Override | 9v Battery Across LED's (+/-) & Master Code | Yes - Using an external plug in power supply | Lock Will Fail Open |
| Compatible with Abloy [®] Cliq™ Remote | No | Yes - Using an RFID Chip inserted in Cliq™ key | Yes - Using a combined 125 kHz / 13,56mHz card or 13.56mHz RFID sticker applied to existing 125kHz cards |

Operation Modes

- 1. Private Use This is the most common function & used where the same code will be repeatedly used 2. Public Use

 - will reopen the lock ready for the next user

Options



a. In this mode the lock operates with a single user code. The user enters their own personal four-digit code. This locks the lock. The same code is entered once only to open the lock before being erased & ready for the next new user b. In this mode the lock operates with a single user card. The presented user card locks the lock. The same user card

Towel Holder

Features

- Self-adhesive towel holder
- Supplied for customer to fit



| Item | Description | Weight |
|------|--------------|--------|
| TH | Towel Holder | 0.1kg |

Removable Plastic Top

Features

- Contoured plastic top
- Raised lip on three sides reduces misplacement of objects
- Clips to the top of the bedside cabinet
- Removable to assist cleaning
- Not compatible with bedside cabinets fitted with wardrobes



Lock Cards

Features

- Compatible with TESA & PS Solo battery powered RFID controlled electronic locks
- Credit card size plastic key cards with an embedded contactless RFID chip
- Available as programming/staff/patient cards
- A programming card is required for each hospital or ward as per hospital procedure
- · A patient card is supplied with the cabinets

| Item | Description |
|------------|------------------|
| RFID020/P | Programming Card |
| RFID020/S | Staff Card |
| RFID020/U | Patient Card |
| RFID020/SF | Staff Fob |
| RFID020/PF | Patient Fob |

Compact Bedside Cabinet / Wardrobe Combination - Single Upper Drawer



BCC1DWR/BE

Alternative Colours Available, Contact Us For Details

-

1-2-2-



Features

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients (with exception of door hinges which are more robust & practical)
- · Constructed from melamine faced chipboard & laminate
- Upper drawer suitable for the storage of patient's medication or personal effects (fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button & RFID locking solutions available for the upper drawer (see page 25)
- Lower storage offered through...
- Cupboard with double doors (adjustable shelf)
- Drawer (recessed shelf)
- · Choice of two colour schemes...
- Two Tone Body grey white. Top, drawer & door fronts available in a range of colours
- Single Colour Beech throughout
- Wardrobe supplied with a hanging rail, available right/left hand opening
- 50mm twin-wheel polypropylene castors

Dimensions (w x d x h)

- Overall Dimensions: 633 x 489 x 1420mm
- Internal Dimensions...
- Upper Drawer: 328 x 381 x 182mm
- Lower Drawer: 328 x 381 x 505mm
- Cupboard: 360 x 443 x 537mm
- Wardrobe: 215 x 443 x 1420mm

Standard Colours

- Two Tone · Body - Grey White
- Top / Drawer / Door Fronts Spearmint
- Single Colour
- Body Beech
- Top / Drawer / Door Fronts Beech

| Grey White | Beech | Description | Weight |
|------------|------------|--|--------|
| Right Hand | | | |
| BCC1CWR/GW | BCC1CWR/BE | Single Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 50.6kg |
| BCC1DWR/GW | BCC1DWR/BE | Single Upper Drawer, Lower Drawer with Recessed Shelf | 56.9kg |
| Left Hand | | | |
| BCC1CWL/GW | BCC1CWL/BE | Single Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 50.6kg |
| BCC1DWL/GW | BCC1DWL/BE | Single Upper Drawer, Lower Drawer with Recessed Shelf | 56.9kg |

Bedside Cabinet / Wardrobe Combination - Single Upper Drawer

Features

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients (with exception of door hinges which are more robust & practical)
- Constructed from melamine faced chipboard & laminate
- Upper drawer suitable for the storage of patient's medication or personal effects (fitted with a CAM lock & supplied with 2 keys)
- Optional personal effects drawer provides secure storage for personal possessions (fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button & RFID locking solutions available for the upper & personal effects drawer (see page 25) Lower storage offered through...
- Cupboard with double doors (adjustable shelf)
- Drawer (recessed shelf model dependant)
- · Choice of two colour schemes...

Dimensions (w x d x h)

Internal Dimensions...

- Two Tone Body grey white. Top, drawer & door fronts available in a range of colours
- Single Colour Beech throughout
- Wardrobe supplied with a hanging rail, available right/left hand opening
- 50mm twin-wheel polypropylene castors

- Overall Dimensions: 718 x 489 x 1420mm
- Upper Drawer: 413 x 381 x 182mm
- Personal Drawer: 413 x 381 x 102mm
- Lower Drawer: 413 x 381 x 505/297mm
- Cupboard: 445 x 443 x 537/394mm
- Wardrobe: 215 x 443 x 1420mm

| Grey White | Beech | Description | Weight |
|-------------|-------------|--|--------|
| Right Hand | | | |
| BC1CWR/GW | BC1CWR/BE | Single Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 54.5kg |
| BC1CPDWR/GW | BC1CPDWR/BE | Single Upper Drawer, Personal Drawer & Lower Cupboard with Double Doors | 56.2kg |
| BC1DWR/GW | BC1DWR/BE | Single Upper Drawer, Lower Drawer with Recessed Shelf | 60.8kg |
| BC1DPDWR/GW | BC1DPDWR/BE | Single Upper Drawer, Personal Drawer & Lower Drawer | 62.5kg |
| Left Hand | | | |
| BC1CWL/GW | BC1CWL/BE | Single Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 54.5kg |
| BC1CPDWL/GW | BC1CPDWL/BE | Single Upper Drawer, Personal Drawer & Lower Cupboard with Double Doors | 56.2kg |
| BC1DWL/GW | BC1DWL/BE | Single Upper Drawer, Lower Drawer with Recessed Shelf | 60.8kg |
| BC1DPDWL/GW | BC1DPDWL/BE | Single Upper Drawer, Personal Drawer & Lower Drawer | 62.5kg |

• PL - PS Solo RFID Electronic Lock (upper drawer)

• TL - TESA RFID Electronic Lock (upper drawer)

CL - Electronic Push Button Lock (upper &

personal effects drawer)

Lock Cards - See page 25

• TH - Towel Holder

Standard Colours

Body - Beech

Body - Grey White

Two Tone

Single Colour

- Top / Drawer / Door Fronts Beech

Top / Drawer / Door Fronts - Spearmint

Options

- PL PS Solo RFID Electronic Lock (upper drawer)
- TL TESA RFID Electronic Lock (upper drawer)
- CL Electronic Push Button Lock (upper &
- personal effects drawer)
- Lock Cards See page 25
- TH Towel Holder





Alternative Colours Available, Contact Us For Details



Alternative Colours Available, Contact Us For Details



Bedside Cabinet / Wardrobe Combination - Double Upper Drawer

Features

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients (with exception of door hinges which are more robust & practical)
- Constructed from melamine faced chipboard & laminate • Upper drawers suitable for the storage of patient's medication or
- personal effects (1 x fitted with a CAM lock & supplied with 2 keys) Optional electronic push button & RFID locking solutions available for
- either drawer (see page 25)
- Lower storage offered through...
- Cupboard with double doors (adjustable shelf)
- Drawer (recessed shelf)
- · Choice of two colour schemes...
- Two Tone Body grey white. Top, drawer & door fronts available in a range of colours
- Single Colour Beech throughout
- Wardrobe supplied with a hanging rail, available right/left hand opening
- 50mm twin-wheel polypropylene castors

Dimensions (w x d x h)

- Overall Dimensions: 718 x 489 x 1420mm
- Internal Dimensions
- Upper Drawer: 180 x 381 x 182mm
- Lower Drawer: 413 x 381 x 505mm
- Cupboard: 445 x 443 x 537mm
- Wardrobe: 215 x 443 x 1420mm

Standard Colours

Two Tone

- Body Grey White
- Top / Drawer / Door Fronts Spearmint
- Single Colour
- Body Beech
- Top / Drawer / Door Fronts Beech

| Grey White | Beech | Description | Weight |
|------------|-----------|--|--------|
| Right Hand | | | |
| BC2CWR/GW | BC2CWR/BE | Double Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 55.0kg |
| BC2DWR/GW | BC2DWR/BE | Double Upper Drawer, Lower Drawer with Recessed Shelf | 61.3kg |
| Left Hand | | | |
| BC2CWL/GW | BC2CWL/BE | Double Upper Drawer, Lower Cupboard with Double Doors & Adjustable Shelf | 55.0kg |
| BC2DWL/GW | BC2DWL/BE | Double Upper Drawer, Lower Drawer with Recessed Shelf | 61.3kg |

Bedside Cabinet / Wardrobe Combination - Upper Section & Rear Side Doors

Features

- Meets BS1765: Part 1: 1990 Specification for general purpose bedside lockers for patients (with exception of door hinges which are more robust & practical)
- Constructed from melamine faced chipboard & laminate
- Upper section suitable for the storage of patient's medication or personal effects (fitted with a CAM lock & supplied with 2 keys)
- Rear side doors allow access to contents
- Optional personal effects drawer provides secure storage for personal possessions (fitted with a CAM lock & supplied with 2 keys)
- Optional electronic push button & RFID locking solutions available for the upper section & personal effects drawer (see page 25) Lower storage offered through...
- Cupboard with double doors (adjustable shelf model dependent)
- Drawer (recessed shelf model dependant)
- · Choice of two colour schemes...
- Two Tone Body grey white. Top, drawer & door fronts available in a range of colours
- Single Colour Beech throughout
- Wardrobe supplied with a hanging rail, available right/left hand opening
- 50mm twin-wheel polypropylene castors

Dimensions (w x d x h)

- Overall Dimensions: 718 x 489 x 1420mm Internal Dimensions...
- Upper Section: 445 x 135 x 249mm
- Rear Side Section: 290 x 445 x 249mm
- Personal Drawer: 413 x 381 x 102mm
- Lower Drawer: 413 x 381 x 505/297mm
- Cupboard: 445 x 443 x 537/394mm
- Wardrobe: 215 x 443 x 1420mm

| Grey White | Beech | Description | Weight |
|-------------|-------------|--|--------|
| Right Hand | | | |
| BC3CWR/GW | BC3CWR/BE | Upper Section, Lower Cupboard with Double Doors & Adjustable Shelf | 52.1kg |
| BC3CPDWR/GW | BC3CPDWR/BE | Upper Section, Personal Drawer & Lower Cupboard with Double Doors | 53.8kg |
| BC3DWR/GW | BC3DWR/BE | Upper Section, Lower Drawer with Recessed Shelf | 58.4kg |
| BC3DPDWR/GW | BC3DPDWR/BE | Upper Section, Personal Drawer & Lower Drawer | 60.1kg |
| Left Hand | | | |
| BC3CWL/GW | BC3CWL/BE | Upper Section, Lower Cupboard with Double Doors & Adjustable Shelf | 52.1kg |
| BC3CPDWL/GW | BC3CPDWL/BE | Upper Section, Personal Drawer & Lower Cupboard with Double Doors | 53.8kg |
| BC3DWL/GW | BC3DWL/BE | Upper Section, Lower Drawer with Recessed Shelf | 58.4kg |
| BC3DPDWL/GW | BC3DPDWL/BE | Upper Section, Personal Drawer & Lower Drawer | 60.1kg |

• PL - PS Solo RFID Electronic Lock (upper drawer)

• TL - TESA RFID Electronic Lock (upper drawer) CL - Electronic Push Button Lock (upper &

BC2CWL/GW

personal effects drawer) Lock Cards - See page 25

TH - Towel Holder

Options

| Grey White | Beech | Description | Weight |
|-------------|-------------|--|--------|
| Right Hand | | | |
| BC3CWR/GW | BC3CWR/BE | Upper Section, Lower Cupboard with Double Doors & Adjustable Shelf | 52.1kg |
| BC3CPDWR/GW | BC3CPDWR/BE | Upper Section, Personal Drawer & Lower Cupboard with Double Doors | 53.8kg |
| BC3DWR/GW | BC3DWR/BE | Upper Section, Lower Drawer with Recessed Shelf | 58.4kg |
| BC3DPDWR/GW | BC3DPDWR/BE | Upper Section, Personal Drawer & Lower Drawer | 60.1kg |
| Left Hand | | | |
| BC3CWL/GW | BC3CWL/BE | Upper Section, Lower Cupboard with Double Doors & Adjustable Shelf | 52.1kg |
| BC3CPDWL/GW | BC3CPDWL/BE | Upper Section, Personal Drawer & Lower Cupboard with Double Doors | 53.8kg |
| BC3DWL/GW | BC3DWL/BE | Upper Section, Lower Drawer with Recessed Shelf | 58.4kg |
| BC3DPDWL/GW | BC3DPDWL/BE | Upper Section, Personal Drawer & Lower Drawer | 60.1kg |

Options

- PL PS Solo RFID Electronic Lock (upper drawer)
- TL TESA RFID Electronic Lock (upper drawer)
- CL Electronic Push Button Lock (upper &
- personal effects drawer)
- Lock Cards See page 25
- TH Towel Holder

Standard Colours

Two Tone Body - Grey White

- Top / Drawer / Door Fronts Spearmint
- Single Colour
- Body Beech
- Top / Drawer / Door Fronts Beech





Alternative Colours Available, Contact Us For Details



Economy Overbed Table

Features

- Cantilever design
- Height/angle adjustable using two hand wheels (650-1080mm)
- Sliding top to allow for central positioning over bed (190mm)
- All round raised plastic lip
- Can be stacked when not in use
- 50mm twin wheel braking castors
- Maximum load 8kg

Dimensions (w x d x h)

• Overall Dimensions - 630 x 375 x 650/1080mm

Standard Colours

- Paint White
- Laminate Top Beech



Top extends to assist positioning

| Item | Description | Weight |
|---------|-----------------------|--------|
| BT/12/A | Economy Overbed Table | 9.0kg |

Overbed & Overchair Table - Polymer, Easy Clean

Features

- Grey White



| Item | Description |
|----------|--|
| OBP/SP/S | Height Adjustable, Polymer, Easy Clean |
| OCP/SP/S | Height Adjustable, Polymer, Easy Clean |

eClean[®]



Overbed & Overchair Table - Laminate

Features

- BS 2483: 1977 Specification for overbed tables (selected tests) • BS EN 438-2: 2016 High-pressure decorative laminates (HPL). Sheets based on thermosetting resins (usually called laminates). Determination of properties
- Height adjustable using a gas damper (735-1045mm)
- Extruded aluminium column provides a smooth linear height adjustment
- Safety feature prevents trapping of the patient between the bed & table
- Cantilever tubular base available as...
- C shaped suitable for use with beds & chairs
- Straight leg suitable for use with beds
- Laminate top with varnished edges available in three configurations
- Flat top
- One varnished beech lip
- Three varnished beech lips
- 50mm twin wheel polypropylene castors
- Maximum load 25kg

Dimensions (w x d x h) Options

- Table Top 880 x 395mm
- Overbed Table 895 x 440 x 735-1045mm
- Overchair Table 918 x 440 x 735-1045mm

Standard Colours

Top - Grey White or Beech Upper Extrusion - Spearmint or Grey White Lower Extrusion/Base - Grey White

- B Braking Castors (x 2)

Top - Grey White / Upright - Spearmint

| Item | Description | Weigh |
|-------------------------|---|--------|
| Overbed Tables | | |
| OBL/GW0/SP/S | Height Adjustable, Laminate, Flat Top | 14.5k |
| OBL/GW1/SP/S | Height Adjustable, Laminate, One Lip | 15.4k |
| OBL/GW3/SP/S | Height Adjustable, Laminate, Three Lips | 15.8kg |
| Overchair Tables | 5 | |
| OCL/GW0/SP/S | Height Adjustable, Laminate, Flat Top | 14.1k |
| OCL/GW1/SP/S | Height Adjustable, Laminate, One Lip | 15.0kg |
| OCL/GW3/SP/S | Height Adjustable, Laminate, Three Lips | 15.4kg |
| | | |

Top - Beech / Upright - Grey White

| Item | Description | Weight |
|-----------------------|---|--------|
| Overbed Tables | | |
| OBL/BE0/GW/S | Height Adjustable, Laminate, Flat Top | 14.5kg |
| OBL/BE1/GW/S | Height Adjustable, Laminate, One Lip | 15.4kg |
| OBL/BE3/GW/S | Height Adjustable, Laminate, Three Lips | 15.8kg |
| Overchair Table | s | |
| OCL/BE0/GW/S | Height Adjustable, Laminate, Flat Top | 14.1kg |
| OCL/BE1/GW/S | Height Adjustable, Laminate, One Lip | 15.0kg |
| OCL/BE3/GW/S | Height Adjustable, Laminate, Three Lips | 15.4kg |
| | | |













Caretray Trolleys & Racks

- Features 36
- Disassembly & Cleaning 38
- Stainless Steel, Low Level Trolleys 40
- Stainless Steel, Standard Level Trolleys 42
 - Stainless Steel, High Level Trolleys 44
 - Stainless Steel, Mobile Racks 46
 - Stainless Steel, Static Racks 48
 - Mild Steel, Low Level Trolleys 50
 - Mild Steel, Standard Level Trolleys 52
 - Mild Steel, High Level Trolleys 54
 - Mild Steel, Mobile Racks 56
 - Mild Steel, Static Racks 58
 - Caretray Trolley Options 60

Caretray Trolley Features

Carerails

• Fitted to both sides of the trolley, compatible with a wide range of accessories

Side Panels

- One piece moulded side panels
- Smooth surfaces, rounded edges & contours assists cleaning (compatible with steam cleaning)
- Runners with positive stops prevent travs from being accidentally withdrawn (symmetrical design allows tray insertion from either side)

Easy Clean

- One piece mouldings used throughout
- Smooth surfaces, rounded edges & contours assists cleaning
- Steam cleaning compatible
- Quick & easy tool-less disassembly assists thorough or deep cleaning programmes



Тор

- One piece moulded top prevents ingress of dirt
- Smooth surfaces, rounded edges & contours assists cleaning (compatible with steam cleaning) Raised edges on three sides provides defined
- grab points to act as push handles, whilst an integral rim at the front contains spillages Maximum load of 30kg

Frame

 Frames are fully welded eliminating any gaps • Available in either 304 grade stainless steel or mild steel that is blast finished & powder coated in a durable epoxy polyester paint



304 grade stainless steel

Ring Buffers

- Ring buffers offer protection to both the trolley & fabric of the building
- Smooth surfaces assist cleaning

Trays - Dividers - Plastic Clip On Label Cover - Labels

- Translucent, impact-resistant trays & dividers
- temperature machine washing)
- Trays feature an integral handle with plastic clip on label cover, (suitable for use with our perforated coloured label strips) • Travs are available in 100 or 150mm depths & can be subdivided in both length & width using simple slot in dividers • Each tray is supplied with 4 removable dividers (2 x long & 2 x short)

- Maximum load of 5kg



Doors & Back Panel

- Solid grade laminate doors & back panel
- Smooth surfaces, rounded edges & contours assists cleaning (compatible with steam cleaning)
- 270 degree opening hinges allow the doors to fold flat against the sides of the trolley
- Three locking options...






CAM with Winged Handle

Security Seal

· Smooth surfaces, rounded edges & contours assists cleaning (compatible with steam cleaning & able to withstand high

Castors

- Trolleys are fitted with 100mm pressed steel castors
- Two braking castors are fitted as standard
- Optional 100mm polymer castors offer reduced rolling resistance making the trolley easier to manoeuvre, lower noise levels & smooth surfaces to assist cleaning



Electronic Push Button Lock



Disassembly & Cleaning

The design of the trolley allows for quick & easy tool-less disassembly providing access for all areas when conducting thorough or deep cleaning programmes. The simple guide below illustrates the disassembly & re-assembly process...



Remove Trays Remove the trays from the unit



Remove Top

The top is held in place with two thumb screws. These are located on the underside of the top. Once these are removed the top can be removed with a firm tug




Remove Back Panel (if fitted)

Once the top has been removed the back panel can be lifted clear of the cross members on the rear of the frame



Remove Side Panels

The side panels can now be removed by lifting off of the positioning inserts. Double height trolleys require thumb screws at the top of the lower panel to be removed.





Clean

The frame, side panels & top can now be cleaned

Reassemble

The trolley can now be re-assembled by repeating the steps in reverse

Caretray Trolleys - Stainless Steel, Low Level, Single & Double Column





CTS1/08/3S1D - Shown with optional CT/WF. CT/AGH & GMBH/15



Easy disassembly - see page 38

Features

- Low level single & double column easy clean procedure trolleys. Height makes it suitable to be stored under bench
- Designed for quick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Seamlessly welded one piece frame from 304 grade stainless steel (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- Carerail fitted to each side of the trolley to support a wide range of accessories (see options)
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- · Each trolley is supplied with perforated coloured label strips (white, blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Single Column 550 x 550 x 807mm
- Double Column 1005 x 550 x 807mm
- Plastic Top
- Single Column 455 x 455 (internal) at height 782mm
- Double Column 915 x 455 (internal) at height 782mm
- Trays
- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)



Options

- CT/AGH Alcohol Gel Holder
- CT/AIR Airway Tube Holder
- CT/BIN Debris Bin & Bracket
- CT/DH/A4 Document Holder, A4
- CT/CR2 & CT/CR5 Catheter Racks
- Single & Double Doors See page 60
- CT/FS Flat Shelf Removable
- CT/IV Infusion Pole & Bracket
- CT/WF Work Flap
- GMBH/15 Glove / Mask Box Holder One Box, 15mm Hook
- SBH/15 Sharps Box Holder One Box, 15mm Hook
- CT/CYD Cylinder Holder C/CD/D Size
 CT/SBB/15 Sharps Box Bracket Universal (Daniels / Frontier / Sharpak) Polymer Castors - See page 69

| Item | Description | Weight |
|--------------|----------------------------|--------|
| CTS1/08/1S2D | 1 x 100mm, 2 x 150mm Trays | 21.2kg |
| CTS1/08/3D | 3 x 150mm Trays | 21.5kg |
| CTS1/08/3S1D | 3 x 100mm, 1 x 150mm Trays | 22.3kg |
| CTS1/08/4S | 4 x 100mm Trays | 21.9kg |

| Item | Description |
|--------------|----------------------------|
| CTS2/08/3S4D | 3 x 100mm, 4 x 150mm Trays |
| CTS2/08/4S3D | 4 x 100mm, 3 x 150mm Trays |
| CTS2/08/6D | 6 x 150mm Trays |
| CTS2/08/6S2D | 6 x 100mm, 2 x 150mm Trays |
| CTS2/08/7S1D | 7 x 100mm, 1 x 150mm Trays |
| CTS2/08/8S | 8 x 100mm Trays |
| | |

Caretray Trolleys - Stainless Steel, Standard Level, Single & Double Column



CTS1/10/3S2D - Shown with optional CT/IV. CT/CR2. GMBH/15 & CT/AGH



Easy disassembly - see page 38

- Options
- · CT/AGH Alcohol Gel Holder
- · CT/AIR Airway Tube Holder
- · CT/BIN Debris Bin & Bracket
- CT/DH/A4 Document Holder, A4
- CT/CR2 & CT/CR5 Catheter Racks
- Single & Double Doors See page 60
- CT/FS Flat Shelf Removable
- CT/IV Infusion Pole & Bracket

- CT/WF Work Flap
- GMBH/15 Glove / Mask Box Holder One Box, 15mm Hook
- SBH/15 Sharps Box Holder One Box, 15mm Hook
- CT/SBB/15 Sharps Box Bracket Universal (Daniels / Frontier / Sharpak)
- CT/LH/S Laptop Holder Swing Arm
- CT/CYD Cylinder Holder C/CD/D Size CT/LH/F Laptop Holder Fixed Arm
 - CT/LH/C5 Motion C5 Holder Swing Arm
 - CT/IV Infusion Pole & Bracket
 - Polymer Castors See page 69

Features

- Standard level single & double column easy clean procedure trolleys
- Designed for guick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Seamlessly welded one piece frame from 304 grade stainless steel (requires no on-site assembly)
- One piece moulded top with raised edges on three sides and integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- Carerails fitted to each side of the trollev to support a wide range of accessories (see options)
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Single Column 550 x 550 x 993mm

ltem

CTS1/10/1S3D

CTS1/10/3S2D

CTS1/10/4D

CTS1/10/6S

CTS1/10/4S1D

- Single Column 455 x 455 (internal) at height 968mm
- Double Column 915 x 455 (internal) at height 968mm
- Travs
- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)



| Item | Description |
|--------------|----------------------------|
| CTS2/10/3S6D | 3 x 100mm, 6 x 150mm Trays |
| CTS2/10/4S5D | 4 x 100mm, 5 x 150mm Trays |
| CTS2/10/5S4D | 5 x 100mm, 4 x 150mm Trays |
| CTS2/10/6S4D | 6 x 100mm, 4 x 150mm Trays |
| CTS2/10/8D | 8 x 150mm Trays |
| CTS2/10/8S2D | 8 x 100mm, 2 x 150mm Trays |

www.bristolmaid.com

- · Each trolley is supplied with perforated coloured label strips (white,

- Double Column 1005 x 550 x 993mm
- Plastic Top

Description

4 x 150mm Trays

6 x 100mm

1 x 100mm, 3 x 150mm Trays

3 x 100mm, 2 x 150mm Trays

4 x 100mm, 1 x 150mm Trays

27.9kg

29.0kg

28.3kg

28.7kg

29.8kg

sales@bristolmaid.com +44 (0)1258 484455

Caretray Trolleys - Stainless Steel, High Level, Single & Double Column





CTS1/14/4S3D - Shown with optional CT/D4/A4



Options

- CT/AGH Alcohol Gel Holder
- CT/AIR Airway Tube Holder
- CT/BIN Debris Bin & Bracket
- CT/DH/A4 Document Holder, A4
- CT/CR2 & CT/CR5 Catheter Racks
- CT/CYD Cylinder Holder C/CD/D Size
 CT/IV Infusion Pole & Bracket CT/FS - Flat Shelf - Removable

Features

- High level single & double column easy clean procedure trolleys
- Designed for guick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Seamlessly welded one piece frame from 304 grade stainless steel (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- Carerail fitted to each side of the trolley to support a wide range of accessories (see options)
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- · Each trolley is supplied with perforated coloured label strips (white, blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Single Column 550 x 550 x 1410mm
- Double Column 1005 x 550 x 1410mm
- Plastic Top

- Single Column 455 x 455 (internal) at height 1385mm
- Double Column 915 x 455 (internal) at height 1385mm Travs
- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)

CTS2/14/10S6D - Shown with optional CT/DH/A4

| Item | Description | Weight |
|--------------|----------------------------|--------|
| CTS1/14/2S5D | 2 x 100mm, 5 x 150mm Trays | 36.2kg |
| CTS1/14/3S4D | 3 x 100mm, 4 x 150mm Trays | 35.8kg |
| CTS1/14/4S3D | 4 x 100mm, 3 x 150mm Trays | 35.4kg |
| CTS1/14/5S3D | 5 x 100mm, 3 x 150mm Trays | 37.8kg |
| CTS1/14/6D | 6 x 150mm Trays | 35.0kg |
| CTS1/14/8S1D | 8 x 100mm, 1 x 150mm Trays | 37.8kg |
| | | |



| Item | Description |
|---------------|----------------------------|
| CTS2/14/4S10D | 2 x 100mm, 10 x 150mm Tra |
| CTS2/14/6S8D | 6 x 100mm, 8 x 150mm Trays |
| CTS2/14/8S6D | 8 x 100mm, 6 x 150mm Trays |
| CTS2/14/10S6D | 10 x 100mm, 6 x 150mm Tra |
| CTS2/14/12D | 12 x 150mm Trays |
| CTS2/14/16S2D | 16 x 100mm, 2 x 150mm Tra |
| | |

- Single & Double Doors See page 60
- CT/WF Work Flap
- GMBH/15 Glove / Mask Box Holder One Box, 15mm Hook
- SBH/15 Sharps Box Holder One Box, 15mm Hook
- CT/SBB/15 Sharps Box Bracket Universal (Daniels / Frontier / Sharpak)
- Polymer Castors See page 69







Caretray Racks - Stainless Steel, Mobile, Single & Double Column



.... 100

Easy disassembly - see page 38

Features

- Single & double column easy clean easy clean rack
- Designed for guick & easy complete disassembly when conducting
- thorough or deep cleaning programmes (including steam cleaning) Seamlessly welded one piece frame from 304 grade stainless steel (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg) • Each trolley is supplied with perforated coloured label strips (white,
- blue, green, yellow & red) • Trays can be subdivided in both length & width using simple slot in
- dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Single Column 550 x 550 x 1782mm
- Double Column 1005 x 550 x 1782mm
- Plastic Top
- Single Column 455 x 455 (internal) at height 1757mm
- Double Column 915 x 455 (internal) at height 1757mm
- Travs
- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)

Options

- CT/FS Flat Shelf Removable
- Polymer Castors See page 69

| ltem | Description | Weight |
|---------------|----------------------------|--------|
| CRSM1/18/2S6D | 2 x 100mm, 6 x 150mm Trays | 41.7kg |
| CRSM1/18/3S6D | 3 x 100mm, 6 x 150mm Trays | 44.1kg |
| CRSM1/18/4S5D | 4 x 100mm, 5 x 150mm Trays | 42.8kg |
| CRSM1/18/5S5D | 5 x 100mm, 5 x 150mm Trays | 44.3kg |
| CRSM1/18/6S4D | 6 x 100mm, 4 x 150mm Trays | 43.9kg |
| CRSM1/18/8D | 8 x 150mm Trays | 43.3kg |
| | | |



| Item | Description |
|-----------------|-------------------------|
| CRSM2/18/4S12D | 4 x 100mm, 12 x 150mm T |
| CRSM2/18/6S12D | 6 x 100mm, 12 x 150mm T |
| CRSM2/18/8S10D | 8 x 100mm, 10 x 150mm T |
| CRSM2/18/10S10D | 10 x 100mm, 10 x 150mm |
| CRSM2/18/12S8D | 12 x 100mm, 8 x 150mm T |
| CRSM2/18/16D | 16 x 150mm Trays |
| | |



- Easy to clean -

See page 69

Caretray Racks - Stainless Steel, Static, Single & Double Column



1 10 10 ile is i

Easy disassembly - see page 38

| Item | Description | Weight |
|---------------|----------------------------|--------|
| CRSS1/16/2S6D | 2 x 100mm, 6 x 150mm Trays | 38.4kg |
| CRSS1/16/3S6D | 3 x 100mm, 6 x 150mm Trays | 40.8kg |
| CRSS1/16/4S5D | 4 x 100mm, 5 x 150mm Trays | 39.5kg |
| CRSS1/16/5S5D | 5 x 100mm, 5 x 150mm Trays | 41.0kg |
| CRSS1/16/6S4D | 6 x 100mm, 4 x 150mm Trays | 40.6kg |
| CRSS1/16/8D | 8 x 150mm Trays | 40.0kg |

Features

- Single & double column easy clean easy clean rack
- Designed for quick & easy complete disassembly when conducting
- thorough or deep cleaning programmes (including steam cleaning) Seamlessly welded one piece frame from 304 grade stainless steel (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
 Each trolley is supplied with perforated coloured label strips (white,
- blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Adjustable feet allow levelling on uneven floors

Dimensions (w x d x h)

- Single Column 536 x 536 x 1670-1685mm
- Double Column 995 x 536 x 1670-1685mm
- Plastic Top
- Single Column 455 x 455 (internal) at height 1645mm
- Double Column 915 x 455 (internal) at height 1645mm
- Trays
- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray -355 x 405 x 150mm (internal)

Options

CT/FS - Flat Shelf - Removable



| CRSS2/1 | 6/8S1 |
|---------|-------|
| | |

| tem | Description | Weight |
|-----------------|------------------------------|--------|
| CRSS2/16/4S12D | 4 x 100mm, 12 x 150mm Trays | 71.1kg |
| CRSS2/16/6S12D | 6 x 100mm, 12 x 150mm Trays | 74.0kg |
| CRSS2/16/8S10D | 8 x 100mm, 10 x 150mm Trays | 72.3kg |
| CRSS2/16/10S10D | 10 x 100mm, 10 x 150mm Trays | 75.3kg |
| CRSS2/16/12S8D | 12 x 100mm, 8 x 150mm Trays | 74.6kg |
| CRSS2/16/16D | 16 x 150mm Trays | 71.5kg |
| | | |



Caretray Trolleys - Mild Steel, Low Level, Single & Double Column





CTM1/08/3S1D - Shown with optional CT/WF, SBH/15 & CT/AGH



Easy disassembly - see page 38

Features

- Low level single & double column easy clean procedure trolleys. Height makes it suitable to be stored underbench
- Designed for quick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Fully welded mild steel frame that is blast finished & powder coated in a durable epoxy paint (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- Carerail fitted to each side of the trolley to support a wide range of accessories (see options)
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- · Each trolley is supplied with perforated coloured label strips (white, blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each trav supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Single Column 550 x 550 x 807mm
- Double Column 1005 x 550 x 807mm
- Plastic Top
- Single Column 455 x 455 (internal) at height 782mm • Double Column - 915 x 455 (internal) at height 782mm
- Travs
- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

· Paint - Spearmint

| Item | Description | Weight |
|--------------|----------------------------|--------|
| CTM1/08/1S2D | 1 x 100mm, 2 x 150mm Trays | 20.9kg |
| CTM1/08/3D | 3 x 150mm Trays | 21.2kg |
| CTM1/08/3S1D | 3 x 100mm, 1 x 150mm Trays | 22.0kg |
| CTM1/08/4S | 4 x 100mm Trays | 21.6kg |



| Item | Description |
|--------------|----------------------------|
| CTM2/08/3S4D | 3 x 100mm, 4 x 150mm Trays |
| CTM2/08/4S3D | 4 x 100mm, 3 x 150mm Trays |
| CTM2/08/6D | 6 x 150mm Trays |
| CTM2/08/6S2D | 6 x 100mm, 2 x 150mm Trays |
| CTM2/08/7S1D | 7 x 100mm, 1 x 150mm Trays |
| CTM2/08/8S | 8 x 100mm Trays |

Options

- CT/AGH Alcohol Gel Holder
- CT/AIR Airway Tube Holder
- CT/BIN Debris Bin & Bracket
- CT/DH/A4 Document Holder, A4
- CT/CR2 & CT/CR5 Catheter Racks
- Single & Double Doors See page 60
- · CT/FS Removable Flat Shelf
- CT/IV Infusion Pole & Bracket
- CT/WF Work Flap
- GMBH/15 Glove / Mask Box Holder One Box, 15mm Hook
- SBH/15 Sharps Box Holder One Box, 15mm Hook
- CT/CYD Cylinder Holder C/CD/D Size CT/SBB/15 Sharps Box Bracket Universal (Daniels / Frontier / Sharpak) Polymer Castors - See page 69

| sales@bristolmaid.com +44 (0)1258 484455 www.bristolmaid.com | sales@bristolmaid.com | +44 (0)1258 484455 | www.bristolmaid.com |
|--|-----------------------|--------------------|---------------------|
|--|-----------------------|--------------------|---------------------|

50

Optional Polymer Castors - Easier to manoeuvre -- Reduced noise levels -- Easy to clean -See page 69





36.8kg

36.4kg

36.0kg

37.6kg

37.2kg

36.8kg

Caretray Trolleys - Mild Steel, Standard Level, Single & Double Column



Options

- · CT/AGH Alcohol Gel Holder
- CT/AIR Airway Tube Holder
- CT/BIN Debris Bin & Bracket
- CT/DH/A4 Document Holder, A4
- CT/CR/2 & CT/CR/5 Catheter Racks
- Single & Double Doors See page 60
- CT/FS Flat Shelf, Removable
- CT/IV Infusion Pole & Bracket

- CT/WF Work Flap
- GMBH/15 Glove / Mask Box Holder One Box, 15mm Hook
- SBH/15 Sharps Box Holder One Box, 15mm Hook
- CT/SBB/15 Sharps Box Bracket Universal (Daniels / Frontier / Sharpak)
- Laptop Holders See page 64
- CT/CYD Cylinder Holder C/CD/D Size CT/LH/C5 Motion C5 Holder, Swing Arm
 - CT/OB/MS Overbridge
 - Tilt Bins See page 70
 - Polymer Castors See page 69

Features

- Standard level single & double column easy clean procedure trolleys
- Designed for guick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Fully welded mild steel frame that is blast finished & powder coated in a durable epoxy paint (requires no on-site assembly)
- One piece moulded top with raised edges on three sides and integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- Carerails fitted to each side of the trolley to support a wide range of accessories (see options)
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- · Each trolley is supplied with perforated coloured label strips (white, blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Single Column 550 x 550 x 993mm
- Double Column 1005 x 550 x 993mm
- Plastic Top
- Single Column 455 x 455 (internal) at height 968mm
- Double Column 915 x 455 (internal) at height 968mm
- Trays
- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

· Paint - Spearmint

| Item | Description | Weight |
|--------------|----------------------------|--------|
| CTM1/10/1S3D | 1 x 100mm, 3 x 150mm Trays | 27.5kg |
| CTM1/10/3S2D | 3 x 100mm, 2 x 150mm Trays | 28.6kg |
| CTM1/10/4D | 4 x 150mm Trays | 27.9kg |
| CTM1/10/4S1D | 4 x 100mm, 1 x 150mm Trays | 28.3kg |
| CTM1/10/6S | 6 x 100mm | 29.4kg |



| Item | Description |
|--------------|----------------------------|
| CTM2/10/3S6D | 3 x 100mm, 6 x 150mm Trays |
| CTM2/10/4S5D | 4 x 100mm, 5 x 150mm Trays |
| CTM2/10/5S4D | 5 x 100mm, 4 x 150mm Trays |
| CTM2/10/6S4D | 6 x 100mm, 4 x 150mm Trays |
| CTM2/10/8D | 8 x 150mm Trays |
| CTM2/10/8S2D | 8 x 100mm, 2 x 150mm Trays |



Options

- CT/AGH Alcohol Gel Holder
- CT/AIR Airway Tube Holder
- CT/BIN Debris Bin & Bracket
- CT/DH/A4 Document Holder, A4
- CT/CR2 & CT/CR5 Catheter Racks
- CT/CYD Cylinder Holder C/CD/D Size
- Single & Double Doors See page 60

Caretray Trolleys - Mild Steel, High Level, Single & Double Column

Features

- High level single & double column easy clean procedure trolleys
- Designed for guick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Fully welded mild steel frame that is blast finished & powder coated in a durable epoxy paint (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- Carerail' fitted to each side of the trolley to support a wide range of accessories (see options)
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- · Each trolley is supplied with perforated coloured label strips (white, blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Single Column 550 x 550 x 1410mm
- Double Column 1005 x 550 x 1410mm
- Plastic Top
- Single Column 455 x 455 (internal) at height 1385mm
- Double Column 915 x 455 (internal) at height 1385mm
- Travs
- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

Paint - Spearmint

| | Item | Description | Weight |
|--|--------------|----------------------------|--------|
| | CTM1/14/2S5D | 2 x 100mm, 5 x 150mm Trays | 36.2kg |
| CT/FS - Flat Shelf - Removable CT/WF - Work Flap | CTM1/14/3S4D | 3 x 100mm, 4 x 150mm Trays | 35.8kg |
| GMBH/15 - Glove / Mask Box Holder - One Box, 15mm Hook | CTM1/14/4S3D | 4 x 100mm, 3 x 150mm Trays | 35.4kg |
| SBH/15 - Sharps Box Holder - One Box, 15mm Hook | CTM1/14/5S3D | 5 x 100mm, 3 x 150mm Trays | 37.3kg |
| CT/SBB/15 - Sharps Box Bracket - Universal (Daniels / Frontier / Sharpak) Polymer Castors - See page 69 | CTM1/14/6D | 6 x 150mm Trays | 35.0kg |
| | CTM1/14/8S1D | 8 x 100mm, 1 x 150mm Trays | 37.2kg |



| ltem | Description |
|---------------|--------------------------|
| CTM2/14/4S10D | 2 x 100mm, 10 x 150mm T |
| CTM2/14/6S8D | 6 x 100mm, 8 x 150mm Tra |
| CTM2/14/8S6D | 8 x 100mm, 6 x 150mm Tr |
| CTM2/14/10S6D | 10 x 100mm, 6 x 150mm T |
| CTM2/14/12D | 12 x 150mm Trays |
| CTM2/14/16S2D | 16 x 100mm, 2 x 150mm T |
| | |

Easy disassembly see page 38

Alternative

Page 410 🔎

Colours



Caretray Racks - Mild Steel, Mobile, Single & Double Column



CRMM1/18/6S4D

Features

- Single & double column easy clean easy clean rack
- Designed for quick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Fully welded mild steel frame that is blast finished & powder coated in a durable epoxy paint (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
 Each trolley is supplied with perforated coloured label strips (white,
- Each trolley is supplied with perforated coloured label strips (white, blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Single Column 550 x 550 x 1782mm
- Double Column 1005 x 550 x 1782mm
- Plastic Top
- Single Column 455 x 455 (internal) at height 1757mm
- Double Column 915 x 455 (internal) at height 1757mm
- Trays
- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)

Standard Colours Options

- Paint Spearmint
- CT/FS Removable Flat shelf
 Polymer Castors See page 69

| ltem | Description | Weight |
|---------------|----------------------------|--------|
| CRMM1/18/2S6D | 2 x 100mm, 6 x 150mm Trays | 41.2kg |
| CRMM1/18/3S6D | 3 x 100mm, 6 x 150mm Trays | 43.6kg |
| CRMM1/18/4S5D | 4 x 100mm, 5 x 150mm Trays | 42.3kg |
| CRMM1/18/5S5D | 5 x 100mm, 5 x 150mm Trays | 43.8kg |
| CRMM1/18/6S4D | 6 x 100mm, 4 x 150mm Trays | 43.4kg |
| CRMM1/18/8D | 8 x 150mm Trays | 42.8kg |



| tem | Description | Weight |
|-----------------|------------------------------|--------|
| CRMM2/18/4S12D | 4 x 100mm, 12 x 150mm Trays | 73.3kg |
| CRMM2/18/6S12D | 6 x 100mm, 12 x 150mm Trays | 76.3kg |
| CRMM2/18/8S10D | 8 x 100mm, 10 x 150mm Trays | 74.6kg |
| CRMM2/18/10S10D | 10 x 100mm, 10 x 150mm Trays | 77.6kg |
| CRMM2/18/12S8D | 12 x 100mm, 8 x 150mm Trays | 76.9kg |
| CRMM2/18/16D | 16 x 150mm Trays | 73.8kg |
| | | |

Optional Polymer Castors - Easier to manoeuvre -- Reduced noise levels -- Easy to clean -See page 69



Caretray Racks - Mild Steel, Static, Single & Double Column



111 444 100

Easy disassembly - see page 38

Features

- Single & double column easy clean easy clean rack
- Designed for quick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning) • Fully welded mild steel frame that is blast finished & powder coated in
- a durable epoxy paint (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
 Each trolley is supplied with perforated coloured label strips (white,
- blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Adjustable feet allow levelling on uneven floors

Dimensions

- Single Column 536 x 536 x 1670-1685mm
- Double Column 995 x 536 x 1670-1685mm
- Plastic Top
- Single Column 455 x 455 (internal) at height 1645mm
- Double Column 915 x 455 (internal) at height 1645mm
- Travs
- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

Paint - Spearmint

Options

CT/FS - Removable Flat Self

| ltem | Description | Weight |
|---------------|----------------------------|--------|
| CRMS1/16/2S6D | 2 x 100mm, 6 x 150mm Trays | 38.4kg |
| CRMS1/16/3S6D | 3 x 100mm, 6 x 150mm Trays | 39.5kg |
| CRMS1/16/4S5D | 4 x 100mm, 5 x 150mm Trays | 39.5kg |
| CRMS1/16/5S5D | 5 x 100mm, 5 x 150mm Trays | 41.0kg |
| CRMS1/16/6S4D | 6 x 100mm, 4 x 150mm Trays | 40.6kg |
| CRMS1/16/8D | 8 x 150mm Trays | 39.6kg |
| | | |



| em | Description | Weight |
|----------------|------------------------------|--------|
| RMS2/16/4S12D | 4 x 100mm, 12 x 150mm Trays | 70.2kg |
| RMS2/16/6S12D | 6 x 100mm, 12 x 150mm Trays | 73.2kg |
| RMS2/16/8S10D | 8 x 100mm, 10 x 150mm Trays | 72.3kg |
| RMS2/16/10S10D | 10 x 100mm, 10 x 150mm Trays | 75.3kg |
| RMS2/16/12S8D | 12 x 100mm, 8 x 150mm Trays | 74.6kg |
| RMS2/16/16D | 16 x 150mm Trays | 71.5kg |
| | | |





SARAN CTM2/14/6S8D/CTD2/14SS

(Double door fitted with Security Seal Lock)

Doors & Back Panels

Features

- Machined solid grade laminate doors allowing for easy cleaning
- Doors fold back against the sides of the trolley when open
- Double door units feature a return to ensure that only one lock is required to secure both doors
- Three locking options
- CAM Lock with Wing Handle
- Supplied with two keys
- Wing handle can operate independently of the key
- Locks keyed to differ
- Keyed alike on request
- Security Seal
 - Accepts standard plastic security seals
- Security seals not included with door
- Electronic Push Button
- Units fitted with a door also supplied with a removable back panel to ensure contents of the trolley are secure Factory fitted

Locking Options





Electronic

Push Button



CTM1/08/3S1D/CTD1/08CA (Single door fitted with CAM Lock with Wing Handle)

Electronic Push Button Lock

- Uses a 10 button keypad and a # button. The # button is hidden beneath the handle when in the locked position
- Easy programming & code change
- Four operating codes with differing levels of authority
 - Master code Sub-master code

 - User code
- Technician code User codes are 4 digits
- Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code


CTM2/10/6S4D/CTD2/10EPB (Double door fitted with Electronic Push Button Lock)

| Description | CAM Lock with Wing Handle | Security Seal | Electronic Push Button | Trolley Compatibility |
|--------------------------|---------------------------|---------------|------------------------|-----------------------|
| Single Door & Back Panel | CTD1/08CA | CTD1/08SS | CTD1/08EPB | CTM1/08, CTS1/08 |
| Double Door & Back Panel | CTD2/08CA | CTD2/08SS | CTD2/08EPB | CTM2/08, CTS2/08 |
| Single Door & Back Panel | CTD1/10CA | CTD1/10SS | CTD1/10EPB | CTM1/10, CTS1/10 |
| Double Door & Back Panel | CTD2/10CA | CTD2/10SS | CTD2/10EPB | CTM2/10, CTS2/10 |
| Single Door & Back Panel | CTD1/14CA | CTD1/14SS | CTD1/14EPB | CTM1/14, CTS1/14 |
| Double Door & Back Panel | CTD2/14CA | CTD2/14SS | CTD2/14EPB | CTM2/14, CTS2/14 |
| | | | | |









Laptop Holder - Removable, Swing Arm

Features

- Three section swing arm offering height, rotation & angle adjustment
 Compatible with models CTM1/10, CTM2/10, CTS1/10, CTS2/10
- Laptop holder fitted to the rear right upright of the trolley
- Adjustable to accept laptops 280-435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Laptop supported in a cradle fitted with adjustable wing arms
- Laptop can be secured using a security clamp (supplied)*
- Design allows easy access to laptop ports
- *Not compatible with the Dell laptops using the N17 wedge lock system

Standard Colours

· Paint - Grey White

| Item | Description | Weight |
|----------|----------------------|--------|
| CT/LHR/S | Removable, Swing Arm | 4.4kg |

Laptop Holder - Non-Removable, Swing Arm

Features

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with models CTM1/10, CTM2/10, CTS1/10, CTS2/10
- Laptop holder fitted to the rear right upright of the trolley
- Adjustable to accept laptops 280-435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Laptop supported in a cradle fitted with a steel mounting bracket that can be fixed using either nylock or wing nuts
- Laptop can be secured using a security clamp (supplied)*
- Design allows easy access to laptop ports
- *Not compatible with the Dell laptops using the N17 wedge lock system

Standard Colours

· Paint - Grey White

| Item | Description | Weight |
|-----------|--------------------------|--------|
| CT/LHNR/S | Non-Removable, Swing Arm | 5.7kg |



+44 (0)1258 484455 sales@bristolmaid.com www.bristolmaid.com

Motion C5 Holder - Swing Arm

Features

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with models CTM1/10, CTM2/10, CTS1/10, CTS2/10
- Motion C5 holder fitted to the rear right upright of the trolley
- Maximum capacity 2.0kg
- Tablet located in a RAM-In-Motion passive dock allowing for easy insertion & removal

Standard Colours

Paint - Grey White



| Item | Description | Weight |
|----------|-------------|--------|
| CT/LH/C5 | Swing Arm | 3.4kg |

Work Flap

Features

- Fold down flap (provides work area 300 x 455mm)
- Machined from solid grade laminate allowing for
- easy cleaningAttaches to upper carerail
- Maximum capacity 2kg

Standard Colours

- Paint Grey White
- Solid Grade Laminate Grey White

Dimensions (w x d x h)



Infusion Pole & Bracket

Features

- Two hook infusion pole & bracket (pole length 1160mm)
- Attaches to upper carerail
- Maximum capacity on hooks 2kg

Standard Colours

Paint - Grey White

Dimensions (w x d x h)

• 238 x 104 x 1157mm

Cylinder Holder - C/CD/D Size

Features

- Compatible with C/CD/D size cylinders (max diameter 100mm)
- Cylinder held in position using a chain, whilst rubber pads
- protect wear points & reduce noise
- Attaches to upper carerail

Standard Colours

Paint - Grey White

Dimensions (w x d x h)



| ltem | Description | Weight |
|--------|-------------------------------|--------|
| CT/CYD | Cylinder Holder - C/CD/D Size | 1.4kg |

Airway Tube Holder

Features

- Internal dimensions of holder 95 x 50 x 510mm
- Outer sleeve removable for cleaning
- Fitted with a 15mm square hook
- Attaches to upper carerail

Standard Colours Dimensions (w x d x h)

Paint - Grey White
 140 x 70 x 555mm



| Item | Description | Weight |
|--------|--------------------|--------|
| CT/AIR | Airway Tube Holder | 1.4kg |

Weigh

1.6kg

Alcohol Gel Holder

Features

- Accepts bottles with maximum 90 x 90mm square base
- Fitted with a 15mm hook
- Attaches to upper/lower carerails

Standard Colours

-)
- Stainless Steel

Dimensions (w x d x h)

• 95 x 110 x 115mm



Alcohol gel not supplied

Glove / Mask Box Holder One Box, 15mm Hook

Features

- · Accommodates glove & mask boxes up to 150mm in width
- One piece construction
- 15mm hook attaches to upper/lower carerails

Standard Colours

· Paint - Grey White

Dimensions (w x d x h)

• 143 x 122 x 177mm







| Item | Description | Weight | Item | Description | Weight |
|--------|--------------------|--------|---------|--|--------|
| CT/AGH | Alcohol Gel Holder | 0.3kg | GMBH/15 | Glove / Mask Box Holder - One Box, 15mm Hook | 0.6kg |

Sharps Box Holder One Box, 15mm Hook

Features

- Accepts sharps boxes up to seven litres capacity
- Sharps box secured using a rubber band fitted with a toggle
- 15mm hook attaches to upper/lower carerails

Standard Colours

· Paint - Grey White

Dimensions (w x d x h)

• 127 x 118 x 230mm



| Item | Description | Weight |
|--------|-------------------|--------|
| SBH/15 | Sharps Box Holder | 0.5kg |

Sharps Box Bracket - Universal

Features

- Compatible with the following brackets (customer to supply the plastic brackets & fit fixings supplied).
- Daniels Healthcare
- POUDS® Wall-Bracket (up to 5L) SHARPSGUARD® mini, 1, 2.5, 3.75, 5 & com plus
- POUDS® Wall-Bracket (7L & above) SHARPSGUARD® 7, 11.5 & 22
- Frontier (Sharpsafe) 2, 3, 4 & 7
- Sharpak Healthcare Sharpak 18, 35, 60, 120 & 220 PLUS
- 15mm hook attaches to upper/lower carerails

Standard Colours

Paint - Grey White

Dimensions (w x d x h)

Debris Bin & Bracket

Features

- Plastic debris bin & bracket
- 15mm hook attaches to upper/lower carerails

Standard Colours

Paint - Grey White

Dimensions (w x d x h)

Overall Dimensions - 210 x 164 x 109mm
Internal Dimensions - 195 x 130 x 95mm



| Item | Description | Weight |
|--------|----------------------|--------|
| CT/BIN | Debris Bin & Bracket | 0.4kg |

Catheter Racks

Features

- Internal dimensions of each section 64 x 75 x 400mm Removable clear perspex panel to aid cleaning
- 15mm hook attaches to upper/lower carerails

Standard Colours

Paint - Grey White

Dimensions (w x d x h)

• CT/CR2 - 130 x 108 x 405mm • CT/CR5 - 325 x 108 x 405mm



| Item | Description |
|--------|--------------|
| CT/CR2 | Two Section |
| CT/CR5 | Five Section |



Plastic brackets not supplied

| Item | Description | Weig |
|-----------|--------------------------------|------|
| CT/SBB/15 | Sharps Box Bracket - Universal | 0.2 |

Document Holder

Features

- A4 size suitable for A4 ring binders or foolscap notes (pocket size 274 x 63mm)
- Open front design allows for easy patient identification
- 15mm hook attaches to upper/lower carerails

Standard Colours

· Paint - Grey White

Dimensions (w x d x h)

• 275 x 82 x 402mm



Features

- Solid grade laminate (max load per shelf 5kg)
- Supported on the runner guides (non pull out)

Standard Colours

Solid Grade Laminate - Grey White

Dimensions (w x d x h)

• 413 x 457 x 4mm







| Veight | |
|--------|--|
| 1.4kg | |
| 3.0kg | |

| ltem | Description | Weight |
|----------|-----------------|--------|
| CT/DH/A4 | Document Holder | 1.7kg |

| Item | Description | Weight |
|-------|------------------------|--------|
| CT/FS | Flat Shelf - Removable | 1.1kg |

Tilt Bins

Features

- Available with 4, 5 or 6 tilting compartments
- Clear tilt bins provide visible storage & assists organisation
- Compartments tilt forward to a 45° angle allowing easy access of contents & can be removed to aid refiling/cleaning
- Label holder for easy identification of the contents
- Retaining bar secures compartments in place, preventing bins from tipping forward during transit
- Supplied with a mild/stainless steel mounting plate
- Attaches to an overbridge (order separately)

Dimensions (w x d x h) Standard Colours

Compartment Dimensions

- TBIN4 143 x 127 x 190mm
- TBIN5 108 x 95 x 146mm
- TBIN6 89 x 67 x 95mm

45° Angle Tilt

Label Holder

- Plastic Cream Paint - Grev White
- Stainless Steel Self Colour









| Mild Steel | Stainless Steel | Description | Dimensions (w x d x h) | Weight |
|------------|-----------------|------------------|------------------------|--------|
| TBIN4/MS | TBIN4/SS | Four Compartment | 600 x 171 x 207mm | 2.3kg |
| TBIN5/MS | TBIN5/SS | Five Compartment | 600 x 135 x 164mm | 1.6kg |
| TBIN6/MS | TBIN6/SS | Six Compartment | 600 x 93 x 113mm | 0.9kg |

Overbridge

Features

- Fixes onto the rear uprights of the trolley (increases overall height by 625mm)
- Available as either painted mild steel or self-colour stainless steel
- Provides a fixing location for tilt bins (order separately)
- Factory fitted

Standard Colours

· Paint - Grev White Stainless Steel - Self Colour



| tem | Description | Weight |
|----------|-----------------------------|--------|
| CT/OB/MS | Overbridge, Mild Steel | 2.0kg |
| CT/OB/SS | Overbridge, Stainless Steel | 2.0kg |

Additional Trays

Features

- 100/150mm deep translucent high impact plastic trays with integral handle (max load per tray 5kg)
- Each tray is supplied with plastic clip on label cover & four removable dividers (2 x long & 2 x short)
- Trays are supplied without coloured label strips (see CT/LAB)

Internal Dimensions (w x d x h)

- CT/TRAY/S 355 x 405 x 100mm
- CT/TRAY/D 355 x 405 x 150mm



| Item | Description |
|-----------|--------------------|
| CT/TRAY/S | 100mm Shallow Tray |
| CT/TRAY/D | 150mm Deep Tray |

Additional Dividers

Features

- 100/150mm deep translucent dividers



| Item | Description | Weight |
|----------|--|--------|
| CT/DIV/S | Set of 4 dividers for 100mm Shallow Tray | 0.6kg |
| CT/DIV/D | Set of 4 dividers for 150mm Deep Tray | 0.4kg |

Additional Label Strips

Coloured Label Strips



Weigh 1.5kg 1.8kg

- Features
- Coloured label strips
- Red
- Yellow
- Green
- Blue

Item

CT/LAB

- White
- Supplied as a set of 6 sheets



Features

- Offers reduced rolling resistance making the trolley easier to manoeuvre, lower noise levels & smooth surfaces to assist cleaning
- Supplied as a set of 4 castors

Polymer Castors

- 2 x braking / 2 x non-braking
- 4 x braking
- 4 x non-braking
- Factory fitted



| Item | Description | Weight |
|------|---|--------|
| C004 | 100mm, Grey White Polymer, 2 x Braking / 2 x Non-Braking (set of 4) | 2.0kg |
| C005 | 100mm, Grey White Polymer, 4 x Braking (set of 4) | 2.0kg |
| C006 | 100mm, Grey White Polymer, 4 x Non-Braking (set of 4) | 2.0kg |

sales@bristolmaid.com +44 (0)1258 484455 www.bristolmaid.com



Chart Trolleys Workstations **Room Dividers** & Screens

- Fixed Height Trolleys 72
- Fixed Height Workstations 74
- Variable Height Workstations 75
 - Chart Boards 76
 - Document Holders 77
 - Room Dividers 78
 - Curtain Screens 80
 - Solid Screens 81



Caretray Trolley - Chart

Features

- Single & double column fixed height chart trolleys
- Angle adjustable chartboard fitted with a spring clip (adjustable 7-45 degrees)
- Usable chart area...
- Single Column 500 x 530mm (w x h)
- Double Column 800 x 530mm (w x h)
- Lower lip on chartboard to support pens etc.
- Designed for guick & easy disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Fully welded mild steel frame that is blast finished & powder coated in a durable epoxy paint (requires no on-site assembly)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- Carerails fitted to each side of the trolley (2 x L/H, 1 x R/H) to support a wide range of accessories (see options)
- Drop down work flap fitted to R/H side
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover (max load per tray 5kg)
- Each trolley is supplied with perforated coloured label strips white, blue, green, yellow & red (not shown on trolley)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions
- Single Column 550 x 665 x 1115-1450mm
- Double Column 1005 665 x 1115-1450mm
- Usable Chart Area
- Single Column 500 x 530mm
- Double Column 800 x 530mm
- Work Flap 420 x 290mm (w x d)
- Trays
- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

- Paint White
- Chartboard/Work Flap White

| em | Description | Weight |
|-----------------|--|--------|
| TCHARM1/10/3S2D | Caretray Trolley - Chart, 3 x 100mm, 2 x 150mm Trays | 33.2kg |
| TCHARM1/10/4S1D | Caretray Trolley - Chart, 4 x 100mm, 1 x 150mm Trays | 32.9kg |

Optional Polymer Castors Easier to manoeuvre Reduced noise levels Easy to clean See page 69



Angle adjustable chartboard fitted with spring clip/s (adjustable 7-45 degrees)

Options

- CT/DH/A4 Document Holder, A4
- CT/AGH Alcohol Gel Holder
- CT/FS Flat Shelf Removable
- DAR/1/15 Apron Roll Holder
- DG/3/15 Three Glove Box Holder
- GMBH/15 Glove / Mask Box Holder One Box, 15mm Hook
- CT/SBB/15 Sharps Box Bracket Universal (Daniels / Frontier / Sharpak)
- SBH/15 Sharps Box Holder One Box, 15mm Hook
- ZP010 100mm Ring Buffers
- C004, C005, C006 Polymer Castors See page 69

Ú,

చ్



| Item | Description | Weight |
|------------------|--|--------|
| CTCHARM2/10/6S4D | Caretray Trolley - Chart, 6 x 100mm, 4 x 150mm Trays | 53.2kg |
| CTCHARM2/10/8S2D | Caretray Trolley - Chart, 8 x 100mm, 2 x 150mm Trays | 51.2kg |



Chart Workstations - Fixed Height

Features

- Angle adjustable chartboard fitted with a spring clip (adjustable 0-60 degrees)
- Usable chart area 800 x 530 (w x h)
- Lower lip on chartboard to support pens etc
- Models available with:
- Chartboard only
- Drawer
- · Rails fitted to each side to support a wide range of accessories (see options)
- 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions 820 x 680 x 1040 1510mm
- Usable Chart Area 800 x 530mm (w x h)
- Drawer 650 x 397 x 100mm (internal)

Standard Colours

- Paint White
- Chartboard/Work Flap White

Options

- CBT900 Work Flap
- CBT905 Lower Tray
- CT/DH/A4 Document Holder, A4
- CT/AGH Alcohol Gel Holder
- DG/3/15 Three Glove Box Holder
- GMBH/15 Glove / Mask Box Holder One Box, 15mm Hook
- CT/SBB/15 Sharps Box Bracket Universal (Daniels / Frontier / Sharpak)
- SBH/15 Sharps Box Holder One Box, 15mm Hook
- ZP010 100mm Ring Buffers

| Item | Description | Weight |
|--------|--|--------|
| CBT050 | Chart Workstation - Fixed Height | 25.3kg |
| CBT055 | Chart Workstation - Fixed Height, Drawer | 35.0kg |

Chart Workstations - Variable Height

Features

- Height adjustable using a gas damper (795-1075mm)
- Angle adjustable chartboard fitted with spring clips (adjustable 0-90 degrees)
- Chartboard available in three sizes (w x h)
- Small 790 x 700mm
- Medium 790 x 800mm
- Large 790 x 900mm
- Lower lip on chartboard to support pens etc
- Slide out shelf at front
- Document storage pocket on side to accept patients notes, ring binders, etc
- 80mm swivel castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions 800 x 700/800/900 x 1095mm
- Document Storage Compartment 75 x 255 x 360mm (internal)

Standard Colours

- Paint White
- Chartboard White

Options

ZP010 - 100mm Ring Buffers

| Item | Description | Weight |
|----------|------------------|--------|
| CB/100/A | 700mm High Board | 40.0kg |
| CB/100/B | 800mm High Board | 41.0kg |
| CB/100/C | 900mm High Board | 42.5kg |





Chart Boards - A4 & A3

Features

- A4 one spring clip (holds one chart 220 x 280mm)
 A3 two spring clips (holds one chart 450 x 280mm or two charts 220 x 280mm)
- Available with a 30 or 50mm mild steel square hook

Standard Colours

White



| Document | Holders | - A4 | & A3 |
|----------|---------|------|------|
|----------|---------|------|------|

Features

- A4 size suitable for A4 ring binders or foolscap notes (pocket size 260 x 65mm)
 A3 size suitable for A3 notes and x-rays (pocket size 398 x 100mm)
- Open front design allows for easy patient identification
 Available for wall fitting or with a 50mm square hook

Standard Colours

Paint - Grey White

| С | B | 0 | 0 | 5 |
|---|---|---|---|---|
| | | | | |

| Item | Description | Overall Dimensions (w x d x h) | Weight |
|--|---------------------------------|--------------------------------|--------|
| A4 Chart BoardCB005A4, 30mm Square Mild Steel HookCB010A4, 50mm Square Mild Steel Hook | | | |
| | | 230 x 80 x 370mm | 0.5kg |
| | | 230 x 100 x 370mm | 0.5kg |
| A3 Chart Board | | | |
| CB015 | A3, 30mm Square Mild Steel Hook | 460 x 80 x 370mm | 0.8kg |
| CB020 | A3, 50mm Square Mild Steel Hook | 460 x 100 x 370mm | 0.8kg |

| ltem | Description | Overall Dimensions (w x d x h) | Weight |
|--------------------|----------------------|--------------------------------|--------|
| A4 Docume | ent Holder | | |
| DH/A4/WF | A4, Wall Fitting | 280 x 70 x 400mm | 1.4kg |
| DH/A4/50 | A4, 50mm Square Hook | 280 x 125 x 405mm | 1.4kg |
| A3 Document Holder | | | |
| DH/A3/WF | A3, Wall Fitting | 405 x 105 x 400mm | 2.9kg |
| DH/A3/50 | A3, 50mm Square Hook | 405 x 160 x 405mm | 3.0kg |
| | | | |

6





Room Dividers

Features

- Fixed height mobile room divider/workstation
- Ideal for use in critical care areas A&E, ICU, ICT, HDU
- Modular design allows the unit to be configured to suit your specific requirements
- Frame, worktop & upstand/panel available as...
- Full width upstand
- Full width upstand & right hand mounted panel
- Full width upstand & left hand mounted panel
- Full width panel
- Four configurations of lower units...
- Three drawer unit, upper drawer fitted with an electronic push button lock (2 x 170, 1 x 285mm drawers)
- Five drawer unit, upper drawer fitted with an electronic push button lock (5 x 125mm drawers)
- Cupboard with double doors & removable shelf
- Cover panel (creates a space within the unit) • Easy clean solid grade laminate work surface & up-
- stand/panel
- Upper panel fitted with
- Triple alove box holder
- Alcohol gel holder (compatible with 500ml bottles, 70 x 70 x 225mm)
- Apron roll dispenser
- 100mm ring buffers to protect the fabric of the building
- 100mm swivel castors, 2 off braking

Standard Colours

- Paint White
- Worksurface / Up-Stand / Panel White


Dimensions (w x d x h)

- RD1 1768 x 698 x 1000mm
- RD2R/RD2L 1768 x 698 x 1588mm
- RD3 1768 x 698 x 1350mm

Options

- CC Colour coded drawers available on request
- DTSF/RD/450/2 Dressing Trolley
- Seamlessly welded 304 grade stainless steel frame
- Fixed shelves with down turned sides
- Shelf Heights 375 & 840mm
- Max load per shelf 40kg

Configure Your Workstation in Four Easy Steps...









| CITERID . | | | |
|-----------|----|---|----|
| | | | 1 |
| | 1. | _ | |
| | | | 1 |
| | | | 11 |
| | | | 1 |



Curtain Screens - Four Section, Folding

Features

- PVC fire retardant curtain (machine washable)
- Available assembled or for self-assembly
- Overall screen height 1708mm
- Distance from floor to lowest point of screening 288mm
- 50mm twin wheel castors

Dimensions (w x d x h)

- Overall Dimensions 2450 x 420 x 1790mm
- Folded 650 x 420 x 1790mm

Standard Colours

- · Paint White
- Curtain Grey

Options

• WS/SPARES/001 - Set of 4 replacement curtains



Features

- Solid plastic screen
- Smooth surface to assist cleaning
- Overall screen height 1680mm
- Distance from floor to lowest point of screening 160mm
- 50mm twin-wheel castors

Dimensions (w x d x h)

- Overall Dimensions
- Three Section 1802 x 438 x 1680mm
- Four Section 2403 x 438 x 1680mm
- Folded 600 x 438 x 1680mm

Standard Colours

- Paint White
- Screen White



| Item | Description | Weight |
|---------|---|--------|
| WSA/4/C | Curtain Screen - Four Section, Folding, Assembled | 14.6kg |
| WSK/4/C | Curtain Screen - Four Section, Folding, Self-Assembly | 14.6kg |

| Item | Description | |
|----------|---------------------------------------|--|
| SS/44/H3 | Solid Screen - Three Section, Folding | |
| SS/44/H4 | Solid Screen - Four Section, Folding | |

Solid Screens - Three & Four Section, Folding

20.6kg 27.0kg









Concealment Trolleys



Variable Height Trolleys 84

- Fixed Height Trolleys 85
- Body Trays & Covers 85

Variable Height - Fixed Body Tray

Features

- Height adjustable, 580 -1280mm to trollev surface
- Scissor action provides vertical lift, ideal for use in confined areas where space is limited Battery powered actuator
- Rechargeable battery is charged through the handset using a UK 3 pin plug
- Audible warning when battery power is at 50%
- Liquid retaining fixed tray, 1900 x 540 x 25mm
- Drop down side-rails, distance between rails 660mm
- Heavy duty reinforced PVC cover (locates over side rails)
- Bumper bar/push handle at either end
- 200mm heavy duty castors, 2 x braking, 2 x directional locking
- Maximum capacity 190kg

Dimensions (w x d x h)

• Overall Dimensions - 2150 x 770 x 580-1280mm

Standard Colours



| Item | Description | Weigh |
|--------|----------------------------------|-------|
| VH/950 | Variable Height, Fixed Body Tray | 82.0k |

| /ariable | Height - | Loose | Body | Tray |
|----------|----------|-------|------|------|
|----------|----------|-------|------|------|

Features

- · Height adjustable, 610 1310mm
- Scissor action provides vertical lift, ideal for use in confined areas where space is limited
- Battery powered actuator
- Rechargeable battery is charged through the handset using a UK 3 pin plug
- Audible warning when battery power is at 50%
- Liquid retaining loose tray, 2100 x 560 x 70mm
- Tray moves on rollers and is positively locked in position by a spring-loaded clamp at each end of the trolley
- Tray guides adjustable 400-450mm to accommodate trays by others (max width 650mm)
- Drop down side-rails, distance between rails 660mm • Heavy duty reinforced PVC cover (locates over side rails)
- Bumper bar/push handle at either end
- Compatible with most other common body trays
- 200mm heavy duty castors, 2 x braking, 2 x directional locking
- Maximum capacity 190kg

Dimensions (w x d x h)

• Overall Dimensions - 2270 x 770 x 610 - 1310mm

Standard Colours

· Paint - White

VH/950/LBT



Variable Height, Loose Body Tray

Weight

92.0kg

Description

Fixed Height - Fixed / Loose Body Tray

Features

- Available with fixed or removable liquid retaining travs Fixed Trav
- Removable Tray Tray rests on top of frame between support lugs
- Drop down side-rails, distance between rails 660mm
- Heavy duty reinforced PVC cover (locates over side rails)
- Bumper bar/push handle at either end
- 200mm heavy duty castors, 2 x braking, 2 x directional locking
- Maximum capacity 190kg



| Item | Description | Overall Dimensions (w x d x h) | Weight |
|---------------|-------------------------------|--------------------------------|--------|
| WT/950/DS | Fixed Height, Fixed Body Tray | 2076 x 751 x 690mm | 55.0kg |
| WT/950/DS/LBT | Fixed Height, Loose Body Tray | 2076 x 751 x 750mm | 61.0kg |

Ĕ

Body Tray

Features

- Heavy duty body tray, mild steel frame & stainless steel tray
- Integral runners
- Fits most standard mortuary racks (check prior to ordering)
- Maximum capacity 190kg

Dimensions (w x d x h)

Overall Dimensions - 2100 x 550 x 70mm

Standard Colours

Paint - White



| Item | Description | Weight |
|--------------|-------------|--------|
| WT/BODY/TRAY | Body Tray | 16.0kg |

Cover

Features

 Heavy duty reinforced PVC cover (locates over side rails)



| Item | Description | Weight |
|----------|-----------------------------------|--------|
| 2FACO024 | Cover - PVC Heavy Duty Reinforced | 5.0kg |



Cylinder Storage & Transportation

- Cylinder Trolleys F/G Size 88
 - Cylinder Trolleys J Size 89
- Bulk Transportation Trolleys 89
- Cylinder Trolleys E, CD/D, AF Size 90
 - Wall Brackets & Cylinder Holders 91

Cylinder Trolley F/G Size

Features

Compatible with F (48 cu ft), HX, LF, VF, ZX & G (120 cu ft) cylinders, max diameter 200mm, height 1350mm

Description

10.0kg

F/G Size

- One piece construction
- Drop over chains secure the cylinder
- 160mm anti-static wheels

Dimensions (w x d x h)

• 510 x 405 x 1160mm

Standard Colours

Paint - Grey

Cylinder Trolley with Stabilisers F/G Size

Features

- Compatible with F (48 cu ft), HX, LF, VF, ZX & G (120 cu ft) cylinders, max diameter 200mm, height 1350mm
- One piece construction
- Drop over chains secure the cylinder
- Stabilising castors at rear of trolley reduce the likelihood of tipping during transportation
- 160mm anti-static wheels, 100mm swivel castors

Dimensions (w x d x h)

520 x 700 x 1160mm

Standard Colours

• Paint - Grey



CYT/FG1/ST



Features

- Compatible with F (48 cu ft), HX, LF, VF, ZX & G (120 cu ft) cylinders, max diameter 200mm, height 1350mm
- One piece construction
- Drop over chains secure the cylinders
- Stabilising castors at rear of trolley reduce the likelihood of tipping during transportation
- 200mm anti-static wheels, 100mm swivel castors

Dimensions (w x d x h)

• 780 x 740 x 1145mm

Standard Colours

Paint - Grey



Cylinder Trolley J Size

Features

- Compatible with J (240 cu ft), AK, L & HL cylinders, max diameter 245mm, height 1550mm
- One piece construction
- Drop over chains secure the cylinder
- 160mm anti-static wheels

Dimensions (w x d x h)

• 510 x 475 x 1280mm

Standard Colours



| Item | Description | Weight |
|--------|-------------|--------|
| CYT/J1 | J Size | 10.5kg |

F/G Size

Weight

14.0kg

CYT/FG1

Cylinder Trolley with Stabilisers J Size

Features

- Compatible with J (240 cu ft), AK, L & HL cylinders, max diameter, max diameter 245mm, height 1550mm
- One piece construction
- · Drop over chains secure the cylinder
- Stabilising castors at rear of trolley reduce the likelihood of tipping during transportation
- 160mm anti-static wheels, 100mm swivel castors

Dimensions (w x d x h)

• 510 x 740 x 1275mm

Standard Colours

Paint - Grey

Bulk Transportation Trolley D & E Cylinders, 12 x Capacity

Features

- Compatible with D, AD, CD, DD, PD, RD, ZD & E (24 cu ft) cylinders (12 x capacity) max diameter 105mm, height 870mm
- One piece construction
- Cylinders locate into individual sections, supported on a non slip rubber base
- 125mm swivel castors, 2 off braking

Dimensions (w x d x h)

Standard Colours

• 460 x 590 x 990mm

Paint - Grey



| Item | Description | Weight | ltem | Description | Weight |
|-----------|-------------|--------|----------|--------------------------------|--------|
| CYT/J1/ST | J Size | 14.8kg | CYT/DE12 | D & E Cylinders, 12 x Capacity | 23kg |

Cylinder Trolley E Size

Features

- Compatible with E (24 cu ft) cylinders. max diameter 115mm, height 900mm
- One piece construction
- Recessed base plate and drop over chain to secure the cylinder
- 40mm hook allows trolley to hang at bed end, a useful feature when transporting patients
- 100mm anti-static wheels

Dimensions (w x d x h)

• 320 x 235 x 980mm

Standard Colours

Paint - Grev

Cylinder Trolley CD/D Size

Features

- · Compatible with CD & D cylinders, max diameter 115mm, height 540mm
- One piece construction
- Recessed base plate and drop over chain to secure the cylinder
- 40mm hook allows trolley to hang at bed end, a useful feature when transporting patients
- 100mm anti-static wheels

Dimensions (w x d x h)

• 320 x 235 x 980mm Standard Colours Paint - Grey

| Item | Description | Weigh |
|---------|-------------|-------|
| CYT/CD1 | CD/D Size | 4.3k |

Cylinder Trolley AF Size

Features

- · Compatible with AF, DF & ZH cylinders, max diameter 178mm, height 900mm
- One piece construction
- Recessed base plate and drop over chain to secure the cylinder 100mm anti-static wheels

Dimensions (w x d x h)

• 390 x 310 x 980mm

Standard Colours

Paint - Grev



| Item | Description | Weight |
|---------|-------------|--------|
| CYT/AF1 | AF Size | 4.9kg |

Wall Brackets Cvlinder

Features

- Max cylinder diameters
- D/E 112mm, F 148mm, G 185mm, J 235mm
- One piece construction
- Drop over chains to secure the cylinder
- Holes for wall mounting (fixings not provided)
- Should be positioned on the wall at approx half cylinder height

Standard Colours

Paint - Grey



| ltem | Description | Dimensions (w x o |
|---------|-------------------|-------------------|
| CYB/DE2 | 2 x D/E Cylinders | 359 x 107 x 88mm |
| CYB/F2 | 2 x F Cylinders | 435 x 126 x 88mm |
| CYB/G2 | 2 x G Cylinders | 511 x 145 x 88mm |
| CYB/J2 | 2 x J Cylinders | 613 x 170 x 88mm |
| | | |



| Item | Description | Weight |
|--------|-------------|--------|
| CYT/E1 | E Size | 4.3kg |

+44 (0)1258 484455 sales@bristolmaid.com www.bristolmaid.com

Features

- Compatible with.
- D Size D, AD, CD, DD, PD, RD & ZD cylinders (max 112mm diameter)
- E Size E (max 112mm diameter)
- One piece construction
- Carry handle and extended hook
- Suitable for beds with 25mm wide side/end members
- Holder can be hung and/or bolted to 25mm section bed-frame

Dimensions (w x d x h)

- CYH/D1 200 x 160 x 290mm
- CYH/E1 200 x 160 x 595mm

Standard Colours

Paint - Grey



Features

- · Painted stainless steel wall plate to support cylinder holders
- Compatible with CYH/D1. CYH/E1. CYHP/D1 & CYHP/E1

Dimensions (w x d x h)

• 240 x 47 x 254mm

Standard Colours

Paint - Grey



| Weight |
|--------|
| 1.1kg |
| 1.4kg |
| 1.7kg |
| 2.1kg |
| |



Weight

1.4kg

2.5kg



Dispensers, Holders & Racks

- Apron Dispensers 94
- Pocket Mask Holders 94
 - Glove Box Holders 95
- Alcohol Gel Holders 96
- Syringe Dispenser 96
- Sharps Box Holder 97
- Urine Bottle Holder 97
- Catheter Bag Holders 97
 - Catheter Racks 98
 - Drainage Racks 99

Apron Dispensers - Roll

Features

Apron roll holders, maximum width 415mm
Holes for wall mounting (fixings not provided)

Dimensions (w x d x h)

- One Roll 420 x 168 x 268mm
- Three Roll 420 x 168 x 615mm

Standard Colours

· Paint - White



| Item | Description | Weight |
|-------|-------------|--------|
| DAR/1 | One Roll | 2.0kg |
| DAR/3 | Three Rolls | 8.2kg |

Apron Dispenser - Flat Pack

Features

- Holder to accept flat packed aprons, maximum width 400mm
- Packs of aprons fold over the centrally positioned divider
- Holes for wall mounting (fixings not provided)

Dimensions (w x d x h)

• 410 x 105 x 290mm

Standard Colours

Paint - White





Features

- Holder to accept Laerdal & Ambu Pocket Masks (hard cases)
- Allows the case to be quickly removed from the holder
- Supplied with self adhesive tape/holes for wall mounting (fixings not provided)

Dimensions (w x d x h)

• 110 x 40 x 110mm

Standard Colours

• Paint - White



PMH005 -Pocket mask not supplied



| Description | Weight | Item | Description | Weight |
|-------------|--------|--------|--------------------|--------|
| Flat Pack | 2.0kg | PMH005 | Pocket Mask Holder | 0.1kg |

Glove Box Holders

Features

- Plastic coated wire frame
- Available for wall mounting or fitted with a 30mm square hook

Standard Colours

Paint - White





| Item | Description | Overall Dimensions | |
|----------|------------------|---------------------------|--|
| One Box | | | |
| DG/1 | Wall Mounting | 290 x 115 x 145mm | |
| DG/1/30 | 30mm Square Hook | 290 x 155 x 145mm | |
| Three Bo | x | | |
| DG/3 | Wall Mounting | 290 x 115 x 405mm | |
| DG/3/30 | 30mm Square Hook | 290 x 155 x 405mm | |

DAF/1

Glove / Mask Box Holder One Box, Wall Mounted

Features

- Accommodates glove and mask boxes up to 140mm in width
- One piece construction
- Holes for wall mounting (fixings not provided)

Dimensions (w x d x h)

• 143 x 100 x 152mm

Standard Colours

Paint - White

Glove Dispenser Three Column, Adjustable

Features

- Packs of sterile gloves separated by two adjustable dividers
- Full width card holder for identification
- Clear front panel assists with stock level inspection
- Hinged lid assists replenishment, whilst top loading ensures stock rotation
- Holes for wall mounting (fixings not provided)

Dimensions (w x d x h)

• 565 x 260 x 315mm

Standard Colours

· Paint - White



Alcohol Gel Holders - Wire Frame

Features

- Compatible with popular 500ml bottles, 70mm x 70mm x 225mm (not supplied)
- Stainless steel construction
- Plastic drip tray
- Available with holes for wall mounting (fixings not provided) or with a 30mm square hook

Standard Colours

Stainless Steel



| Item | Description | Overall Dimensions (w x d x h) | Weight |
|-------|--------------------------------------|--------------------------------|--------|
| GH005 | Lever Operated, Wall Mounted, Angled | 77 x 300 x 278mm | 0.4kg |
| GH015 | 30mm Hook | 77 x 117 x 180mm | 0.3kg |

Syringe Dispenser - Five Column, Adjustable

Features

- Syringes separated by four adjustable dividers
- Full width card holder for identification
- Clear front panel assists with stock level inspection
- Hinged lid assists replenishment, whilst top loading ensures stock rotation
- Holes for wall mounting (fixings not provided)

Dimensions (w x d x h)

• 565 x 317 x 320mm

Standard Colours

· Paint - White

www.bristolmaid.com

+44 (0)1258 484455

sales@bristolmaid.com



| Item | Description | Weight |
|-------|-------------------------|--------|
| SS/36 | Five Column, Adjustable | 9.2kg |

Sharps Box Holder One Box, Wall Mounted

Features

- · Accepts sharps boxes up to seven litres capacity
- Sharps box secured using a rubber band fitted with a toggle
- · Holes for wall mounting (fixings not provided

Dimensions (w x d x h)

• 127 x 100 x 190mm

Standard Colours

Paint - White



| Item | Description | Weight |
|--------|-----------------------|--------|
| SBH/WM | One Box, Wall Mounted | 0.4kg |

Urine Bottle Holder Hanging

Features

- Compatible with bottles up to 115mm Diameter
- Plastic coated wire frame
- 30mm square hook

Dimensions (w x d x h)

• 150 x 182 x 282mm

Standard Colours

· Plastic - White



Catheter Bag Holders Free Standing/Hanging

Features

- Plastic coated wire frame
- Accepts most widely used bags
- Available freestanding or hanging (30mm square hook)
- Upper clip to hold the inlet tube
- Lower clip angles the drainage tube upwards to prevent drips should the tap leak

Standard Colours

Plastic - White



| Item | Description | Weight |
|------|--------------|--------|
| UB/9 | Urine Bottle | 0.3kg |

| ltem | Description | Overall Dimensions (w x d x h) | Weight |
|-------|--------------|--------------------------------|--------|
| UB005 | Freestanding | 250 x 250 x 550mm | 0.6kg |
| UB010 | Hanging | 250 x 77 x 410mm | 0.4kg |

Catheter Racks - Two & Five Section

Features

- Internal dimensions of each section 65 x 75 x 405mm
- Removable clear perspex panel to aid cleaning
- Available with holes for wall mounting (fixings not provided) or with a 15mm square hook

Standard Colours



| Item Description Overall Dimensions (w x d x h) | | Weight | |
|---|--------------|-------------------|-------|
| Two Section | | | |
| RC/2 | Wall Fitting | 135 x 85 x 410mm | 1.3kg |
| RC/2/15 | 15mm Hook | 135 x 115 x 410mm | 1.4kg |
| Five Section | | | |
| RC/5 | Wall Fitting | 320 x 85 x 410mm | 3.3kg |
| RC/5/15 | 15mm Hook | 320 x 115 x 410mm | 3.4kg |

Catheter Racks - One Section

Features

- Outer sleeve removable for cleaning
- 20mm square hook

Dimensions (w x d x h) Standard Colours

• 185 x 55 x 555mm Paint - White



Drainage Racks - Urinal

Features

- Holds five urine bottles
- Stainless steel construction
- Supplied with a stainless steel drip tray
- Holes for wall mounting (fixings not provided)

Dimensions (w x d x h)

• 535 x 205 x 415mm



| ltem | Description | Weight |
|----------|-------------|--------|
| SS/40/MS | One Section | 2.0kg |

| Item | Description |
|----------|-------------------------|
| RD/U5/SS | Urinal, Stainless Steel |

Drainage Racks - Bedpan

Features

- Holds five bedpans
- Stainless steel construction
- Supplied with a stainless steel drip tray
- Holes for wall mounting (fixings not provided)

Dimensions (w x d x h)

• 530 x 295 x 375mm



| Weight | Item | Description | Weight |
|--------|----------|-------------------------|--------|
| 1.7kg | RD/B5/SS | Bedpan, Stainless Steel | 3.3kg |



Dressing Instrument Trolleys



- Removable Shelves, Stainless Steel 102
 - Removable Shelves, Mild Steel 103
 - Fixed Shelves, Sides Up 104
 - Fixed Shelves, Sides Down 105
 - Drawer Units 107
 - Options 108



Removable Shelves, Stainless Steel, 450mm Deep

Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Flat or reversible 304 grade stainless steel shelves, cushioned on anti-static spiral buffers (removable for cleaning)
- Reversible shelves allow the end user to choose between a flat shelf or a shelf with three upturned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x nonanti-static)

Options

- C001 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- C013 75mm, Pressed Steel Castors (4 x braking)
- Drawer Units Range of sizes / configurations (see page 107)
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers

Removable Shelves, Mild Steel, 450mm Deep

Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded mild steel frame
- Flat or reversible 304 grade stainless steel shelves, cushioned on anti-static spiral buffers (removable for cleaning)
- Reversible shelves allow the end user to choose between a flat shelf or a shelf with three upturned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x non-anti-static)

Standard Colours

Paint - White

Options

- C001 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- C013 75mm, Pressed Steel Castors (4 x braking)
- Drawer Units Range of sizes / configurations (see page 107)
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers

| Flat Shelves | Reversible Shelves | Size | Shelf Size (w x d) | Overall Dimensions (w x d x h) | Weight |
|--------------|---------------------------|---------|--------------------|--------------------------------|--------|
| DTSR/450 | DTSR/450/SR | Small | 450 x 450mm | 470 x 470 x 890mm | 11.5kg |
| DTSR/600 | DTSR/600/SR | Medium | 600 x 450mm | 620 x 470 x 890mm | 16.0kg |
| DTSR/750 | DTSR/750/SR | Large | 750 x 450mm | 770 x 470 x 890mm | 18.0kg |
| DTSR/900 | DTSR/900/SR | X-Large | 900 x 450mm | 920 x 470 x 890mm | 20.0kg |
| | | | | | |

| Flat Shelves | Reversible Shelves | Size | Shelf Size (w x d) | Overall Dimensions (w x d x h) | Weight |
|--------------|---------------------------|--------|--------------------|--------------------------------|--------|
| DTMR/450 | DTMR/450/SR | Small | 450 x 450mm | 470 x 470 x 890mm | 11.5kg |
| DTMR/600 | DTMR/600/SR | Medium | 600 x 450mm | 620 x 470 x 890mm | 16.0kg |





Fixed Shelves, Sides Up, Stainless Steel, 450mm Deep

Features

 Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed shelves with three sides up & one down
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x nonanti-static)

Options

- C001 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- C013 75mm, Pressed Steel Castors (4 x braking)
- Drawer Units Range of sizes / configurations (see page 107)
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers

Fixed Shelves, Sides Down, Stainless Steel, 450mm Deep

Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings
 & for anaesthetists' use
- · Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed shelves with down turned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x non-anti-static)

Options

- C001 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- C013 75mm, Pressed Steel Castors (4 x braking)
- Drawer Units Range of sizes / configurations (see page 107)
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers

| Item | Size | Shelf Size (w x d) | Overall Dimensions (w x d x h) | Weight |
|----------|---------|--------------------|--------------------------------|--------|
| DTSF/450 | Small | 450 x 450mm | 470 x 470 x 890mm | 11.5kg |
| DTSF/600 | Medium | 600 x 450mm | 620 x 470 x 890mm | 16.0kg |
| DTSF/750 | Large | 750 x 450mm | 770 x 470 x 890mm | 18.0kg |
| DTSF/900 | X-Large | 900 x 450mm | 920 x 470 x 890mm | 20.0kg |

| Size | Shelf Size (w x d) | Overa |
|---------|--------------------|--|
| | | |
| Small | 450 x 450mm | 470 x 4 |
| Medium | 600 x 450mm | 620 x 4 |
| Large | 750 x 450mm | 770 x 4 |
| X-Large | 900 x 450mm | 920 x 4 |
| | Large | Small 450 x 450mm Medium 600 x 450mm |

 $\overline{\mathbf{O}}$



Fixed Shelves, Sides Down, Stainless Steel, 600mm Deep

• Fixed shelves with down turned sides

Maximum capacity per shelf 40kg

DP/C - Infusion Pole & Bracket

ZP005 - Ring Buffers

Meets BS4068:1977 Specification for hospital trolleys for

Seamlessly welded from 304 grade stainless steel for strength

C002 - 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
 C003 - 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking)

instruments, dressings & for anaesthetists' use

100mm pressed steel anti-static non-marking castors

C015 - 100mm, Pressed Steel Castors (4 x braking)

Features

Options

& ease of cleaning

Shelf heights 405 & 870mm



| tem | Size | Shelf Size (w x d) | Overall Dimensions (w x d x h) | Weight |
|---------------|--------|--------------------|--------------------------------|--------|
| DTSF/66/2/SD | Small | 600 x 600mm | 620 x 620 x 870mm | 17.5kg |
| DTSF/96/2/SD | Medium | 900 x 600mm | 920 x 620 x 870mm | 24.0kg |
| DTSF/126/2/SD | Large | 1200 x 600mm | 1220 x 620 x 870mm | 29.6kg |

Drawer Units

Features

- Compatible with removable & fixed shelf ranges, with exception of 600mm deep trolleys
- Available in three widths...
 DRW260 Set of two drawers
- DRW260 Set of two draw
 DRW350 Single drawer
- DRW350 Single drawer
 DRW650 Single drawer
- DRVV650 Single drawe
 Maximum capacity 10.0kg
- Factory fitted & must be ordered at the same time as the trolley

Removable Shelf Trolleys

| | One Tier Combination Code | | | Two Tier Combination Code | | |
|------------------|---------------------------|--------------|--------------|---------------------------|--------------------|--------------------|
| | DRW260/R | DRW350/R | DRW650/R | DRW260/R DRW260 | DRW350/R DRW350 | DRW650/R DRW650 |
| Trollow | | | | | | |
| Trolley Width | 1 1 | | 1 1 | 1 1 | 1 1 | |
| 450mm | | \checkmark | | | \checkmark | |
| 600mm | \checkmark | \checkmark | | \checkmark | \checkmark | |
| 750mm | \checkmark | \checkmark | \checkmark | \checkmark | \checkmark | \checkmark |
| 900mm | \checkmark | \checkmark | \checkmark | \checkmark | \checkmark | \checkmark |

Fixed Shelf Trolleys

| | DRW260/F | DRV |
|-----------------|--------------|-----|
| rolley Width | | |
| 450mm | | |
| 600mm | \checkmark | |
| 750mm | \checkmark | |
| 900mm | \checkmark | |

ō

Internal Dimensions (w x d x h)

DRW260 - 230 x 410 x 100mm
DRW350 - 350 x 410 x 100mm
DRW650 - 650 x 410 x 100mm



DTSR/750 with DRW260/R/DRW260



Castors

Features

- Polymer & pressed steel castors
- Anti-static & non-marking
- Supplied as a set of 4 castors
- 80mm polymer castors increase the height of the trolley by approx. 5mm
- 100mm polymer castors increase the height of the trolley by approx. 30mm
- 100mm steel castors increase the height of the trolley by approx. 30mm
- Factory fitted



Item Description

- C001 Castor 80mm, Grey Polymer, 2 x Braking / 2 x Non-Braking (set of 4) C002 Castor - 100mm, Grey Polymer, 2 x Braking / 2 x Non-Braking (set of 4)
- C003 Castor 100mm, Pressed Steel, 2 x Braking / 2 x Non-Braking (set of 4)
- C013 Castor 75mm, Pressed Steel (4 x braking)
- C015 Castor 100mm, Pressed Steel (4 x braking)

Reversible Shelves

Features

- Allows the end user to choose between a flat shelf or one with three upturned sides
- Removable for cleaning
- Compatible with DTMR & DTSR ranges only



| Item | Description | Weight |
|------|--------------------|--------|
| SR | Reversible Shelves | 5-10kg |

Clamping Infusion Stands - Trolley

Features

- Clamp secures to 25mm vertical round tube (designed to fit our range of dressing trolleys)
- When fitting to other trolleys, we strongly recommend that a risk assessment is carried out to assess stability
- When fitted protrudes 45mm
- Supplied with a two pigtail hook infusion pole
- Maximum capacity on hooks 2kg

Dimensions (w x d x h)

• 55 x 220 x 1200mm





| Item | Description |
|------|------------------------------------|
| DP/C | Clamping Infusion Stands - Trolley |

Dre

Ring Buffers

Features

- Set of four ring buffers
- Offers protection to both the trolley & fabric of the building
- Increases the length & width of trolley by approx. 80mm
- Factory fitted





Weight





Drug & Medicine Dispensing

- Dispensing Trolleys 112
- CLIQ[®] & RFID Solutions 115
- Patient Administration Trolleys 116
 - Drawer Trolleys 118
 - Tray Trolleys 119
 - Blister Pack (MDS) Trolleys 120
 - Nomad[®] Cassette Trolleys 122
- Biodose[®] & MultiMeds Trolleys 124
 - Medinoxx Trolleys 125
 - Laptop Holders 126
 - Security Clamp 128
 - Distribution Trolleys 129



Dispensing Trolleys - Wooden, Standard Capacity

Features

- Trollevs suitable for the storage & dispensing of drugs & medicines
- Lockable hinged lids
 - Upper Assisted using a gas strut, providing document support
- Lower Provides a work area, at height 980mm (max capacity 4.0kg)
- Available with either an electronic push button or high security bolt lock
- Electronic Push Button Uses a four digit code (available with master/ submaster/user & technicians codes)
- High Security Bolt Lock Supplied with two keys
- Hinged writing flap provides a work area at height 930mm (max capacity 2.0kg)
- Lower shelf
- Five tier plastic rack (395mm)
- Two slide out storage / waste bins
- 100mm swivel pressed steel castors
- When not in use, trolley should be secured to the building - see optional MT/1/6F

Dimensions (w x d x h)

- Overall...
- Small 750 x 490 x 1130mm
- Medium 890 x 490 x 1130mm
- Large 1030 x 525 x 1130mm
- Internal...
 - Small 530 x 427 x 263mm sloping to 112mm at front
- Medium 683 x 427 x 263mm sloping to 112mm at front
- Large 805 x 461 x 263mm sloping to 112mm at front

Standard Colours

Options

- Paint Grev White
 - MT/1/6F Security Wall Clamp
- Cabinet/Shelf/Flap Spearmint/White ZP005 - Ring Buffers

| ltem | Description | Weight | | | | |
|-----------------------------|-------------------------------------|--------|--|--|--|--|
| High Security I | High Security Bolt Lock | | | | | |
| MT100 | Small, High Security Bolt Lock | 26.0kg | | | | |
| MT110 | Medium, High Security Bolt Lock | 35.0kg | | | | |
| MT120 | Large, High Security Bolt Lock | 42.5kg | | | | |
| Electronic Push Button Lock | | | | | | |
| MT100CL | Small, Electronic Push Button Lock | 26.0kg | | | | |
| MT110CL | Medium, Electronic Push Button Lock | 35.0kg | | | | |
| MT120CL | Large, Electronic Push Button Lock | 42.5kg | | | | |

Dispensing Trolleys - Wooden, High Capacity

Features

- Trollevs suitable for the storage & dispensing of drugs & medicines Lockable hinged lids
- Upper Assisted using a gas strut, providing document support
- Lower Provides a work area, at height 921mm (max capacity 4.0kg)
- Available with either an electronic push button or high security bolt lock Electronic Push Button - Uses a four digit code (available with master/ submaster/user & technicians codes)
- High Security Bolt Lock Supplied with two keys
- Hinged writing flap provides a work area at height 850mm (max capacity 2.0kg)
- Lower shelf
- Five tier plastic rack (395mm)
- Three angled shelves, complete with six adjustable U shaped dividers
- Two slide out storage / waste bins
- Ring buffers
- 125mm swivel pressed steel castors
- When not in use, trolley should be secured to the building see optional MT/1/6F

Dimensions (w x d x h)

- Overall...
- Medium 760 x 560 x 1290mm
- Large 905 x 560 x 1290mm
- Internal...
- Medium 530 x 461 x 513mm sloping to 137mm at the front
- Large 683 x 461 x 513mm
- sloping to 137mm at the front

Standard Colours

- Paint Grey White
- · Cabinet/Shelf/Flap Spearmint/White

Options

MT/1/6F - Security Wall Clamp

| Item | Description |
|------------|-------------------------------------|
| High Secur | ity Bolt Lock |
| MT220 | Medium, High Security Bolt Lock |
| MT230 | Large, High Security Bolt Lock |
| Electronic | Push Button Lock |
| MT220CL | Medium, Electronic Push Button Lock |
| MT230CL | Large, Electronic Push Button Lock |





Dispensing Trolleys - Metal, Standard Capacity

Features

- Trolleys suitable for the storage & dispensing of drugs & medicines
- Lockable hinged lid provides document support
- Fitted with a high security bolt lock (supplied with two keys)
- Hinged writing flap provides a work area at height 907mm (max capacity 2.0kg)
- Lower shelf
- Five tier plastic rack (395mm)
- Two slide out storage / waste bins
- 100mm swivel pressed steel castors
- When not in use, trolley should be secured to the building see optional MT/1/6F

Dimensions (w x d x h)

- Overall...
- Small 750 x 459 x 1112mm
- Medium 891 x 459 x 1112mm
- Large 1030 x 492 x 1112mm
- Internal.
- Small 489 x 457 x 292mm sloping to 127mm at the front
- Medium 706 x 457 x 292mm sloping to 127mm at the front
- Large 832 x 490 x 292mm sloping to 127mm at the front

Standard Colours

- Paint Grey White
- Shelf/Flap Spearmint
- Options
- MT/1/6F Security Wall Clamp • ZP005 - Ring Buffers



| ltem | Description | Weight |
|-------|---------------------------------|--------|
| MT130 | Small, High Security Bolt Lock | 22.0kg |
| MT140 | Medium, High Security Bolt Lock | 25.0kg |
| MT150 | Large, High Security Bolt Lock | 31.0kg |

CLIQ[®] & RFID Solutions

The use of CLIQ[®] coupled with unique RFID tag identifiers allows a single key to be used with Medicines Management and for Accessing Patient Records.

- Controlled Drug Cabinets (audit trail available)
- Drug & Medicine Cabinets (audit trail available)
- Self-Administration Cabinets
- Bedside Cabinets
- Drug & Medicine Dispensing Trolleys
- Medical Records Trollevs


The use of RFID cards allows a single card to be used with Self-Administration Schemes.

- Self-Administration Cabinets
- Bedside Cabinets





ď



D



Pharmacy Trolleys - Patient Administration

Features

- Trolleys suitable for the storage & dispensing of drugs & medicines using individual compartments dedicated for each patient
- Available with either an electronic push button or high security bolt lock
- Electronic Push Button Uses a four digit code (available with master/ submaster/user & technicians codes)
- High Security Bolt Lock Supplied with two keys
- Trolleys supplied with 4 / 8 trays subdivided into equal compartments
- 24 Compartment Trolley supplied with 4 trays subdivided into 6 equal compartments, each compartment - 190mm x 132mm x 220mm (5.5 Litres).
 Label holders provided for each compartment
- 32 Compartment Trolley supplied with 8 trays subdivided into 4 equal compartments, each compartment - 192mm x 143mm x 220mm (6.1 Litres). Label holders provided for each compartment
- Upper work surface
- Push handles incorporate integral buffers to protect the fabric of the building
- Plastic debris bin & bracket
- 125mm polymer castors offer reduced rolling resistance making the trolley easier to manoeuvre, lowers noise levels, whilst smooth surfaces assist cleaning
- When not in use, trolley should be secured to the building see optional MT/1/6F

Dimensions (w x d x h)

Standard Colours

- Single Door 627 x 579 x 1165mm
 Double Door 820 x 579 x 1165mm
- Paint Grey
 Work Surface Spearmint

Options

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm
- PT/LH/C5 Motion C5 Holder



PTS/EPB/PA24

Single Door

| High Security Bolt Lock | Electronic Push Button Lock | Description | Weight |
|-------------------------|-----------------------------|---|--------|
| PTS/HSB/PA24 | PTS/EPB/PA24 | Patient Administration, 24 Compartments | 51.4kg |



Double Door

| High Security Bolt Lock | Electronic Push Butto |
|-------------------------|-----------------------|
| PTD/HSB/PA32 | PTD/EPB/PA32 |

PTD/HSB/PA32



1000





PTD/HSB/PA32



Pharmacy Trolleys - Drawers

Features

- Trollevs suitable for the storage & dispensing of drugs & medicines using drawers
- Available with either an electronic push button or high security bolt lock • Electronic Push Button - Uses a four digit code (available with master/ submaster/user & technicians codes)
- High Security Bolt Lock Supplied with two keys
- Internal drawer dimensions 335 x 295 x 113mm (w x d x h)
- Upper work surface
- Push handles incorporate integral buffers to protect the fabric of the building
- Adjustable door shelves
- Single Door 3 door shelves
- Plastic debris bin & bracket
- 125mm polymer castors offer reduced rolling resistance making the trolley easier to manoeuvre, lowers noise levels, whilst smooth surfaces assist cleaning
- When not in use, trolley should be secured to the building see optional MT/1/6F

Standard Colours Dimensions (w x d x h)

- 511 x 542 x 1072mm
- Paint Grey Work Surface - Spearmint

Options

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm
- PT/LH/C5 Motion C5 Holder Swing Arm



PTS/EPB/DR6

| Item | Description | Weight |
|-------------|--|--------|
| PTS/HSB/DR6 | High Security Bolt Lock, Six Drawers | 53.0kg |
| PTS/EPB/DR6 | Electronic Push Button Lock, Six Drawers | 53.0kg |

Pharmacy Trolleys - Trays

Features

- Trollevs suitable for the storage & dispensing of drugs & medicines using individual patient travs supported on adjustable shelves
- Available with either an electronic push button or high security bolt lock Electronic Push Button - Uses a four digit code (available with master/submaster/user & technicians codes)
- High Security Bolt Lock Supplied with two keys
- Travs supported on adjustable shelves
- Tray Sizes (w x d x h)
- LA 86 x 268 x 57mm (aluminium)
- A 85 x 181 x 54mm (plastic)
- C 129 x 181 x 54mm (plastic)
- LP 115 x 266 x 95mm (plastic)
- Upper work surface
- Push handles incorporate integral buffers to protect the fabric of the buildina
- Adjustable door shelves
- Single Door 3 door shelves
- Double Door 6 door shelves
- Plastic debris bin & bracket
- 125mm polymer castors offer reduced rolling resistance making the trolley easier to manoeuvre, lowers noise levels, whilst smooth surfaces assist cleaning
- When not in use, trolley should be secured to the building - see optional MT/1/6F

Dimensions (w x d x h)

- Single Door 511 x 542 x 1072mm
- Double Door 711 x 568 x 1072mm

Options

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm
- PT/LH/C5 Motion C5 Holder Swing Arm

| High Security Bolt Lock | Electronic Push Button |
|-------------------------|-------------------------------|
| Single Door | |
| PTS/HSB/LA16 | PTS/EPB/LA16 |
| Double Door | |
| PTD/HSB/A18C12 | PTD/EPB/A18C12 |
| PTD/HSB/LP24 | PTD/EPB/LP24 |
| | |





Pharmacy Trolleys - Blister Pack (MDS)

Features

- Trolleys suitable for the storage & dispensing of drugs & medicines using the blister packed monitored dosage system (MDS)
- Available with either an electronic push button or high security bolt lock Electronic Push Button - Uses a four digit code (available with master/submaster/user & technicians codes)
- · High Security Bolt Lock Supplied with two keys
- Frames supported on either a fixed panel mounted to the inside rear of the trolley or on three adjustable slide out shelves (frames not supplied)
- Upper work surface
- Push handles incorporate integral buffers to protect the fabric of the buildina
- Adjustable door shelves
- Single Door 3 door shelves
- Double Door 6 door shelves
- Fold over MDS dispensing surface
- Plastic debris bin & bracket
- Lower drawer (refer to table for models)
- Single Door 335 x 295 x 66mm
- Double Door 535 x 324 x 66mm
- 125mm polymer castors offer reduced rolling resistance making the trolley easier to manoeuvre, lowers noise levels, whilst smooth surfaces assist cleaning
- When not in use, trolley should be secured to the building - see optional MT/1/6F

Dimensions (w x d x h)

- Single Door 511 x 542 x 1072mm
- Paint Grey Double Door - 711 x 568 x 1072mm
 Work Surface - Spearmint

Standard Colours

Options

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm
- PT/LH/C5 Motion C5 Holder Swing Arm

Single Door

PTS/HSB/BP6

(frames not included)

PTS/EPB/BP4

| High Security Bolt Lock | Electronic Push Button Lock | Description | Weight |
|-------------------------|-----------------------------|----------------------------------|--------|
| PTS/HSB/BP4 | PTS/EPB/BP4 | Blister Pack, Four Frames | 31.0kg |
| PTS/HSB/BP6 | PTS/EPB/BP6 | Blister Pack, Six Frames | 38.0kg |
| PTS/HSB/DR/BP6 | PTS/EPB/DR/BP6 | Drawer, Blister Pack, Six Frames | 39.0kg |



PTS/HSB/DR/BP6



Double Door

| High Security Bolt Lock | Electronic Push Butt |
|-------------------------|----------------------|
| PTD/HSB/BP6 | PTD/EPB/BP6 |
| PTD/HSB/BP9 | PTD/EPB/BP9 |
| PTD/HSB/DR/BP9 | PTD/EPB/DR/BP9 |
| | |

120

PTD/HSB/DR/BP9

| _ock | Description | Weight |
|------|-----------------------------------|--------|
| | Blister Pack, Six Frames | 38.0kg |
| | Blister Pack, Nine Frames | 42.0kg |
| | Drawer, Blister Pack, Nine Frames | 46.0kg |

Tab prevents the shelf from being removed



Pharmacy Trolleys - Nomad[®] Cassettes

Features

- Trolleys suitable for the storage & dispensing of drugs & medicines using the Nomad[®] cassette monitored dosage system
- Available with either an electronic push button or high security bolt lock
- Electronic Push Button Uses a four digit code (available with master/ submaster/user & technicians codes)
- High Security Bolt Lock Supplied with two keys
- Cassettes supported on three adjustable shelves
- Upper work surface
- Push handles incorporate integral buffers to protect the fabric of the building
- Adjustable door shelves
- Single Door 3 door shelves
- Double Door 6 door shelves
- Plastic debris bin & bracket
- 125mm polymer castors offer reduced rolling resistance making the trolley easier to manoeuvre, lowers noise levels, whilst smooth surfaces assist cleaning
- When not in use, trolley should be secured to the building see optional MT/1/6F

Dimensions (w x d x h)

Standard Colours

Nomad[®], Thirty Cassettes

Weight

38.0kg

- Single Door 511 x 542 x 1072mm • Double Door - 711 x 568 x 1072mm
 - Paint Grey Work Surface - Spearmint

- Options
- MT/1/6F Security Wall Clamp

PTS/EPB/NC30

- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm
- PT/LH/C5 Motion C5 Holder Swing Arm



PTS/EPB/NC30

Double Door

| High Security Bolt Lock | Electronic Push Button Lock | Description | Weight |
|-------------------------|-----------------------------|--|---------------|
| PTD/HSB/NC56 | PTD/EPB/NC56 | Nomad [®] , Fifty Six Cassettes | 48.0kg |
| | | | 0)4050 404455 |

+44 (0)1258 484455 www.bristolmaid.com sales@bristolmaid.com

PTS/HSB/NC30

ö

ž ంర

b





PTD/HSB/NC56 (cassettes not included)



sales@bristolmaid.com +44 (0)1258 484455 www.bristolmaid.com



Pharmacy Trolleys - Biodose[®] & MultiMeds

Features

- Trollevs suitable for the storage & dispensing of drugs & medicines using the Biodose[®] & MultiMeds monitored dosage system
- Available with either an electronic push button or high security bolt lock Electronic Push Button - Uses a four-digit code (available with master/submaster/user & technicians codes)
- High Security Bolt Lock Supplied with two keys
- Travs supported on fixed runners (travs not supplied)
- Upper work surface
- Push handles incorporate integral buffers to protect the fabric of the buildina
- Adjustable door shelves
- Single Door 3 door shelves
- Double Door 6 door shelves
- Plastic debris bin & bracket
- 125mm polymer castors offer reduced rolling resistance making the trollev easier to manoeuvre. lowers noise levels, whilst smooth surfaces assist cleaning
- When not in use, trolley should be secured to the building see optional MT/1/6F

Dimensions (w x d x h)

- Single Door 581 x 568 x 1072mm
- Double Door 812 x 568 x 1072mm

Standard Colours

- Paint Grev
- · Work Surface Spearmint

Options

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm
- PT/LH/C5 Motion C5 Holder Swing Arm

| | Pharmacy | Trolleys - | Medinoxy |
|--|----------|-------------------|----------|
|--|----------|-------------------|----------|

Features

- Trollevs suitable for the storage & dispensing of drugs & medicines using the Medinoxx monitored dosage system (compatible with 7 x 4 medi-travs)
- Available with either an electronic push button or high security bolt lock • Electronic Push Button - Uses a four-digit code (available with master/ submaster/user & technicians codes)
- High Security Bolt Lock Supplied with two keys
- Travs supported on fixed runners (travs not supplied)
- Upper work surface
- Push handles incorporate integral buffers to protect the fabric of the buildina
- Adjustable door shelves
- Single Door 3 door shelves
- Double Door 6 door shelves
- Plastic debris bin & bracket
- 125mm polymer castors offer reduced rolling resistance making the trolley easier to manoeuvre. lowers noise levels, whilst smooth surfaces assist cleaning
- When not in use, trolley should be secured to the building - see optional MT/1/6F

Dimensions (w x d x h)

- Single Door 581 x 568 x 1072mm
- Double Door 812 x 568 x 1072mm

Standard Colours

- Paint Grev
- Work Surface Spearmint

Options

- MT/1/6F Security Wall Clamp
- PT/LHR/S Laptop Holder Removable, Swing Arm
- PT/LHNR/S Laptop Holder Non-Removable, Swing Arm
- PT/LH/C5 Motion C5 Holder Swing Arm

| High Security Bolt Lock | Electronic Push Butto |
|-------------------------|-----------------------|
| Single Door | |
| PTS/HSB/MX40 | PTS/EPB/MX40 |
| Double Door | |
| PTD/HSB/MX60 | PTD/EPB/MX60 |


| High Security Bolt Lock | Electronic Push Button Lock | Description | Weight |
|-------------------------|-----------------------------|--|--------|
| Single Door | | | |
| PTS/HSB/BD46 | PTS/EPB/BD46 | Biodose [®] & MultiMeds, 46 Trays | 37.3kg |
| Double Door | | | |
| PTD/HSB/BD69 | PTD/EPB/BD69 | Biodose® & MultiMeds, 69 Trays | 52.0kg |



n Lock







Laptop Holder - Removable, Swing Arm

Features

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the PT range of single & double door pharmacy trolleys
- Laptop holder fitted to an upright bolted to the rear left hand corner of the trolley (upright & holder supplied for customer to fit)
- Adjustable to accept laptops 280-435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Laptop supported in a cradle fitted with adjustable wing arms
- Laptop can be secured using a security clamp (supplied)*
- Design allows easy access to laptop ports

*Not compatible with the Dell laptops using the N17 wedge lock system

Standard Colours

Paint - Grey

| Item | Description | Weight |
|----------|----------------------|--------|
| PT/LHR/S | Removable, Swing Arm | 3.7kg |

Laptop Holder - Non-Removable, Swing Arm

Features

- Three section swing arm offering height, rotation & angle adjustment
- · Compatible with the PT range of single & double door pharmacy trolleys
- Laptop holder fitted to an upright bolted to the rear left hand corner of the trolley (upright & holder supplied for customer to fit)
- Adjustable to accept laptops 280-435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Laptop supported in a cradle fitted with a steel mounting bracket that can be fixed using either nylock or wing nuts
- Laptop can be secured using a security clamp (supplied)*
- Design allows easy access to laptop ports

*Not compatible with the Dell laptops using the N17 wedge lock system

Standard Colours

Paint - Grey

Description Weight Item PT/LHNR/S Non-Removable, Swing Arm 5.0kg



ò





Motion C5 Holder - Swing Arm

Features

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the PT range of single & double door pharmacy trolleys
- Motion C5 holder fitted to an upright bolted to the rear left hand corner of the trolley (upright & holder supplied for customer to fit)
- Maximum capacity 2.0kg
- Tablet located in a RAM-In-Motion passive dock allowing for easy insertion & removal

Standard Colours

Paint - Grev



| Item | Description | Weight |
|----------|-------------|--------|
| PT/LH/C5 | Swing Arm | 2.7kg |

Fitting Your Security Clamp



Securely fasten one end of the clamp to the wall



The 'hook shaped' end is trapped by the door when closed & locked



Security Clamp

Features

- Allows the trolley to be secured to the 'fabric of the building'
- · Requires no additional keys, utilises the trolleys own locking system
- One end is secured to the wall, the 'hook shaped' end is trapped by the door/lid of the trolley when closed & locked
- Recommended wall fitting height approx 1000mm



Important Note - Security Clamps It is a requirement that when the trolley is not in use it should be secured to the 'fabric of the building'. Not only does this prevent casual theft, i.e. trolleys being pushed to an isolated area & broken into, but also encourages good nursing practice, giving the trolley a "permanent home".

"Safe & Secure Handling of Medicines - Royal Pharmaceutical Society, GB-March 2005"

| Item | Description | Cable Length | Weight |
|---------|----------------|--------------|--------|
| MT/1/6F | Security Clamp | 950mm | 0.75kg |

Distribution Trolleys - Lockable Drawers

Features

- Trolleys suitable for the distribution of drugs & medicines from pharmacy to the ward & also with self-administration/patient own drug schemes
- Provides secure drug delivery using drawers which when locked are securely fastened to the trolley
- Boxes are fitted with a CAM lock, supplied with two keys Supplied as standard with locks to differ, allowing each ward to have its
- own key & be responsible for removing the drawer from the trolley • The remaining key should be kept at the pharmacy
- Max capacity per drawer 3kg
- Master keys are available (additional cost)
- Upper work surface
- Push handles
- Document pocket fitted to one end of the trolley
- Low level corner bumpers to protect the fabric of the building
- 125mm swivel castors, 2 off braking

Dimensions (w x d x h)

- SA005 1026 x 523 x 1002mm
- SA025 1180 x 523 x 1002mm
- Drawer 145 x 423 x 230mm (Internal)

Standard Colours

- Paint Grey
- Work Surface Spearmint

Options

Master Keys

Ö

& Me

Drug







Drug & Medicine Storage

- Drug & Medicine Cabinets 132
- Controlled Drug Cabinets 137
 - CLIQ[®] 140
- Drug & Medicine Cupboards 142
 - Options 144
- Patients/Residents Own Medication Cabinets 146
 - Tablet Counter & Storage Boxes 150
 - Pharmacy & Vaccine Refrigerators 151



Single Door Cabinets

Drug & Medicine Cabinets - Shelves & Door Shelves

Features

- Independently tested to BS 2881: 1989 (Security Level One) Specification for Cupboards for the Storage of Medicines in Health Care Premises
- BS 3621 anti-pick & drill, key retaining five pin cylinder lock (supplied with 3 keys)
- Adjustable & removable shelves (190mm deep) & door shelves (75mm deep)
- Key profile unique to Hospital Metalcraft Ltd, offering 250,000 differs
- Compatible with CLIQ[®], integrating electronics into high security mechanical disc cylinders & keys, combining electronic security & remote management in real time
- Optional Digital Push Button Lock (does not conform to BS 2881: 1989)
- 12v DC LED lighting/warning system (model dependent) with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug. CE marked & conforms to the low voltage directives 73/23/EEC & 89/336/EC
- Lighting transformer housing adds 105mm to the overall height
- 0.9mm (20 gauge) mild steel welded construction
- Mounting holes provided in back (fixings not included)
- Modular design allows for easy interchange with HTM 63 storage systems
- Compatible with our range of CD cabinets for the storage of Controlled Drugs (see page 138)

Standard Colours

Paint - White

Options

- DGL Digital Push Button Lock
- SL Suited Locks (contact us for more information)
- PC/500/ST Sloping Top (500mm Wide Cabinets)
- PC/600/ST Sloping Top (600mm Wide Cabinets)
- PC/800/ST Sloping Top (800mm Wide Cabinets)
- PC/1000/ST Sloping Top (1000mm Wide Cabinets)
- PC/RLBS Remote Light Box (Surface Mounted)
- PC/RLBF Remote Light Plate (Flush Mounted)
- PCB Buzzer Alarm (Factory Fit To Cabinets With Lights Only)



Double Door Cabinets

| Single Lock (Se | ecures Both Doors) | Each Door F | itted With A Lock | | ***CD015 not | compatible with PC839/IND & PC | 839/IND/LX |
|-----------------|--------------------|--------------|-------------------|----------------------------|--------------------------------|--------------------------------|------------|
| With Light* | Without Light | With Light* | Without Light | Shelves/Door Shelves | Overall Dimensions (w x d x h) | Compatible CD Cabinets | Weight |
| PC836 | PC836/LX | PC836/IND** | PC836/IND/LX** | 6 Shelves / 5 Door Shelves | 800 x 300 x 600mm | CD005, CD010** | 40.9kg |
| PC839 | PC839/LX | PC839/IND*** | PC839/IND/LX*** | 8 Shelves / 8 Door Shelves | 800 x 300 x 900mm | CD005, CD010, CD015*** | 45.9kg |
| PC1036 | PC1036/LX | PC1036/IND | PC1036/IND/LX | 6 Shelves / 6 Door Shelves | 1000 x 300 x 600mm | CD005, CD010 | 48.9kg |
| PC1039 | PC1039/LX | PC1039/IND | PC1039/IND/LX | 8 Shelves / 8 Door Shelves | 1000 x 300 x 900mm | CD005, CD010, CD015 | 53.9kg |

| Right Hand Opening | | Left Hand Opening | | * Transformer casing adds 105mm to overall | | | erall height |
|--------------------|---------------|-------------------|---------------|--|--------------------------------|------------------------|--------------|
| With Light* | Without Light | With Light* | Without Light | Shelves/Door Shelves | Overall Dimensions (w x d x h) | Compatible CD Cabinets | |
| PC536 | PC536/LX | PC536/LH | PC536/LH/LX | 3 Shelves / 3 Door Shelves | 500 x 300 x 600mm | CD005, CD010 | 29.9kg |
| PC539 | PC539/LX | PC539/LH | PC539/LH/LX | 4 Shelves / 4 Door Shelves | 500 x 300 x 900mm | CD005, CD010, CD015 | 34.0kg |
| PC636 | PC636/LX | PC636/LH | PC636/LH/LX | 3 Shelves / 3 Door Shelves | 600 x 300 x 600mm | CD005, CD010 | 33.9kg |
| PC639 | PC639/LX | PC639/LH | PC639/LH/LX | 4 Shelves / 4 Door Shelves | 600 x 300 x 900mm | CD005, CD010, CD015 | 38.0kg |

PC836/LX

*Transformer casing adds 105mm to overall height **CD010 not compatible with PC836/IND & PC836/IND/LX ***CD015 not compatible with PC839/IND & PC839/IND/LX



Diskt Hand Onemian Laft Hand Onemian

Single Door Cabinets

Drug & Medicine Cabinets - Deep Shelves

Features

- Independently tested to BS 2881: 1989 (Security Level One) Specification for Cupboards for the Storage of Medicines in Health Care Premises
- BS 3621 anti-pick & drill, key retaining five pin cylinder lock (supplied with 3 keys)
- Adjustable & removable deep shelves (250mm deep)
- Key profile unique to Hospital Metalcraft Ltd, offering 250,000 differs
- Compatible with CLIQ[®], integrating electronics into high security mechanical disc cylinders & keys, combining electronic security & remote management in real time
- Optional Digital Push Button Lock (does not conform to BS 2881: 1989)
- 12v DC LED lighting/warning system (model dependent) with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug. CE marked & conforms to the low voltage directives 73/23/EEC & 89/336/EC
- Lighting transformer housing adds 105mm to the overall height
- 0.9mm (20 gauge) mild steel welded construction
- Mounting holes provided in back (fixings not included)
- Modular design allows for easy interchange with HTM 63 storage systems
- Compatible with our range of CD cabinets for the storage of Controlled Drugs (see page 138)

Standard Colours

Paint - White

Options

- DGL Digital Push Button Lock
- SL Suited Locks (contact us for more information)
- PC/300/ST Sloping Top (300mm Wide Cabinets)
- PC/500/ST Sloping Top (500mm Wide Cabinets)
- PC/600/ST Sloping Top (600mm Wide Cabinets)
- PC/800/ST Sloping Top (800mm Wide Cabinets)
- PC/1000/ST Sloping Top (1000mm Wide Cabinets)
- PC/RLBS Remote Light Box (Surface Mounted)
- PC/RLBF Remote Light Plate (Flush Mounted)
- PCB Buzzer Alarm (Factory Fit To Cabinets With Lights Only)

* Transformer casing adds 105mm to overall height



Double Door Cabinets

| Single Loci | k (Secures Doors) | Each Door Fi | itted With A Lock | *** CD015 not compatible with PC839/FS/IND & PC839/FS | | | /FS/IND/LX | |
|-------------|-------------------|-----------------|--------------------|---|-------------------------|--------------------------------|------------------------|--------|
| With Light* | Without Light | With Light* | Without Light | | Nomad Cassette Capacity | Overall Dimensions (w x d x h) | Compatible CD Cabinets | Weight |
| PC836/FS | PC836/FS/LX | PC836/FS/IND** | PC836/FS/IND/LX** | 6 | - | 800 x 300 x 600mm | CD005, CD010** | 40.9kg |
| PC839/FS | PC839/FS/LX | PC839/FS/IND*** | PC839/FS/IND/LX*** | 8 | - | 800 x 300 x 900mm | CD005, CD010, CD015*** | 45.9kg |
| PC1036/FS | PC1036/FS/LX | PC1036/FS/IND | PC1036/FS/IND/LX | 6 | 58 Cassettes | 1000 x 300 x 600mm | CD005, CD010 | 48.9kg |
| PC1039/FS | PC1039/FS/LX | PC1039/FS/IND | PC1039/FS/IND/LX | 8 | 80 Cassettes | 1000 x 300 x 900mm | CD005, CD010, CD015 | 53.9kg |

| Right Ha | na Opening | Lett Hai | nd Opening | | | | 5 | |
|-------------|---------------|-------------|----------------|---|-------------------------|--------------------------------|------------------------|--------|
| With Light* | Without Light | With Light* | Without Light | | Nomad Cassette Capacity | Overall Dimensions (w x d x h) | Compatible CD Cabinets | Weight |
| PC336/FS | PC336/FS/LX | PC336/FS/LH | PC336/FS/LH/LX | 3 | 13 Cassettes | 300 x 300 x 600mm | CD005 | 16.9kg |
| PC536/FS | PC536/FS/LX | PC536/FS/LH | PC536/FS/LH/LX | 3 | 24 Cassettes | 500 x 300 x 600mm | CD005, CD010 | 29.9kg |
| PC636/FS | PC636/FS/LX | PC636/FS/LH | PC636/FS/LH/LX | 3 | 32 Cassettes | 600 x 300 x 600mm | CD005, CD010 | 33.9kg |
| PC539/FS | PC539/FS/LX | PC539/FS/LH | PC539/FS/LH/LX | 4 | 35 Cassettes | 500 x 300 x 900mm | CD005, CD010, CD015 | 34.9kg |
| PC639/FS | PC639/FS/LX | PC639/FS/LH | PC639/FS/LH/LX | 4 | 44 Cassettes | 600 x 300 x 900mm | CD005, CD010, CD015 | 38.9kg |



Drug & Medicine Storaç

*Transformer casing adds 105mm to overall height **CD010 not compatible with PC836/FS/IND & PC836/FS/IND/LX *** CD015 not compatible with PC839/FS/IND & PC839/FS/IND/LX



PC536/MDS (Frames not included)

Drug & Medicine Cabinets - Blister Pack (MDS)

Features

- Independently tested to BS 2881: 1989 (Security Level One) Specification for Cupboards for the Storage of Medicines in Health Care Premises • BS 3621 anti-pick & drill, key retaining five pin cylinder lock (supplied with 3 keys)
- Hooks at rear of cabinet support frames for the popular blister pack Monitored Dosage System (MDS)
- Key profile unique to Hospital Metalcraft Ltd, offering 250,000 differs
- Compatible with CLIQ[®], integrating electronics into high security mechanical disc cylinders & keys, combining electronic security & remote management in real time Optional Digital Push Button Lock (does not conform to BS 2881: 1989)
- 12v DC LED lighting/warning system (model dependent) with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug. CE marked & conforms to the low voltage directives 73/23/EEC & 89/336/EC
- Lighting transformer housing adds 105mm to the overall height
- 0.9mm (20 gauge) mild steel welded construction
- Mounting holes provided in back (fixings not included)
- Modular design allows for easy interchange with HTM 63 storage systems

Standard Colours

Paint - White

Options

- DGL Digital Push Button Lock
- SL Suited Locks (contact us for more information)
- PC/500/ST Sloping Top (500mm Wide Cabinets)
- PC/600/ST Sloping Top (600mm Wide Cabinets)
- PC/1000/ST Sloping Top (1000mm Wide Cabinets)
- PC/RLBS Remote Light Box (Surface Mounted)
- PC/RLBF Remote Light Plate (Flush Mounted)
- PCB Buzzer Alarm (Factory Fit To Cabinets With Lights Only)

Single Door Cabinets

| Right Hand Opening | | Left Hand Opening | | | | | |
|--------------------|-------------|-------------------|--------------|-----------------|--------------------|--------------------------------|--------|
| | With Light* | Without Light | With Light* | Without Light | MDS Frame Capacity | Overall Dimensions (w x d x h) | Weight |
| | PC536/MDS | PC536/MDS/LX | PC536/LH/MDS | PC536/LH/MDS/LX | 2 Frames | 500 x 300 x 600mm | 26.9kg |
| | PC539/MDS | PC539/MDS/LX | PC539/LH/MDS | PC539/LH/MDS/LX | 4 Frames | 500 x 300 x 900mm | 31.4kg |
| | PC639/MDS | PC639/MDS/LX | PC639/LH/MDS | PC639/LH/MDS/LX | 6 Frames | 600 x 300 x 900mm | 35.2kg |

Double Door Cabinets

| Single Lock (Secures Both Doors) | | Each Door Fitted With A Lock | | * Transformer casing adds 105mm to overall heigh | | |
|----------------------------------|---------------|------------------------------|-------------------|--|--------------------------------|--------|
| With Light* | Without Light | With Light* | Without Light | MDS Frame Capacity | Overall Dimensions (w x d x h) | Weight |
| PC1039/MDS | PC1039/MDS/LX | PC1039/IND/MDS | PC1039/IND/MDS/LX | 8 Frames | 1000 x 300 x 900mm | 26.9kg |

Controlled Drug Cabinets - Multi-Point Locking

Features

- Exceeds the statutory requirement for the storage of controlled drugs, "The Misuse of Drugs (Safe Custody) Regulations 1973"
- Independently security tested by Sold Secure attaining "Silver Accreditation" based on a realistic evaluation of the attack resistance of Controlled Drug Cabinets • BS 3621 anti-pick & drill, key retaining six pin cylinder lock (supplied with 2 keys & identification card)
- Multi-point locking, operated using a handle mounted on the door (independently tested)
- Key profile unique to Hospital Metalcraft Ltd. offering 250.000 differs Compatible with CLIQ[®], integrating electronics into high security mechanical disc
- cylinders & keys, combining electronic security & remote management in real time Lock protected by hardened steel plates & stainless steel escutcheon
- 3.0mm (10 gauge) mild steel, seamlessly welded construction
- Mounting holes provided in back (fixings not included)
- Fixed shelf

Standard Colours

Paint - White

Options

- MS910 12v DC LED Warning Light System (compatible with MS005 & MS010)
- MS911 12v DC LED Warning Light System (compatible with MS015)

| ltem | Description | Shelves | Lighting / Warnin |
|-------|-------------|---------|-------------------|
| MS005 | Small | 0 | MS910 |
| MS010 | Medium | 1 | MS910 |
| MS015 | Large | 1 | MS911 |
| | | | |



MS005



ng System (Optional) Overall Dimensions (w x d x h) Weight 335 x 265 x 300mm 17.0kg 335 x 265 x 550mm 26.0kg 37.0kg 500 x 325 x 550mm



Controlled Drug Cabinets - Single Point Locking

Features

- Exceeds the statutory requirement for the storage of controlled drugs "The Misuse of Drugs (Safe Custody) Regulations 1973" (with the exception of CD040, CD045, CD140 & CD145)
- BS 3621 anti-pick & drill, key retaining six pin cylinder lock (supplied with 2 kevs & identification card)
- Single-point locking
- Key profile unique to Hospital Metalcraft Ltd, offering 250,000 differs per profile
- Compatible with CLIQ[®], integrating electronics into high security mechanical disc cylinders & keys, combining electronic security & remote management in real time
- Can be supplied with a micro-switch to allow connection to remote lighting/nurse call system by others (CD905), or supplied with a 12v DC LED lighting/warning system with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug. CE marked & conforms to the low voltage directives 73/23/EEC & 89/336/EC
- Lighting transformer casing adds 105mm to the overall height
- 1.5mm (16 gauge) mild steel, seamless welded construction
- Mounting holes provided in back (fixings not included)
- Adjustable & removable shelves
- Security of models CD005, CD010 & CD015 can be enhanced by concealing inside one of our 'Drug & Medicine Cabinets'
- Modular design allows for easy interchange with HTM 63 storage systems
- Supplied right hand hinging as standard, left hand available on request

Standard Colours

· Paint - White

Options

- CD905 Micro-switch to connect to remote light box / alarm
- CD910 / CD911 / CD912 / CD913 12v DC LED Lighting/ Warning System (for compatibility see tables below)
- PCB Buzzer Alarm (Factory Fit To Cabinets With Lights Only)

| R/H Door | L/H Door | Shelves | Lighting / Warning System (Optional) | Overall Dimensions (w x d x h) | Weight |
|----------|----------|---------|--------------------------------------|--------------------------------|--------|
| CD005 | - | 1 | CD910 | 210 x 270 x 300mm | 8.0kg |
| CD010 | CD110 | 1 | CD910 | 335 x 270 x 300mm | 9.0kg |
| CD015 | CD115 | 2 | CD910 | 335 x 270 x 550mm | 17.0kg |
| CD020 | CD120 | 2 | CD910 | 500 x 300 x 550mm | 20.0kg |
| CD025 | CD125 | 3 | CD910 | 500 x 300 x 850mm | 37.0kg |

- additional light/nurse call system
- Optional buzzer



| L/H Door | Shelves | Lighting / Warni |
|----------|-------------------------|--------------------|
| CD130 | 3 | CD911 |
| CD135 | 4 | CD911 |
| CD140 | 6 | CD912 |
| CD145 | 6 | CD913 |
| | CD130 CD135 CD140 | CD135 4 CD140 6 |

CD125

| System (Optional) | Overall Dimensions (w x d x h) | Weight |
|-------------------|--------------------------------|--------|
| | 500 x 450 x 850mm | 42.0kg |
| | 500 x 450 x 1250mm | 59.0kg |
| | 500 x 315 x 1900mm | 73.0kg |
| | 500 x 465 x 1900mm | 91.0kg |
| | | |





CD145

Auditable Access Control for Medicines Management



What is CI IQ[®]?

CLIQ® is the very latest in user-friendly, web-based security management. Part of the Protec² CLIQ® from Ablov range, CLIQ® is designed for remote sites, enabling the commissioning, programming and full administration of security management centrally via web-based software.

How Does It Work?

CLIQ[®] integrates electronics into mechanical cylinders and keys, providing a security solution that can be managed remotely. In short, this works like an electronic handshake: when the key is inserted it engages with the cylinder, sending a random number back to the key. The key replies with an encrypted response which is decrypted by the cylinder and if valid the cylinder then allows access. It is important to remember that the mechanical combination must match as well as the digital combination. CLIQ® allows a single key to be used with multiple cabinets or trolleys.be used with multiple cabinets or trolleys.

System Features

- Secure Provides dual security through the patented Protec² CLIQ[®] from Ablov mechanical disc cylinders and keys along with DES encryption, a recognised electronic standard
- · Easily Programmable Keys can be added guickly and programmed to specific time and date schedules
- Fully Auditable- Audit trails can be accessed from both the cylinder and key
- · Reduces Keys One key can be programmed to open multiple cabinets or trolleys
- The CLIQ[®] environment works in an Amazon cloud hosted server
- · maintained and supported by Abloy allowing remote user access via a secure
- · connection from anywhere in the world
- · Cost Effective Low installation & maintenance costs along with a long lifecycle (requires no wires or system updates, power for the cylinder is provided by the key)
- Easily Retrofitable CLIQ[®] can be used with a variety of locking solutions including euro lock profile cases, padlocks, camlocks & furniture locks

System Benefits

- Save valuable nursing time which can be used to deliver patient care
- Reduces delay in administering medication
- Improves drug & medicine security
- No more unlocked drug & medicine cabinets
- Eliminates the need to replace locks as a result of lost keys
- · One key can be programmed to open multiple cabinets or trolleys
- · The key is made personal to that member of staff, resulting in increased ownership and a reduction in the number of keys lost
- Easy to apply remote management, saves money and has the ability to cover several hospitals
- Time and date based access rights, with a full audit trail allow for full visibility of when visits were made
- · New keys can be added to the system quickly and easily, removal of lost ones is hassle-free
- Reduction in time spent looking for the person holding the key



in the state



During The Shift

Operation

The key will only open cabinets where access permissions have been granted. Each time the key is used in a cabinet fitted with an PROTEC² cylinder an audit trail is created on both the key and cylinder



System Components



CLIQ® Kev

· Provides dual security through the patented PROTEC² mechanical disc cylinders and keys along with DES encryption, a recognised electronic standard Thin, light but incredibly strong key

PROTEC² Lock Barrel

- A bump proof mechanism
- Exceptionally smooth key entry and operation
- Can be master keyed and keyed alike with cylinders, cam and furniture locks as well as padlocks



Wall Programming Device

- Updates keys remotely (linking the server and users keys)
- Offers high level data encryption
- Allows the use of recurring validation, forcing users to validate their keys at the start of their shift whilst automatically downloading audit trails which are relayed to the server

Traka21

- Sophisticated stand-alone key management system to control the issuing of keys
- · Keys are accessed using PIN number technology
- Holds up to 21 keys or keysets



End of Shift

Start of Shift

Kev Validation & Authorisation

authorisations are uploaded to the key

Recurring validation forces users to

place their key into the Wall PD. User

Automated Data Transfer

Key is placed into the Wall PD. Previous key transactions are automatically downloaded from the key and stored in the audit trail archive

Questions & Answers Answers to common guestions about CLIQ[®]



O How secure is it?

- A CLIQ[®] uses a DES encryption (Data Encryption Standard) so it is just as safe as online banking
- Q Can I use this with an existing Abloy master key system?
- A Yes: CLIQ[®] can be integrated with an existing PROTEC² master key system
- What is the benefit?
- A CLIQ[®] will give you some of the benefits of access control, such as adding, deleting and amending access rights, take audit trails from keys and cylinders as well as setting the key to expire on a specific date at a specific time. Also, lost or stolen key can be easily be deleted electronically

Q Is the system expensive to install?

- A The system has a low installation & maintenance cost along with a long life cycle. CLIQ[®] cylinders are designed to fit into Euro Profile mechanical lock cases. No wiring is required to the cabinet. The wall programming device requires an Ethernet connection and power through a 12-24V adapter or alternatively power over Ethernet
- Q Who maintains the system?
- A System maintenance is normally assigned to a member of the hospital staff who will be responsible for programming/issuing of keys along with interrogation of audit trails
- Q How many audit trails can be collected from a single key?
- A The programmable user key stores the last 1.800 transactions and the last 20 foreign transactions
- Q How large is the cylinder memory?
- A 200 key groups; each key group holding 65,000 individual unique key combinations. This provides up to 13 million unique possible key combinations in one cylinder and can retain a black list of 3.000 keys
- Q Can I take an audit trail from a cylinder?
- A Yes, the cylinder stores an audit trail of the last 1,000 key insertions
- Q What happens if a key is lost?
- A We recommend that time schedules are set to each key, so if the key is not returned it will expire. Keys can be removed from the cylinder with the programming key or directly from the system. If subsequently found keys can be added back and re-issued. In the event that an audit trail is required for a lost key. lock cylinders can be interrogated using a programming key (each cylinder stores an audit trail of the last 1,000 key insertions)

Q How long do the batteries last?

- A The battery is located in the key, not the cylinder. Under normal usage a battery will operate up to 20.000 cvcles or 2 vears
- O How is the information sent?
- A Via the internet to a wall mounted PD
- Can the information be taken off the Wall mounted PD by anyone?
- No, there is no information stored on the PD. it is stored on a server until the correct key is inserted and only then is the information securely transferred

Q Is it possible to access the system remotely?

- A The CLIQ[®] environment works in an Amazon cloud hosted server environment ISO27001 Accredited. maintained and supported by Ablov allowing remote user access via a secure connection from anywhere in the world
- Where is the server?
- A Amazon Republic Of Ireland
- Q Can I host this myself?
- A Yes, but there would be additional costs
- What if there is no internet connection on a particular site?
- A mobile PD is available which will work on the GSM network
- Q How can I better control the issuing of keys?
- A Through the use of Traka21 (key safe) keys can be automatically assigned to specific users. Keys are issued using either the staff members own security card or via pin number technology.
- Q Is it possible to integrate the CLIQ[®] key with other systems?
- A Unique RFID tag identifiers can be inserted into the key allowing it to be used with other MIFARE card based RFID systems
 - Self-Administration Schemes
 - Medicines Management
- Patient Records



Drug & Medicine Cupboards - Single Door

Features

- Independently tested to BS2881: 1989 (Security Level One) Specification for Cupboards for the Storage of Medicines in Health Care Premises BS3621 anti-pick & drill, key retaining five pin cylinder lock (supplied with 3 keys)
- Single door cupboards available as either 870 / 1870mm high
- Available with either...
- Plastic Travs
- Shelves
- Shelves & Dividers
- Key profile unique to Hospital Metalcraft Ltd, offering 250,000 differs
- Compatible with CLIQ[®], integrating electronics into high security mechanical disc cylinders & keys, combining electronic security & remote management in real time
- Optional Digital Push Button Lock (does not conform to BS 2881: 1989)
- Worksurface supplied with 870mm cupboards (removable to assist with positioning under existing worktops)
- 0.9mm (20 Gauge) mild steel welded construction
- Modular design allows for easy interchange with HTM 63 storage systems

Dimensions (w x d x h)

- Overall Dimensions
- 870mm 500 x 536 x 870mm
- 1870mm 500 x 536 x 1870mm
- Travs
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal

Standard Colours

Paint - White

Options

- SL Suited Locks (contact us for more information)
- DGL Digital Push Button Lock
- MC/JK/870 Joining Kit & Full Width Worktop (required when connecting 870mm cupboards)
- MC/JK/1870 Joining Kit (required when connecting 1870mm cupboards)
- MC/ST Sloping Top (1870mm cupboards)
- MC/DIV Set of 16 Dividers

| R/H Hinged | L/H Hinged | Description | Weight |
|---------------|---------------|--|--------|
| MC/R/870/3S2D | MC/L/870/3S2D | Single Door, Five Pin Cylinder Lock, 870mm High, Plastic Trays | 45.6kg |
| MC/R/870/FS | MC/L/870/FS | Single Door, Five Pin Cylinder Lock, 870mm High, Flat Shelves | 37.0kg |



| R/H Hinged | L/H Hinged | Description | Weight |
|----------------|----------------|---|--------|
| MC/R/1870/7S5D | MC/L/1870/7S5D | Single Door, Five Pin Cylinder Lock, 1870mm High, Plastic Trays | 88.6kg |
| MC/R/1870/FS | MC/L/1870/FS | Single Door, Five Pin Cylinder Lock, 1870mm High, Flat Shelves | 69.0kg |
| MC/R/1870/FSD | MC/L/1870/FSD | Single Door, Five Pin Cylinder Lock, 1870mm High, Flat Shelves & Dividers | 79.0kg |

Plastic Trays

- Interchangeable with our popular CT range (post 2016)
- 870mm Cupboards 3 x 100mm, 2 x 150mm trays
- 1870mm Cupboards 7 x 100mm. 5 x 150mm travs
- Translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- Each cupboard is supplied with perforated coloured label strips (white, blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)

Shelves

- Flat shelves adjustable at 50mm increments
- 870mm Cupboards 3 x shelves
- 1870mm Cupboards 8 x shelves
- Suitable for the storage of MDS frames (if required)
- 870mm Cupboards 6 x frames
- 1870mm Cupboards 16 x frames

Shelves & Dividers

- Flat shelves adjustable at 50mm increments
- Removable dividers 100mm high adjustable at 25mm increments
- 870mm Cupboards N/A
- 1870mm Cupboards 8 x shelves & 16 x dividers

Sloping Top

Features

- Compatible with 300/500/600/800 & 1000mm wide PC cabinets & 500mm wide MC cupboards
- Prevents items being stored on the top of the cabinet & assists cleaning
- Increases the overall height of the cabinet by 193mm
- Void within the top provides a space to locate the transformer for the lighting/warning system
- Access to the transformer provided via a removable plate

Standard Colours

· Paint - White

Warning Light System

Features

- Suitable for use with Bristol Maid MS & CD Controlled Drug cabinets
- Provides a visual indication when the cabinet is open
- LED illuminates...
 - · Green when the door is closed
 - Red when the door is open
 - Flashes red when open for 2 ¹/₂ minutes
- Operates using a door activated micro-switch connected to a 12v DC LED lighting/warning system (model dependent) with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug, CE marked & conforms to the low voltage directives 73/23/EEC & 89/336/EC
- Lighting transformer housing adds 105mm to the overall height



| Item | Description |
|------------|--|
| MS Range C | Compatible |
| MS910 | 12v DC LED Lighting/Warning System - Compatible with MS005 & MS010 |
| MS911 | 12v DC LED Lighting/Warning System - Compatible with MS015 |
| CD Range C | ompatible |
| CD910 | 12v DC LED Lighting/Warning System - Refer To Page 138 For Compatibility |
| CD911 | 12v DC LED Lighting/Warning System - Refer To Page 139 For Compatibility |
| CD912 | 12v DC LED Lighting/Warning System - Refer To Page 139 For Compatibility |
| CD913 | 12v DC LED Lighting/Warning System - Refer To Page 139 For Compatibility |

| | 12v DC LED Lighting/Warning System - Compatible with MS015 |
|-------|--|
| je Co | |
| | 12v DC LED Lighting/Warning System - Refer To Page 138 For Compatibility |
| | 12v DC LED Lighting/Warning System - Refer To Page 139 For Compatibility |
| | 12v DC LED Lighting/Warning System - Refer To Page 139 For Compatibility |
| | 12v DC LED Lighting/Warning System - Refer To Page 139 For Compatibility |
| | |

Joining Kits - MC Range of Cabinets

Features

- Allows two units to be joined together
- 870mm joining kit supplied with a full width worktop (2 units)



| Item | Description | |
|------------|----------------------------------|----|
| MC/JK/870 | Joining Kit & Full Width Worktop | 8 |
| MC/JK/1870 | Joining Kit | 18 |


Remote Light Box

Features

- Connects to all Bristol Maid PC, MS & CD cabinets fitted with a light)
- Provides a remote visual indication when the cabinet is open
- Models available for surface or flush mounting
- Units fitted without a light require a separate micro-switch & a 12v DC supply



| Item | Description | Ove |
|---------|------------------------------------|-----|
| PC/RLBS | Remote Light Box - Surface Mounted | 1 |
| PC/RLBF | Remote Light Plate - Flush Mounted | 1 |



| Ite | | Description | Weight |
|-----|----------|-------------------------|--------|
| PC | /300/ST | 300mm Wide PC Cabinets | 1.3kg |
| PC | /500/ST | 500mm Wide PC Cabinets | 1.5kg |
| PC | /600/ST | 600mm Wide PC Cabinets | 1.6kg |
| PC | /800/ST | 800mm Wide PC Cabinets | 1.7kg |
| PC | /1000/ST | 1000mm Wide PC Cabinets | 1.8kg |
| MC | C/ST | 500mm Wide MC Cupboards | 2.9kg |

Digital Push Button Lock

Features

- Requires a five digit code to operate (set by customer)
- Cabinets fitted with this lock will not conform to BS2881: 1989 (Security Level One) Specification for Cupboards for the Storage of Medicines in Health Care Premises



Micro Switch

Features

 Micro-switch to allow connection to remote lighting/nurse call system by others



370mm High Units 370mm Hiah Units Digital Push Button Lock

Dividers - MC Range of Cabinets

Features

DGL

- Spare / replacement dividers for MC cupboards supplied with flat shelves & dividers
- Removable 100mm high dividers locate into slots on the shelves (adjustable at 25mm increments)
- Compatible with 870mm & 1870mm high units

erall Dim

100 x 50 x 100 x 50 x



| lensions | | | |
|----------|--------|--------------------|--|
| x 150mm | Item | Description | |
| x 150mm | MC/DIV | Set of 16 Dividers | |

| em | Description |
|-------|----------------------|
| 1S905 | Micro-switch (MS Ran |

nge Compatible) CD905 Micro-switch (CD Range Compatible) PC905 Micro-switch (PC Range Compatible)

Buzzer Alarm

Features

- Replacement circuit board with buzzer
- Compatible for PC, MS & CD cabinets fitted with a light unit
- Buzzer pulses on/off after the door has been open for 2 1/2 minutes
- Factory fit only to cabinets fitted with lights

| Item | Description |
|------|--------------|
| PCB | Buzzer Alarm |

Suited Locks

Features

- · Suited locks (allows multiple locks / cabinets to operate using the same key)
- · Not all cabinets can have suited locks contact customer service

| Item | Description |
|------|--------------|
| SL | Suited Locks |

Patients / Residents Own Medication Cabinets - Self Administration





SAC/313/S/VP

SAC/FS/2 - Shelf, 200mm Wide

SAC/FS/3 - Shelf, 300mm Wide

SAC/FS/4 - Shelf, 400mm Wide

SAC/WP - Mounting Plate, Wall

SAC/LT - Mounting Plate, Locker Top

Options

Features

Suitable for use with...

- Drugs retained in original packs
- Individual MDS packs with dividers & name cards
- Nomad cassettes
- Roll type blister packs
- Available as side or bottom hinged fitted with either a CAM or electronic push button lock
- CAM Supplied with two keys & uses 40 different master suites (for security we recommend that each ward/wing operates a different series, master keys available at additional cost)
- Electronic Push Button Battery operated
- Available with a solid door or with an inset perspex visibility panel
- 0.9mm (20 gauge) mild steel, spot welded construction
 Mounting holes in back/bottom (fixings not included)
- Range of mounting plates (optional)
- Shelves not included (optional)
- Supplied with a self-adhesive name card holder (length 60mm)

Electronic Push Button Lock

- Uses a 10 button keypad and a # button. The # button is hidden beneath the handle when in the locked position
- Easy programming & code change
- Four operating codes with differing levels of authority Master code
 - Sub-master code
 - User code
 - Technician code
- User codes are 4 digits
- Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code

Side Hinged

| | CAM Lock | Electric | Push Button Lock | | | |
|------------|--------------------------|--------------|--------------------------|--------------------------------|---------------------|----|
| Solid Door | Perspex Visibility Panel | Solid Door | Perspex Visibility Panel | Overall Dimensions (w x d x h) | Shelf Compatibility | |
| SAC/213/S | SAC/213/S/VP | SAC/213/S/CL | SAC/213/S/CL/VP | 210 x 155 x 315mm | SAC/FS/2 | 2. |
| SAC/313/S | SAC/313/S/VP | SAC/313/S/CL | SAC/313/S/CL/VP | 310 x 155 x 315mm | SAC/FS/3 | 3. |
| SAC/413/S | SAC/413/S/VP | SAC/413/S/CL | SAC/413/S/CL/VP | 410 x 155 x 315mm | SAC/FS/4 | 4. |



Bottom Hinged

| | CAM Lock | Electric | Push Button Lock | | | |
|------------|--------------------------|--------------|--------------------------|--------------------------------|---------------------|--------|
| Solid Door | Perspex Visibility Panel | Solid Door | Perspex Visibility Panel | Overall Dimensions (w x d x h) | Shelf Compatibility | Weight |
| SAC/213/B | SAC/213/B/VP | SAC/213/B/CL | SAC/213/B/CL/VP | 210 x 155 x 315mm | SAC/FS/2 | 2.6kg |
| SAC/313/B | SAC/313/B/VP | SAC/313/B/CL | SAC/313/B/CL/VP | 310 x 155 x 315mm | SAC/FS/3 | 3.4kg |
| SAC/413/B | SAC/413/B/VP | SAC/413/B/CL | SAC/413/B/CL/VP | 410 x 155 x 315mm | SAC/FS/4 | 4.1kg |

Shelves

Features

- Flat shelves (100mm deep)
- Removable to assist cleaning, allow the storage of larger items

Standard Colours

Paint - White





| ltem | Description | Wei |
|----------|---|-----|
| SAC/FS/2 | Shelf - 200mm Wide Cabinets (SAC/213 Range) | 0. |
| SAC/FS/3 | Shelf - 300mm Wide Cabinets (SAC/313 Range) | 0. |
| SAC/FS/4 | Shelf - 400mm Wide Cabinets (SAC/413 Range) | 0. |

Mounting Plates

Features

- Compatible with all SAC cabinets using existing holes (fixings not provided)
- · SAC/WP Cabinet secured to wall plate using a small thumbscrew. Allows the cabinet to be repositioned
- SAC/LT Mounts to upper surface of traditional bedside cabinets

Standard Colours



Patients / Residents Own Medication Cabinets - Monitored Dosage, MDS Frames

Features

- Suitable for the storage of MDS frames with dividers & name cards (max 100mm depth)
- Models available to hold one or two frames (not included)
- Small One frame
- Large Two frames
- Available with either a CAM or electronic push button lock
- · CAM Supplied with two keys, master keys available at additional cost Electronic Push Button - Battery operated
- Mounting holes provided in back (fixings not included)
- Supplied with an adjustable shelf (SA045 & SA045/CL only)
- Constructed from melamine faced chipboard offering an extremely strong & durable, yet aesthetically pleasing appearance
- MDS frames not supplied

Standard Colours

Body/Door - Beech Woodgrain

Electronic Push Button Lock

- Uses a 10 button keypad and a # button. The # button is hidden beneath the handle when in the locked position
- Easy programming & code change
- Four operating codes with differing levels of authority
- Master code
- Sub-master code
- User code
- Technician code
- User codes are 4 digits
- Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code

| m | Description | Weight |
|--------|---|--------|
| C/FS/2 | Shelf - 200mm Wide Cabinets (SAC/213 Range) | 0.2kg |
| C/FS/3 | Shelf - 300mm Wide Cabinets (SAC/313 Range) | 0.2kg |
| C/FS/4 | Shelf - 400mm Wide Cabinets (SAC/413 Range) | 0.3ka |

| Item | Description | Overall Dimensions (w x d x h) | Weight |
|--------|-------------|--------------------------------|--------|
| SAC/WP | Wall | 190 x 12 x 162mm | 0.5kg |
| SAC/LT | Locker Top | 190 x 150 x 150mm | 0.9kg |

| CAM Lock | Electronic Push Button Lock | Capaci |
|----------|-----------------------------|---------|
| SA040 | SA040/CL | 1 x Fra |
| SA045 | SA045/CL | 2 x Fra |

www.bristolmaid.com







SA040 - frames not included

Overall Dimensions (w x d x h) Weight ame 325 x 245 x 420mm 7.5kg 12.1kg

ames 325 x 245 x 775mm

Tablet Counter

Features

- Triangular stainless steel tablet counter
- Etched legend to interpret rows (max 22) into number of tablets



| Storage | Boxes | - | Lig | htwe | eigł | nt |
|---------|-------|---|-----|------|------|----|
|---------|-------|---|-----|------|------|----|

Features

- Lightweight aluminium riveted construction
- Fitted with handle
- Single-point locking using a CAM lock (supplied with 2 keys)
- · Cannot be wall mounted

Standard Colours

· Paint - White



| tem | Description | Overall Dimensions (w x d x h) | Neight |
|--------|-------------|--------------------------------|--------|
| SB/221 | Small | 200 x 200 x 150mm | 0.9kg |
| SB/321 | Medium | 300 x 200 x 150mm | 1.1kg |
| SB/421 | Large | 400 x 200 x 150mm | 1.4kg |

Pharmacy & Vaccine Refrigerators - IntelliCold[™], Solid & Glass Doors

Features

- Compliant with the Department of Health Green Book
- Recommendations in the General Pharmaceutical Council Standards
- Royal College of Veterinary Surgeons Guidelines
- Designed specifically for the safe storage of vaccines and pharmaceuticals
- Available as bench top, under bench, wall fitting or free-standing
- Temperature range +2°C to +8°C
- Second wandering temperature probe (can be placed anywhere within the refrigerator, ideal for pharmacists needing to comply with MHRA regulations)
- IntelliCold[®] touch screen controller
- Automatic 24/7 temperature & fridge operation data recording (micro SD card required)
- Min/max temperature control, display & recording
- Displays temperature from both probes
- High/low temperature alarm
- Door open alarm
- Mains failure alarm (requires batteries)
- Alarm battery back-up
- Internal LED light
- Access port & remote alarm contacts
- Fan assisted cooling
- Off cycle auto defrost
- Low GWP environmentally friendly refrigerant
- Self-closing doors with door open alarm
- Available with either a CAM or digital push button lock CAM - Supplied with two keys
- Digital Push Button Mains operated (uses a 4-digit pin)
- Adjustable shelves
- Micro SD card (optional)
- UKAS traceable calibration service (optional)
- Repeat alarm (optional)

Standard Colours

- Exterior White/Blue
- Interior White

| CAM Lock | Digital Push Button | Description | Shelves | Overall Dimensions (w x d x h) | Capacity (Litres) | Weight |
|------------|---------------------|----------------------|---------|--------------------------------|-------------------|--------|
| Solid Door | | | | | | |
| 5RLDF0119 | 5RLDF0119DIGILOCK | Bench/Wall Mounting* | 2 | 450 x 510 x 538mm | 36 Litres | 30.0kg |
| 5RLDF0219 | 5RLDF0219DIGILOCK | Bench/Wall Mounting* | 3 | 450 x 510 x 740mm | 66 Litres | 34.0kg |
| 5RLDF0519 | 5RLDF0519DIGILOCK | Under Bench | 3 | 600 x 600 x 835mm | 150 Litres | 51.0kg |
| Glass Door | | | | | | |
| 5RLDG0519 | 5RLDG0519DIGILOCK | Under Bench | 3 | 600 x 600 x 835mm | 150 Litres | 54.0kg |

*Integral Wall Mounting Brackets

| · | |
|----|--|
| b | |
| | |
| د | |
| 5 | |
| 2 | |
| - | |
| 5 | |
| 20 | |
| - | |
| د | |
| | |
| | |

| | / | |
|---|---|--|
| / | | |
| | | |
| | | |
| | | |

| Item | Description | Dimensions (w x d x h) |
|-------|----------------|------------------------|
| PE/30 | Tablet Counter | 205 x 195 x 10mm |





5RLDF0219












Examination Treatment Couches & Chairs

- Examination & Treatment Couches 154
 - Leg & Arm Rests 162
 - Couch Steps 163
- Bariatric Examination & Treatment & Couches 164
 - Therapy Couches 170
 - Examination & Treatment Chairs 174
 - Bariatric Examination & Treatment Chairs 176
 - Phlebotomy Chairs 178
 - Options 181







Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Hydraulic in Aquamarine - EC1205/VAQ



Examination & Treatment Couch Two Section

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric or hydraulic height adjustment
- Height adjustment from 490-950mm
- Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Anti-trap gas strut fitted with a release lever used on the backrest section (optional electric)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless upholstery using flame retardant anti-microbial vinyl
- Frames prepared to accept optional accessories
- Maximum capacity 225kg

Dimensions (w x d x h)

- 630 x 1790 x 490-950mm
- Paint White Vinyl - Bristol Blue

Standard Colours

Options

- E02 Paper Roll Holder
- E03 Side Support Rails (Pair)
- E04 Side Support Rails C/W Cushions (Pair)
- E05 Breather Hole & Plug
- E06 Battery Back Up
- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole
- E11 Lithotomy Stirrups (Pair)
- E12 Leg Supports (Pair)

| ltem | Description | Weight |
|------------|---|--------|
| EC1205/VBB | Hydraulic | 68.8kg |
| EC1210/VBB | Electric, Hand Switch | 68.8kg |
| EC1215/VBB | Electric, Foot Switch | 68.8kg |
| EC1220/VBB | Electric with Electric Backrest - Hand Switch | 71.1kg |
| EC1225/VBB | Electric with Electric Backrest - Foot Switch | 71.1kg |

Examination & Treatment Couch Two Section, Mobile

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric or hydraulic height adjustment
- Height adjustment from 490-950mm
- Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Anti-trap gas strut fitted with a release lever used on the backrest section (optional electric)
- 125mm swivel castors featuring a central locking & steering facility
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless upholstery using flame retardant anti-microbial vinvl
- Swing down side support rails fitted to either side
- Frames prepared to accept optional accessories
- Maximum capacity 225kg

Standard Colours Dimensions (w x d x h)

- 630 x 1790 x 490-950mm
- Paint White Vinvl - Bristol Blue

Options

- E02 Paper Roll Holder
- E04/U Side Support Rails C/W Cushions Upgrade (Pair)
- E05 Breather Hole & Plug
- E06 Battery Back Up
- E11 Lithotomy Stirrups (Pair)
- E12 Leg Supports (Pair)

| Description |
|---|
| Hydraulic |
| Electric, Hand Switch |
| Electric, Foot Switch |
| Electric with Electric Backrest - Hand Switch |
| Electric with Electric Backrest - Foot Switch |
| |

154



Weight

75.0kg

75.0kg

75.0kg

77.3ka

77.3kg

sales@bristolmaid.com +44 (0)1258 484455





Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric, Foot Switch in Onyx Black - EC1315/VOB



Examination & Treatment Couch **Three Section**

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric or hydraulic height adjustment
- Height adjustment from 490-950mm
- Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Adjustable angle footrest from -25 to +75 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless upholstery using flame retardant anti-microbial vinyl
- · Frames prepared to accept optional accessories
- Maximum capacity 225kg

Dimensions (w x d x h)

- Standard Colours
- 630 x 1790 x 490-950mm
- Paint White Vinvl - Bristol Blue

Options

- E02 Paper Roll Holder
- E03 Side Support Rails (Pair)
- E04 Side Support Rails C/W Cushions (Pair)
- E05 Breather Hole & Plug
- E06 Battery Back Up
- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole

| Item | Description | Weight |
|------------|---|--------|
| EC1305/VBB | Hydraulic | 70.2kg |
| EC1310/VBB | Electric, Hand Switch | 70.2kg |
| EC1315/VBB | Electric, Foot Switch | 70.2kg |
| EC1320/VBB | Electric with Electric Backrest - Hand Switch | 72.5kg |
| EC1325/VBB | Electric with Electric Backrest - Foot Switch | 72.5kg |

Examination & Treatment Couch Three Section, Mobile

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric or hydraulic height adjustment
- Height adjustment from 490-950mm
- Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Adjustable angle footrest from -25 to +75 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- 125mm swivel castors featuring a central locking & steering facility
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless & contoured upholstery using flame retardant & anti-microbial vinyl
- Swing down side support rails fitted to either side
- Frames prepared to accept optional accessories
- Maximum capacity 225kg

Dimensions (w x d x h)

- 630 x 1790 x 490-950mm
- Paint White Vinvl - Bristol Blue

Options

- E02 Paper Roll Holder
- E04/U Side Support Rails C/W Cushions Upgrade (Pair)
- E05 Breather Hole & Plug
- E06 Battery Back Up

| Item | Description |
|------------|--|
| EC1330/VBB | Hydraulic |
| EC1335/VBB | Electric, Hand Switch |
| EC1340/VBB | Electric, Foot Switch |
| EC1345/VBB | Electric with Electric Backrest - Hand |
| EC1350/VBB | Electric with Electric Backrest - Foot S |
| | |

156



- Service Contract ontact us fo



Standard Colours

Weight 76.4kg 76.4kg 76.4kg Switch 78.7ka Switch 78.7kg



sales@bristolmaid.com +44 (0)1258 484455



Examination & Treatment Couch **Four Section**

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric or hydraulic height adjustment
- Height adjustment from 490-950mm
- Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Adjustable angle footrest from -25 to +75 degrees
- · Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- · Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless upholstery using flame retardant anti-microbial vinyl
- Frames prepared to accept optional accessories
- Maximum capacity 225kg

Dimensions (w x d x h)

- Standard Colours
- 630 x 1790 x 490-950mm Paint - White Vinvl - Bristol Blue

Options

- E02 Paper Roll Holder
- E03 Side Support Rails (Pair)

• E04 - Side Support Rails C/W Cushions (Pair)

- E05 Breather Hole & Plug
- E06 Battery Back Up
- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole

| tem | Description | Weight |
|-----------|---|--------|
| C1405/VBB | Hydraulic | 70.3kg |
| C1410/VBB | Electric, Hand Switch | 70.3kg |
| C1415/VBB | Electric, Foot Switch | 70.3kg |
| C1420/VBB | Electric with Electric Backrest - Hand Switch | 72.6kg |
| C1425/VBB | Electric with Electric Backrest - Foot Switch | 72.6kg |
| | | |

Examination & Treatment Couch Four Section, Mobile

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric or hydraulic height adjustment
- Height adjustment from 490-950mm
- Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Adjustable angle footrest from -25 to +75 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- 125mm swivel castors featuring a central locking & steering facility
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless & contoured upholstery using flame retardant & anti-microbial vinvl
- Swing down side support rails fitted to either side
- Frames prepared to accept optional accessories
- Maximum capacity 225kg

Dimensions (w x d x h)

- 630 x 1790 x 490-950mm
- Paint White Vinyl - Bristol Blue

Options

- E02 Paper Roll Holder
- E04/U Side Support Rails C/W Cushions Upgrade (Pair)
- E05 Breather Hole & Plug

| Item | Description |
|------------|--|
| EC1430/VBB | Hydraulic |
| EC1435/VBB | Electric, Hand Switch |
| EC1440/VBB | Electric, Foot Switch |
| EC1445/VBB | Electric with Electric Backrest - Hand |
| EC1450/VBB | Electric with Electric Backrest - Foot S |
| | |

Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Hydraulic in Sand Dune - EC1405/VSD



sales@bristolmaid.com +44 (0)1258 484455

- E06 Battery Back Up

| Description | Weight |
|--|--------|
| 3B Hydraulic | 70.3kg |
| 3B Electric, Hand Switch | 70.3kg |
| 3B Electric, Foot Switch | 70.3kg |
| 3B Electric with Electric Backrest - Hand Switch | 72.6kg |
| BB Electric with Electric Backrest - Foot Switch | 72.6kg |
| | |

Service Contract ontact us fo details

Standard Colours



76.6kg 76.6kg 76.6ka Switch 78.9kg Switch 78.9kg

Weight

VAQ sales@bristolmaid.com +44 (0)1258 484455

VBB

www.bristolmaid.com

VSD

VOB

VL

VGH

VFR

VCB

VCH

159



Examination & Treatment Couch Two Section, Vertical Lift

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Scissor action provides vertical lift, ideal for use in confined areas where space is limited
- Height adjustment from 560-1000mm
- Adjustable angle backrest 0 to +90 degrees
- Anti-trap gas strut fitted with a release lever used on the backrest section
- 125mm swivel braking castors
- Clearance beneath frame allows access for mobile hoists
- 230V ~. 50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- Seamless upholstery using flame retardant anti-microbial vinyl Paper roll holder (suitable for rolls max 510mm wide x 200mm) diameter)
- Frames prepared to accept optional accessories
- Maximum capacity 190kg

Dimensions (w x d x h)

• 610 x 1860 x 560-1000mm

Standard Colours

 Paint - White Vinyl - Bristol Blue

Options

- EC905 Side Support Rails (Pair)
- EC920 Breather Hole & Plug

Examination & Treatment Couch **Two Section, Fixed Height**

Features

- Fixed height approx. 786mm
- Adjustable angle backrest 0 to +90 degrees using a ratchet
- Adjustable feet to allow levelling on uneven floors
- Seamless upholstery using flame retardant anti-microbial vinyl
- Paper roll holder fitted to foot end (suitable for rolls max 510mm wide x 200mm diameter)
- Self-assembly model flat packed making it easier to ship
- Frames prepared to accept optional accessories
- Maximum capacity 125kg

Dimensions (w x d x h)

• 610 x 1860 x 795mm

Standard Colours

- Paint White
- Vinyl Bristol Blue

Options

- EC900 Drawer
- EC935 Drop Down Side Rails (Pair)



| Item | Description |
|-----------|-----------------------------------|
| EC005/VBB | Examination Couch - Assembled |
| EC010/VBB | Examination Couch - Self Assembly |

Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric, Foot Switch in Chocolate - EC020/VCH



| Item | Description | Weight |
|-----------|-----------------------|--------|
| EC015/VBB | Electric, Hand Switch | 54.0kg |
| EC020/VBB | Electric, Foot Switch | 54.0kg |

Service Contract Contact us fo details



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Assembled Couch in Chocolate - EC005/VCH



www.bristolmaid.com

161

Leg Rest - Adjustable

Features

- Independently adjustable legs
- Height & angle adjustment from 355-510mm
- Seamless upholstery using flame retardant anti-microbial vinyl
- Maximum capacity 25kg

Dimensions (w x d x h)

Standard Colours

FIRE

- Overall Dimensions 460 x 305 x 350-510mm
- Upholstered Pad 460 x 305mm

- Paint White
- Vinyl Bristol Blue



Features

- Four leg stand
- Height & angle adjustment from 495-750mm / 17 degrees
- Seamless upholstery using flame retardant anti-microbial vinyl
- Maximum capacity 15kg

Dimensions (w x d x h)

- Overall Dimensions 420 x 490 x 495-750mm
- Upholstered Pad 380 x 130mm





| Item | Description | Overall Di |
|-------|--------------------------------|-------------------|
| CS005 | Couch Step - One Tread | 519 x 284 |
| CS010 | Couch Step - One Tread, Handle | 519 x 284 |
| CS015 | Couch Step - Two Tread | 519 x 599 |
| | | |



VL

VOB

VCH

VCB

VPE

VBB

Item

G/16/TILT/VBB

VAQ

VFR

VGH

VSD

Description

VOB

Arm / Leg Rest - Adjustable

VL

VCH

VPE

Weight

3.5kg

VCB

VBB

Item

CS/28/VBB

VAQ

VFR

VGH

VSD





Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric with Electric Backrest - Hand Switch in Peach - EC2215/VPE



Examination & Treatment Couch **Two Section, Bariatric**

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric height adjustment
- Height adjustment from 520-950mm
- Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- · Anti-trap gas strut fitted with a release lever used on the backrest section (optional electric)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- 230V ~, 50/60H actuator fitted with a UK 3 pin plug & available foot switch
- Seamless upholstery using flame retardant antimicrobial viny
- Frames prepared to accept optional accessories
- Maximum capacity 325kg

Dimensions (w x d x h) Standard Colours

- 850 x 1790 x 520-950mm
- Paint White Vinvl - Bristol Blue

Options

- E03 Side Support Rails (Pair)
- E04 Side Support Rails C/W Cushions (Pair)
- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole
- E18 Battery Back Up
- E19 Paper Roll Holder
- E21 Breather Hole & Plug



| able with a hand or | 230V ~, 50/ foot switch |
|---------------------|---|
| | Seamless u |
| ıyl | Swing dowr |
| | |

- wing down side support rails fitted to either side Frames prepared to accept optional accessories
- Maximum capacity 325kg

Dimensions (w x d x h)

• 850 x 1790 x 520-950mm

Standard Colours

Paint - White

Features

& safety specifications

Electric height adjustment

 Adjustable angle backrest Manual -25 to +90 degrees

(optional electric)

Height adjustment from 520-950mm

Electric -22 to +90 degrees

Vinyl - Bristol Blue

Options

- E04/U Side Support Rails C/W Cushions Upgrade (Pair)
- E18 Battery Back Up
- E19 Paper Roll Holder
- E21 Breather Hole & Plug

| ltem | Description |
|------------|--|
| EC2225/VBB | Electric, Hand Switch |
| EC2230/VBB | Electric, Foot Switch |
| EC2235/VBB | Electric with Electric Backrest - Hand |
| EC2240/VBB | Electric with Electric Backrest - Foot |
| | |

Description Item EC2205/VBB Electric, Hand Switch EC2210/VBB Electric, Foot Switch EC2215/VBB Electric with Electric Backrest - Hand Switch EC2220/VBB Electric with Electric Backrest - Foot Switch

| - | |
|---|--|
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| 5 | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |

Service

- details

Contract ontact us fo





79.0kg


81.3kg

81.3kg











Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric with Electric Backrest - Hand Switch in Peach - EC2315/VPE



Examination & Treatment Couch Three Section, Bariatric

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric height adjustment
- Height adjustment from 520-950mm
- · Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Adjustable angle footrest from -25 to +75 degrees
- Anti-trap gas strut fitted with a release lever used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- 230V ~, 50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- Seamless upholstery using flame retardant antimicrobial vinyl
- · Frames prepared to accept optional accessories
- Maximum capacity 325kg

Standard Colours Dimensions (w x d x h)

- 850 x 1790 x 490-950mm
- Paint White Vinvl - Bristol Blue

Options

Item

EC2305/VBB

EC2310/VBB

EC2315/VBB

EC2320/VBB

- E03 Side Support Rails (Pair)
- E04 Side Support Rails C/W Cushions (Pair)
- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole
- E18 Battery Back Up
- E19 Paper Roll Holder
- E21 Breather Hole & Plug

Service Contract details



Examination & Treatment Couch Three Section, Mobile, Bariatric

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric height adjustment
- Height adjustment from 520-950mm
- Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Adjustable angle footrest from -25 to +75 degrees
- Anti-trap gas strut fitted with a release lever used on the back & footrest sections (optional electric backrest)
- 125mm swivel castors featuring a central locking & steering facility
- Clearance beneath frame allows access for mobile hoists
- 230V ~, 50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- Seamless upholstery using flame retardant antimicrobial vinvl
- Swing down side support rails fitted to either side
- Frames prepared to accept optional accessories
- Maximum capacity 325kg

Dimensions (w x d x h)

850 x 1790 x 520-950mm

Standard Colours

- Paint White
- Vinyl Bristol Blue

Options

- E04/U Side Support Rails C/W Cushions Upgrade (Pair)
- E18 Battery Back Up
- E19 Paper Roll Holder
- E21 Breather Hole & Plug

| ltem | Description |
|------------|--|
| EC2325/VBB | Electric, Hand Switch |
| EC2330/VBB | Electric, Foot Switch |
| EC2335/VBB | Electric with Electric Backrest - Hand |
| EC2340/VBB | Electric with Electric Backrest - Foot |
| | |

166

- Contact us fo



d Switch Switch

Weight 88.0kg 88.0kg 90.3kg 90.3kg VBB VAQ VFR VGH VSD VOB VL VCH

www.bristolmaid.com

167

VCB





Foot Switch

Hand Switch

Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric with Electric Backrest - Hand Switch in Graphite - EC2415/VGH



Examination & Treatment Couch Four Section, Bariatric

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric height adjustment
- Height adjustment from 520-950mm
- Adjustable angle backrest Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Adjustable angle footrest from -25 to +75 degrees
- Anti-trap gas strut fitted with a release lever used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- 230V ~, 50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- Seamless upholstery using flame retardant antimicrobial vinyl
- Frames prepared to accept optional accessories

Dimensions (w x d x h)

Standard Colours

- 850 x 1790 x 520-950mm
- Paint White

Options

- E03 Side Support Rails (Pair)
- E04 Side Support Rails C/W Cushions (Pair)
- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole
- E18 Battery Back Up
- E19 Paper Roll Holder
- E21 Breather Hole & Plug

Examination & Treatment Couch Four Section, Mobile, Bariatric

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric height adjustment
- Height adjustment from 520-950mm
- Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Adjustable angle footrest from -25 to +75 degrees
- Anti-trap gas strut fitted with a release lever used on the back & footrest sections (optional electric backrest)
- 125mm swivel castors featuring a central locking & steering facility
- Clearance beneath frame allows access for mobile hoists
- 230V ~, 50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- Seamless upholstery using flame retardant antimicrobial vinvl
- Swing down side support rails fitted to either side
- Frames prepared to accept optional accessories
- Maximum capacity 325kg

Dimensions (w x d x h)

850 x 1790 x 520-950mm

Standard Colours

- Paint White
- Vinyl Bristol Blue

Options

- E04/U Side Support Rails C/W Cushions Upgrade (Pair)
- E18 Battery Back Up
- E21 Breather Hole & Plug

| ltem | Description |
|------------|--|
| EC2425/VBB | Electric, Hand Switch |
| EC2430/VBB | Electric, Foot Switch |
| EC2435/VBB | Electric with Electric Backrest - Hand |
| EC2440/VBB | Electric with Electric Backrest - Foot |
| | |

Description Item EC2405/VBB Electric, Hand Switch EC2410/VBB Electric, Foot Switch EC2415/VBB Electric with Electric Backrest - Hand Switch Electric with Electric Backrest - Foot Switch

168

- Vinyl Bristol Blue

Weight

82.0kg

82.0kg

84.3ka

84.3kg

- E19 Paper Roll Holder

EC2420/VBB

- Maximum capacity 325kg

Service Contract ontact us fo details



Weight 88.1kg 88.1kg Switch 90.4ka Switch 90.4kg

sales@bristolmaid.com +44 (0)1258 484455

www.bristolmaid.com


169







Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Hydraulic in Sand Dune - EC3205/VSD



Examination & Treatment Couch **Bobath Therapy**

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric or hydraulic height adjustment
- Height adjustment from 520-950mm
- Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Anti-trap gas strut fitted with a release lever used on the backrest section (optional electric)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless upholstery using flame retardant anti-microbial vinyl
- Frames prepared to accept optional accessories
- Maximum capacity 225kg

Dimensions (w x d x h) Standard Colours

- 1000 x 1790 x 520-950mm
- · Paint White • Vinyl - Bristol Blue

Options

- E16 Side Support Rails (Pair)
- E17 Side Support Rails C/W Cushions (Pair)
- E18 Battery Back Up
- E19 Paper Roll Holder
- E22 Breather Hole & Plug

| Item | Description | Weight |
|------------|---|--------|
| EC3205/VBB | Hydraulic | 95.2kg |
| EC3210/VBB | Electric - Hand Switch | 95.2kg |
| EC3215/VBB | Electric - Foot Switch | 95.2kg |
| EC3220/VBB | Electric with Electric Backrest - Hand Switch | 97.5kg |
| EC3225/VBB | Electric with Electric Backrest - Foot Switch | 97.5kg |

Examination & Treatment Couch **Bobath Therapy, Mobile**

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric or hydraulic height adjustment
- · Height adjustment from 520-950mm
- Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Anti-trap gas strut fitted with a release lever used on the backrest section (optional electric)
- 125mm swivel castors featuring a central locking & steering facility
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless upholstery using flame retardant anti-microbial vinvl
- Swing down side support rails fitted to either side
- Frames prepared to accept optional accessories
- Maximum capacity 225kg

Dimensions (w x d x h)

• 1000 x 1790 x 520-950mm

Standard Colours

- Paint White
- Vinvl Bristol Blue

Options

- E17/U Side Support Rails C/W Cushions Upgrade (Pair)
- E18 Battery Back Up
- E19 Paper Roll Holder
- E22 Breather Hole & Plug

| Item | Description |
|------------|---|
| EC3230/VBB | Hydraulic |
| EC3235/VBB | Electric - Hand Switch |
| EC3240/VBB | Electric - Foot Switch |
| EC3245/VBB | Electric with Electric Backrest - Hand Switch |
| EC3250/VBB | Electric with Electric Backrest - Foot Switch |

sales@bristolmaid.com +44 (0)1258 484455



Weight

103.2kg

103.2kg

103.2kg

105.5ka

105.5kg



1750mm Hand Switch



Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Electric - Foot Switch in Chocolate - EC3115/VCH



Electric - Hand Switch

Electric - Foot Switch

Description

Hydraulic

Weight

96.9kg

96.9kg

96.9kg

Examination & Treatment Couch Mat Therapy

Features

- BS 2838-2: 1991 Examination & general treatment couches. build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric or hydraulic height adjustment
- Height adjustment from 520-950mm
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless upholstery using flame retardant anti-microbial vinyl
- Frames prepared to accept optional accessories
- Maximum capacity 225kg

Dimensions (w x d x h)

1000 x 1790 x 520-950mm

Standard Colours

- Paint White
- Vinyl Bristol Blue

Options

ltem

EC3105/VBB

EC3110/VBB

EC3115/VBB

- E16 Side Support Rails (Pair)
- E17 Side Support Rails C/W Cushions (Pair)
- E18 Battery Back Up
- E19 Paper Roll Holder

Examination & Treatment Couch Mat Therapy, Mobile

Features

- BS 2838-2: 1991 Examination & general treatment couches, build. performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric or hydraulic height adjustment
- Height adjustment from 520-950mm
- 125mm swivel castors featuring a central locking & steering facility
- Clearance beneath frame allows access for mobile hoists Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug &
- available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless upholstery using flame retardant anti-microbial vinyl
- Swing down side support rails fitted to either side
- Frames prepared to accept optional accessories
- Maximum capacity 225kg

Dimensions (w x d x h)

• 1000 x 1790 x 520-950mm

Standard Colours

- Paint White
- Vinvl Bristol Blue

Options

- E17/U Side Support Rails C/W Cushions Upgrade (Pair)
- E18 Battery Back Up
- E19 Paper Roll Holder

| Item | Description |
|------------|------------------------|
| EC3120/VBB | Hydraulic |
| EC3125/VBB | Electric - Hand Switch |
| EC3130/VBB | Electric - Foot Switch |
| | |


172

Foot Switch

Service Contract ontact us fo details



Weigh 104.9kg 104.9kg 104.9kg



sales@bristolmaid.com +44 (0)1258 484455



Examination & Treatment Chair **Drop End, Three Section**

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric or hydraulic height adjustment
- Height adjustment from 490-950mm
- · Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- · Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless upholstery using flame retardant anti-microbial vinyl
- Frames prepared to accept optional accessories
- Maximum capacity 225kg

Dimensions (w x d x h)

- Options
 - E02 Paper Roll Holder
 - E05 Breather Hole & Plug
 - E06 Battery Back Up
 - E07 Height Adjustable Armrests (Pair)
 - E08 Infusion Pole
 - E09 Bowl
 - E10 Bowl Support
 - E13 Phlebotomy Armrests (Pair) • E14 - Side Support Loops (Pair)

Weight

66.5kg

66.5kg

66.5kg

68.8kg

68.8kg

| Item | Description |
|------------|---|
| EC4305/VBB | Hydraulic |
| EC4310/VBB | Electric - Hand Switch |
| EC4315/VBB | Electric - Foot Switch |
| EC4320/VBB | Electric with Electric Backrest - Hand Switch |
| EC4325/VBB | Electric with Electric Backrest - Foot Switch |

Examination & Treatment Chair Drop End, Four Section

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric or hydraulic height adjustment
- Height adjustment from 490-950mm
- Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a
- hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless upholstery using flame retardant anti-microbial vinyl Frames prepared to accept optional accessories
- Maximum capacity 225kg

Dimensions (w x d x h)

- 630 x 1710 x 490-950mm · E02 - Paper Roll Holder
 - E06 Battery Back Up

Options

- E08 Infusion Pole
- E09 Bowl
- E10 Bowl Support
- E14 Side Support Loops (Pair)

| Item | Description |
|------------|--|
| EC4405/VBB | Hydraulic |
| EC4410/VBB | Electric - Hand Switch |
| EC4415/VBB | Electric - Foot Switch |
| EC4420/VBB | Electric with Electric Backrest - Hand |
| EC4425/VBB | Electric with Electric Backrest - Foot S |
| | |

- Vinvl Bristol Blue
- · Paint White
- 630 x 1710 x 490-950mm
- Standard Colours

- E13 Phlebotomy Armrests (Pair)

Standard Colours

Paint - White

Vinyl - Bristol Blue

 E05 - Breather Hole & Plug • E07 - Height Adjustable Armrests (Pair)











Examination & Treatment Chair **Drop End, Three Section, Bariatric**

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric height adjustment
- Height adjustment from 520-950mm
- Adjustable angle backrest Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas strut fitted with a release lever used on the back & footrest sections
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- 230V ~, 50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- · Seamless upholstery using flame retardant antimicrobial vinyl
- · Frames prepared to accept optional accessories
- Maximum capacity 325kg

Dimensions (w x d x h)

- **Standard Colours**
- 850 x 1710 x 520-950mm
- Paint White Vinvl - Bristol Blue

Options

- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole
- E09 Bowl
- E10 Bowl Support
- E13 Phlebotomy Armrests (Pair)
- E14 Side Support Loops (Pair)
- E18 Battery Back Up
- E19 Paper Roll Holder
- E21 Breather Hole & Plug

| Item | Description | Weight |
|------------|---|--------|
| EC4330/VBB | Electric - Hand Switch | 76.7kg |
| EC4335/VBB | Electric - Foot Switch | 76.7kg |
| EC4340/VBB | Electric with Electric Backrest - Hand Switch | 79.0kg |
| EC4345/VBB | Electric with Electric Backrest - Foot Switch | 79.0kg |

Examination & Treatment Chair Drop End, Four Section, Bariatric

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric height adjustment
- Height adjustment from 520-950mm
- Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas strut fitted with a release lever used on the back & footrest sections
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- 230V ~, 50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- Seamless upholstery using flame retardant antimicrobial vinvl
- Frames prepared to accept optional accessories
- Maximum capacity 325kg

Dimensions (w x d x h)

- 850 x 1710 x 520-950mm
- Paint White Vinvl - Bristol Blue

Options

- E07 Height Adjustable Armrests (Pair)
- E08 Infusion Pole
- E09 Bowl
- E10 Bowl Support
- E13 Phlebotomy Armrests (Pair)
- E14 Side Support Loops (Pair)
- E18 Battery Back Up
- · E19 Paper Roll Holder
- E21 Breather Hole & Plug

| Description |
|---|
| Electric - Hand Switch |
| Electric - Foot Switch |
| Electric with Electric Backrest - Hand Switch |
| Electric with Electric Backrest - Foot Switch |
| |

360mm

www.bristolmaid.com

176

- Service Contract Contact us fo details

Standard Colours



sales@bristolmaid.com

Weight

77.6kg

77.6kg

79.9kg

79.9kg

+44 (0)1258 484455





Phlebotomy Chair Variable Height

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric or hydraulic height adjustment
- Height adjustment from 490-950mm
- Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless upholstery using flame retardant anti-microbial vinyl
- Adjustable polyurethane phlebotomy armrests
- Side support loops assist with patient mount & dismount (self-assembly)
- Paper roll holder (suitable for rolls max 510mm wide x 200mm diameter)
- Frames prepared to accept optional accessories
- Maximum capacity 225kg

Dimensions (w x d x h) Options

• 770 x 1710 x 490-950mm

Standard Colours

Paint - White

Vinyl - Bristol Blue

| Item | Description | Weight |
|------------|---|--------|
| EC5305/VBB | Hydraulic | 67.0kg |
| EC5310/VBB | Electric - Hand Switch | 67.0kg |
| EC5315/VBB | Electric - Foot Switch | 67.0kg |
| EC5320/VBB | Electric with Electric Backrest - Hand Switch | 69.3kg |
| EC5325/VBB | Electric with Electric Backrest - Foot Switch | 69.3kg |

Phlebotomy Chair Variable Height, Bariatric

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric height adjustment
- Height adjustment from 490-950mm
- Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- 230V ~, 50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- Seamless upholstery using flame retardant anti-microbial vinyl
- Adjustable polyurethane phlebotomy armrests
- Side support loops assist with patient mount & dismount (self-assembly)
- Paper roll holder (suitable for rolls max 510mm wide x 200mm diameter)
- Frames prepared to accept optional accessories
- Maximum capacity 325kg

Dimensions (w x d x h)

• 990 x 1710 x 520-950mm

Standard Colours

- Paint White
- Vinyl Bristol Blue

Options

· E18 - Battery Back Up

| Item | Description |
|------------|---|
| EC5330/VBB | Electric - Hand Switch |
| EC5335/VBB | Electric - Foot Switch |
| EC5340/VBB | Electric with Electric Backrest - Hand Switch |
| EC5345/VBB | Electric with Electric Backrest - Foot Switch |

178

sales@bristolmaid.com +44 (0)1258 484455

- · E06 Battery Back Up









Phlebotomy Chair **Fixed Height**

Features

- BS EN 16139: 2013 Level 2. Furniture. Strength, durability & safety. Requirements for non-domestic seating
- BS 4875: Part 1: 2001 Test Level 4 (strength & durability)
- BS 7945: 1999 (stability)
- Certified for use up to 32 stone (203kg)
- Stable four leg design
- Black finish tubular frame
- Seamless upholstery using flame retardant anti-microbial vinyl
- Adjustable polyurethane phlebotomy armrest
- Armrest can be positioned on either side of the chair (additional armrest available if required)
- Maximum capacity 203kg

Dimensions (w x d x h)

• 620 x 580 x 800mm

Standard Colours

- Paint Black
- Vinyl Bristol Blue

Options

• E15 - Phlebotomy Armrest (Single)







E07

Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Fixed Height Phlebotomy Chair in Aquamarine - EC5000/VAQ



| Item | Description | Seat Height | Weight |
|------------|--------------|-------------|--------|
| EC5000/VBB | Fixed Height | 485mm | 15.0kg |

sales@bristolmaid.com +44 (0)1258 484455 www.bristolmaid.com

Paper Roll Holder E02 - E19

 Fitted to the head section Suitable for rolls max 510mm wide x 200mm diameter



Side Support Rails (Pair)

E03 - EC905 - E16

- Swing down side support rails
- Can be added to either one or both sides of the couch to provide additional safety for patients



E05 - E21 - E22 - EC920

· Breather hole located in backrest section to improve patient comfort when in the prone position Supplied with a plug/insert



Cushioned Side Support Rails (Pair)

E04 - E04/U - E17 - E17/U

- · Vinyl upholstered swing down side support rails
- Can be added to either one or both sides of the couch to provide additional safety for patients
- Upgrade option for factory fitted side support rails (E04/U & E17/U)

Height Adjustable Armrests (Pair)

 Height adjustable easy clean moulded polyurethane armrests



Battery Back Up

E06 - E18

- Allows the couch to be lowered in a controlled manner in the event of a power failure
- Provides enough charge for 1-3 cycles



Infusion Pole

E08

- Height adjustable two hook infusion pole
 Cannot be fitted in conjunction with side support rails E03, E04, E04/U, E16, E17, E17/U



Bowl

E09

 360mm diameter 304 Grade stainless steel bowl



E14



E15



E12

 Fully adjustable easy clean moulded polyurethane leg supports



Bowl Support

E10

- Support to hold 360mm diameter 304 Grade stainless steel bowl
- Order bowl separately



Lithotomy Stirrups (Pair) E11

Traditional lithotomy stirrups



Phlebotomy Armrests (Pair)

E13

 Fully adjustable easy clean moulded polyurethane armrests

Side Support Loops (Pair)

 Side support loops to assist with mount/ dismount & provide patient security



Drawer EC900

- Mounts beneath the couch frame offering additional storage space • Compatible with EC005 & EC010

Phlebotomy Armrest (Single)

• Fully adjustable easy clean moulded Provide the set of the s



Drop Down Side Rails (Pair) EC935

- Drop down side support rails
- Can be added to either one or both sides of the couch to provide additional safety for patients
- Compatible with EC005 & EC010









General Purpose Trolleys



Stainless Steel Trolleys 186

Polymer Trolleys 188

Mild Steel Trolleys 189



General Purpose Trolley Stainless Steel, Heavy Duty

Features

- Seamlessly welded from 304 grade stainless steel for strength & ease
 of cleaning
- Fixed shelves with three sides up & one down (shelf size 920 x 535mm)
- Available with two or three shelves, shelf heights
- Two shelf 280 & 860mm
- Three shelf 220, 540 & 860mm
- Maximum capacity per shelf 75kg
- 125mm pressed steel castors

Dimensions (w x d x h)

- No Ring Buffers 970 x 540 x 1020mm
- Ring Buffers 1045 x 615 x 1020mm

Standard Colours

Stainless Steel

Options

• C008 - 125mm Castors, Heavy Duty Pressed Steel, 2 x Braking, 2 x Non-Braking (set of 4)

General Purpose Trolley Stainless Steel, Light Duty

Features

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed liquid retaining sunken shelves (shelf size 540 x 430mm)
- Available with two or three shelves, shelf heights
- Two shelf 250 & 870mm
- Three shelf 250, 560 & 870mm
- Maximum capacity per shelf 50kg
- 100mm pressed steel castors

Dimensions (w x d x h)

- No Ring Buffers 990 x 510 x 1035mm
- Ring Buffers 1065 x 585 x 1035mm

Standard Colours

Stainless Steel

| Item | Description | Weight |
|---------------|---|--------|
| WT/37/2 | Stainless Steel - Heavy Duty, Two Shelves | 28.0kg |
| WT/37/3 | Stainless Steel - Heavy Duty, Three Shelves | 37.5kg |
| WT/37/2/ZP010 | Stainless Steel - Heavy Duty, Two Shelves, Ring Buffers | 28.2kg |
| WT/37/3/ZP010 | Stainless Steel - Heavy Duty, Three Shelves, Ring Buffers | 37.7kg |

| Item | Description | Weight |
|------------------|---|--------|
| WT/33/SS/2 | Stainless Steel - Light Duty, Two Shelves | 16.0kg |
| WT/33/SS/3 | Stainless Steel - Light Duty, Three Shelves | 22.0kg |
| WT/33/SS/2/ZP005 | Stainless Steel - Light Duty, Two Shelves, Ring Buffers | 16.2kg |
| WT/33/SS/3/ZP005 | Stainless Steel - Light Duty, Three Shelves, Ring Buffers | 22.2kg |



General Purpose Trolley Polymer, Shallow Shelves

Features

- Easy to clean corrosion resistant high density polyethylene shelves (not suitable for use with automated washers)
- · 11mm raised edge retains spills & prevents items from sliding off during transport
- Microban® product protection inhibits growth of bacteria, mould, mildew & funai
- Chrome plated posts offers shelf adjustment in 25mm increments
- Available with two or three shelves in three sizes
- Small Shelf size 620 x 345mm
- Easy grip handle & utility tray
- Maximum capacity per shelf 60kg
- 100mm pressed steel castors
- Medium Shelf size 690 x 475mm
- · Handle has built in recessed support channel to manage odd shaped items
- Maximum capacity per shelf 60kg
- 100mm pressed steel castors
- Large Shelf size 840 x 585mm
- Handle has built in slot to hold scanners
- Maximum capacity per shelf 75kg
- 125mm pressed steel castors
- Self-assembly (supplied boxed)



Standard Colours

METRO

- · Uprights Chrome
- Shelves Blue

| Item | Shelves | Weight |
|-------------------|------------|--------|
| Small - 800 x 465 | x 900mm | |
| 8MY1627-24BU | 2 | 12.9kg |
| 8MY1627-34BU | 3 | 15.4kg |
| Medium - 873 x 5 | 95 x 900mr | n |
| 8MY2030-24BU | 2 | 14.9kg |
| 8MY2030-34BU | 3 | 18.8kg |
| Large - 1022 x 70 |)3 x 936mm | |
| 8MY2636-25BU | 2 | 18.1kg |
| 8MY2636-35BU | 3 | 23.1kg |
| | | |

General Purpose Trolley Polymer, Deep Shelves

Features

- Easy to clean corrosion resistant high density polyethylene shelves (not suitable for use with automated washers)
- 70mm deep shelves retains spills & prevents items from sliding off during transport
- Microban® product protection inhibits growth of bacteria, mould, mildew & funai
- · Chrome plated posts offers shelf adjustment in 25mm increments Available with two or three shelves in three sizes
- Small 800 x 546mm
- Large 948 x 686mm
- Max capacity shelf 68kg
- 100mm pressed steel castors
- Self-assembly (supplied boxed)

Standard Colours

- Uprights Chrome
- Shelves Blue

| Item | Shelves | Weight |
|-------------------|----------|--------|
| Small - 832 x 546 | x 1041mm | |
| 8BC2030-2DMB | 2 | 13.0kg |
| 8BC2030-3DMB | 3 | 17.0kg |
| Large - 984 x 685 | x 1041mm | |
| 8BC2636-2DMB | 2 | 15.9kg |
| 8BC2636-3DMB | 3 | 21.1kg |
| | | |



General Purpose Trolley Mild Steel, Laminate Shelves

Features

- Laminate shelves
- Available with two or three shelves in two sizes
- Shelf heights
- Two shelf 210 & 845mm
- Three shelf 210, 525 & 845mm
- Shelf size
- Small 705 x 445mm
- Large 855 x 445mm
- Maximum capacity per shelf 40kg
- 100mm pressed steel castors



| ltem | Description |
|-------------------|-------------------------------------|
| Small - 972 x 533 | 3 x 1008mm |
| GP005 | Mild Steel - Two Laminate Shelves, |
| GP015 | Mild Steel - Three Laminate Shelves |
| Large - 1122 x 53 | 33 x 1008mm |
| GP025 | Mild Steel - Two Laminate Shelves, |
| GP035 | Mild Steel - Three Laminate Shelves |
| | |

General Purpose Trolley Mild Steel, Aluminium Trays

Features

- Liquid retaining painted aluminium travs (depth 50mm)
- Available with two or three shelves in two sizes
- Shelf heights
- Two shelf 210 & 845mm
- Three shelf 210, 525 & 845mm
- Shelf size
- Small 705 x 445mm
- Large 855 x 445mm
- Maximum capacity per tray 40kg
- 100mm pressed steel castors

Standard Colours

Paint - Grev

Options

ZP005 - Ring Buffer

| arge | | |
|------|--|--|

Standard Colours

Paint - Grev

Options

Shelf - Spearmint

ZP005 - Ring Buffer



| ltem | Description | Weight |
|-----------------|---|--------|
| Small - 972 x 5 | 33 x 1008mm | |
| GP010 | Mild Steel - Two Aluminium Trays, Small | 16.5kg |
| GP020 | Mild Steel - Three Aluminium Trays, Small | 20.6kg |
| Large - 1122 x | 533 x 1008mm | |
| GP030 | Mild Steel - Two Aluminium Trays, Large | 18.5kg |
| GP040 | Mild Steel - Three Aluminium Trays, Large | 23.7kg |
| | | |













HTM71 Storage Solutions

- ScanCell[®] Carts 192
- Distribution Trolleys 197
- Open Transfer/Transport Trolleys 199
 - HTM71 Racking 200
- Trays, Baskets, Dividers & Accessories 202



Type A - 600mm Insertion





ScanCell[®] Carts - Open Fronts



Features

- Mobile open front storage carts for use with 600 x 400mm modular travs & baskets Lightweight easyclean design
- · Injection moulded ABS top/chassis with aluminium corner profiles, laminate sides & back

See page 202 for

Trays Baskets

Accessories

- · Available as single, double and treble columns in a range of heights
- Fitted with angled/horizontal side panels to suit 400 or 600mm insertions
- Type A 600mm Insertions
- Type B 400mm Insertions
- 125mm castors, 2 off braking

Standard Colours

· Body - Grey

Options

- 4UR1A Utility Rail, Type A (Cart Heights 905mm, 1055mm, 1205mm)
- 4UR1B Utility Rail, Type B (Cart Heights 905mm, 1055mm, 1205mm)
- 4PR1A Push Rail, Type A (Cart Heights 1503mm, 1803mm)
- 4PR1B Push Rail, Type B (Cart Heights 1503mm, 1803mm)



ScanCell[®] Carts - Roller Shutters

Features

- Mobile storage carts fitted with roller shutters for use with 600 x 400mm modular travs & baskets
- Lightweight easyclean design
- Injection moulded ABS top/chassis with aluminium corner profiles, laminate sides & back, polypropylene roller shutter
- Available as single, double and treble columns in a range of heights
- Fitted with angled/horizontal side panels to suit 400 or 600mm insertions
- Type A 600mm Insertions
- Type B 400mm Insertions
- 125mm castors, 2 off braking

Standard Colours

See page 202 for **Travs Baskets** Accessories

· Body - Grey Options

- 4UR1A Utility Rail, Type A (Cart Heights 905mm, 1055mm, 1205mm)
- 4UR1B Utility Rail, Type B (Cart Heights 905mm, 1055mm, 1205mm)
- 4PR1A Push Rail, Type A (Cart Heights 1503mm, 1803mm)
- 4PR1B Push Rail, Type B (Cart Heights 1503mm, 1803mm)





(trays & baskets not included)











Type A - 600mm Insertion



ScanCell[®] Carts - Doors

Features

 Mobile storage carts fitted with doors for use with 600 x 400mm modular travs & baskets

- Lightweight easy clean design
- Injection moulded ABS top/chassis with aluminium corner profiles, laminate sides & back, post form/plastic finish doors
- · Available as single, double and treble columns in a range of heights
- Fitted with angled/horizontal side panels to suit 400 or 600mm insertions
- Type A 600mm Insertion Single door for each column Single Column - R/H Hinged
- Double Column R/H & L/H Hinged
- Treble Column 1 x R/H & 2 x L/H Hinged
- Type B 400mm Insertion Double doors for each column
- 125mm castors, 2 off braking
- For travs/baskets/stops & accessories please refer to page 202

Standard Colours

Body - Grey

Options

- 4UR1A Utility Rail, Type A (Cart Heights 905mm, 1055mm, 1205mm)
- 4UR1B Utility Rail, Type B (Cart Heights 905mm, 1055mm, 1205mm)
- 4PR1A Push Rail, Type A (Cart Heights 1503mm, 1803mm)
- 4PR1B Push Rail, Type B (Cart Heights 1503mm, 1803mm)



ScanCell[®] Carts - Vario Drawers

Features

- Mobile storage carts, lightweight easyclean design, fitted with vario drawers and supplied with travs
- Injection moulded ABS top/chassis with aluminium corner profiles, laminate sides & back
- Type B 400mm Insertions
- Available as single, double and treble columns in a range of heights
- 50mm/100mm drawers can be subdivided (see page 203 for dividers & label holders)
- 125mm castors, 2 off braking

Dimensions (w x d x h)

- Sinale Column 662 x 494 x Height
- Double Column 1324 x 494 x Height
- Treble Column 1986 x 494 x Height

Options

• 4UR1B - Utility Rail (Cart Height - 905mm, 1055mm, 1205mm) 4PR1B - Push Rail (Cart Height - 1503mm)


| Item | 50mm Tray | 100mm Tray | 20 | | |
|--------------------------------|-----------|------------|----|--|--|
| Height - 905mm | | | | | |
| 4B1VD09/A | 1 | 2 | | | |
| 4B1VD09/B | | 4 | | | |
| 4B1VD09/C | 2 | 3 | | | |
| 4B2VD09/AA* | 2 | 4 | | | |
| 4B3VD09/AAA** | 3 | 6 | | | |
| Height - 1055mm | | | | | |
| 4B1VD10/D | 1 | 3 | | | |
| 4B1VD10/E | | 5 | | | |
| 4B1VD10/F | 2 | 4 | | | |
| 4B2VD10/DD* | 2 | 6 | | | |
| 4B3VD10/DDD** | 3 | 9 | | | |
| Height - 1205mm | | | | | |
| 4B1VD12/G | 1 | 4 | | | |
| 4B1VD12/H | 2 | 2 | | | |
| Height - 1503mm | | | | | |
| 4B1VD15/I | 1 | 6 | | | |
| 4B1VD15/J | 1 | 4 | | | |
| *double column **treble column | | | | | |


HTM71



Standard Colours

 Body/Drawers - Grey Drawer Handle - Blue











ScanCell[®] Carts - Telescopic Drawers

Features

- Mobile storage carts with telescopic drawers supplied with 600 x 400mm modular travs & baskets
- Lightweight easy clean design
- Injection moulded ABS top/chassis with aluminium corner profiles, laminate sides & back
- Type B 400mm Insertions
- Available as single, double and treble columns in a range of heights
- Fitted with horizontal side panels to suit 400 insertions
- Drawers can be subdivided (see page 203 for dividers & label holders)
- 125mm castors, 2 off braking

Dimensions (w x d x h)

- Single Column 662 x 494 x Height
- Double Column 1324 x 494 x Height

Standard Colours

· Body/Drawers - Grev

Options

• 4UR1B - Utility Rail (Cart Height - 905mm, 1055mm, 1205mm)

| Item | 50mm Tray | 100mm Tray | 200mm Basket |
|-----------------|-----------|------------|--------------|
| Height - 905mm | | | |
| 4B1TD09/A | 1 | 2 | 1 |
| 4B1TD09/B | | 4 | |
| 4B1TD09/C | 2 | 3 | |
| 4B2TD09/AA * | 2 | 4 | 2 |
| Height - 1055mm | | | |
| 4B1TD10/D | 1 | 3 | 1 |
| 4B1TD10/E | | 5 | |
| 4B1TD10/F | 2 | 4 | |
| 4B1TD10/G | 10 | | |
| Height - 1205mm | | | |
| 4B1TD12/H | 1 | 4 | 1 |
| 4B1TD12/I | 2 | 5 | |
| 4B1TD12/J | 6 | 3 | |
| double column | | | |

ScanCell[®] Carts - Compact

Features

- Single column mobile compact carts for use with 300 x 400mm modular travs & baskets
- Lightweight easyclean design
- Injection moulded ABS top/chassis with aluminium corner profiles, laminate sides & back
- Available with drawers with telescopic runners
- Type D 400mm Insertions
- 125mm castors, 2 off braking
- Travs/baskets/drawers can be subdivided

Dimensions (w x d x h) Options

- 462 x 494 x 905mm
- 4UR1B Utility Rail (Type B)

Standard Colours

· Body/Drawers - Grey

This is a small election of a large

configurable range Please contact our

customer services team to discuss your requirements



| Item | 50mm Tray | 100mm Tra |
|-----------|-----------|-----------|
| 4D1TD09/A | 1 | 2 |
| 4D1TD09/B | | 4 |



Distribution Trolleys - Compact

Features

- Single column compact distribution trollevs for use with 300 x 400mm modular travs. baskets & stops
- Available in mild/stainless steel
- Mild Steel Welded & powder coated
- · 304 Grade Stainless Steel Seamlessly welded for ease of cleaning
- Travs & baskets can be supported on the upper & lower frame (optional solid grade laminate worktop)
- 4 stops required per tray / basket
- Runners allow trays to be accessed from either side (4 stops required per tray/basket) Type D - 400mm insertion
- 75mm castors
- For travs/baskets/stops & accessories please refer to page 202

Dimensions (w x d x h)

• 360 x 460 x 850mm

Standard Colours

Paint - Grev

Options

- ZP010 Ring Buffers
- 4D1SLWT- Solid Grade Laminate Worktop





1

4D1DSTCM 4D1DSTCS

Description

Compact, Single Column, Mild Steel Compact, Single Column, Stainless Steel



Distribution Trolleys

Features

- Single & double column distribution trolleys for use with 600 x 400mm modular travs & baskets Available in mild/stainless steel
- Mild Steel Welded & powder coated
- 304 Grade Stainless Steel Seamlessly welded for ease of cleaning
- Travs & baskets can be supported on the upper & lower frame (optional solid grade laminate worktop)
- Runners allow trays to be accessed from either side (4 stops required per tray/basket)
- Type A 600mm insertion
- Type B 400mm insertion
- 75mm castors
- · For trays/baskets/stops & accessories please refer to page 202

Standard Colours

· Paint - Grey

Options

- ZP010 Ring Buffers
- 4A1SLWT Solid Grade Laminate Worktop, Single Column (Type A)
- 4A2SLWT Solid Grade Laminate Worktop, Double Column (Type A)
- 4B1SLWT Solid Grade Laminate Worktop, Single Column (Type B)

TYPE A - 600mm Insertion

TYPE B - 400mm Insertion



| Item | Description | Overall Dimensions (w x d x h) |
|-----------------|--|--------------------------------|
| Mild Steel | | |
| 4A1DSTM | Single Column, Type A, Mild Steel | 460 x 655 x 850mm |
| 4A2DSTM | Double Column, Type A, Mild Steel | 895 x 655 x 850mm |
| 4B1DSTM | Single Column, Type B, Mild Steel | 660 x 455 x 850mm |
| Stainless Steel | | |
| 4A1DSTS | Single Column, Type A, Stainless Steel | 460 x 655 x 850mm |
| 4A2DSTS | Double Column, Type A, Stainless Steel | 895 x 655 x 850mm |
| 4B1DSTS | Single Column, Type B, Stainless Steel | 660 x 455 x 850mm |

Open Transfer / Transport Trolleys

Features

- Single, double & treble column transfer/transport trollevs for use with 600 x 400mm modular travs & baskets
- Available in mild/stainless steel
- · Mild Steel Welded & powder coatedx · 304 Grade Stainless Steel - Seamlessly welded for ease of
- cleaning
- Runners allow travs to be accessed from either side (4 stops required per tray/basket)
- Type A 600mm insertion
- Push handles at either end
- · Ring buffers to protect the trolley/fabric of the building
- Castors
- Single/Double Column 125mm castors, 2 off braking
- Treble Column 200mm heavy duty cushioned castors, 2 off fixed. 2 off swivel braking
- For travs/baskets/stops and accessories please refer to page 202

Standard Colours

Paint - Grev

Options

- 4A1OTTNC Nylon Cover (Single Column)
- 4A2OTTNC Nylon Cover (Double Column)
- 4A3OTTNC Nylon Cover (Treble Column)

| ltem | Description |
|-----------------|--------------------------------------|
| Mild Steel | |
| 4A1OTTM | Single Column, Type A, Mild Steel |
| 4A2OTTM | Double Column, Type A, Mild Steel |
| 4A3OTTM | Treble Column, Type A, Mild Steel |
| Stainless Steel | |
| 4A10TTS | Single Column, Type A, Stainless Ste |
| 4A2OTTS | Double Column, Type A, Stainless St |
| 4A3OTTS | Treble Column, Type A, Stainless Ste |
| | |

198



Overall Dimensions (w x d x h)

665 x 725 x 1555mm 1100 x 725 x 1555mm 1535 x 725 x 1635mm

teel Steel

665 x 725 x 1555mm 1100 x 725 x 1555mm 1535 x 725 x 1635mm This is a small selection of a large configurable range. Please contact our customer services team to discuss your requirements





Room Design & Installation Service Available



2HH-7532293 Wall Mounting Bracket



2HH-7532277 10mm Plastic Cover Cap

4HH-HTM71/S1/400 Runners, trays, baskets, dividers & angle brackets not included (order separately)

Components

- 2HH-7532301-HH Base Unit, Type A, 600mm Insertion (x 1)
- 2HH-7532311 Extension Unit, Type A, 600mm Insertion (x 1)
- 2HH-7532356 Runners, Type A, 600mm Insertion (pair)
- 2HH-7532329 Angle Adaptor Plate, Type A, 600mm Insertion (pair)
- 2HH-7532304-HH Base Unit, Type B, 400mm Insertion (x 1)
- 2HH-7532314 Extension Unit, Type B, 400mm Insertion (x 1)
- 2HH-7532355 Runners, Type B, 400mm Insertion (pair)
- 2HH-7532328 Angle Adaptor Plate, Type B, 400mm Insertion (pair)
- 2HH-7532320 Adjustable Foot Set (x 4)
- 2HH-7532321 Adjustable Foot Set (x 2)
- 2HH-7532322 50mm Pressed Steel Castors (2 x Braking / 2 x Non-Braking)
- 2HH-7532323 50mm Pressed Steel Castors (1 x Braking / 1 x Non-Braking)
- 2HH-7532324 80mm Pressed Steel Castors (2 x Braking / 2 x Non-Braking)
- 2HH-7532325 80mm Pressed Steel Castors (1 x Braking / 1 x Non-Braking)

HTM71 Racking

Features

- Fully configurable freestanding modular racking for use with 600 x 400mm modular trays & baskets
 Available to suit 600 or 400mm insertions
- Type A 600mm insertion
- Type B 400mm insertion
- · Trays & baskets supported on adjustable moulded plastic/glass fibre runners
- Runners can be positioned horizontal or angled (approx. 25 degrees requires angle adaptor brackets)
- · Integral stops on the runners allow trays/baskets to be accessible in two positions
- Straight (approx. 5 degrees)
- Angled (approx. 25 degrees)
- Runners offer pass through facility allowing access on either side
- Anodised extruded aluminium uprights & cross members
- Bolted construction using high quality moulded fixings provides rigidity & allows the system to be altered easily
- Available as either static or mobile units
- Static Adjustable feet allow levelling on uneven floors
- Mobile 50/80mm pressed steel castors
- Optional wall mounting brackets allow static units to be fixed securely to the wall, (when units are loaded over 100kg

Mobile

Type A - 600mm Insertion

Type B - 400mm Insertion

Options

Single Bay - 494 x 642 x 1996mm

• Double Bay - 937 x 642 x 1996mm

Single Bay - 694 x 442 x 1996mm

• Double Bay - 1337 x 442 x 1996mm

2HH-7532293 - Wall Mounting Bracket

• 2HH-7532277 - 10mm Plastic Cover Cap

Supplied for self-assembly

Dimensions (w x d x h)

Static

- Type A 600mm Insertion
- Single Bay 494 x 642 x 1940mm
 Double Bay 937 x 642 x 1940mm
- Type B 400mm Insertion
- Single Bay 694 x 442 x 1940mm
 Double Bay 1337 x 442 x 1940mm

Standard Colours

- Uprights & Cross Members Anodised Aluminium
- Plastic Mouldings Green/Grey

Static - Adjustable Feet

| Item | Description | Weight |
|--------------------------|-------------------------------------|--------|
| Type A - 600mm Insertion | | |
| 4HH-HTM71/S1/600 | Static, Single Bay, 600mm Insertion | 10.3kg |
| 4HH-HTM71/S2/600 | Static, Double Bay, 600mm Insertion | 16.1kg |
| Type B - 400mm Insertion | | |
| 4HH-HTM71/S1/400 | Static, Single Bay, 400mm Insertion | 10.2kg |
| 4HH-HTM71/S2/400 | Static, Double Bay, 400mm Insertion | 16.3kg |



Mobile - 50mm Castors

| Item | Description | Weight |
|--------------------------|-------------------------------------|--------|
| Type A - 600mm Insertion | | |
| 4HH-HTM71/M51/600 | Mobile, Single Bay, 600mm Insertion | 11.1kg |
| 4HH-HTM71/M52/600 | Mobile, Double Bay, 600mm Insertion | 17.3kg |
| Type B - 400mm Insertion | | |
| 4HH-HTM71/M51/400 | Mobile, Single Bay, 400mm Insertion | 11.0kg |
| 4HH-HTM71/M52/400 | Mobile, Double Bay, 400mm Insertion | 17.5kg |





Trays & Baskets

| Type D - Size 300 x 400mm | | Type A & B - Size 600 x 400m | m | Type D - Size 300 x 400mm | | | | |
|---------------------------|--|------------------------------|--|---|----------------------------|--------------------------|--|----------------------------|
| Item | Description | Item | Description | Item Description | Tray/Basket Compatibility | Item | Description | Fray/Basket Compatibility |
| 2HH-7014050 | Tray - ABS, 50mm Dimensions (w x d x h): 300 x 400 x 50mm Material: ABS Grey | 2HH-7011050 | Tray - ABS, 50mm Dimensions (w x d x h): 600 x 400 x 50mm Material: ABS Grey | Divider - PC Dimensions (w x h): 300 x 50mm Material: Polycarbonate, Transparent 2HH-7810103 | 2HH-7014050 | | Divider - ABS Dimensions (w x h): 400 x 200mm Material: ABS Grey | 2HH-7014200 2HH-7012200 |
| | Basket - ABS, 100mm Dimensions (w x d x h): 300 x 400 x 100mm Material: ABS Grey | | Tray - ABS, 100mm Dimensions (w x d x h): 600 x 400 x 100mm Material: ABS Grey | Divider - PC Dimensions (w x h): 400 x 50mm Material: Polycarbonate, Transparent 2HH-7810102 | 2HH-7014050 2HH-7011050 | 2HH-7811305 | Sliding Divider - SAN Dimensions (w x h): 170 x 45mm Material: Plastic, Transparent | 2HH-7014050 2HH-7011050 |
| 2HH-7014100 | | 2HH-7011100 | | Samuel | | 2HH-241302 | | |
| | Basket - ABS, 200mm Dimensions (w x d x h): 300 x 400 x 200mm Material: ABS Grey | 2НН-7012100 | Basket - ABS, 100mm Dimensions (w x d x h): 600 x 400 x 100mm Material: ABS Grey | Divider - ABS Dimensions (w x h): 300 x 100mm Material: ABS Grey 2HH-7811206N | 2HH-7014100 | | Sliding Divider - SAN (including cardholder) Dimensions (w x h): 170 x 45mm Material: Plastic, Transparent | 2HH-7014050 2HH-7011050 |
| 2HH-7014200 | | | Basket - ABS, 200mm Dimensions (w x d x h): 600 x 400 x 200mm Material: ABS Grey | Divider - ABS Dimensions (w x h): 400 x 100mm Material: ABS Grey 2HH-7811205N | 2HH-7014100 2HH-7012100 | 2HH-241301 | Sliding Divider - SAN Dimensions (w x h): 170 x 100mm Material: Plastic, Transparent | 2HH-7014100 |
| | | 2НН-7012200 | Tray/Lid - ABS, 15mm Dimensions (w x d x h): 600 x 400 x 15mm Material: Plastic Grey | Divider - ABS Dimensions (w x h): 300 x 200mm Material: ABS Grey 2HH-7811306 | 2HH-7014200 | 2HH-240820 2HH-240821 | Sliding Divider - SAN (including cardholder) Dimensions (w x h): 170 x 100mm Material: Plastic, Transparent | 2HH-7014100 |

Dividers

HTM71 Storage Soluti

Dividers

Accessories

Type A & B - Size - 600 x 400mm



MTH

Tray/Basket Compatibilit

All Trays & Baskets

(Scan Modul Compatible)

All Trays & Baskets

(H+H System Compatible)

All Trays & Baskets

Credit Card Size: (85 x 54mm)

All Trays & Baskets

Credit Card Size: (85 x 54mm) Material: Plastic, White

All Trays & Baskets



Infusion Stands & Accessories

- Mobile Pump Stands 208
- Mobile Infusion Stands 210
- Clamping Infusion Stands 215
- Wall Mounted Infusion Holders 215



Mobile Pump Stands - Heavy Duty, Easy Clean

Features

- Range of stands designed specifically to hold blood warmers, pumps & syringe drivers
- Fully welded heavy duty five leg base with low centre of gravity to assist stability
- 32mm diameter 304 grade stainless steel upright & 25mm diameter pump support bows
- Polyurethane overmoulded base with smooth edges & contours provides...
- Water-resistant finish, allowing for easy & effective cleaning (including steam cleaning)
- Shield for the twin-wheel castors
- Impact protection to patients & staff, reducing the risk of skin abrasions & bruising
 Protection to the fabric of the building
- Height adjustable upright using a handwheel (1300 2000mm)
- Available with four/six U shaped hooks
- Maximum recommended load...
- Hooks 12kg per stand plus
- Fixed Upright/s 22kg per stand
- 75mm polypropylene/thermoplastic rubber twin-wheel castors (1 x anti-static)

Dimensions (w x d x h)

• 620 x 620 x 1300-2000mm

Vertical Clamping Area

- EPS005 975mm
- EPS010 & EPS015 771mm (upright) 2 x 646mm (support bows)

Standard Colours

- Base Grey
- Upright & Support Bows Stainless Steel



Smooth edges & contours, allows easy & effective cleaning

| Item | Description | Weight |
|--------|--|--------|
| EPS005 | Mobile Pump Stand - Heavy Duty, Easy Clean, Standard Capacity, Four Hook | 11.0kg |
| EPS010 | Mobile Pump Stand - Heavy Duty, Easy Clean, High Capacity, Four Hook | 13.0kg |
| EPS015 | Mobile Pump Stand - Heavy Duty, Easy Clean, High Capacity, Six Hook | 13.1kg |





Four or six U shaped hooks



Height adjustable upright using a heavy duty handwheel





Mobile Infusion Stands - Polyurethane, Easy Clean

Features

- Meets appendix A (Stability Requirements) of BS 3619:1976 Mobile Infusion Stands
- Five leg base with low centre of gravity to assist stability
- 25mm diameter 304 grade stainless steel upright
- · Polyurethane overmoulded base with smooth edges & contours provides...
- Water-resistant finish, allowing for easy & effective cleaning (including steam cleaning)
- Shield for the twin-wheel castors
- Impact protection to patients & staff, reducing the risk of skin abrasions & bruising
- Protection to the fabric of the building
- Coloured indicator caps provide department identification & segregation, reducing the risk of cross contamination
- Height adjustable upright using a handwheel (1260 2100mm)
- Available with two/four pigtail hooks
- Maximum recommended load...
- Hooks 4kg per stand
- Fixed Upright 4kg per stand
- Optional handle keeps feet clear of the base, whilst also assisting manoeuvrability & improving patient safety
- 50mm polypropylene/thermoplastic rubber twin-wheel non-marking castors (1 x anti-static)

Dimensions (w x d x h) Vertical Clamping Area

• 620 x 620 x 1260 - 2100mm • IS/P/2H & IS/P/4H - 970mm • IS/P/2H/H & IS/P/4H/H - 840mm


Standard Colours

- Base Grey
- Upright & Handle Stainless Steel
- Coloured Indicator Cap Spearmint

Options

· Coloured Indicator Cap - Blue / Green / Red / Yellow







4 hook pole supplied with models with '4H' product code suffix





Mobile Infusion Stands - Stainless Steel

Features

Meets appendix A (Stability Requirements) of BS 3619:1976 Mobile Infusion Stands

- Seamless welded 304 grade stainless steel construction
- Height adjustable upright using a handwheel (1270 2000mm)
- Available with two or four pigtail hooks
- 9kg plastic coated weighted base available for additional stability
- Maximum recommended load...
- Hooks 4kg per stand plus
- Fixed Upright
 - Non-Weighted Base 2kg per stand
 - Weighted Base 4kg per stand
- 50mm polypropylene/thermoplastic rubber twin wheel castors (1 x anti-static)

Dimensions (w x d x h)

• 650 x 650 x 1270-2000mm

Standard Colours

 Base/Upright - Stainless Steel Weighted Base - Grey



4 hook pole supplied with models with '4H' product code suffix

| Item | Description | Weight |
|--------------|--------------------------|--------|
| DPWS/1 | Two Hook | 3.8kg |
| DPWS/1/4H | Four Hook | 3.9kg |
| DPWS/1/W9 | Two Hook, Weighted Base | 12.8kg |
| DPWS/1/W9/4H | Four Hook, Weighted Base | 12.9kg |

Mobile Infusion Stands - Mild Steel

Features

Meets appendix A (Stability Requirements) of BS 3619:1976 Mobile Infusion Stands

- Welded mild steel construction
- Height adjustable upright using a handwheel (1270 2000mm)
- Available with two or four pigtail hooks
- 9kg plastic coated weighted base available for additional stability
- Maximum recommended load...
- Hooks 4kg per stand plus
- Fixed Upright
 - Non-Weighted Base 2kg per stand
- Weighted Base 4kg per stand
- 50mm polypropylene/thermoplastic rubber twin wheel castors (1 x anti-static)

Dimensions (w x d x h)

• 650 x 650 x 1270-2000mm

Standard Colours

- Paint White
- Weighted Base Grey



4 hook pole supplied with models with '4H' product code suffix

| Item | Description |
|--------------|--------------------------|
| DPWM/1 | Two Hook |
| DPWM/1/4H | Four Hook |
| DPWM/1/W9 | Two Hook, Weighted Base |
| DPWM/1/W9/4H | Four Hook, Weighted Base |

0000 (ind HEIGHT ADJUSTABLE Weight DPWM/1/W9/4H 3.8kg 3.9kg 0 12.8kg 12.9kg DPWM/1

5

Mobile Infusion Stands - Knockdown



| Item | Description | Weight |
|-----------|--|--------|
| DPKS/1 | Mobile Infusion Stand - Knockdown, Two Hook | 3.8kg |
| DPKS/1/4H | Mobile Infusion Stand - Knockdown, Four Hook | 3.9kg |

Infusion Poles - Replacement

Features

- Two & four pigtail hook stainless steel infusion poles
- Compatible with our drip stands (19mm diameter upright)
- Pole length 1200mm

Standard Colours

Stainless Steel

Item

DP/2H

DP/4H



Description

Two Hook Infusion Pole

Four Hook Infusion Pole

0.6kg

0.8kg

Features

- Clamp secures to 25mm vertical round tube (designed to fit our range of dressing trolleys)
- When fitting to other trollevs, we strongly recommend that a risk assessment is carried out to assess stability
- When fitted protrudes 45mm
- Supplied with a two pigtail hook infusion pole
- Maximum load on hooks 2kg

Dimensions (w x d x h)

• 55 x 220 x 1200mm





| Item | Description |
|------|------------------------------------|
| DP/C | Clamping Infusion Stands - Trolley |

+44 (0)1258 484455 sales@bristolmaid.com www.bristolmaid.com

Clamping Infusion Stands - Trolley

Wall Mounted Infusion Holders

Features

- · Wall mounted infusion bag holder
- Available with one/two U shaped hooks
- 25mm diameter "D" shaped tube to accept pumps (max load 10kg)
- Maximum load per hook 2kg
- Fixings not included

Dimensions (w x d x h)

- IV001 344 x 75 x 500mm
- IV002 344 x 75 x 570mm




IV001

Description Item IV001 One Hook, Wall Mounted Infusion Holder 1.2kg 1.6kg IV002 Two Hook, Wall Mounted Infusion Holder

+44 (0)1258 484455 sales@bristolmaid.com www.bristolmaid.com

1.8kg











Linen Trolleys & Bags

- Soiled Linen Trolleys 218
 - Linen Trolleys 222
 - Hamper Trolleys 223
 - Soiled Linen Bags 223
- Clean Linen Trolleys 224

Soiled Linen Trolleys - Stainless Steel, Cantilever Frame



Features

- Seamlessly welded from from 304 grade stainless steel for strength & ease
- of cleaning
- Suitable for use with fabric linen bags using a tubular C shaped ring*
- Bags supported on an easy clean solid base
- Available with a plastic, pedal operated lid, when depressed prevents the trolley from moving
- 75mm pressed steel castors

*Recommended minimum bag circumference 1360mm, height 750mm

Dimensions (w x d x h)

- Single 440 x 460 x 820mm
- Double 835 x 460 x 820mm
- Single, Pedal Operated 440 x 560 x 830mm
- Double, Pedal Operated 835 x 560 x 830mm

Options

Linen Bags:

- R/BAG/WHITE White
- R/BAG/BLUE Blue
- R/BAG/RED Red
- R/BAG/GREEN Green
- R/BAG/YELLOW Yellow
- R/BAG/GREY Grey
- R/BAG/ORANGE Orange

See Page 223 For Linen Bag



| Item | Description | Weight |
|--------------|--|--------|
| WT/6/SS | Stainless Steel Cantilever, Single | 6.5kg |
| WT/6/2/SS | Stainless Steel Cantilever, Double | 10.0kg |
| WT/6/PL/SS | Stainless Steel Cantilever, Single, Pedal Operated | 8.5kg |
| WT/6/2/PL/SS | Stainless Steel Cantilever, Double, Pedal Operated | 14.5kg |

Soiled Linen Trolleys - Mild Steel, Cantilever Frame

See Page 223 For Linen Bags

Features

- Welded mild steel frame
- Suitable for use with fabric linen bags using a tubular C shaped ring*
- Bags supported on an easy clean solid base
- Available with a plastic, pedal operated lid, when depressed prevents the trolley from moving
- 75mm pressed steel castors
- *Recommended minimum bag circumference 1360mm, height 750mm

Dimensions (w x d x h)

- Single 440 x 460 x 820mm
- Double 835 x 460 x 820mm
- Single, Pedal Operated 440 x 560 x 830mm
- Double, Pedal Operated 835 x 560 x 830mm

Standard Colours

Paint - Spearmint

Options

Linen Bags:

- R/BAG/WHITE White
- R/BAG/BLUE Blue
- R/BAG/RED Red
- R/BAG/GREEN Green
- R/BAG/YELLOW Yellow
- R/BAG/GREY Grey
- R/BAG/ORANGE Orange





WT/9/PL/MS

Soiled Linen Trolleys - Mild Steel, Sprung Frame

Features

- · Welded mild steel frame
- Suitable for use with either fabric or disposable polythene bags using a sprung rod frame*
- Bags supported on an easy clean solid base
- Available with a plastic, pedal operated lid, when depressed prevents the trolley from moving
- 75mm pressed steel castors

*Factory set for standard 1360mm circumference bags. Can be manually adjusted to suit other bags. Recommended minimum bag height 750mm

Dimensions (w x d x h)

- Single 440 x 460 x 825mm
- Double 835 x 460 x 825mm
- Single, Pedal Operated 440 x 560 x 830mm
- Double, Pedal Operated 835 x 560 x 830mm

Standard Colours

Paint - Spearmint

See Page 223 For Linen Bags

Options

- Linen Bags:
- R/BAG/WHITE White
- R/BAG/BLUE Blue
- R/BAG/RED Red
- R/BAG/GREEN Green
- R/BAG/YELLOW Yellow
- R/BAG/GREY Grev R/BAG/ORANGE - Orange

Description Weight Item WT/9/MS Mild Steel, Sprung Frame, Single 6.5kg WT/9/PL/MS Mild Steel, Sprung Frame, Single, Pedal Operated 8.5kg WT/9/2/MS Mild Steel, Sprung Frame, Double 10.0kg WT/9/2/PL/MS Mild Steel, Sprung Frame, Double, Pedal Operated 14.5kg

Soiled Linen Trolleys Stainless & Mild Steel. Three Castor

Features

- Stainless or mild steel construction
- Stainless Steel (304 Grade) Seamlessly welded for strength & ease of cleaning Mild Steel - Welded & painted
- Suitable for use with fabric linen bags using a tubular C shaped ring*
- Bag supported on a base frame
- 100mm pressed steel castors
- *Recommended minimum bag circumference 1360mm, height 750mm

Dimensions (w x d x h)

• 500 x 425 x 890mm

Standard Colours

- Stainless Steel
- Paint Spearmint

Options

Linen Bags:

- R/BAG/WHITE White
- R/BAG/BLUE Blue
- R/BAG/RED Red
- R/BAG/GREEN Green
- R/BAG/YELLOW Yellow
- R/BAG/GREY Grev
- R/BAG/ORANGE Orange



| Item | Description | Weight |
|----------------|---|--------|
| WT/10/SS | Stainless Steel, Three Castor | 5.5kg |
| WT/10/SS/ZP005 | Stainless Steel, Three Castor, Ring Buffers | 6.0kg |
| WT/10/MS | Mild Steel, Three Castor | 5.5kg |
| WT/10/MS/ZP005 | Mild Steel, Three Castor, Ring Buffers | 6.0kg |
| | | |

WT/9/2/PL/MS



Soiled Linen Trolleys Stainless & Mild Steel, Four Castor

Features

- · Stainless or mild steel construction
- Stainless Steel (304 Grade) Seamlessly welded for strength & ease of cleaning
- Mild Steel Welded & painted
- Suitable for use with fabric linen bags using a tubular C shaped ring*
- Bag supported on a base frame
- 100mm pressed steel castors
- *Recommended minimum bag circumference 1360mm, height 750mm

Dimensions (w x d x h)

- Single 460 x 460 x 890mm
- Double 840 x 460 x 890mm

Standard Colours

- Stainless Steel
- Paint Spearmint

Options

- ZP005 100mm Ring Buffers
- Linen Bags:
- R/BAG/WHITE White
- R/BAG/BLUE Blue
- R/BAG/RED Red
- R/BAG/GREEN Green
- R/BAG/YELLOW Yellow
- R/BAG/GREY Grey
- R/BAG/ORANGE Orange



| Item | Description | Weight |
|----------|--------------------------------------|--------|
| Single | | |
| WT/14/SS | Stainless Steel, Four Castor, Single | 7.0kg |
| WT/14/MS | Mild Steel, Four Castor, Single | 7.0kg |
| Double | | |
| WT/15/SS | Stainless Steel, Four Castor, Double | 8.0kg |
| WT/15/MS | Mild Steel, Four Castor, Double | 8.0kg |

Soiled Linen Trolley - Mild Steel, Folding

Features

- Folds to minimize space
- Suitable for use with fabric linen bags
- 50mm polypropylene/thermoplastic rubber twin-wheel castors
- *Recommended minimum bag circumference 1360mm, height 750mm

WT/7/MS

0.

Dimensions (w x d x h)

- Overall Dimensions 430 x 430 x 1030mm
- Folded 430 x 200 x 1040mm

Standard Colours

Paint - Grev

Options

- Linen Bags:
- R/BAG/WHITE White
- R/BAG/BLUE Blue
- R/BAG/RED Red
- R/BAG/GREEN Green
- R/BAG/YELLOW Yellow
- R/BAG/GREY Grev

See Page 223 For Linen Bags

R/BAG/ORANGE - Orange

Linen Trolley - Bed Changing

Features

- Welded mild steel frame
- Three removable, sunken, liquid retaining trays
- Tray heights 230, 590 & 960mm
- Dividing panel separates clean & soiled linen
- Suitable for use with fabric linen bags using a tubular C shaped ring*
- Bag supported on a base frame
- 100mm pressed steel castors

*Recommended minimum bag circumference 1360mm, height 750mm

Dimensions (w x d x h)

Standard Colours

- Overall Dimensions 830 x 490 x 1040mm
 Paint Grey
- Tray 390 x 400 x 25mm

Options

ZP005 - 100mm Ring Buffers

- Linen Bags:
- R/BAG/WHITE White
- R/BAG/BLUE Blue
- R/BAG/RED Red
- R/BAG/GREEN Green
- R/BAG/YELLOW Yellow
- R/BAG/GREY Grev
- R/BAG/ORANGE Orange

See Page 223





WT/13/MS

| Item | Description | Weight | Item | Description | Weight |
|---------|---------------------|--------|----------|--------------|--------|
| WT/7/MS | Folding, Mild Steel | 3.25kg | WT/13/MS | Bed Changing | 15.5kg |

Hamper Trolleys - Stainless & Mild Steel

Features

- Stainless or mild steel construction
- Stainless Steel (304 Grade) Seamlessly welded for strength & ease of cleaning Mild Steel - Welded & painted
- Suitable for use with fabric hamper bags supported on 4 hooks Bag supported on an easy clean solid base
- 100mm pressed steel castors

Dimensions (w x d x h)

• 815 x 514 x 1000mm

Standard Colours

- Stainless Steel
- Paint Spearmint

Options

- ZP005 100mm Ring Buffers • 2FACO013 - Hamper Bag (190 litre)
- WT/11/SS

| Item | Description | Weight |
|----------|---------------------------------|--------|
| WT/11/SS | Hamper Trolley, Stainless Steel | 12.5kg |
| WT/11/MS | Hamper Trolley, Mild Steel | 12.5kg |

48



Soiled Linen Bags

Features

- Washable 7oz polvester soiled linen bag
- Circumference 1360mm, height 750mm, 120 litres
- Fitted with eyelets and draw cord with plastic toggle
- Available in a range of colours to allow compliance with local, national & international guidelines







R/BAG/RED

R/BAG/GREEN

Weigh

0.3kg

0.3kg

0.3kg

0.3kg

0.3kg

0.3kg 0.3kg



R/BAG/YELLOW



R/BAG/GREY





Description tem R/BAG/WHITE White Linen Bag R/BAG/BLUE Blue Linen Bag R/BAG/RED Red Linen Bag R/BAG/GREEN Green Linen Bag R/BAG/YELLOW Yellow Linen Bag R/BAG/GREY Grey Linen Bag R/BAG/ORANGE Orange Linen Bag



Clean Linen Trolleys - Small (390 Litre)

Features

- Suitable for ward clean linen storage
- Mild steel frame
- Lightweight aluminium cladding & shelves
- · Powder coated for durability
- Two adjustable shelves plus base (max 45kg per shelf) • Shelf size - 1000 x 450mm (w x d)
- Push handles at either end of trollev
- Available with a flexible cover or fold back doors, preventing dust ingress Flexible Cover
 - Fire retardant cover
 - Secured to the top of the trolley & drops down to enclose the front
 - Weighted bottom hem
 - Velcro style side fastenings
- Fold Back Doors
- Fold back doors encloses the front
- . When opened, doors fold against the push handles at either end of the trolley
- Secured with a drop over latch, complete with clasp to allow locking with a padlock (not supplied)
- Drop down bolt secures base of doors
- 125mm pressed steel castors
- Maximum capacity 135kg per trolley

Dimensions (w x d x h)

- Overall Dimensions 1175 x 530 x 1065mm
- Internal Dimensions 1000 x 500 x 810mm

Standard Colours

Paint - Grey

Options

• C008 - 125mm Castors, Heavy Duty Pressed Steel, 2 x Braking, 2 x Non-Braking (set of 4) • ZP010 - 100mm Ring Buffers

| Item | Description | Weight |
|-------------|--|--------|
| WT/42 | Clean Linen Trolley - Small | 39.7kg |
| WT/42/COVER | Clean Linen Trolley - Small, Flexible Cover | 42.4kg |
| WT/42/FBD | Clean Linen Trolley - Small, Fold Back Doors | 46.0kg |

Clean Linen Trolleys - Medium (690 Litre)

Features

- Widely used for bulk linen delivery, part of a linen exchange system, or as an alternative to dedicated ward linen rooms
- Mild steel frame
- Lightweight aluminium cladding & shelves
- Powder coated for durability
- Three adjustable shelves plus base (max 45kg per shelf) • Shelf size - 1180 x 460mm (w x d)
- Inset push handles at either end of trolley
- Low level bumper rail offers protection to both the trolley & fabric of the building
- Available with a flexible cover or fold back doors, preventing dust ingress Flexible Cover
 - Fire retardant cover
 - Secured to the top of the trolley & drops down to enclose the front
 - Weighted bottom hem
 - Velcro style side fastenings
- Fold Back Doors
- · Fold back doors encloses the front
- When opened, doors fold back flat at either end of the trolley*
- Secured with a drop over latch, complete with clasp to allow locking with a padlock (not supplied)
- Drop down bolt secures base of doors
- Optional tow bar fitting
- Tow bar fitted to the front (swivel castor)
- Antiluce fastener fitted to rear (fixed castor)
- Tow bar position approx. 360-380mm above floor
- Recommended max speed 7kph for towed trolleys
- 200mm heavy duty castors, 2 x fixed, 2 x swivel
- Maximum capacity 180kg per trolley

Dimensions (w x d x h)

*Doors will not fold back flat when a tow bar is fitted

Standard Colours

- Overall Dimensions 1430 x 600 x 1495 Internal Dimensions - 1180 x 500 x 1180

| Item | Description | Weight |
|------------------|--|---------|
| Non-Towing | | |
| WT/50 | Clean Linen Trolley - Medium | 87.0kg |
| WT/50/COVER | Clean Linen Trolley - Medium, Flexible Cover | 90.0kg |
| WT/50/FBD | Clean Linen Trolley - Medium, Fold Back Doors | 97.3kg |
| Towing | | |
| WT/50/TBAR | Clean Linen Trolley - Medium, Tow Bar | 89.5kg |
| WT/50/COVER/TBAR | Clean Linen Trolley - Medium, Flexible Cover, Tow Bar | 92.5kg |
| WT/50/FBD/TBAR | Clean Linen Trolley - Medium, Fold Back Doors, Tow Bar | 100.8kg |







Paint - Grey



AGV Compatible Trolleys & Carts Bristol Maid™ is your ideal partner fo transport trolleys & carts for Automated Guided Vehicle (AGV) systems. We can provide solutions that are compatible with Sterile Supplies, Laundry, Medical Records & Pharmacy distribution



AGV Compatible Trolleys & Carts

Bristol Maid[™] is your ideal partner for transport trolleys & carts for Automated Guided Vehicle (AGV) systems. We can provide solutions that are compatible with Sterile Supplies, Laundry, Medical Records & Pharmacy distribution

Clean Linen Trolleys - Large (1100 Litre)

Features

- Widely used for bulk linen delivery, part of a linen exchange system, or as an alternative to dedicated ward linen rooms
- Mild steel frame
- Lightweight aluminium cladding & shelves
- · Powder coated for durability
- Three adjustable shelves plus base (max 45kg per shelf)
- Shelf size 1430 x 555mm (w x d)
- · Inset push handles at either end of trolley
- Low level bumper rail offers protection to both the trolley & fabric of the building
- Available with a flexible cover or fold back doors, preventing dust ingress
- Flexible Cover
- Fold Back Doors
- Optional tow bar fitting
- 200mm heavy duty castors, 2 x fixed, 2 x swivel
- Maximum capacity 180kg per trolley

Dimensions (w x d x h)

- Overall Dimensions 1710 x 700 x1620mm
- Internal Dimensions 1420 x 610 x 1295mm

Standard Colours

Paint - Grey

| Item | Description | Weight |
|------------------|---|---------|
| Non-Towing | | |
| WT/48 | Clean Linen Trolley - Large | 99.6kg |
| WT/48/COVER | Clean Linen Trolley - Large, Flexible Cover | 102.9kg |
| WT/48/FBD | Clean Linen Trolley - Large, Fold Back Doors | 111.5kg |
| Towing | | |
| WT/48/TBAR | Clean Linen Trolley - Large, Tow Bar | 103.1kg |
| WT/48/COVER/TBAR | Clean Linen Trolley - Large, Flexible Cover, Tow Bar | 106.4kg |
| WT/48/FBD/TBAR | Clean Linen Trolley - Large, Fold Back Doors, Tow Bar | 115.0kg |



Flexible Cover

- Fire retardant cover
- Secured to the top of the trolley & drops down to enclose the front
- Weighted bottom hem
- Velcro style side fastenings

Tow Bar Fitting

- Tow bar fitted to the front (swivel castor)
- Antiluce fastener fitted to rear (fixed castor)
- Tow bar position approx. 360-380mm above floor
- Recommended max speed 7kph for towed trolleys

Fold Back Doors

- Fold back doors, encloses the front
- When opened, doors fold back flat at either end of the trolley*
- Secured with a drop over latch, complete with clasp to allow locking with a padlock (not supplied)
- Drop down bolt secures base of doors

*Doors will not fold back flat when a tow bar is fitted









Maternity Equipment



Variable Height Baby Cribs 230

- Fixed Height Baby Cribs 232
- Variable Height Baby Bath 233



Baby Cribs - Variable Height, Battery Operation

Features

- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety and essential performance
- Rechargeable battery powered height adjustment
- Height adjustment from 700-910mm
- Crib can be tilted at either end, horizontal to approx 6 degrees
- Two push handles
- Lower cupboard with fixed shelf
- Two lift up work flaps (max capacity 2.0kg)
- Supplied with a clear plastic crib (PETG) & flame retardant mattress (max capacity 6.0kg)
- Supplied with one rechargeable battery & handset*
- 100mm pressed steel braking castors

*Battery has to be removed from the unit for charging (see options for additional batteries)

Dimensions (w x d x h)

- Overall Dimensions 940 x 510 x 700-910mm
- Mattress 600 x 330 x 60mm
- Cupboard 210 x 265 x 315mm (internal, either side of the column)
- Work Flap 415 x 305mm

Standard Colours

- Paint Spearmint
- Cabinet/Lower Shelf/Work Flaps- Spearmint/Grev

Options

- 2LIDE033 Additional Rechargeable Battery
- 2LIDE037 Charger for Rechargable Battery
- 2PLCO014 Clear Plastic Crib, Solid Side (Spare/Replacement)
- 2FACO052 Mattress, Flame Retardant, Washable (Spare/Replacement)

| Item | Description | Weight |
|----------|-----------------------------------|--------|
| BA/VH | Battery Operated, With Charger | 49.5kg |
| BA/VH/NC | Battery Operated, Without Charger | 49.5kg |

Baby Crib - Variable Height, Gas Strut Operation

Features

- Gas damper powered height adjustment
- Height adjustment from 710-1025mm
- Crib can be tilted at either end, horizontal to approx 6 degrees
- Two push handles
- Supplied with a clear plastic crib (PETG) & flame retardant mattress (max capacity 6.0kg)
- 75mm twin-wheel braking castors

Dimensions (w x d x h)

- Overall Dimensions 930 x 476 x 775-1025mm
- Mattress 600 x 330 x 60mm

Standard Colours

Paint - Spearmint/Grev

Options

- BA902 Slide Out Debris Bin
- BA903 Wire Basket
- BA904 Lift Up Flap
- 2PLCO014 Clear Plastic Crib, Solid Side (Spare/Replacement)
- 2FACO052 Mattress, Flame Retardant, Washable (Spare/Replacement)



| BA904 - |
|---------|

| Item | Description |
|-------|--------------------------------------|
| BA030 | Variable Height, Gas Strut Operation |

sales@bristolmaid.com





20.2kg

- Lift Up Flap





Baby Cribs - Fixed Height

Features

- Constructed from melamine faced chipboard & laminate
- Available with...
- Lower shelf
- Lower shelf & small cupboard (with adjustable shelf)
- Large cupboard (with adjustable shelf)
- Crib can be tilted at either end, horizontal to approx 6 degrees
- Two lift up work flaps (max capacity 2.0kg)
- Supplied with a clear plastic crib (PETG) & flame retardant mattress (max capacity 6.0kg)
- 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions 952 x 500 x 929mm
- Mattress 600 x 330 x 60mm
- Large Base Shelf 810 x 400mm at height 270mm (BA/3/B)
- Small Base Shelf 460 x 355mm at height 270mm (BA/6)
- Small Cupboard 355 x 420 x 345mm (internal) at height 270mm (BA/6)
- Large Cupboard 710 x 420 x 345mm (internal) at height 270mm (BA/6/A)
- Work Flap 415 x 305mm

Standard Colours

- Paint Grey
- Cabinet/Lower Shelf/Work Flaps Spearmint/Grey

Options

- 2PLCO014 Clear Plastic Crib, Solid Side (Spare/Replacement)
- 2FACO052 Mattress, Flame Retardant, Washable (Spare/Replacement)

| Item | Description | Weigh |
|--------|------------------------------|-------|
| BA/3/B | Fixed Height, Lower Shelf | 23.5k |
| BA/6 | Fixed Height, Small Cupboard | 34.0k |
| BA/6/A | Fixed Height, Large Cupboard | 42.0k |

Baby Bath - Variable Height, Gas Strut Operation

Features

- Gas damper powered height adjustment
- Height adjustment from 720-970mm
- Supplied with a plastic bath (max capacity 20.0kg)
- 75mm twin-wheel braking castors

Dimensions (w x d x h)

• 930 x 476 x 720-970mm

Standard Colours

- Paint Spearmint/Grey
- Laminate Spearmint

Options

- BA900 Wire Basket
- BA901 Lift Up Flap
- 2PLCO013 Baby Bath (Spare/Replacement)



BA900 - Wire Basket

| Item | Description |
|-------|--------------------------------------|
| BA015 | Variable Height, Gas Strut Operation |















Medical Records

- Medical Records Trolleys 236
- Medical Record Workstations 242
 - Laptop/Motion C5 Holders 244
 - Carefile Filing Pockets 246
 - CLIQ[®] & RFID Solutions 247
- Medical Record Filing Trolleys 248
- Medical Records Transfer Trolleys 249
 - Medical Records Cabinet 251



Medical Records Trolleys - Small

Features

- Suitable for use with our popular narrow & wide filing pockets (order separately)
- Available as open or with a locking lid...
- Open sides, lower shelf & work flap
- Locking lid & lower shelf
- Locking lid fitted with either a...
- CAM lock (supplied with 2 keys)
- Electronic push button lock
- Maximum load on lid 2.0kg
- Document pocket fitted to end of trolley
- Optional ring buffers to protect the fabric of the building (increases length & width of trolley by approx. 80mm)
- 100mm pressed steel castors

Standard Colours

- Paint Spearmint/Grey
- · Lids/Flaps Spearmint

Options

- ZP010 Ring Buffers
- CF/F/PK10 Narrow Foolscap Filing Pockets (x 10)
- CF/X/PK10 Narrow X-Ray Filing Pockets (x 10)
- CF/XF/PK10 Narrow X-Ray/Foolscap Filing Pockets (x 10)
- CF/WF/PK10 Wide Foolscap Filing Pockets (x 10)
- R/C/T Paper Tabs (x 500)
- R/C/TH Tab Holders (x 25)
- C004 100mm Polymer Castor, 2 x Braking / 2 x Non-Braking (set of 4)
- C005 100mm Polymer Castor, 4 x Braking (set of 4)

Filing Capacity

| Item | Standard |
|----------------|----------|
| Narrow Pockets | 470mm |
| Wide Pockets | 380mm |

*When fully open trolley length is approximately 1190mm

| Item | Description | Overall Dimensions (w x d x h) | Weight |
|----------|--|--------------------------------|--------|
| MR200 | Open Sides, Lower Shelf & Work Flap | 627 x 522 x 997mm | 18.0kg |
| MR210* | Locking Lid & Lower Shelf, CAM Lock | 671 x 522 x 1056mm | 27.2kg |
| MR210CL* | Locking Lid & Lower Shelf, Electronic Push Button Lock | 671 x 522 x 1056mm | 27.2kg |

************ - 2

MR210CL





Carefile Pocket Compatibility

| | | Narrow | | Wide |
|---------|--------------|--------------|--------------|--------------|
| Item | CF/F/PK10 | CF/X/PK10 | CF/XF/PK10 | CF/WF/PK10 |
| MR200 | \checkmark | \checkmark | \checkmark | \checkmark |
| MR210 | \checkmark | \checkmark | \checkmark | \checkmark |
| MR210CL | \checkmark | \checkmark | \checkmark | \checkmark |



Electronic Push Button Lock

- Uses a 10-button keypad & a # button. The # button is hidden beneath the handle when in the locked position
- Easy programming & code change
- Four operating codes with differing levels of authority -Master / Sub-master / User / Technician
- User codes are 4 digits
- Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code



Medical Records Trolleys - Medium

Features

- Suitable for use with our popular narrow filing pockets (order separately)
- Available as open or with locking lids...
- Open sides, lower shelf & work flap
- Locking lids & lower shelf
- · Locking lids fitted with either a...
- CAM lock (supplied with 2 keys)
- Electronic push button lock
- Maximum load on lid 2.0kg
- Document pocket fitted to end of trollev
- Optional ring buffers to protect the fabric of the building (increases length & width of trolley by approx. 80mm)
- 100mm pressed steel castors

Standard Colours

- Paint Spearmint/Grev
- Lids/Flaps Spearmint

Options

- ZP010 Ring Buffers
- CF/F/PK10 Narrow Foolscap Filing Pockets (x 10)
- CF/X/PK10 Narrow X-Ray Filing Pockets (x 10)
- CF/XF/PK10 Narrow X-Ray/Foolscap Filing Pockets (x 10)
- R/C/T Paper Tabs (x 500)
- R/C/TH Tab Holders (x 25)
- C004 100mm Polymer Castor, 2 x Braking / 2 x Non-Braking (set of 4)
- C005 100mm Polymer Castor, 4 x Braking (set of 4)

Filing Capacity



*When fully open trolley length is approximately 1420mm

| Item | Description | Overall Dimensions (w x d x h) | Weight |
|----------|---|--------------------------------|--------|
| MR300 | Open Sides, Lower Shelf & Work Flap | 851 x 522 x 997mm | 20.8kg |
| MR310* | Locking Lids & Lower Shelf, CAM Lock | 920 x 522 x 1056mm | 34.4kg |
| MR310CL* | Locking Lids & Lower Shelf, Electronic Push Button Lock | 920 x 522 x 1056mm | 34.4kg |



Wide

x

x

X

- Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code

 \checkmark

 \checkmark

 \checkmark



Medical Records Trolleys - Large

Features

- Suitable for use with our popular narrow & wide filing pockets (order separately)
- · Available as open or with locking lids...
- Open sides, lower shelf & work flap
- Locking lids & lower shelf
- · Locking lids fitted with either a ...
- CAM lock (supplied with 2 keys)
- Electronic push button lock
- Maximum load on lid 2.0kg
- Document pocket fitted to end/side of trollev
- Optional ring buffers to protect the fabric of the building (increases length & width of trolley by approx. 80mm)
- 100mm pressed steel castors

Standard Colours

- Paint Spearmint/Grev
- · Lids/Flaps Spearmint

Options

- ZP010 Ring Buffers
- CF/F/PK10 Narrow Foolscap Filing Pockets (x 10)
- CF/X/PK10 Narrow X-Ray Filing Pockets (x 10)
- CF/XF/PK10 Narrow X-Ray/Foolscap Filing Pockets (x 10)
- CF/WF/PK10 Wide Foolscap Filing Pockets (x 10)
- R/C/T Paper Tabs (x 500)
- R/C/TH Tab Holders (x 25)
- C004 100mm Polymer Castor, 2 x Braking / 2 x Non-Braking (set of 4)
- C005 100mm Polymer Castor, 4 x Braking (set of 4)

Filing Capacity


| Item | Standard |
|----------------|----------|
| Narrow Pockets | 940mm |
| Wide Pockets | 760mm |



*When fully open trolley length is approximately 2085mm

| Item | Description | Overall Dimensions (w x d x h) | Weight |
|----------|---|--------------------------------|--------|
| MR400 | Open Sides, Lower Shelf & Work Flap | 1158 x 584 x 997mm | 28.0kg |
| MR410* | Locking Lids & Lower Shelf, CAM Lock | 1240 x 522 x 1056mm | 46.4kg |
| MR410CL* | Locking Lids & Lower Shelf, Electronic Push Button Lock | 1240 x 522 x 1056mm | 46.4kg |





Carefile Pocket Compatibility

| | | Narrow | | Wide |
|---------|--------------|--------------|--------------|--------------|
| Item | CF/F/PK10 | CF/X/PK10 | CF/XF/PK10 | CF/WF/PK10 |
| MR400 | \checkmark | \checkmark | \checkmark | \checkmark |
| MR410 | \checkmark | \checkmark | \checkmark | \checkmark |
| MR410CL | \checkmark | \checkmark | \checkmark | \checkmark |



Electronic Push Button Lock

- Uses a 10-button keypad & a # button. The # button is hidden beneath the handle when in the locked position
- Easy programming & code change
- · Four operating codes with differing levels of authority -Master / Sub-master / User / Technician
- User codes are 4 digits
- Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code



Electronic Push Button Lock

- Uses a 10-button keypad & a # button. The # button is hidden beneath the handle when in the locked position
- Easy programming & code change
- Four operating codes with differing levels of authority -Master / Sub-master / User / Technician
- User codes are 4 digits
- Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code

Medical Records Workstations

Features

- Constructed from melamine faced chipboard / laminate offering an extremely strong & durable, yet aesthetically pleasing appearance Available as single & double column units with...
- One / two drawers & x-ray pigeon hole (5 / 10 sections)
- One / two drawers & open cupboard
- Two / four drawers
- Two locking options
- CAM lock (supplied with 2 keys)
- Electronic push button lock
- Where relevant one lock secures all drawers
- Double drawer units fitted with an anti-tilt mechanism preventing two drawers from being opened at the same time
- Max load per drawer 16.0kg
- Push handle & document pocket fitted to end of trolley
- Suitable for use with our filing pockets (order separately)
- Provides 375mm filing capacity per drawer
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Single Column 604 x 538 x 1028mm
- Double Column 1056 x 538 x 1034mm

Standard Colours

- Paint Grey White
- Base / Sides / Back / Centre Panels Grey White
- Top / Drawer Fronts Spearmint
- Edging Grey White



Double Column

| CAM Lock | Electronic Push Button Lock | Description | Weight |
|---------------|-----------------------------|--------------------------------|--------|
| MRW/D2DXP/CAM | MRW/D2DXP/CL | Two Drawers, X-Ray Pigeon Hole | 62.0kg |
| MRW/D2DOC/CAM | MRW/D2DOC/CL | Two Drawers, Open Cupboard | 58.0kg |
| MRW/D4D/CAM | MRW/D4D/CL | Four Drawers | 72.0kg |


Single Column

MRW/S1DOC/CAM

| CAM Lock | Electronic Push Button Lock | Description | Weight |
|---------------|-----------------------------|-------------------------------|--------|
| MRW/S1DXP/CAM | MRW/S1DXP/CL | One Drawer, X-Ray Pigeon Hole | 34.2kg |
| MRW/S1DOC/CAM | MRW/S1DOC/CL | One Drawer, Open Cupboard | 34.0kg |
| MRW/S2D/CAM | MRW/S2D/CL | Two Drawers | 40.6kg |







Laptop Holder - Removable, Swing Arm

Features

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the MRW range of medical record workstations
- Laptop holder fitted to an upright bolted to the rear left hand corner of the trolley (upright & holder supplied for customer to fit)
- Adjustable to accept laptops 280-435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Laptop supported in a cradle fitted with adjustable wing arms
- Laptop can be secured using a security clamp (supplied)*
- Design allows easy access to laptop ports

*Not compatible with the Dell laptops using the N17 wedge lock system

Standard Colours

· Paint - Grey White

| Item | Description | Weight |
|-----------|----------------------|--------|
| MRW/LHR/S | Removable, Swing Arm | 3.7kg |

Laptop Holder - Non-Removable, Swing Arm

Features

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the MRW range of medical record workstations
- Laptop holder fitted to an upright bolted to the rear left hand corner of the trolley (upright & holder supplied for customer to fit)
- Adjustable to accept laptops 280-435 x 22mm (w x h)
- Maximum capacity 2.0kg
- Laptop supported in a cradle fitted with a steel mounting bracket that can be fixed using either nyloc or wing nuts
- Laptop can be secured using a security clamp (supplied)*
- Design allows easy access to laptop ports

*Not compatible with the Dell laptops using the N17 wedge lock system

Standard Colours

· Paint - Grey White

| Item | Description | Weight |
|------------|--------------------------|--------|
| MRW/LHNR/S | Non-Removable, Swing Arm | 5.0kg |



Laptop Holder - Motion C5 Holder

Features

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the MRW range of medical record workstations
- Motion C5 holder fitted to an upright bolted to the rear left hand corner of the trolley (upright & holder supplied for customer to fit)
- Maximum capacity 2.0kg
- Tablet located in a RAM-In-Motion passive dock allowing for easy insertion & removal

Standard Colours

· Paint - Grey White





MRW/S2D/CAM fitted with MRW/LHR/S (laptop not included)

| tem | Description | Weight |
|------------|------------------------------|--------|
| /IRW/LH/C5 | Motion C5 Holder - Swing Arm | 2.7kg |



Carefile Filing Pockets

Features

- Filing pockets suitable with Bristol Maid[™] trolleys and most standard filing suspension systems
- Provides a durable alternative to paper or card types in narrow and wide versions Narrow - 380mm spacings/centres
- Wide 470mm spacings/centres
- Fire retardant PVC coated nylon fabric with stainless steel runners
- Unique 'stay-put' tab holders and pink / blue reversible card tabs
- Supplied in packs of ten

Standard Colours

Vinyl - Grey

CLIQ[®] & RFID Solutions

The use of CLIQ[®] coupled with unique RFID tag identifiers allows a single key to be used with Medicines Management and for Accessing Patient Records.

- Controlled Drug Cabinets (audit trail available)
- Drug & Medicine Cabinets (audit trail available)
- Self-Administration Cabinets
- Bedside Cabinets
- Drug & Medicine Dispensing Trolleys
- Medical Records Trollevs

The use of RFID cards allows a single card to be used with Self-Administration Schemes.


- Self-Administration Cabinets
- Bedside Cabinets





| Item | Description | Internal Pocket Dimensions (w x h) |
|--------------|-----------------------|---|
| Narrow Pocke | ts | |
| CF/F/PK10 | Foolscap (x 10) | 360 x 260mm |
| CF/X/PK10 | X-Ray (x 10) | 360 x 460mm |
| CF/XF/PK10 | X-Ray/Foolscap (x 10) | X-ray - 360 x 460mm Foolscap - 360 x 260mm |
| Wide Pockets | | |
| CF/WF/PK10 | Wide Foolscap (x 10) | 450 x 260mm |
| Paper Tabs & | Tab Holders | |
| R/C/T | Paper Tabs (x 500) | |

Paper Tabs (x 500) R/C/TH Tab Holders (x 25)





| | Case Note Files | Ring Binders | Lever Arch Files |
|----------------|------------------|------------------|------------------|
| | 235 x 50 x 320mm | 250 x 40 x 315mm | 290 x 75 x 315mm |
| Small - 900mm | 18 | 22 | 12 |
| Large - 1800mm | 36 | 44 | 24 |

Electronic Push Button Lock

- Uses a 10-button keypad & a # button. The # button is hidden beneath the handle when in the locked position • Easy programming & code change
- · Four operating codes with differing levels of authority -Master / Sub-master / User / Technician
- User codes are 4 digits
- Master & sub-master codes are 8 digits
- 15.000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code

Medical Record Filing Trolleys - Tiered

Features

- Suitable for use with ring binders, case note & lever arch files (max file height 350mm)
- Available in two sizes with either an open front or locking door/s
- Small 900mm filing capacity
- Large 1800mm filing capacity
- Trolleys supplied with adjustable dividers (small 4 / large 8)
- Doors fold against the sides of the trolley (not flat)
- Large trolleys use one lock to secure both doors
- Two locking options
- CAM lock (supplied with 2 keys)
- Electronic push button lock
- Optional polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning, (increases length & width of trolley by approx. 80mm)

Standard Colours

Paint - Spearmint/Grey

 Top - Spearmint Edging - Grev

Door/s - Grev

100mm pressed steel castors

Dimensions (w x d x h)

- Single Column 604 x 538 x 1028mm
- Double Column 1056 x 538 x 1034mm

Options

- C004 100mm Polymer Castor, 2 x Braking / 2 x Non-Braking (set of 4)
- C005 100mm Polymer Castor, 4 x Braking (set of 4)
- C006 100mm Polymer Castor, 4 x Non-Braking (set of 4)
- ZP010 Ring Buffers (set of 4)
- MRF/LHR/S Laptop Holder Removable, Swing Arm
- MRF/LHNR/S Laptop Holder Non-Removable, Swing Arm
- MRF/LH/C5 Motion C5 Holder Swing Arm

| Item | Description | Weight |
|--------------|--|--------|
| Small | | |
| MRF/S/OC | Small - Open Cupboard | 21.5kg |
| MRF/S/LD/CAM | Small - Locking Door, CAM Lock with Wing Handle | 26.5kg |
| MRF/S/LD/CL | Small - Locking Door, Electronic Push Button Lock | 26.5kg |
| Large | | |
| MRF/L/OC | Large - Open Cupboard | 30.0kg |
| MRF/L/LD/CAM | Large - Locking Doors, CAM Lock with Wing Handle | 41.0kg |
| MRF/L/LD/CL | Large - Locking Doors, Electronic Push Button Lock | 41.0kg |

Medical Records Transfer Trolleys - Single Box

Features

- Suitable for use the transportation of bulk case notes & X-ravs
- Records stored in an upper box internal dimensions 600 x 450 x 250mm (bottom of box positioned 520mm to floor)
- Lower trav
- Available with a locking lid fitted with either a CAM or electronic push button lock providing secure storage patient records
- Trollev max capacity 50kg
- Optional ring buffers to protect the fabric of the building (increases length & width of trolley by approx. 80mm)
- 100mm pressed steel castors

Dimensions (w x d x h)

817 x 523 x 975mm

Standard Colours

Paint - Spearmint/Grev

Options

ZP005 - Ring Buffers

| Item | Description |
|------------|--|
| Open | |
| MR130 | Open Top |
| Locking Li | id |
| MR140 | Locking Lid, CAM Lock |
| MR140CL | Locking Lid, Electronic Push Button Lock |

🛋 Alternative 🦱

Colours

Page 410

00





19.0kg 19.0kg

Electronic Push Button Lock

- Uses a 10-button keypad & a # button. The # button is hidden beneath the handle when in the locked position
- Easy programming & code change
- · Four operating codes with differing levels of authority -Master / Sub-master / User / Technician
- User codes are 4 digits
- Master & sub-master codes are 8 digits
- 15.000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code



Medical Records Transfer Trolley **Double Box**

Features

- Trolley for the transportation of bulk case notes & x-rays
- Records stored in two removable boxes internal dimensions 535 x 344 x 290mm
- Lower Top of box 440mm to floor
- Upper Top of box 975mm to floor
- Clearance between upper & lower boxes, 235mm
- Max capacity 25kg (per box)
- Optional ring buffers to protect the fabric of the building (increases length & width of trolley by approx. 80mm)
- 100mm pressed steel castors

Dimensions (w x d x h)

• 940 x 425 x 1040mm

Standard Colours

Paint - Spearmint

Options

- ZP005 Ring Buffers
- 2PLCO327 Removable Box





Description Weigh tem MR105 Double Box 24.0kg

Medical Records Transfer Trolley **Double Box, Extending**

Features

- Trollev for the transportation of bulk case notes & x-rays
- Records stored in two removable boxes internal dimensions 568 x 368 x 316mm
- Lower Top of box 440mm to floor
- Upper Top of box 975mm to floor
- Clearance between upper & lower boxes, 235mm
- Max capacity 25kg (per box)
- 100mm pressed steel castors

Dimensions (w x d x h) Options

- 875 x 598 x 1040mm
- 2PLCO271 Removable Box

Standard Colours





MR110 Double Box, Extending 31.6kg

Medical Records Transfer Trolley **Double Box, Duo**

Features

- Unique 'two' trolley system for the bulk transportation of case notes & x-rays
- Records stored in fixed boxes internal dimensions 600 x 420 x 255mm
- Allows trolleys to be loaded independently, before being linked together & transported
- Trolleys link together using a simple clevis pin
- Max capacity 25kg (per box)
- 100mm pressed steel castors (2 off braking) & 200mm wheels (braking)

Dimensions (w x d x h) Standard Colours

- 866 x 522 x 970mm
- Paint Spearmint/Grey



| Item | Description | Weight |
|-------|-----------------|--------|
| MR120 | Double Box, Duo | 34.6kg |

tem

Medical Records Cabinet X-Ray Filing

Features

- Lockable x-ray filing cabinet (complete with dividers)
- Three drawers fitted with heavy duty ball bearing slides
- Anti-tilt device prevents more than one drawer being opened at the same time
- Offers total filing capacity of 1740mm, 580mm per drawer* *Suspended filing is not compatible with this unit

Dimensions (w x d x h)

Standard Colours

- 590 x 660 x 1320mm
- Paint Grey



| Item | Description | Weight |
|----------|----------------------|--------|
| 5AM3XRAY | X-Ray Filing Cabinet | 75.0kg |



Phlebotomy



Caretray Trolley - Phlebotomy 254

Phlebotomy Chairs 256



Caretray Trolley - Phlebotomy

Features

- Designed for guick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Available in mild or stainless steel (requires no on-site assembly)
- Fully welded mild steel frame that is blast finished & powder coated in a durable epoxy paint Seamlessly welded one piece frame from 304 grade stainless steel
- Ergonomic push handle fitted to one side of the trolley ensures correct posture & keeps feet away from the frame, ideal for mobile phlebotomy or in cases where the trolley is wheeled for long periods of time
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- Dre piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- Carerails fitted to each side of the trolley to support a wide range of accessories (see options)
- Trolley pre-configured with the following accessories...
- Single Sharps Box Holder Accepts sharps boxes up to seven litres capacity. Sharps box secured using a rubber band fitted with a toggle
- Alcohol Gel Holder Accepts bottles with maximum 90 x 90mm square base
- Single Glove / Mask Box Holder Accommodates glove & mask boxes up to 150mm in width
- Debris Bin & Bracket Plastic debris bin & bracket (internal dims 200 x 135 x 100mm)
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- Each trolley is supplied with perforated coloured label strips (white, blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm polymer castors offering smooth surfaces to assist cleaning, 2 off braking

Description

6 x 100mm Trays

3 x 100mm, 2 x 150mm Trays

4 x 100mm, 1 x 150mm Trays

Dimensions (w x d x h)

Standard Colours

Weight

28.1kg

27.8kg

28.9kg

• Overall Dimensions - 610 x 550 x 993mm • Shallow Tray - 355 x 405 x 100mm (internal)

• Deep Tray - 355 x 405 x 150mm (internal)

- Paint Spearmint
- Plastic Top Single Column 455 x 455 (internal) at height 968mm



Stainless Steel

| Item | Description |
|------------------|----------------------|
| CTPHLES1/10/3S2D | 3 x 100mm, 2 x 150mm |
| CTPHLES1/10/4S1D | 4 x 100mm, 1 x 150mm |
| CTPHLES1/10/6S | 6 x 100mm Trays |
| | |

Mild Steel

CTPHLEM1/10/3S2D

CTPHLEM1/10/4S1D

CTPHLEM1/10/6S

Item

Weight m Travs 28.1kg m Trays 27.8kg 28.9kg

Options



DG/3/15 Triple Glove Box Holder



CT/WF Work Flap



CTD1/10EPB Single Door & Back Panel (Electronic Push Button Lock)



CT/PARM Phlebotomy Armrest



CTD1/10CA Single Door & Back Panel (CAM Lock with Wing Handle)



Phlebotomy Chair Variable Height

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- · Electric or hydraulic height adjustment
- Height adjustment from 490-950mm
- Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections (optional electric backrest)
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- Electric models (230V ~, 50/60H) fitted with a UK 3 pin plug & available with a hand or foot switch
- Hydraulic models use a foot pump pedal positioned on both sides of the couch
- Seamless upholstery using flame retardant anti-microbial vinyl
- Adjustable polyurethane phlebotomy armrests
- Side support loops assist with patient mount & dismount (self-assembly)
- Paper roll holder (suitable for rolls max 510mm wide x 200mm diameter)
- Frames prepared to accept optional accessories
- Maximum capacity 225kg

Dimensions (w x d x h)

• 770 x 1710 x 490-950mm

Options

• E06 - Battery Back Up

| Item | Description | Weight |
|------------|---|--------|
| EC5305/VBB | Hydraulic | 67.0kg |
| EC5310/VBB | Electric - Hand Switch | 67.0kg |
| EC5315/VBB | Electric - Foot Switch | 67.0kg |
| EC5320/VBB | Electric with Electric Backrest - Hand Switch | 69.3kg |
| EC5325/VBB | Electric with Electric Backrest - Foot Switch | 69.3kg |

Phlebotomy Chair Variable Height, Bariatric

Features

- BS 2838-2: 1991 Examination & general treatment couches, build, performance & safety specifications
- EN 60601-1:2006+A12:2014 Medical electrical equipment. General requirements for basic safety & essential performance
- Electric height adjustment
- Height adjustment from 490-950mm
- Adjustable angle backrest
- Manual -25 to +90 degrees
- Electric -22 to +90 degrees
- Adjustable angle footrest from -90 to +10 degrees
- Anti-trap gas struts fitted with release levers used on the back & footrest sections
- Retractable castors fitted to both ends to assist with movement & positioning
- One adjustable foot fitted to allow levelling on uneven floors
- Clearance beneath frame allows access for mobile hoists
- 230V ~, 50/60H actuator fitted with a UK 3 pin plug & available with a hand or foot switch
- Seamless upholstery using flame retardant anti-microbial vinyl
- Adjustable polyurethane phlebotomy armrests
- Side support loops assist with patient mount & dismount (self-assembly)
- Paper roll holder (suitable for rolls max 510mm wide x 200mm diameter)
- Frames prepared to accept optional accessories
- Maximum capacity 325kg

Dimensions (w x d x h)

• 990 x 1710 x 520-950mm

Standard Colours

- Paint White
- Vinyl Bristol Blue

Options

E18 - Battery Back Up

| Item | Description |
|------------|---|
| EC5330/VBB | Electric - Hand Switch |
| EC5335/VBB | Electric - Foot Switch |
| EC5340/VBB | Electric with Electric Backrest - Hand Switch |
| EC5345/VBB | Electric with Electric Backrest - Foot Switch |

VBB VAQ VFR

VGH

VSD VOB VL

VCH

sales@bristolmaid.com +44 (0)1258 484455

- EC5310/VE EC5315/VE
- EC5320/VE EC5325/VE

VPF

VCB

Standard Colours Paint - White Vinvl - Bristol Blue







Weight

+44 (0)1258 484455 www.bristolmaid.com

Alternative colours available at extra cost. Change suffix 'VBB' on the product code with any of the following i.e. Fixed Height Phlebotomy Chair in Aguamarine - EC5000/VAQ



Phlebotomy Chair Fixed Height

Features

- BS EN 16139: 2013 Level 2. Furniture. Strength, durability &
- BS 7945: 1999 (stability)
- Certified for use up to 32 stone (203kg)
- Stable four leg design
- Black finish tubular frame
- Seamless upholstery using flame retardant anti-microbial vinyl
- Adjustable polyurethane phlebotomy armrest
- Armrest can be positioned on either side of the chair (additional armrest available if required)
- Maximum capacity 203kg

Dimensions (w x d x h)

• 620 x 580 x 800mm

Standard Colours

- Paint Black
- Vinyl Bristol Blue

Options

• E15 - Phlebotomy Armrest (Single)



| Item | Description | Seat Heigl |
|------------|--------------|------------|
| EC5000/VBB | Fixed Height | 485mm |



Portering Chairs



- Portering Chair Rear Steer 262
- Portering Chair Front Steer 263
- Portering Chair Bariatric, Rear Steer 265
 - Options 266



Alternative colours available at additional cost. Change suffix 'VBB' on the product code with any of the following i.e. Rear Steer, Nesting, Sliding Foot Rest in Bristol Blue - G/200/RS/VBB



Portering Chair Rear Steer, Nesting, Sliding Foot Rest

Features

- Nests to occupy reduced floor space and encourage tidy storage
- · Seamless upholstery using flame retardant antimicrobial vinyl
- Moulded drop down arm rests allows side access for patients and helpers
- Sliding rubber foot rest aids access as well as supporting feet during transit
- Design ensures chair will not tip should patient stand on the step
- X-ray/notes pocket on rear of seat
- Two hooks for patients belongings
- · Fittings to accept IV poles/cylinder holders mounted to rear uprights
- 200mm fixed wheels/150mm total stop braking castors (cushioned)
- Maximum capacity 190.0kg

Dimensions (w x d x h)

Overall Dimensions - 660 x 840 x 1040mm Paint - Grey

- Seat 540 x 420mm at Height 530mm Upholstery - Bristol Blue
- Back 490 x 455mm

Options

- CYHP/D1 Cylinder Holder, CD/D
- CYHP/E1 Cvlinder Holder. E
- G/CLS Calf/Leg Supports (pair)
- G/DP Infusion Pole
- G/SB Lap Seatbelt, Factory Fitted



X-Ray/Notes Pocket

| Item | Description | Weight |
|--------------|--|--------|
| G/200/RS/VBB | Rear Steer, Nesting, Sliding Foot Rest | 29.0kg |

Portering Chair Front Steer, Nesting, Sliding Foot Rest

Features

- Nests to occupy reduced floor space and encourage tidy storage
- · Seamless upholstery using flame retardant antimicrobial vinyl
- Moulded drop down arm rests allows side access for patients and helpers
- Sliding rubber foot rest aids access as well as supporting feet during transit
- Design ensures chair will not tip should patient stand on the step
- X-rav/notes pocket on rear of seat
- Two hooks for patients belongings
- Fittings to accept IV poles/cylinder holders mounted to rear uprights
- 125mm swivel castors/200mm braking fixed wheels (cushioned)
- Maximum capacity 190.0kg

Dimensions (w x d x h)

- Overall Dimensions 660 x 1005 x 1090mm
- Seat 540 x 420mm at Height 600mm Back - 490 x 455mm

Options

- CYHP/D1 Cylinder Holder, CD/D
- CYHP/E1 Cylinder Holder, E
- G/CLS Calf/Leg Supports (pair)
- G/DP Infusion Pole
- G/SB Lap Seatbelt, Factory Fitted
- G/CARD/LOCK Card/Coin Lock, Factory Fitted
- G/CARD/STATION Head of Station
- 2LOKE011 Release Card




X-Ray/Notes Pocket

| Item | Description |
|--------------|----------------------------------|
| G/305/FS/VBB | Front Steer, Nesting, Sliding Fo |

www.bristolmaid.com

Standard Colours

- G/CARD/LOCK Card/Coin Lock, Factory Fitted
- · G/CARD/STATION Head of Station
- 2LOKE011 Release Card



Standard Colours

 Paint - Grey Upholstery - Bristol Blue



Alternative colours available at additional cost. Change suffix 'VBB' on the product code with any of the following i.e. Front Steer, Nesting, Sliding Foot Rest in Lime - G/305/FS/VLI



www.bristolmaid.com

+44 (0)1258 484455

sales@bristolmaid.com



Alternative colours available at additional cost. Change suffix 'VBB' on the product code with any of the following i.e. Front Steer, Nesting, Hinged Foot Rests in Flame Red - G/205/FS/VFR



Portering Chair Front Steer, Nesting, Hinged Foot Rests

Features

- Nests to occupy reduced floor space and encourage tidy storage
- Seamless upholstery using flame retardant antimicrobial vinyl
- Moulded drop down arm rests allows side access for patients and helpers
- Folding/hinging foot rests aids access as well as supporting feet during transit

Standard Colours

Upholstery - Bristol Blue

Paint - Grev

- X-ray/notes pocket on rear of seat
- Two hooks for patients belongings
- Fittings to accept IV poles/cylinder holders mounted to rear uprights
- 125mm swivel castors/200mm braking fixed wheels (cushioned)
- Maximum capacity 190.0kg

Dimensions (w x d x h)

- Overall Dimensions 660 x 710 x 1055mm
- Seat 540 x 420mm at Height 530mm
- Back 490 x 455mm

Options

- CYHP/D1 Cylinder Holder, CD/D
- CYHP/E1 Cylinder Holder, E
- G/CLS Calf/Leg Supports (pair)
- G/DP Infusion Pole
- G/SB Lap Seatbelt, Factory Fitted
- G/CARD/LOCK Card/Coin Lock, Factory Fitted
- G/CARD/STATION Head of Station
- 2LOKE011 Release Card



X-Rav/Notes Pocket

| tem | Description | Weight |
|--------------|---|--------|
| G/205/FS/VBB | Front Steer, Nesting, Hinged Foot Rests | 27.5kg |

Portering Chair Bariatric, Rear Steer, Sliding Foot Rest

Features

- Heavy duty bariatric portering chair
- Seamless upholstery using flame retardant antimicrobial vinyl
- Moulded drop down arm rests allows side access for patients and helpers
- Sliding rubber foot rest aids access as well as supporting feet during transit
- Design ensures chair will not tip should patient stand on the step
- X-rav/notes pocket on rear of seat
- Two hooks for patients belongings
- Fittings to accept IV poles mounted to rear uprights
- Retractable diagonal seatbelt
- 200mm fixed wheels/200mm total stop braking castors (cushioned)
- · We recommend that two members of staff are in attendance when used with a patient in excess of 190kg
- Maximum capacity 380.0kg

Dimensions (w x d x h)

- Overall Dimensions 930 x 1000 x 1090mm
- Seat 740 x 670mm at Height 530mm
- Back 700 x 480mm

Options

- G/500/RS/CYL Cylinder Holder, CD/D
- G/CLS/500 Calf/Leg Supports (pair)
- G/DP Infusion Pole
- G/CARD/LOCK Card/Coin Lock, Factory Fitted
- G/CARD/STATION Head of Station
- 2LOKE011 Release Card



X-Ray/Notes Pocket

| Item | Description |
|--------------|-----------------------------------|
| G/500/RS/VBB | Bariatric, Rear Steer, Sliding Fo |

264



Standard Colours

 Paint - Grey Upholstery - Bristol Blue



Alternative colours available at additional cost. Change suffix 'VBB' on the product code with any of the following i.e. Bariatric, Rear Steer, Sliding Foot Rest in Coastal Blue - G/500/RS/VCB



+44 (0)1258 484455

www.bristolmaid.com

Cylinder Holders

Features

- Compatible with 100mm diameter cylinders (C/CD/D/E)
- Mounts to IV support tube at rear of chair using a handwheel
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models

Standard Colours

· Paint - Grey

Calf / Leg Support - Pair

Features

- Moulded calf/leg support on captive retractable support tube
- Retracts beneath the seat when not in use
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models
- Supplied as a pair

Standard Colours

Infusion Pole

Features

- Stainless steel infusion pole with a U shaped hook (length of pole 980mm)
- Mounts to IV support tube at rear of chair using a handwheel
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models



Features

- Lap seat belt with snap buckle
- Compatible with G/200/RS, G/205/FS & G/305/FS models
- Factory fitted

Standard Colours





| Item | Description | Weight |
|------|--------------|--------|
| G/SB | Lap Seatbelt | 0.4kg |



Card/Coin Lock

Features

- Compatible with old & new one pound coin (euro compatible)
- Coin/card remains in the lock until the chair is connected to another chair or the head of station
- Plastic release card available for use by porters
- Factory fitted
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models
 Replacement card/coin lock cassette available (requires cassette removal tool)



ItemDescriptionWeightG/CARD/LOCKCard/Coin Lock, Factory Fitted1.0kg2LOKE011Plastic Release Card0.01kg2LOCK217Card/Coin Lock Cassette, New £1 Coin Compatible0.2kg2LOCK218Cassette Removal Tool0.01kg

| Item | Description | Weight |
|----------------|-----------------|--------|
| G/CARD/STATION | Head of Station | 1.6kg |

Card/Coin Lock Wall Bracket

Features

- Head of station for chairs fitted with card/coin locks
- Requires fitting to the fabric of the building (wall fixings not supplied)
- Chain length 1000mm
- One required for each row of chairs, (G/200/RS, G/205/FS & G/305/FS)
- Compatible with G/200/RS, G/205/FS, G/305/FS & G/500/RS models





Procedure Trolleys & Carts



Plaster Trolleys 270

- Mobile Workstations 271
 - Procedure Carts 272

Plaster Trolley - Stainless Steel

Features

- Seamlessly welded easy clean stainless steel plaster trolley (304 grade)
- Supplied with three drawers (max load 10kg per drawer)
- Integral bucket ring (supplied with a 15 litre stainless steel bucket)
- Lower shelf (max load 40kg)
- 100mm pressed steel anti-static non marking castors

Dimensions (w x d x h)

Overall Dimensions - 1150 x 540 x 915mm
 Stainless

Standard Colours

- Drawers 650 x 410 x 100mm (internal)
- Work Surface 745 x 515 at height 915mm
- Lower Shelf 746 x 515 at height 175mm

Options

- WT/96/FLAP Lift Up Flap
- WT/96/SS/LP Louver Panel
- 2PLCO024 Plastic Lin Bin
- ZP010 Ring Buffers
- 2BK 2 x Braking Castors / 2 x Non-Braking Castors
- 4BK 4 x Braking Castors



2PLCO024 optional lin bin

Low level mobile workstation constructed from mild steel

- Recessed liquid retaining top
- Single point locking

Features

- To maximise safety, the unit is fitted with an anti-tilt mechanism, allowing only one drawer to be opened at any one time
- Full extension drawers, heights 85, 170 & 255mm (max load 30kg per drawer)
- Two tone colour keyway, drawer fronts available in a range of colours
- Low level corner bumpers to protect the fabric of the building
- 80mm heavy duty castors, 2 off braking

Dimensions (w x d x h)

- Overall 615 x 500 x 810mm
- Drawers 480 x 375 x 85/170/255mm (internal)

Standard Colours

- Body Light Grey
- Drawers Dark Grey



| Item | Description | Weight |
|------------|------------------------|--------|
| WT/96/3/SS | Three Drawers & Bucket | 51.0kg |

Mobile Workstation - Mild Steel, Low Level



Drawers, Cam Lock



Procedure Carts - Polymer, Metro Flexline[™]

Features

- Range of pre-configured carts to suit a wide range of applications
- Lightweight polymer construction with smooth rounded corners to assist cleaning
- Upper work surface provides a writing surface/space for preparation work
- Three locking options (CAM/electronic locking & security seal compatible)
- Full extension lightweight polymer drawers, depths 76, 152, 229mm (max load 7-20kg per drawer (dependent on drawer size)
- Incorporates Microban offering antimicrobial protection
- Two tone colour keyway, drawer handles available in a range of colours
- 127mm castors, 2 x Directional Castors / 2 x Braking Castors

Dimensions (w x d x h)

• Drawers - 514 x 432 x 76/152/229mm (internal)

Standard Colours

- Body Light Taupe
- Drawers Handles Model Dependant

Options

• 8FL190 - Drawer Label Holders (x 10)



Drawer Label Holders

| Item | Description | Overall Dimensions (w x d x h) | Weight |
|-----------|--|--------------------------------|--------|
| 8FLCRP4 | Emergency Cart - 3 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Code Blue Handles, Security Seal Lock | 1056 x 607 x 1273mm | 91.0kg |
| 8FLPED1 | Paediatric Emergency Cart - 8 x 76mm, 1 x 152mm Drawers, Coloured Handles, Security Seal Lock | 1056 x 607 x 1273mm | 91.0kg |
| 8FLTMENT2 | Treatment Cart - 2 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Slate Blue Handles, CAM Lock | 819 x 568 x 1064mm | 51.0kg |
| 8FLCOM-LT | Laptop Cart - 1 x Keyboard Tray, 1 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Slate Blue Handles, Keyless Entry Touch Pad | 819 x 568 x 1064mm | 83.0kg |
| 8FLCLINE | Central Line Cart - 3 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Dark Taupe Handles, Keyless Entry Touch Pad | 819 x 568 x 1064mm | 91.0kg |
| 8FLCCU1 | Critical Care Cart - 3 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Red Handles, Keyless Entry Touch Pad | 840 x 569 x 1807mm | 90.0kg |





Fitted To Trolley

Drawer Seal Lock Bar



| mergency 3FLCRP4 | Paediatric Emer- gency 8FLPED1 | Treatment 8FLTMENT2 | Laptop 8FLCOM-LT | Central Line 8FLCLINE | Critical Care 8FLCCU1 |
|---------------------|-----------------------------------|------------------------|---------------------|--------------------------|--------------------------|
| | | 2 | | 2 | |
| | | 2 | | 1 | |
| 3 | | | | | |
| 2 | | | | | |
| 4 | | | | 2 | 4 |
| | | | 3 | 2 | 4 |
| | | | | \checkmark | |
| \checkmark | \checkmark | | | | |
| | | | | √ | |
| \checkmark | | | | | |
| \checkmark | \checkmark | | | | |
| \checkmark | | | | | |
| \checkmark | \checkmark | | | | |
| \checkmark | \checkmark | | | | |
| | | | \checkmark | | |
| \checkmark | \checkmark | | | | |
| \checkmark | \checkmark | | | | |
| \checkmark | \checkmark | | | | |
| | \checkmark | | | | |
| | | | √ | √ | \checkmark |
| | | | | | \checkmark |
| | | | | | \checkmark |



Procedure Carts - Polymer, Metro Starsys[™]

Features

- Range of pre-configured carts to suit a wide range of applications
- Lightweight polymer construction with smooth rounded corners to assist cleaning
- Upper work surface provides a writing surface/space for preparation work
- Five locking options, (electronic touchpad/pushbutton/keyless, CAM & security seal compatible)
- Full extension lightweight polymer drawers, depths 76, 152, 229mm (max load 7-20kg per drawer (dependent on drawer size)
- Incorporates Microban offering antimicrobial protection
- Two tone colour keyway, drawer handles available in a range of colours
- 127mm castors, 2 x Swivel Castors / 2 x Braking Castors

Standard Colours

- Body Light Taupe
- Drawers Handles Model Dependant

Options

• 8SXLABKIT - Drawer/Door Label Holders, (x 10)



Drawer / Door Label Holders

| Item | Description | Overall Dimensions (w x d x h) | Weight |
|------------|--|--------------------------------|---------|
| 8SXRSANES | Anaesthesia Cart - 4 x 152mm, 1 x 229mm Drawers, Dark Taupe Handles, CAM Lock | 882 x 633 x 1136mm | 104.0kg |
| 8SXRSDRS | Dressing Cart - 3 x 76mm, 2 x 152mm, 1 x 229mm Drawers, Dark Taupe Handles, CAM Lock | 577 x 633 x 1060mm | 83.0kg |
| 8SXRSTRAMA | Trauma Cart - 3 x 76mm, 4 x 152mm Drawers, Slate Blue Handles, CAM Lock | 882 x 633 x 1136mm | 116.0kg |
| 8SXRSCCU | Critical Care Cart - 1 x 76mm, 1 x 152mm, 3 x 229mm Drawers, Slate Blue Handles | 882 x 633 x 1212mm | 107.0kg |
| 8SXRSARTH | Arthroscopic Cart - 6 x 152mm, 1 x 229mm Drawers, Dark Taupe Handles, CAM Lock | 882 x 633 x 1136mm | 110.0kg |
| 8SXRSIV | IV Therapy Cart - 5 x 76mm, 2 x 229mm Drawers, Green Handles, CAM Lock | 882 x 633 x 1136mm | 106.0kg |
| 8SXRSCAST | Cast Cart - 2 x 76mm, 2 x 152mm, 2 x 229mm Drawers, Slate Blue Handles | 882 x 633 x 1212mm | 109.0kg |



Fitted To Trolley

Pod Tray With Dividers Full Size Utility Bin (Metal)



| | Anaesthesia 8SXRSANES | Trauma 8SXRSTRAMA | Critical Care 8SXRSCCU | Arthroscopic 8SXRSARTH | IV Therapy 8SXRSIV | Cast 8SXRSCAST |
|-------|--------------------------|----------------------|---------------------------|---------------------------|-----------------------|-------------------|
| | | 2 | 1 | | 4 | 1 |
| | | | \checkmark | | | √ |
| ock) | | | | | \checkmark | |
| () | | \checkmark | | \checkmark | | |
| | | \checkmark | \checkmark | \checkmark | \checkmark | √ |
| | | \checkmark | \checkmark | \checkmark | | √ |
| | | | \checkmark | | | √ |
| | | | | \checkmark | | |
| Lock) | | | | | \checkmark | |
| :k) | | \checkmark | | | | |
| | | | | \checkmark | | |
| | | \checkmark | | | | |
| | | \checkmark | | \checkmark | | |
| ock) | \checkmark | | | | | |
| Lock) | \checkmark | | | | | |
| | | \checkmark | \checkmark | | | |
| | \checkmark | | | \checkmark | | √ |
| | | \checkmark | | | \checkmark | |
| | | \checkmark | \checkmark | \checkmark | | \checkmark |
| | \checkmark | | | | | |
| | \checkmark | | | | | |
| | \checkmark | | | | | |
| | \checkmark | | | | | |
| | \checkmark | | | | | |
| | \checkmark | | \checkmark | \checkmark | | |



\checkmark 1 \checkmark \checkmark \checkmark \checkmark \checkmark \checkmark \checkmark \checkmark \checkmark

Procedure Carts - Metro Basix[™]

Features

- Lightweight metal procedure carts
- Recessed liquid retaining plastic work surface
- Three locking options
- CAM (Compact)
- Push button (Trauma)
- Security seal compatible (Anaesthesia)
- Full extension drawers, depths 76, 152 & 230mm
- Compact Max load 7kg per drawer
- Trauma/Anaesthesia Max load 24kg per drawer
- Drawer fronts available in a range of colours
- 125mm castors, total lock at front & directional lock at rear

Dimensions (w x d x h)

- Standard Colours
- Compact 489 x 470 x 876mm Body - Light Taupe Drawers - Dark Blue
- Drawers (internal)
- (see below for alternative colours) • 76mm - 381 x 432 x 76mm
- 152mm 381 x 432 x 152mm
- 230mm 381 x 432 x 230mm
- Trauma/Anaesthesia Cart 762 x 613 x 1440mm
- Drawers (internal)
- 76mm 432 x 565 x 76mm
- 152mm 432 x 565 x 152mm
- 230mm 432 x 565 x 230mm

Options

8MBA172 - Label Holders (x 12)

Alternative drawer colours are available. Change suffix 'DB' on the product code with any of the following i.e. - Compact Basix Cart with red drawers 8MBC1210RE



| | Description | |
|------------------|--|--------|
| 8MBC1210DB | Compact - 1 x 76mm, 2 x 152mm, 1 x 230mm Drawers | 63.0kg |
| 8MBX3210TRAUMADB | Trauma - 3 x 76mm, 2 x 152mm, 1 x 230mm Drawers | 78.0kg |
| 8MBX1310ANES4DB | Anaesthesia - 1 x 76mm, 3 x 152mm, 1 x 230mm Drawers | 81.0kg |













Resuscitation Trolleys & Carts

Resuscitation Trolleys 280

- Caretray Trolley Resuscitation 282
- Resuscitation Carts Polymer 284

Airways / Difficult Intubation Trolley 287



Resuscitation Trolleys - Mild Steel

Features

- · Popular resuscitation trollev
- Drawer contents are secured when the lower cupboard is closed (security seal compatible)
- · Upper drawer supplied with a plastic drawer liner, with dividers
- Recessed liquid retaining stainless steel work surface (usable 630 x 445)
- Lift up work flap (usable 305 x 455mm)
- Two hook infusion pole mounted to rear of trolley
- Low level bumper bar
- Push handles one end
- 125mm pressed steel castors

Dimensions (w x d x h)

- Drawers (internal) 510
 - 70mm 600 x 430 x 70mm
 - 96mm 600 x 430 x 96mm
 - 155mm 600 x 430 x 155mm
 - Lower Cupboard (internal) • 625 x 475 x 295mm

Standard Colours

- Paint Red
- Lift Up Flap White

Drop down cupboard door prevents drawers from being opened. Can be secured using security seals.

shown with optional cylinder holder WT/CYL/D/E

| ltem | Description | Overall Dimensions (w x d x h) | Weight |
|-------|---|--------------------------------|--------|
| WT200 | 2 x 70mm & 1 x 155mm Deep Drawers, Lower Cupboard | 825 x 575 x 999mm | 74.0kg |
| WT210 | 3 x 96mm Deep Drawers, Lower Cupboard | 825 x 575 x 999mm | 75.1kg |
| WT250 | 4 x 70mm & 1 x 155mm Deep Drawers, Lower Cupboard | 825 x 575 x 1165mm | 86.0kg |

Cardiac Board

Features

- Laminate cardiac board approx 500 x 400 x 17mm, complete with fixing brackets
- Positioned at the rear of the trolley for
- easy access
- Factory fitted

| Item | Description | Weight |
|-----------|---------------|--------|
| WT/CBOARD | Cardiac Board | 1.5kg |

Cylinder Holders

Features

- WT/CYL/D/E compatible with AD, CD, D, DD, PD, RD, ZD & E cylinders, max diameter 100mm
- WT/CYL/F compatible with F, LF, VF, HX & ZX cylinders, max diameter 145mm
- Maximum two cylinder holders per trolley

| Item | Description | Weight |
|------------|------------------------------|--------|
| WT/CYL/D/E | D/E Size, Max Diameter 100mm | 3.0kg |
| WT/CYL/F | F Size, Max Diameter 145mm | 3.4kg |

Drawer Tidy

Features

 Plastic drawer liner with removable dividers One tidy supplied as standard with each trolley

| Item | Description | Weight |
|-------|-------------|--------|
| WT/DT | Drawer Tidy | 1.6kg |



Four Gang Socket

Features

- Switched & fused four gang socket
- Supplied with 3 metres of beige cable & fitted with a UK 3 pin plug
- Fitted to the rear of the trolley, complete with a cleat to secure the cable when not in use
- · Factory fitted



| Item | Description | Weight |
|------------|------------------|--------|
| E/SOCKET/4 | Four Gang Socket | 1.4kg |

Guard Rail

Features

- Stainless steel quard rail
- Increases the height of the trolley by 60mm
- Compatible with WT200, WT210 and WT250 models
- Factory fitted




| Item | Description | Weight |
|----------|-------------|--------|
| G/RAIL-C | Guard Rail | 1.5kg |

Security Seals

Features

Snap on plastic security



Description

Security Seals

Castors Features

- Pressed steel castors
- Supplied as a set of 4 castors
- Factory fitted

| Item | Description | Weight |
|------|---|--------|
| C012 | 2 x Braking, 2 x Non-Braking (set of 4) | 5.0kg |



WT/CYL/D/E

WT/CYL/F



sales@bristolmaid.com +44 (0)1258 484455 www.bristolmaid.com

Item 2PLCO026



Options (refer to page 62 for additional information)

- · CT/AGH Alcohol Gel Holder
- CT/AIR Airway Tube Holder
- CT/BIN Debris Bin & Bracket
- CT/DH/A4 Document Holder, A4
- CT/CR5 Catheter Rack, Five Section
- CT/FS Flat Shelf Removable
- CT/IV Infusion Pole & Bracket
- CT/OB/MS Overbridge

- TBIN4/MS Tilt Bin, 4 Compartment (overbridge required)
- TBIN5/MS Tilt Bin, 5 Compartment (overbridge required)
- TBIN6/MS Tilt Bin, 6 Compartment (overbridge required)
- CT/WF Work Flap
- CT/CR2 Catheter Rack, Two Section
 GMBH/15 Glove / Mask Box Holder One Box, 15mm Hook
 - SBH/15 Sharps Box Holder One Box, 15mm Hook
 - CT/SBB/15 Sharps Box Bracket Universal (Daniels / Frontier / Sharpak)
 - 2PLCO026 Security Seals, Pack of 100

Caretray Trolley - Resuscitation

Features

- Single column standard level resuscitation trollevs
- Designed for quick & easy complete disassembly when conducting thorough or deep cleaning programmes (including steam cleaning)
- Fully welded mild steel frame that is blast finished & powder coated in a durable epoxy paint (requires no on-site assembly)
- One piece moulded top with raised edges on three sides & integral rim at the front to contain any spillages (max load 30kg)
- One piece moulded side panels incorporate runners with positive stops, preventing trays from being accidentally withdrawn
- Carerails fitted to each side of the trollev to support a wide range of accessories (see options)
- Fitted with a back panel & tamper evident door providing security for the contents of the trays (security seal compatible)
- Trolley pre-configured with the following...
- Cylinder Holder Compatible with C/CD/D size cylinders
- Single Glove / Mask Box Holder Accommodates glove & mask boxes up to 140 mm x 95 mm
- 100/150mm deep translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- Each trollev is supplied with six sheets of perforated coloured label strips - white, blue, green, yellow & red (not shown on trolley)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building; smooth surfaces assist cleaning
- 100mm polymer castors offering smooth surfaces to assist cleaning, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions 550 x 550 x 993mm
- Overall Dimensions (with overbridge) 600 x 550 x 1620mm
- Plastic Top
- Single Column 455 x 455 (internal) at height 968mm
- Travs
- Shallow 355 x 405 x 100mm (internal)
- Deep 355 x 405 x 150mm (internal)

Standard Colours

· Paint - Red

| Item | Description | Weight |
|----------------------|----------------------------|--------|
| CTRESUSM1/10/3S2D/SS | 3 x 100mm, 2 x 150mm Trays | 33.3kg |
| CTRESUSM1/10/4S1D/SS | 4 x 100mm, 1 x 150mm Trays | 32.9kg |
| CTRESUSM1/10/6S/SS | 6 x 100mm Trays | 34.0kg |

Overbridge

Features

- Fixes onto the rear uprights of the trolley.
- Provides a fixing location for tilt bins (order separately) Factory fitted

Standard Colours

Paint - Grev White





Overbridge, Mild Steel

Ö ంర

Tilt Bins

Features

- Available with 4, 5 or 6 tilting compartments
- Clear tilt bins provide visible storage & assists organisation
- · Compartments tilt forward to a 45° angle allowing easy access of contents & can be removed to aid refiling/cleaning
- Label holder for easy identification of the contents
- Retaining bar secures compartments in place, preventing bins from tipping forward during transit

Plastic - Cream

· Paint - Grey White

Attaches to an overbridge (order separately)

Dimensions (w x d x h)

Standard Colours

Compartment Dimensions

- TBIN4/MS 143 x 127 x 190mm
- TBIN5/MS 108 x 95 x 146mm
- TBIN6/MS 89 x 67 x 95mm



Label Holder

45° Anale Tilt

Top: TBIN6/MS Middle: TBIN5/MS Bottom: TBIN4/MS

| Item | Description | Overall Dimensions (w x d x h) | Weight |
|----------|------------------|--------------------------------|--------|
| TBIN4/MS | Four Compartment | 600 x 171 x 207mm | 2.3kg |
| TBIN5/MS | Five Compartment | 600 x 135 x 164mm | 1.6kg |
| TBIN6/MS | Six Compartment | 600 x 93 x 113mm | 0.9kg |

sales@bristolmaid.com

Weight

2.0kg



Resuscitation Carts - Polymer, Metro Lifeline[™]

Features

- Range of pre-configured emergency carts
- Incorporates Microban offering antimicrobial protection
- Lightweight polymer construction with smooth rounded corners to assist cleaning
- Single seal mechanism secures the top compartment, drawers and tilt out side bins (security seal compatible)
- Top storage compartment with clear removable cover provides instant access to first line equipment
- Strapless defibrillator platform on swing arm adjusts to accommodate various models
- Full extension lightweight polymer drawers, depths 76, 152, 229 & 305mm, (max load 7-20kg per drawer (dependent on drawer size)
- Ergonomic handle with fifth wheel steering provides maximum control when attending an emergency situation
- Two tone colour keyway
- 127mm castors, 2 x Braking

Dimensions (w x d x h)

• Drawers - 514 x 432 x 76 / 152 / 229 / 305mm (internal)

Standard Colours

Body - Code BlueDrawers - Light Taupe/Code Blue

Options

 8FL190 - Drawer Label Holders (x 10)



Drawer Label Holders

| Item | Description | Overall Dimensions (w x d x h) | Weight |
|----------|--|--------------------------------|--------|
| 8LECCRP2 | 1 x 76mm, 3 x 152mm Drawers | 965 x 575 x 984mm | 84.0kg |
| 8LECCRP3 | 2 x 76mm, 1 x 152mm, 1 x 229mm Drawers | 965 x 575 x 984mm | 88.0kg |
| 8LECCRP4 | 2 x 76mm, 1 x 152mm, 1 x 305mm Drawers | 965 x 575 x 1070mm | 87.0kg |



ANTIMICROBIAL



Top Cavity Tray Locking Side Bin Waste Basket (28 Litre) Lockable Sharps Container Adjustable Defibrillator Tray Suction Pump Shelf Backboard (Front Fitting) Backboard (Rear Fitting) Oxygen Cylinder Holder Two Hook Infusion Pole Storage Bin Plastic Security Seals (x 100)

g



| 8LECCRP2 | 8LECCRP3 | 8LECCRP4 |
|--------------|----------|--------------|
| | √ | |
| | 3 | 3 |
| \checkmark | | |
| \checkmark | | |
| \checkmark | ✓ | \checkmark |
| | ✓ | |
| | ✓ | \checkmark |
| \checkmark | | |
| \checkmark | ✓ | \checkmark |



Resuscitation Carts - Polymer, Metro Flexline[™]

Features

- Lightweight polymer construction with smooth rounded corners to assist cleaning
- Incorporates Microban® offering antimicrobial protection
- Upper work surface provides a writing surface/space for preparation work
- Single seal mechanism secures the drawers (security seal compatible)
- Full extension lightweight polymer drawers, depths 76, 152, 229mm (max load 7-20kg per drawer (dependent on drawer size)
- Two tone colour keyway
- 127mm castors, 2 x Directional / 2 x Braking

Dimensions (w x d x h)

- Overall Dimensions 819 x 568 x 914mm
- Drawers 514 x 432 x 76 / 152 / 229mm (internal)

Standard Colours

- Body Light Taupe
- Drawers Light Taupe/Red

Options

- 8FL141 76mm Egg Crate Style Drawer Divider Kit (10 x Vertical, 5 x Horizontal)
- 8FL142 152mm Egg Crate Style Drawer Divider Kit (10 x Vertical, 5 x Horizontal)

8FL21P3369 2 x 76mm, 1 x 152mm, 1 x 229mm Drawers, Red Drawer Pulls, Security Seal Lock 68.0kg

- 8FL315 Two Hook Infusion Pole
- 8FL234 100mm Diameter Oxygen Cylinder Holder
- 8FL190 Drawer Label Holders (x 10)
- 8LEC320 Plastic Security Seals (x 100)
- 8LEC304 Defribillator Platform

Airways / Difficult Intubation Trolley - Mild Steel

Features

- Recessed liquid retaining top
- contents of the drawers (security seal compatible)
- 135mm deep drawers, fitted on full extension glides (max load 30kg per drawer)
- Two accessory rails fitted to L/H side
- Low level corner bumpers to protect the fabric of the building
- 80mm heavy duty castors, 2 off braking

Dimensions (w x d x h)

- Overall 615 x 500 x 915mm
- Drawers 480 x 375 x 135mm (internal)

Standard Colours

- Body White
- Drawers White

Options

Weight



| Item | Description |
|-------|---|
| MW010 | 5 x 125mm Drawers, Non Locking |
| MW015 | 5 x 125mm Drawers, Tamper Evident Panel |

8FL21P3369 (fitted with optional two hook infusion pole, oxygen cylinder holder & defibrillator platform. Defibrillator not included)

Description

 Designed specifically for airways / difficult intubation in conjunction with resuscitation officers • Available as non-locking or fitted with a side mounted tamper evident panel to secure the MW015 Tamper Evident Panel CT/AGH GMBH/15 CT/CR5 & CT/CR2 Weight

55.0kg 57.0kg





Sackholders & Waste Bins

- Removable Body Bins 290
- Free Standing Sackholders 295
- Fixed Body Plastic Sackholders 296
- Removable Body Plastic Sackholders 297



Removable Body Bins - 20 / 50 Litre

Features

- Premium range of removable body bins, fire retardant to HTM83
- Radiused corners to assist cleaning
- Captive drip tray to collect spills
- Positive bag retention through the use of a drop over clamp
- 'Hands free' design with recessed lid discourages opening by hand & encourage use of the pedal
- Stainless steel pedal & linkages
- Plastic feet to keep the unit clear of the floor
- Wheels fitted to rear of bin assist moving when cleaning
- Available with 'Colour Coding' to assist refuse identification, (HTM07)
- Choice of construction...
- Mild steel tubular frame & sheet metal body, painted for extra durability
 'Rust Free' All stainless steel construction, painted for extra durability
- Silent closing versions are fitted with a gas damper to minimise noise when the lid closes

Dimensions (w x d x h)

20 Litre - 267 x 346 x 552mm
50 Litre - 444 x 397 x 551mm

Standard Colours

Paint - Grey Frame/White Body




'Rust Free' versions are constructed from painted stainless steel (frame, body & lid)

20 Litre

| tem | Description | Weight |
|-----------|--|--------|
| 3R015/WHL | Hands Free | 8.5kg |
| 3R020/WHL | Hands Free & Silent Closing | 8.5kg |
| 3R025/WHL | Hands Free & Rust Free | 8.5kg |
| 3R030/WHL | Hands Free, Rust Free & Silent Closing | 8.5kg |
| | | |



Positive bag retention through the use of a drop over clamp



Raised flange & recessed lid used on hands free bins

50 Litre

| ltem | Description |
|-----------|---------------------------------------|
| BR045/WHL | Hands Free |
| BR050/WHL | Hands Free & Silent Closing |
| BR055/WHL | Hands Free & Rust Free |
| BR060/WHL | Hands Free, Rust Free & Silent Closin |

sales@bristolmaid.com +44 (0)1258 484455 www.bristolmaid.com



Gas strut used on silent closing bins





sales@bristolmaid.com +44 (0)1258 484455 www.bristolmaid.com



Removable Body Bins - 75 / 90 Litre

Features

- Premium range of removable body bins, fire retardant to HTM83
- Radiused corners to assist cleaning
- Captive drip tray to collect spills
- Positive bag retention through the use of a drop over clamp
- 'Hands free' design with recessed lid discourages opening by hand & encourage use of the pedal
- Stainless steel pedal & linkages
- Plastic feet to keep the unit clear of the floor
- Wheels fitted to rear of bin assist moving when cleaning
- Available with 'Colour Coding' to assist refuse identification, (HTM07)
 Choice of construction...
- Mild steel tubular frame & sheet metal body, painted for extra durability
 'Rust Free' All stainless steel construction, painted for extra durability
- Silent closing versions are fitted with a gas damper to minimise noise
 when the lid closes

Dimensions (w x d x h)

- 75 Litre 444 x 397 x 781mm
- 90 Litre 444 x 397 x 910mm

Standard Colours

• Paint - Grey Frame/White Body




'Rust Free' versions are constructed from painted stainless steel (frame, body & lid)

75 Litre

| Item | Description | Weight |
|-----------|--|--------|
| BR075/WHL | Hands Free | 15.0kg |
| BR080/WHL | Hands Free & Silent Closing | 15.0kg |
| BR085/WHL | Hands Free & Rust Free | 15.0kg |
| BR090/WHL | Hands Free, Rust Free & Silent Closing | 15.0kg |



Positive bag retention through the use of a drop over clamp



Raised flange & recessed lid used on hands free bins

90 Litre

| Item | Description |
|-----------|---------------------------------------|
| BR105/WHL | Hands Free |
| BR110/WHL | Hands Free & Silent Closing |
| BR115/WHL | Hands Free & Rust Free |
| BR120/WHL | Hands Free, Rust Free & Silent Closin |
| | |

sales@bristolmaid.com +44 (0)1258 484455 www.bristolmaid.com

m



Gas strut used on silent closing bins





sales@bristolmaid.com +44 (0)1258 484455 www.bristolmaid.com

Alternative lid colours are available. Change suffix 'WHL' on the product code with any of the following i.e. Removable Body Bin - 75ltr, Hands free with a yellow lid - BR075/YEL

BLL YTL Blue Yellow & Black

GNL Green

BNL

Brown







Positive bag retention through the use of a drop over clamp



Gas strut used on silent closing bins



Raised flange & recessed lid used on hands free bins

Removable Body Bins - 75 Litre, Front Opening

Features

- Premium range of removable body bins, fire retardant to HTM83
- Radiused corners to assist cleaning
- Captive drip trav to collect spills
- Positive bag retention through the use of a drop over clamp
- 'Hands free' design with recessed lid discourages opening by hand & encourage use of the pedal
- Stainless steel pedal & linkages
- · Plastic feet to keep the unit clear of the floor
- · Wheels fitted to rear of bin assist moving when cleaning
- Available with 'Colour Coding' to assist refuse identification, (HTM07) Choice of construction...
- Mild steel tubular frame & sheet metal body, painted for extra durability
- 'Rust Free' All stainless steel construction, painted for extra durability
- Silent closing versions are fitted with a gas damper to minimise noise when the lid closes

Dimensions (w x d x h)

• 75 Litre - 444 x 397 x 781mm

Standard Colours

Paint - Grey Frame/White Body



'Rust Free' versions are constructed from painted stainless steel (frame, body & lid)

Alternative lid colours are available. Change suffix 'WHL' on the product code with any of the following i.e. Removable Body Bin - 75ltr, Hands Free & Silent Closing with a grey lid - BR140/GYL



| | Description | Weight |
|----|--|--------|
| ΗL | Hands Free | 15.0kg |
| ΗL | Hands Free & Silent Closing | 15.0kg |
| ΗL | Hands Free & Rust Free | 15.0kg |
| ΗL | Hands Free, Rust Free & Silent Closing | 15.0kg |

Features

- Range of freestanding sackholders
- Mild steel tubular/sheet metal frame, painted for extra durability
- Captive drip tray to collect spills
- Positive bag retention through the use of an elasticated cord
- Plastic feet to keep the unit clear of the floor
- Stainless steel pedal and linkages
- Wheels fitted to rear of bin assist moving when cleaning
- Available with 'Colour Coding' to assist refuse identification, (HTM07)

Dimensions (w x d x h)

- 50 Litre 442 x 397 x 551mm
- 75 Litre 442 x 397 x 780mm

Standard Colours

Paint - Grev Frame

| Item | Description |
|-----------|----------------------------|
| BH015/WHL | Elasticated Cord, 50 Litre |
| BH030/WHL | Elasticated Cord, 75 Litre |

Free Standing Sackholders - 50 / 75 Litre



Alternative lid colours are available. Change suffix 'WHL' on the product code with any of the following i.e. Free Standing Bag Holder - 50 Litre, Elastic Cord with brown lid - BH015/BNL



+44 (0)1258 484455

www.bristolmaid.com



Fixed Body Sackholders - Plastic, 20 Litre

Features

- Meets NHS England 'Managing Healthcare Fire Safety' (HTM-05-01)
- Fixed body
- Plastic construction with no metal parts
- Smooth surfaces to assist cleaning
- One piece plastic lid, available in a range of colours to assist refuse identification (HTM07)
- Soft closing lid
- Bag positively secured using a drop over sack retention ring

Dimensions (w x d x h)

• 350 x 220 x 660mm

Standard Colours

 Body - White Lid - White

Options

- Alternative Lid Colours Black/Orange/Yellow/Yellow & Black
- Moulded In Graphics Available on request

Alternative lid colours available. Change suffix 'WH' on the item no. with any of the following, i.e. Fixed Body Sackholder - Plastic, Silent Closing, 20 Litre with a black lid - 5FSH20BL

| WH | BL | OR | YE | YB |
|-------|-------|--------|-----------------|----------------|
| White | Black | Orange | Ye ll ow | Yellow & Black |

Removable Body Sackholders - Plastic, 42 & 70 Litre

Features

- Meets NHS England 'Managing Healthcare Fire Safety' (HTM-05-03)
- One piece removable body
- Plastic construction with no metal parts
- Smooth surfaces to assist cleaning
- One piece plastic lid, available in a range of colours to assist refuse identification (HTM07)
- Soft operating lid on both opening & closing (moving parts enclosed)
- Recessed drip tray to collect unwanted spills
- Bag positively secured using a drop over sack retention ring
- Non-marking rubber feet prevents the bin from sliding
- Wheels to assist when cleaning

Dimensions (w x d x h)

- 42 Litre 460 x 420 x 580mm
- 70 Litre 460 x 450 x 790mm

Standard Colours

- Body White
- Lid White

Options

- Alternative Lid Colours Black/Orange/Yellow/Yellow & Black
- Moulded In Graphics Available on request

Alternative lid colours available. Change suffix 'WH' on the item no. with any of the following, i.e. Removable Body Sackholder - Plastic, Silent Closing, 42 Litre with a black lid - 5RSH42BL



| ltem | Description | Weight |
|----------|---|--------|
| 5RSH42WH | Removable Body Sackholder - Plastic, Silent Closing, 42 Litre | 8.0kg |
| 5RSH70WH | Removable Body Sackholder - Plastic, Silent Closing, 70 Litre | 9.0kg |

| Item | Description | Weight |
|----------|---|--------|
| 5FSH20WH | Fixed Body Sackholder - Plastic, Silent Closing, 20 Litre | 3.0kg |







Seating Solutions

- Office, Task & Operators Chairs 300
- 24 Hour Task & Operator Chairs 302
 - Bariatric Chairs 304
 - Saddle Stools 306
 - Low Back Patient Chairs 308
- Medium & High Back Patient Armchairs 309
 - Reclining Armchairs 310
 - Bariatric Reclining Armchairs 320
 - Reclining Porters Armchairs 322
 - Rise & Recline Porters Armchairs 324
- Seat, Back Cushion & Headrest Options 326
 - Overnight Chair 328
 - Soft Seating 330
 - Visitor & Waiting Room Seating 332
 - Waiting Room Seating 334
 - Tech Chairs & Stools 340
 - PU Tech Chairs & Stools 342
 - Static Safe Tech Chairs & Stools 344
- Static Safe & Sterile Tech Chairs & Stools 346
- Medical Grade Clean Room Tech Chairs 348



Office, Task & Operators Chairs - Century

Features

- Available as...
- Medium back
- High back
- Draughtsman
- Medium & high back available with fixed ring arms & non-braking castors
 Draughtsman has a medium back, foot ring & glides on a 27" base
- Vinyl or fabric upholstery
 Vinyl Fire retardant & antimicrobial
- Fabric Fire retardant
- Two lever mechanism
- Black polypropylene five star base
- Pneumatic seat & back rest height adjustment
 Fully upholstered seat & back with deep CMHR foam
- Maximum capacity 120kg

Replace suffix on the product code with any of the following to select your colour i.e. Medium Back Chair in Aztec - 5160/FAC

Vinyl - Ambla New Richmond in Colour





| Item | Description | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|-------|---------------------------------|--------------------------------|------------------------|--------|
| 5160 | Century, Medium Back | 670 x 670 x 800-1040mm | 440-570mm | 10.0kg |
| 5161 | Century, Medium Back, Ring Arms | 670 x 670 x 800-1040mm | 440-570mm | 12.0kg |
| 5162 | Century, High Back | 670 x 670 x 940-1070mm | 440-570mm | 11.0kg |
| 5163 | Century, High Back, Ring Arms | 670 x 670 x 980-1110mm | 440-570mm | 13.0kg |
| 5160D | Century, Draughtsman | 460 x 605 x 1010-1200mm | 536-686mm | 12.0kg |





| | 5160D |
|--|-------|
| | |
| | |

| 6 | |
|---|---|
| | 4 |
| | |
| | |
| | |
| | - |
| | |
| | |
| Ū | - |
| 6 | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |



24 Hour Task & Operator Chairs - Eon

Features

- BS 5459: Part 2: 2000. Specification for performance requirements & tests for office furniture. Office pedestal seating for use by persons weighing up to 150kg & for use up to 24 hours a day S1 level tested (severe) for heavy duty use
- Certified for use up to 150kg
- Available as...
- Medium back
- High back
- 24/7 Flax fire retardant heavy duty fabric upholstery
- Black polypropylene five-star base
- Easy lift back, ratchet height adjustable
- Independent seat & back mechanism with tension control
- Fully upholstered seat & back with deep CMHR foam
- Optional
- Height adjustable, fold back or ring arms
- 12" or 14" draughtsman's kit
- Maximum capacity 150kg

Replace suffix on the product code with any of the following to select your colour i.e. Eon, Medium Back in Zone - 5260/FYR09

Fabric - Camira 24/7 Flax

| 5 FYR03 FYR21 al Equinox While |
|-----------------------------------|
| |

13.0kg 14.0kg

24 Hour Task & Operator Chairs - Zenith

Features

- BS 5459: Part 2: 2000, Specification for performance requirements & tests for office furniture. Office pedestal seating for use by persons weighing up to 150kg & for use up to 24 hours a day
- S1 level tested (severe) for heavy duty use
- Certified for use up to 150kg
- 24/7 Flax fire retardant heavy duty fabric upholstery
- Black polypropylene five-star base
- Easy lift back, ratchet height adjustable
- Independent seat & back mechanism with tension control & seat slider
- Fully upholstered seat & back with deep CMHR foam
- Chrome swivel height adjustable arms & seat slider
- Maximum capacity 150kg

Options

- 5/Z/HR Headrest
- 5/Z/PLS Pump-Up Lumbar Support

Replace suffix on the product code with any of the following to select your colour i.e. Zenith, High Back in Vernal - 59003/FYR02

Fabric - Camira 24/7 Flax



| 5260 Eon, Medium Back | 710 x 710 x 835-925mm | 470-600mm |
|-----------------------|------------------------|-----------|
| 5262 Eon, High Back | 710 x 710 x 925-1030mm | 470-600mm |

| Item | Description | Overall Dimensions |
|-------|-------------------|-----------------------|
| 59003 | Zenith, High Back | 710 x 710 x 1160-1240 |





While

59003/FYR02/5/Z/HR/5/Z/PLS shown with optional headrest & pump-up lumbar support

21.0kg 0mm 450-530mm



Bariatric High Back Chairs

Features

FCT

FGE

- Deep CMHR foam
- Vinyl or fabric upholstery
- Vinyl Fire retardant & antimicrobial
- Fabric Fire retardant
- Black heavy duty steel five star base
- Black heavy duty gas lift
- Black height adjustable soft touch armrests
- Easy lift ratchet height adjustable back
- Independent seat & back mechanism with side tension control
- Maximum capacity 315kg

Replace suffix on the product code with any of the following to select your colour i.e. Bariatric High Back in Lime - 5660XL/VLI

Vinyl - Ambla New Richmond in Colour



FBK

Bariatric Chairs

Features

- Deep CMHR foam
- Vinyl or fabric upholstery
- Vinyl Fire retardant & antimicrobial
- Fabric Fire retardant
- Steel or beech four leg frames
- Steel Powder coated
- Beech Varnished
- Maximum capacity 315kg

Replace suffix on the product code with any of the following to select your colour i.e. Bariatric, Steel Frame in Lime - 5670/VLI





| Item | Description | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|--------|-------------------|--------------------------------|------------------------|--------|
| 5660XL | High Back, Swivel | 750 x 750 x 1080-1190mm | 425-535mm | 29.0kg |

| Item | Description | Overall Dimensions |
|------|------------------------|---------------------------|
| 5670 | Bariatric, Steel Frame | 920 x 600 x 885mm |
| 5671 | Bariatric, Beech Frame | 844 x 580 x 890mm |







Saddle Stools

Features

- Vinyl or fabric upholstery
 Vinyl Fire retardant & antimicrobial
 - Fabric Fire retardant
- Black polypropylene five star base
- Black gas lift with pneumatic seat height adjustment
- Free wheeling castors
- Maximum capacity 120kg

Replace suffix on the product code with any of the following to select your colour i.e. Saddle Stool in Lime - 5202/VLI

Vinyl - Ambla New Richmond in Colour



Fabric - Camira Advantage



PU Saddle & Perch Stools

Features

- Soft touch black polyurethane seat
 Black polypropylene five star base
 Black gas lift with pneumatic seat height adjustment
 Free wheeling castors
- Maximum capacity 120kg

| tem | Description | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|-----|--------------|--------------------------------|------------------------|--------|
| 202 | Saddle Stool | 670 x 670 x 555-730mm | 555-730mm | 8.0kg |

| ltem | Description | Overall Dimensions (w x d > |
|------|-----------------|-----------------------------|
| 5201 | PU Saddle Stool | 670 x 670 x 590-780mm |
| 5205 | PU Perch Stool | 670 x 670 x 650-730mm |





Alternative colours available at additional cost. Change suffix 'VVDEL' on the product code with any of the following i.e. Alderbury, Low Back Chair in Baltic - 5ND/LC/ALDE/VVBAL

Vinyl - Spradling Valencia C5



5ND/LC/ALDE/VVDEL Alderbury, Low Back Chair, Fixed Seat

Low Back Patient Chairs - Alderbury

BS EN 16139: 2013 Level 2. Furniture. Strength. durability & safety.

• Hygiene gap prevents the trapping of dirt & bacteria between the seat &

Handgrips assist positioning from standing & seated positions

Requirements for non-domestic seating

Solid beech frame (PEFC accredited)

· Plastic glides fitted to chair legs

Maximum capacity 160kg

Dimensions (w x d x h)

• Overall - 580 x 540 x 890mm

Width Between Arms - 470mm

Arm Height - 655mm

Seat - 470 x 430 at Height 470mm

Standard Colours

Vinyl - Delft

back, provides access for cleaning • Fire retardant foam (BS 5852; Part 2: 1982)

Frame finished with anti-bacterial varnish

Contoured seat & back, seat section fixed

Ergonomic design to ensure correct sit & posture

Mortice & tenon joint construction, glued & pinned

Upholstered with fire retardant & antimicrobial vinvl

Features

8.2kg

Medium & High Back Patient Armchairs - Alderbury

Features

- BS EN 16139; 2013 Level 2, Furniture, Strength, durability & safety, Requirements for non-domestic seating
- Eraonomic design to ensure correct sit & posture
- Handgrips assist positioning from standing & seated positions
- Solid beech frame (PEFC accredited)
- Mortice & tenon joint construction, glued & pinned
- Frame finished with anti-bacterial varnish
- Plastic glides fitted to chair legs
- Contoured seat & back
- Fixed Seat is fixed to the frame
- Removable Seat can be removed to assist cleaning (smooth viny) upholstered underside) & is reversible to allow the seat section to be repositioned to provide additional clearance at the front of the chair
- Hygiene gap prevents the trapping of dirt & bacteria between the seat & back, provides access for cleaning
- Fire retardant foam (BS 5852; Part 2: 1982)
- Upholstered with fire retardant & antimicrobial vinyl
- Maximum capacity 160kg

Dimensions (w x d x h)

- Medium Back 660 x 650 x 1000mm
- High Back 660 x 650 x 1195mm
- Seat 450 x 435/475* at Height 500mm
- Width Between Arms 530mm
- Arm Height 655mm
- *Seat depth can be altered by repositioning the removable seat pad

Standard Colours

Vinvl - Delft

Options

HW - Housekeeping Wheels

| Item | Description | Weight |
|--------------------|---|--------|
| Medium Back | | |
| 5ND/MAF/ALDE/VVDEL | Alderbury, Medium Back Armchair, Fixed Seat | 13.6kg |
| 5ND/MAR/ALDE/VVDEL | Alderbury, Medium Back Armchair, Removable Seat | 13.6kg |
| High Back | | |
| 5ND/HAF/ALDE/VVDEL | Alderbury, High Back Armchair, Fixed Seat | 14.4kg |
| 5ND/HAR/ALDE/VVDEL | Alderbury, High Back Armchair, Removable Seat | 14.4kg |

sales@bristolmaid.com +44 (0)1258 484455





Alternative colours available at additional cost. Change suffix 'VVDEL' on the product code with any of the following i.e. Alderbury, High Back Armchair, Removable & Reversible Seat in Olive - 5ND/HAR/ALDE/VVOLI



www.bristolmaid.com

Reclining Armchair - Shapwick

Features

- Available as either manual or electric operation
- Manual Tilt in space recliner
- Electric Single tilt in space reclining only or single tilt in space reclining with rise function
- Manual chairs use a smooth easy to use push operation
- Electric chairs are fitted with a mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- Upholstered with either
- Flame retardant antimicrobial vinyl throughout • A combination of flame retardant antimicrobial vinyl & vapour permeable
- (VP) multiway stretch fabric. VP fabric applied to seat & back sections
- Easiform reflex foam provides comfort when sitting (fixed seat)
- Foam leg rest
- 4 x 50mm swivel braking castors
- Maximum capacity 127kg

Dimensions (w x d x h)

Standard Colours

- Overall Dimensions 737 x 902 x 1130mm Vinyl - Delft Seat - 510 x 510mm at Height 460mm
 - VP Fabric Navy Blue



Alternative colours available at additional cost. Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VPNB' i.e. Shapwick, Reclining Armchair, Electric in Olive Vinyl & Green Vapour Permeable (VP) Multi-Way Stretch Fabric - 5RE/RAE/SHAP/VVOLI/VPGR

Vinyl - Spradling Valencia C5



| Item | Descriptio |
|--------------------------|-------------|
| 5RE/RAM/SHAP/VVDEL/VPNB | Reclining A |
| 5RE/RAE/SHAP/VVDEL/VPNB | Reclining A |
| 5RE/RRAE/SHAP/VVDEL/VPNB | Rise & Rec |

| Item | Description | Weight |
|---------------------|-------------------------------------|--------|
| 5RE/RAM/SHAP/VVDEL | Reclining Armchair, Manual | 41.0kg |
| 5RE/RAE/SHAP/VVDEL | Reclining Armchair, Electric | 55.0kg |
| 5RE/RRAE/SHAP/VVDEL | Rise & Reclining Armchair, Electric | 57.0kg |



CE



Delivery to Chosen Location

- Unwrapping & Removal of Packaging
- Product Positioning & Set-Up
- Full Demonstration of Use to Staff

Fabric - VP Multi-Way Stretch



| | Weight |
|-------------------------|--------|
| nchair, Manual | 41.0kg |
| nchair, Electric | 55.0kg |
| ning Armchair, Electric | 57.0kg |



FIRE

CE

5RE/RAE/ALWE/VVDEL (fitted with optional pocket)

Reclining Armchair - Alweston, Wide Arm

Features

- Ideal for use in oncology, dialysis or haematology for intravenous cannulation, where a patient may be seated for a prolonged period of time, or alternatively in recovery following minor treatment
- Wide supportive cushioned arm (150mm)
- Available as either manual or electric operation Manual - Tilt in space recliner
- Electric Single tilt in space reclining only, or single tilt in space reclining with rise function
- Manual chairs use a smooth easy to use push operation
- Electric chairs are fitted with a mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- Upholstered with either
- Flame retardant antimicrobial vinyl throughout
- A combination of flame retardant antimicrobial vinyl & vapour permeable (VP) multiway stretch fabric. VP fabric applied to seat & back sections
- Easiform reflex foam provides comfort when sitting (fixed seat)
- Foam leg rest
- 4 x 50mm swivel braking castors
- Maximum capacity 127kg

Dimensions (w x d x h)

Standard Colours

- Overall Dimensions 737 x 902 x 1130mm Seat - 510 x 510mm at Height 460mm
 - Vinyl Delft VP Fabric - Navy Blue

Options

- Seat Pressure Management Coolform, Viscoform, Liquiform, Celliform
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- Table with Wireless Charging See page 327
- Full Length Arm Caps Positioned over the arm rests
- AccuPack System Allows mains free chair movements
- USB Charging Socket Located on outside of arm
- · Pocket Positioned on the right-hand side of the chair

| Item | Description | Weight |
|---------------------|---|--------|
| 5RE/RAM/ALWE/VVDEL | Reclining Armchair, Wide Arm, Manual | 41.0kg |
| 5RE/RAE/ALWE/VVDEL | Reclining Armchair, Wide Arm, Electric | 55.0kg |
| 5RE/RRAE/ALWE/VVDEL | Rise & Reclining Armchair, Wide Arm, Electric | 57.0kg |



(fitted with optional pocket)

Alternative colours available at additional cost. Change vinyl suffix 'VVDEL' then change VP fabric suffix VPNB' i.e. Ashington, Reclining Armchair, Waterfall Back, Wide Arm, Electric in Olive Vinyl & Green VP Multi-Way Stretch Fabric - 5RE/RAE/ASHI/VVOLI/VPGR

Vinyl - Spradling Valencia C5



| Item | Description |
|--------------------------|------------------|
| 5RE/RAM/ALWE/VVDEL/VPNB | Reclining Armcha |
| 5RE/RAE/ALWE/VVDEL/VPNB | Reclining Armcha |
| 5RE/RRAE/ALWE/VVDEL/VPNB | Rise & Reclining |

sales@bristolmaid.com +44 (0)1258 484455 www.bristolmaid.com

Fabric - VP Multi-Way Stretch



| | Weight |
|-----------------------------|--------|
| r, Wide Arm, Manual | 41.0kg |
| r, Wide Arm, Electric | 55.0kg |
| rmchair, Wide Arm, Electric | 57.0kg |





- Delivery to Chosen Location
- Unwrapping & Removal of Packaging
- Product Positioning & Set-Up
- Full Demonstration of Use to Staff



Reclining Armchair - Twyford, Waterfall Back

Features

- Available as either manual or electric operation
- Manual Tilt in space recliner
- Electric Single tilt in space reclining only or single tilt in space reclining with rise function
- Manual chairs use a smooth easy to use push operation
- Electric chairs are fitted with a mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- · Upholstered with either
- Flame retardant antimicrobial vinyl throughout
- A combination of flame retardant antimicrobial vinyl & vapour permeable (VP) multiway stretch fabric. VP fabric applied to seat & back sections
- Easiform reflex foam provides comfort when sitting (fixed seat)
- Foam leg rest
- Fixed deep padded waterfall back
- 2 x 50mm swivel braking castors
- Maximum capacity 127kg

Dimensions (w x d x h)

Standard Colours

- Overall Dimensions 765 x 865 x 1145mm
 Vinyl Delft Seat - 520 x 510mm at Height 482mm
 - VP Fabric Navy Blue

Options

- Seat Pressure Management Coolform, Viscoform, Liquiform
- · Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- Full Length Arm Caps Positioned over the arm rests
- AccuPack System Allows mains free chair movements



| IVERY mainland only | Unwrapping & Removal of Packaging Product Positioning & Set-Up Full Demonstration of Use to Staff |
|------------------------|---|
| | |

| Item | Description | Weight |
|---------------------|---|--------|
| 5RE/RAM/TWYF/VVDEL | Reclining Armchair, Waterfall Back, Manual | 43.0kg |
| 5RE/RAE/TWYF/VVDEL | Reclining Armchair, Waterfall Back, Electric | 57.0kg |
| 5RE/RRAE/TWYF/VVDEL | Rise & Recline Armchair, Waterfall Back, Electric | 59.0kg |



5RE/RAM/TWYF/VVMEL/VPGR

Alternative colours available at additional cost. Change vinyl suffix "VVDEL' then change VP fabric suffix 'VPNB' i.e. Twyford, Reclining Armchair, Electric in Olive Vinyl & Green Vapour Permeable (VP) Multi-Way Stretch Fabric - 5RE/RAE/TWYF/VVOLI/VPGR

Vinyl - Spradling Valencia C5



| Item | Description | |
|--------------------------|---|--------|
| 5RE/RAM/TWYF/VVDEL/VPNB | Reclining Armchair, Waterfall Back, Manual | 43.0kg |
| 5RE/RAE/TWYF/VVDEL/VPNB | Reclining Armchair, Waterfall Back, Electric | 57.0kg |
| 5RE/RRAE/TWYF/VVDEL/VPNB | Rise & Recline Armchair, Waterfall Back, Electric | 59.0kg |

Fabric - VP Multi-Way Stretch





Reclining Armchair - Ashington, Waterfall Back, Wide Arm

Features

- · Ideal for use in oncology, dialysis or haematology for intravenous cannulation, where a patient may be seated for a prolonged period of time, or alternatively in recovery following minor treatment
- Wide supportive cushioned arm (150mm)
- Available as either manual or electric operation
- Manual Tilt in space recliner
- Electric Single tilt in space reclining only, or single tilt in space reclining with rise function
- Manual chairs use a smooth easy to use push operation
- Electric chairs are fitted with a mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- Upholstered with either
- Flame retardant antimicrobial vinyl throughout
- A combination of flame retardant antimicrobial vinyl & vapour permeable (VP) multiway stretch fabric. VP fabric applied to seat & back sections

Standard Colours

VP Fabric - Navy Blue

Vinvl - Delft

- Easiform reflex foam provides comfort when sitting (fixed seat)
- Foam leg rest
- Fixed deep padded waterfall back
- 2 x 50mm swivel braking castors
- Maximum capacity 127kg

Dimensions (w x d x h)

- Overall Dimensions 765 x 865 x 1145mm
- Seat 520 x 510mm at Height 482mm

Options

- · Seat Pressure Management Coolform, Viscoform, Liquiform
- · Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- Table with Wireless Charging See page 327
- Full Length Arm Caps Positioned over the arm rests
- AccuPack System Allows mains free chair movements
- USB Charging Socket Located on outside of arm
- · Pocket Positioned on the right-hand side of the chair

| tem | Description | Weight |
|--------------------|---|--------|
| RE/RAM/ASHI/VVDEL | Reclining Armchair, Waterfall Back, Wide Arm, Manual | 43.0kg |
| RE/RAE/ASHI/VVDEL | Reclining Armchair, Waterfall Back, Wide Arm, Electric | 57.0kg |
| RE/RRAE/ASHI/VVDEL | Rise & Recline Armchair, Waterfall Back, Wide Arm, Electric | 59.0kg |



(fitted with optional pocket)

Alternative colours available at additional cost. Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VPNB' i.e. Alweston, Reclining Armchair, Wide Arm, Electric in Olive Vinyl & Green Vapour Permeable VP Multi-Way Stretch Fabric - 5RE/RAE/ALWE/VVOLI/VPGR

Vinyl - Spradling Valencia C5



| | Description | Weight |
|--------------------------|---|--------|
| 5RE/RAM/ASHI/VVDEL/VPNB | Reclining Armchair, Waterfall Back, Wide Arm, Manual | 43.0kg |
| 5RE/RAE/ASHI/VVDEL/VPNB | Reclining Armchair, Waterfall Back, Wide Arm, Electric | 57.0kg |
| 5RE/RRAE/ASHI/VVDEL/VPNB | Rise & Recline Armchair, Waterfall Back, Wide Arm, Electric | 59.0kg |



FIRE

Fabric - VP Multi-Way Stretch





 Product Positioning & Set-Up • Full Demonstration of Use to Staff



Reclining Armchair - Beckington, Waterfall Back

Features

- Electric reclining only, or dual tilt in space reclining with rise function
- Mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Upholstered with flame retardant vinyl
- · Easiform reflex foam provides comfort when sitting
- Foam leg rest
- Fixed waterfall back
- Side pocket & loop for handset
- 2 x 50mm swivel braking castors
- Maximum capacity 127kg

Dimensions (w x d x h)

- Overall Dimensions 840 x 930 x 1160mm
- Seat 510 x 510mm at Height 483mm
- Back Height 710mm
- Arm Height 180mm

Standard Colours

Vinyl - Ultra



 Delivery to Chosen Location Unwrapping & Removal of Packaging
 Product Positioning & Set-Up Full Demonstration of Use to Staff



Alternative colours available at additional cost. Change suffix 'VMULT' on the product code with any of the following i.e. Beckington, Reclining Armchair, Waterfall Back, Electric in Garden - 5RE/RAE/BECK/VMGAR

Memphis Faux Leather



| Item | Description | Weight |
|---------------------|---|--------|
| 5RE/RAE/BECK/VMULT | Reclining Armchair, Waterfall Back, Electric | 58.0kg |
| 5RE/RRAE/BECK/VMULT | Rise & Recline Armchair, Waterfall Back, Electric | 62.0kg |

Reclining Armchair - Hadleigh, Waterfall Back, Bariatric



Features

- Specifically designed for comfort & pressure management control
- Dual standard rise & recliner, provides a full range of posture positive recline positions as well as support to a standing position
- Mains powered actuator (230V ~. 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset
- Upholstered with a combination of flame retardant antimicrobial vinyl & vapour permeable (VP) multiway stretch fabric. VP fabric applied to seat & back sections
- · Easiform reflex foam provides comfort when sitting
- Seat & leg sections generously padded with memory foam & fitted with a removable cover
- Separate waterfall back cushions can be individually repositioned & adjusted for comfort & support
- Back blockers provided for additional lumbar support
- Flat wide arms for deep comfort & support
- Loop for handset & deep side pocket
- Heavy duty glides
- Maximum capacity 254kg

Dimensions (w x d x h)

- Chair 1000 / 1100 x 940 x 1145mm
- Seat 610 / 710mm x 510 x 460mm

Standard Colours

- Vinvl Delft
- VP Fabric Navy Blue

Options

- Seat Pressure Management Coolform, Viscoform
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- USB Charging Socket Located on outside of arm
- Full Length Arm Caps Positioned over the arm rests





Alternative colours available at additional cost. Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VPNB' i.e. Hadleigh, Rise & Reclining Armchair, Waterfall Back, Bariatric, Electric, 710mm in Olive Vinyl & Green VP Multi-Way Stretch Fabric - 5RE/RRAE/HADL710/VVOLI/VPGR

Vinyl - Spradling Valencia C5



| Item | Description | Weight |
|-----------------------------|---|--------|
| 5RE/RRAE/HADL610/VVDEL/VPNB | Rise & Reclining Armchair, Waterfall Back, Bariatric, Electric, 610mm | 88.0kg |
| 5RE/RRAE/HADL710/VVDEL/VPNB | Rise & Reclining Armchair, Waterfall Back, Bariatric, Electric, 710mm | 88.0kg |



VPBL

VPRE

VPGR

Reclining Porters Armchair - Bathampton, Waterfall Back



Features

- Single tilt in space recliner, provides a full range of posture positive recline positions • Mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied
- with a handset (positioned at rear of chair/optional side of chair)
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- Upholstered with either
- Flame retardant antimicrobial vinyl throughout
- A combination of flame retardant antimicrobial vinyl & vapour permeable (VP) multi-way stretch fabric. VP fabric applied to seat & back sections
- · Upholstered seat & back sections are fitted with removable covers, allows seat section to be easily upgraded to aid pressure care
- Removable internal Easiform reflex seat cushion
- Seat & leg sections generously padded with memory foam & fitted with a removable cover
- · Separate waterfall back cushions can be individually repositioned & adjusted for comfort & support
- · Push handle at rear of the chair to assist when positioning/portering
- Slide out footplate to support feet
- 4 x 75mm swivel braking castors
- Maximum capacity 127kg

Dimensions (w x d x h)

- Overall Dimensions 685 x 990 x 1220mm
- Seat 520 x 510mm at Height 470mm

Standard Colours

- Vinyl Delft
- VP Fabric Navy Blue

Options

- Seat Pressure Management Coolform, Viscoform, Liquiform, Celliform, Airform, Dvnaform
- Back Cushions Comfort Lateral, Support Lateral, Large Profile with Waterfall, Large Profile Headrest with Comfort Lateral, Large Profile Headrest with Support Lateral
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- Full Length Arm Caps Positioned over the arm rests
- AccuPack System Allows mains free chair movements

| Item | Description | Weight |
|---------------------|--|--------|
| 5RE/RPAH/BATH/VVDEL | Reclining Porters Armchair, Waterfall Back, Electric | 75.0kg |



Alternative colours available at additional cost. Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VPNB' i.e. Bathampton, Waterfall Back, Handset Operated in Olive Vinyl & Green VP Multi-Way Stretch Fabric - 5RE/RPAH/BATH/VVOLI/VPGR

Vinyl - Spradling Valencia C5



| 5RE/RPAH/BATH/VVDEL/VPNB | Reclining Porters Ar |
|--------------------------|----------------------|

- Delivery to Chosen Location
- Unwrapping & Removal of Packaging
- Product Positioning & Set-Up
- Full Demonstration of Use to Staff

Fabric - VP Multi-Way Stretch



Armchair, Waterfall Back, Electric 75.0kg



Rise & Recline Porters Armchair - Dorchester, Waterfall Back



Features

- Dual tilt in space rise & recliner, provides a full range of posture positive recline positions as well as support to a standing position
- Mains powered actuator (230V ~, 50/60Hz) fitted with a UK 3 pin plug & supplied with a handset (positioned at side of chair)
- Optional AccuPack system allows the chair to be used without being connected to an electrical supply (average of 100 cycles based on a 24hr charge)
- Upholstered with either
- Flame retardant antimicrobial vinyl throughout
- A combination of flame retardant antimicrobial vinvl & vapour permeable (VP) multi-way stretch fabric. VP fabric applied to seat & back sections
- Upholstered seat & back sections are fitted with removable covers, allows seat section to be easily upgraded to aid pressure care
- Removable internal Easiform reflex seat cushion
- Seat & leg sections generously padded with memory foam & fitted with a removable cover
- Separate waterfall back cushions can be individually repositioned & adjusted for comfort & support
- Push handle at rear of the chair to assist when positioning/portering
- Optional sliding footplate to support feet
- 4 x 100mm swivel braking castors
- Maximum capacity 127kg

Dimensions (w x d x h)

Standard Colours

- Overall Dimensions 775 x 1055 x 1245mm Vinvl - Delft Seat - 510 x 510mm at Height 480mm
 - VP Fabric Navy Blue

Options

- Sliding Footplate
- Seat Pressure Management Coolform, Viscoform, Liquiform, Celliform, Airform, Dynaform
- Back Cushions Comfort Lateral, Support Lateral, Large Profile with Waterfall, Large Profile Headrest with Comfort Lateral, Large Profile Headrest with Support Lateral
- Headrests Head Pillow, Head Roll, Small Profile, Large Profile, Midline
- Full Length Arm Caps Positioned over the arm rests
- AccuPack System Allows mains free chair movements



Alternative colours available at additional cost.

Vinyl - Spradling Valencia C5

ÉLIVER

| Description |
|------------------|
| Rise & Recline I |

| tem | Description | Weigh |
|---------------------|---|-------|
| RE/RRPAH/DORC/VVDEL | Rise & Recline Porters Armchair, Waterfall Back, Electric | 76.5k |

5RE/RRPAH/DORC/VVBAL/VPNB (with optional footplate)

Change vinyl suffix 'VVDEL' then change VP fabric suffix 'VPNB' i.e. Dorchester, Waterfall Back, Handset Operated in Leinen Vinyl & Grev VP Multi-Way Stretch Fabric - 5RE/RPAH/DORC/VVOLI/VPGY

Fabric - VP Multi-Way Stretch




5RE/RRPAH/DORC/VVDEL/VPNB Rise & Recline Porters Armchair, Waterfall Back, Electric 76.5kg




Seat Options For Pressure Management

Cushions are covered in multi-way stretch, vapour permeable fabric which is wipeable, removable, washable at 95°C, incontinence proof & breathable.



For those sitting for long periods with specific comfort issues, such as a high risk of pressure sores. Uses the latest medical grade silicone gel technology allowing pressure to be evenly distributed.

The gel conforms to the shape of the body, equalises to body temperature, allows air to flow easily & is totally odourless & hypoallergenic.



Celliform




Pressure risk levels detailed are for guidance only. Each individual user has different posture, medical needs & health issues which contribute

to the performance of any pressure cushion. If in doubt, please consult a GP or Health Care Professional.


For those sitting for long periods of time with specific comfort issues, such as a high risk of pressure sores & sheer forces.

Gelacell's unique grid design eliminates pressure points, contouring in response to the transfer of weight, evenly distributes pressure & reduces heat absorbtion, keeping the seat cooler & fresher.



For those who have good posture & mobility, but are seated for periods longer than one hour, require pressure relief & supported comfort.

High density foam provides support, whilst a top layer of memory foam allows submersion into the cushion, providing a large surface area to reduce pressure levels.



Pressure Risk HIGH


For those who have good posture, but are unable to alter their position easily & require high level pressure care.

This cushion provides the benefits of an air cushion, however without the sensation of alternating changes in air pressure.

Suitable for users 40 - 127kg

For those who have good posture & mobility, but are seated for long periods of time, require pressure relief & supported comfort. Excellent for people with little fatty tissue.

High density foam provides support, topped with layers of fluid liquid gel & memory foam to provide high levels of pressure relief. The liquid gel helps to dissipate heat, whilst constantly adapting to movement.



Dynaform Pressure Risk VERY HIGH

For those at very high risk of developing pressure sores, or have a history of sores, where the priority is to prevent further pressure damage.

Dynaform is a dynamic alternating air cushion, which constantly changes air levels within its cells to recreate the effect of pressure relieving & prevent pressure build up.

Air cell features

- 6 dynamic cells 10cm wide
- 3 cell cycle repeated every 9 minutes

Suitable for users 40 - 121kg.

Please note the Dynaform cushion does need an electric supply at all times to operate.

Back Cushion Options



Waterfall Back Cushion Individual waterfall cushions, fully adjustable to accommodate maximum postural support



Large Profile Headrest with Waterfall Back

For those with limited head control & requiring a soft adjustable lumbar support

Headrest Options






Head Pillow

Head Roll



Comfort Lateral Back Cushion Soft fibre filled removable lateral cushion for lower back & side lumber support, without feeling restricted



Large Profile Headrest with Comfort Lateral Back

For those with limited head control & requiring soft lateral support without feeling restricted



Support Lateral Back Cushion A foam filled back providing firmer lateral support with a soft fibre filled headrest



Large profile headrest with Support Lateral Back

For those with limited head control, requiring firm foam lateral lumbar & midline support

Table with Wireless Charging

360-degree rotation & angle adjustable table with a wipe clean surface. Supports QI wireless charging & fitted with a USB socket. Removable magnetic bar provides support & ensures correct positioning on the wireless receiver coil, wipe clean surface (w 345 x h 245mm)





Small Profile Headrest



Large Profile Headrest



Midline Headrest

Overnight Chair - Portland, Three Position



Features

- Multi-functional compact chair that provides comfortable seating by day, transforming into a bed by night
- Suitable for use in maternity units, children's wards, relative's rooms, hospices & care homes
- Provides visitors with a comfortable place to sleep during longer visits
- Chair is pulled out in three motions to a bed
 - Chair
 - Lounge Chair
 - Bed
- Upholstered using flame retardant antimicrobial vinyl
- · Easiform reflex foam provides comfort when sitting or sleeping
- 50mm swivel braking castors
- Maximum capacity 160kg

Dimensions (w x d x h)

- Overall Dimensions 870 x 840 x 870mm (when open overall bed length is 2060mm)
- Seat 690 x 540mm at Height 410mm

Standard Colours

- Vinyl Delft
- Vinyl & Patterned Fabric Vinyl Lavendel & Savanna Lavender



Delivery to Chosen Location
Unwrapping & Removal of Packaging
Product Positioning & Set-Up
Full Demonstration of Use to Staff





Alternative colours available at additional cost. Change suffix 'VVDEL' on the product code with any of the following i.e. Portland, Three Position in Olive - 5RE/OC/PORT/VVOLI

Vinyl - Spradling Valencia C5



| Item | Description | Weight |
|-----------------------|---|--------|
| 5RE/OC/PORT/VVDEL | Three Position, Vinyl Delft | 61.5kg |
| 5RE/OC/PORT/VVLAVSAVL | Three Position, Vinyl Lavendel & Savanna Lavender | 61.5kg |





Replace suffix on the product code with any of the following to select your colour i.e. Single Unit in Cobalt - 5ZETA1/FCT

Vinyl - Ambla New Richmond in Colour

FAC

FGE

FRD

FCT



FKR

FBK



Soft Seating - Zeta, Modular

Features

- BS EN 16139: 2013 Level 2. Furniture. Strength, durability & safety. Requirements for non-domestic seating
- Zeta is a stylish, fully modular range of soft seating based around eight standard elements suitable for use in seating environments, ranging from reception areas, breakout spaces to waiting areas & cafeterias
- FSC® certified wood core, constructed from 18mm boards
- Traditional webbed seat for additional comfort & durability
- High density CMHR foam
- Vinyl or fabric upholstery
- Vinyl Fire retardant & antimicrobial
- Fabric Fire retardant
- Chrome finish legs
- Arms available on single & double units with backs
- · Backs are available with a raised feature to allow fluid run off (overall height 810mm)
- Maximum capacity
- Single 146kg

Options

- RA Right Arm
- LA Left Arm
- 2A Left & Right Arms

| | Description | Overall Dimensions (w x d x h) | Seat Height | Seat Depth | |
|-------------|---|--------------------------------|-------------|------------|--------|
| 5ZETA1 | Single Unit | 650 x 650 x 415mm | 415mm | 650mm | 12.0kg |
| 5ZETA2 | Double Unit | 1300 x 650 x 415mm | 415mm | 650mm | 22.0kg |
| 5ZETA1/B | Single Unit with Back | 650 x 650 x 765mm | 415mm | 485mm | 16.5kg |
| 5ZETA2/2B | Double Unit with Back | 1300 x 650 x 765mm | 415mm | 485mm | 31.0kg |
| 5ZETA1/CO | Single Unit with Corner | 650 x 650 x 765mm | 415mm | 485mm | 16.5kg |
| 5ZETA2/B/RC | Double Unit with Right Corner & Left Back | 1300 x 650 x 765mm | 415mm | 485mm | 35.5kg |
| 5ZETA2/B/LC | Double Unit with Left Corner & Right Back | 1300 x 650 x 765mm | 415mm | 485mm | 35.5kg |
| 5ZETAQ | Quarter Segment | 650 x 650 x 415mm | 415mm | 650mm | 12.0kg |
| 5ZETAH | Half Segment | 1300 x 650 x 415mm | 415mm | 650mm | 23.0kg |
| 5ZETAH/B | Half Segment with Back | 1300 x 650 x 765mm | 415mm | 485mm | 25.0kg |
| 5ZETAC45 | 45 Degree Curve | 950 x 650 x 415mm | 415mm | 650mm | 11.5kg |
| 5ZETAC45/B | 45 Degree Curve with Back | 950 x 650 x 765mm | 415mm | 485mm | 16.5kg |
| 5ZETAC90 | 90 Degree Curve | 1625 x 650 x 415mm | 415mm | 650mm | 22.0kg |
| 5ZETAC90/B | 90 Degree Curve with Back | 1625 x 650 x 765mm | 415mm | 485mm | 32.0kg |
| 5ZETA930C | Single Corner Unit | 930 x 650 x 765mm | 415mm | 650mm | 28.0kg |
| | | | | | |





Visitor & Waiting Room Seating - Mata Chairs

Features

- BS EN 1729-1: Part 1 & 2: 2015 (functional dimensions & safety requirements)
- BS 4875: Part 1: 2001 Test Level 4 (strength & durability)
- BS 7945: 1999 (stability)
- Available with
- Four leas
- Sled base
- Moulded polypropylene shell available in seven colours
- Silver finish four lea
- Chrome finish sled base
- Four leg chair stackable five high
- Maximum capacity 120kg

Replace suffix on the product code with any of the following to select your colour i.e. Four Leg Chair in Olive - 5SM6L/SOE

Plastic Shell



Visitor & Waiting Room Seating - Prima Chairs

Features

- BS EN 16139: 2013 Level 2 (strength, durability & safety)
- BS 4875: Part 1: 2001 Test Level 4 (strength & durability)
- BS 7945: 1999 (stability)
- Certified for use up to 203kg
- Choice of vinyl, fabric upholstery or polypropylene seat & back Vinyl - Fire retardant & antimicrobial
- Fabric Fire retardant
- Stable four leg design
- Black finish tubular frame
- Stacks up to ten high (transporter available)
- Built-in linking system (optional extension links)
- Maximum capacity 203kg

Replace suffix on the product code with any of the following to select your colour i.e. Four Leg Chair in Cobalt - 5301/FCT

Vinyl - Ambla New Richmond in Colour



| ltem | Description | Overall Dimensions (w x d x h) | Seat Height | Weight |
|-------|-------------|--------------------------------|-------------|--------|
| 5SM6L | Four Leg | 500 x 490 x 805mm | 460mm | 5.0kg |
| 5SL6 | Sled Base | 510 x 570 x 805mm | 460mm | 7.0kg |

| ltem | Description | Overall Dimensions (|
|--------|---------------------|----------------------|
| 5301 | Four Leg | 520 x 580 x 800mm |
| 5301/A | Four Leg, Arm Rests | 620 x 580 x 800mm |

485mm

485mm





Waiting Room Seating - Nova, Tub

Features

- Nova single & double tub chairs are suitable for a wide range of seating environments, ranging from cafeterias to waiting areas
- FSC[®] certified wood core with moulded deep CMHR seat cushions
- Vinyl or fabric upholstery
- Vinvl Fire retardant & antimicrobial
- Fabric Fire retardant
- Chrome finish leas
- Maximum capacity
- Single 146kg
- Double 292kg

Replace suffix on the product code with any of the following to select your colour i.e. Double Seat in Cobalt - 5NOVA/3/FCT

Vinyl - Ambla New Richmond in Colour





| Item | Description | Overall Dimensions (w x d x h) | Seat Height | Weight |
|---------|-------------------------|--------------------------------|-------------|--------|
| 5NOVA/1 | Nova - Tub, Single Seat | 680 x 680 x 800mm | 450mm | 16.0kg |
| 5NOVA/2 | Nova - Tub, Double Seat | 1360 x 680 x 800mm | 450mm | 31.0kg |

Waiting Room Seating - Horizon, Beech Frame

Features

- Horizon chairs, sofas & tables are designed to suit a variety of waiting & occasional seating areas
- Beech frame with plywood seat shells
- Arms assist positioning from standing & seated positions
- Vinyl or fabric upholstery
- Vinyl Fire retardant & antimicrobial
- · Fabric Fire retardant
- Fixed seat, with a loose back cushion to assist cleaning
- Back cushion has a zipped cover, allowing for easy replacement if damaged
- Maximum capacity

FGE

FCT

- Single chair 146kg
- Two seater sofa 120kg per seat
- Three seater sofa 120kg per seat

Replace suffix on the product code with any of the following to select your colour i.e. Two Seater Sofa with Arms in Cobalt - 5HORI/2/A/FCT





FBK

| ltem | Description | Overall Dimensions (w x d x h) | Seat Height | Weight |
|-----------|-----------------------------------|--------------------------------|-------------|--------|
| 5HORI/1/A | Horizon - Easy Chair, Arms | 750 x 650 x 800mm | 450mm | 19.0kg |
| 5HORI/2/A | Horizon - Two Seater Sofa, Arms | 1350 x 650 x 800mm | 450mm | 33.0kg |
| 5HORI/3/A | Horizon - Three Seater Sofa, Arms | 1950 x 650 x 800mm | 450mm | 47.0kg |





5HORI/2/A/FCT

FIRE RETARDAN



Waiting Room Seating - Prima Beam

Features

- Prima beam is suitable for all environments requiring durable & robust waiting room seating solutions
- Certified for use up to 32 stone (203kg) per seat
- Available as...
- Fixed Seat
- Tip Up Seat (must be fixed to the floor)
- Choice of vinyl, fabric upholstery or polypropylene seat & back
 Vinyl Fire retardant & antimicrobial

Replace suffix on the product code with any of the following to select your colour i.e. Fixed Beam Seat, Two Seat in Cobalt - 5PF332/FCT

- Fabric Fire retardant
- Black finish beam & legs
- Drilled for floor fixing (fixings not supplied)
- Maximum capacity 203kg per seat

Options

• 510065 - Armrest (x1)

5PT338/FRD Shown with 510065 (optional arms sold individually) - contact us for details



Fixed Seat Beam

| | Description | Overall Dimensions (w x d x h) | Seat Height | Weight |
|--------|-------------|--------------------------------|-------------|--------|
| 5PF332 | Two Seat | 1200 x 560 x 790mm | 480mm | 20.0kg |
| 5PF333 | Three Seat | 1800 x 560 x 790mm | 480mm | 25.0kg |
| 5PF334 | Four Seat | 2400 x 560 x 790mm | 480mm | 30.0kg |
| 5PF335 | Five Seat | 3000 x 560 x 790mm | 480mm | 35.0kg |





Tip-Up Seat Beam

| Item | Description | Overall Dimensions (w x d x h) | Seat Height | Weight |
|--------|-------------|--------------------------------|-------------|--------|
| 5PT336 | Two Seat | 1200 x 560 x 480mm | 490mm | 20.0kg |
| 5PT337 | Three Seat | 1800 x 560 x 480mm | 490mm | 25.0kg |
| 5PT338 | Four Seat | 2400 x 560 x 480mm | 490mm | 30.0kg |
| 5PT339 | Five Seat | 3000 x 560 x 480mm | 480mm | 35.0kg |



Waiting Room Seating - Voyage, Beam

Features

- BS EN 16139: 2013 Level 2. Furniture. Strength, durability & safety. Requirements for non-domestic seating
- Voyage beam is suitable in heavy usage waiting room environments
- Cast aluminium polished leg frames
- One piece metal seat face with cast aluminium polished edges
- Optional...
- Floor fixing plates
- Upholstered seat & back pads
- Arms
- Tables Maple finish, rectangular / round with swivel linking joint
- Maximum capacity 203kg per seat

Options

- 510046 Floor Fixing Plates Mild Steel, Painted
- 510068 Seat Pad Upholstered, Vinyl
- 510069 Back Pad Upholstered, Vinyl
- 5VOYAGE ARM Arm Cast Aluminium, Polished
- 5VOYAGECT Table Round, Maple Finish with Swivel Linking Joint
- 5VOYAGERT Table Rectangle, Maple Finish

Seat & Back Pad Vinyl

Replace suffix on the product code with any of the following to select your colour i.e. Beam, Two Seat, Seat Pad in Bristol Blue - 5VOYAGE2/510068/VBB

Vinyl - Ambla New Richmond in Colour



| Item | Description | Overall Dimensions (w x d x h) | Seat Height | Weight |
|----------|---------------------------|--------------------------------|-------------|--------|
| 5VOYAGE2 | Voyage - Beam, Two Seat | 1300 x 690 x 770mm | 345mm | 27.5kg |
| 5VOYAGE3 | Voyage - Beam, Three Seat | 1900 x 690 x 770mm | 345mm | 38.5kg |
| 5VOYAGE4 | Voyage - Beam, Four Seat | 2473 x 690 x 770mm | 345mm | 49.5kg |





5VOYAGERT Table - Rectangle, Maple Finish

5VOYAGECT Table - Round, Maple Finish with Swivel Linking Joint



Tech Chairs

Features

- Vinyl or fabric upholstery
- Vinyl Fire retardant & antimicrobial
- Fabric Fire retardant
- Black polypropylene five star base
- Models
- Medium/High Foot ring & glides
- Low Free wheeling castors
- Black gas lift with pneumatic seat height adjustment
- Two lever mechanism
- Maximum load 120kg

Options

• 5BLSTC - Brake loaded castors

Replace suffix on the product code with any of the following to select your colour i.e. High Chair in Peach - 5TC101/VPE

Vinyl - Ambla New Richmond in Colour



Fabric - Camira Advantage



| Item | Description | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|--------|-------------|--------------------------------|------------------------|--------|
| 5TC103 | Low | 670 x 670 x 890-1090mm | 440-580mm | 13.0kg |
| 5TC102 | Medium | 710 x 710 x 950-1210mm | 520-710mm | 14.0kg |
| 5TC101 | High | 710 x 710 x 980-1320mm | 570-820mm | 14.0kg |

Tech Stools

Features

- Vinyl or fabric upholstery
- Vinyl Fire retardant & antimicrobial
 Fabric Fire retardant
- Black polypropylene five star base
- Models
 - Medium/High Foot ring & glides
 - Low Non-braking castors
- Black gas lift with pneumatic seat height adjustment
- Maximum load 120kg

Options

5BLSTC - Brake loaded castors

Replace suffix on the product code with any of the following to select your colour i.e. Low Stool in Kingfisher - 5TC110/FKR

Vinyl - Ambla New Richmond in Colour



| Item | Description | Overall Dimensions (w x d x l |
|--------|-------------|-------------------------------|
| 5TC110 | Low | 670 x 670 x 440-570mm |
| 5TC109 | Medium | 670 x 670 x 570-690mm |
| 5TC108 | High | 670 x 670 x 650-810mm |





| Item | Description | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|------|----------------------------|--------------------------------|------------------------|--------|
| 5PU6 | Low, Castors | 670 x 670 x 750-980mm | 430-570mm | 13.0kg |
| 5PU5 | Medium, Foot Ring & Glides | 710 x 710 x 830-1120mm | 510-700mm | 14.0kg |
| 5PU4 | High, Foot Ring & Glides | 710 x 710 x 910-1270mm | 580-840mm | 14.0kg |

| ltem | Description | Overall Dimer |
|--------|----------------------------|----------------|
| 5PU153 | Low, Castors | 670 x 670 x 41 |
| 5PU152 | Medium, Foot Ring & Glides | 670 x 670 x 47 |
| 5PU151 | High, Foot Ring & Glides | 670 x 670 x 53 |





| | Seat Height Adjustment | Weight |
|----------|------------------------|--------|
| 10-540mm | 410-540mm | 8.0kg |
| 70-660mm | 470-660mm | 10.0kg |
| 30-780mm | 530-780mm | 10.0kg |



| Item | Description | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|--------|----------------------------|--------------------------------|------------------------|--------|
| 5TC203 | Low, Castors | 670 x 670 x 890-1090mm | 440-580mm | 13.0kg |
| 5TC202 | Medium, Foot Ring & Glides | 710 x 710 x 950-1210mm | 520-710mm | 14.0kg |
| 5TC201 | High, Foot Ring & Glides | 710 x 710 x 980-1320mm | 570-820mm | 14.0kg |

Static Safe Tech Stools

Features

- Hi-Tech ESD (anti-static) fire retardant black fabric upholstery
- Aluminium five star base
- Models
- Medium/High Foot ring & glides
- Low ESD free wheeling castors
- Chrome gas lift with pneumatic seat height adjustment
 Supplied with ESD (anti-static) castors or glides
- Maximum capacity 120kg

Options

Brake loaded ESD castors



| Item | Description | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|--------|----------------------------|--------------------------------|------------------------|--------|
| 5TC210 | Low, Castors | 670 x 670 x 386-556mm | 386-556mm | 10.0kg |
| 5TC209 | Medium, Foot Ring & Glides | 670 x 670 x 490-680mm | 490-680mm | 11.0kg |
| 5TC208 | High, Foot Ring & Glides | 670 x 670 x 540-790mm | 540-790mm | 11.0kg |





| Item | Description | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|--------|----------------------------|--------------------------------|------------------------|--------|
| 5TC306 | Low, Castors | 670 x 670 x 890-1090mm | 440-580mm | 13.0kg |
| 5TC305 | Medium, Foot Ring & Glides | 710 x 710 x 950-1210mm | 520-710mm | 14.0kg |
| 5TC304 | High, Foot Ring & Glides | 710 x 710 x 980-1320mm | 570-820mm | 14.0kg |

Static Safe & Sterile Tech Stools

- black vinyl upholstery

- adjustment



| Item | Description | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|--------|----------------------------|--------------------------------|------------------------|--------|
| 5TC310 | Low, Castors | 670 x 670 x 386-556mm | 386-556mm | 10.0kg |
| 5TC309 | Medium, Foot Ring & Glides | 670 x 670 x 490-680mm | 490-680mm | 11.0kg |
| 5TC308 | High, Foot Ring & Glides | 670 x 670 x 540-790mm | 540-790mm | 11.0kg |



| ltem | Description | Overall Dimensions (w x d x h) | Seat Height Adjustment | Weight |
|--------|----------------------------|--------------------------------|------------------------|--------|
| 5TC313 | Low, Castors | 670 x 670 x 890-1090mm | 440-580mm | 13.0kg |
| 5TC312 | Medium, Foot Ring & Glides | 710 x 710 x 950-1210mm | 520-710mm | 14.0kg |
| 5TC311 | High, Foot Ring & Glides | 710 x 710 x 980-1320mm | 570-820mm | 14.0kg |

Medical Records Trolleys





Sepsis Trolleys

Sepsis Trolleys 352 Options 353



Sepsis Trolleys

Features

. .

- Trolleys designed for use when treating Sepsis
- Supports delivery of 'The Sepsis Six' with clearly defined locations for equipment used in each of the stages
- · Door/s fitted with either:
- Electronic (battery operated) Push Button Lock Customer can set their own separate user, master, submaster & technician access codes
- · High Security Bolt Lock Supplied with two keys
- Trolleys supplied with door shelves (useful storage for antibiotics etc) & an A4 document holder
- Single Door 2 x door shelves / 1 x A4 document holder
- Double Door 5 x door shelves / 1 x A4 document holder Six labelled drawers (Oxygen / Blood Cultures / Intravenous Antibiotics /
- Intravenous Fluid /Haemoglobin & Lactates / Urine Measurement) Optional 'Sharps Box Holder' or 'Universal Sharps Box Bracket Mounting
- Plate' (fits to the I/h side of the trolley)
- Sharps Box Holder Accepts up to seven litre capacity sharps boxes (customer to fit - fixings supplied)
- Universal Sharps Box Bracket Mounting Plate Compatible with
- (customer to supply the plastic brackets & fit fixings supplied) Daniels Healthcare
 - POUDS® Wall-Bracket (up to 5L) SHARPSGUARD® mini, 1, 2.5, 3.75, 5 & com plus
 - POUDS® Wall-Bracket (7L & above) SHARPSGUARD® 7, 11.5 & 22
- Frontier (Sharpsafe) 2, 3, 4 & 7
- Sharpak Healthcare Sharpak 18, 35, 60, 120 & 220 PLUS
- Upper work surface
- Push handles incorporate integral buffers
- 125mm polymer swivel castors
- When not in use, trolleys storing medicines should be secured to the wall. See option MT/1/6F Security Clamp

Dimensions (w x d x h)

- Single Trolley: 512 x 542 x 1072mm
- Double Trolley: 712 x 542 x 1072mm
 - · Work Surface Orange

Standard Colours

Paint - Grev

- Drawers (internal) Single Width Drawer: 338 x 297 x 115mm
- Double Width Drawer: 538 x 297 x 115mm

| High Security Bolt Lock | Electronic Push Button Lock | Description | Weight |
|-------------------------|-----------------------------|-------------------------|---------|
| SEPS/6/HSB | SEPS/6/EPB | Single Door, Six Drawer | 53.4kg |
| SEPD/6/HSB | SEPD/6/EPB | Double Door, Six Drawer | 62.2 kg |

Sharps Box Holder

Features

- Accepts up to seven litre capacity sharps boxes
- Sharps box secured using a rubber band fitted with a toggle
- Attaches to the left hand side of the unit (customer to fit fixings supplied)





| Item | Description |
|---------|-------------------|
| SEP/SBH | Sharps Box Holder |

121 0

High Security

Bolt Lock

352

Sharps Box Bracket Mounting Plate - Universal

Features

- Compatible with the following brackets (customer to supply the plastic brackets & fit fixings supplied)
- Daniels Healthcare
- POUDS[®] Wall-Bracket (up to 5L) SHARPSGUARD[®] mini, 1, 2.5, 3.75, 5 & com plus
- POUDS[®] Wall-Bracket (7L & above) SHARPSGUARD[®] 7, 11.5 & 22
- Frontier (Sharpsafe) 2, 3, 4 & 7
- Sharpak Healthcare Sharpak 18, 35, 60, 120 & 220 PLUS



Security Clamp

Features

- · Allows the trolley to be secured to the 'fabric' of the building
- Requires no additional keys, utilises the trolleys own locking system
- One end is secured to the wall, the 'hook shaped' end is trapped by the door of the trolley when closed and locked
- Recommended wall fitting height 1000mm



Important Note - Security Clamps

It is a UK requirement that when the trolley is not in use it should be secured to the 'fabric of the building'. Not only does this prevent casual theft, i.e. trolleys being pushed to an isolated area and broken into, but also encourages good nursing practice, giving the trolley a 'permanent home'.


'Safe & Secure Handling of Medicines - Royal Pharmaceutical Society, GB- March 2005'

| Item | Description | Cable Length | Weight |
|---------|----------------|--------------|--------|
| MT/1/6F | Security Clamp | 950mm | 0.75kg |

Drawer Inserts

Features

- Plastic drawer inserts, height 100mm
- Assists segregation of drawer contents
- Suitable for single & double door trolleys
- Insert supplied with removable & repositionable dividers
- SEP/DIS 2 x dividers providing 4 separate areas
- SEP/DID 4 x dividers providing 6 separate areas
- Supplied for self-assembly
- Clip on label holders & plastic labels available to assist identification




| Item | Description | Weight |
|---------|--|--------|
| SEP/DIS | Drawer Insert - Set - Single Door Trolleys | 0.49kg |
| SEP/DID | Drawer Insert - Set - Double Door Trolleys | 0.55kg |

Label Holders & Labels

Features

- Plastic label holders & labels
- Assists identification of drawer contents
- Label holders snap onto drawer inserts/dividers
- Accepts ISO labels 40 x 60mm (w x h)
- Sold in packs of ten




| Item | Description | Weight |
|--------|--------------|---------|
| SEP/LH | Label Holder | 0.006kg |
| SEP/L | Label | 0.001kg |

Label - Sepsis Tiger

Features


- Sammy the sepsis tiger label developed for use on Paediatric trolleys
- Compatible with both single & double door trolleys














Shelving & Racking Solutions

Polymer Racking - MetroMax I[®] 360 Polymer Racking - MetroMax Q[®] 362

- Zinc Chrome Racking Olympic[®] 364
- High Density Storage Top-Track[®] 366

Engineered for a Demanding Enviroment





Microban[®] antimicrobial product protection helps keep shelves "cleaner between cleanings" by inhibiting the growth of mold, mildew, fungus and bacteria that cause odors and stains on the shelf surface.



METRO



Stationary Units

| 8MX546SG 8MX546SS 85MX546SG 85MX546SS 1066 x 610mm 8MX546MG 8MX556SG 8MX556SG 85MX556SS 1219 x 610mm 8MX556MG | 4 Shelves - 1590mm High | | 5 Shelves - 1 | 5 Shelves - 1895mm High | |
|---|-------------------------|---------------|---------------|-------------------------|--------------|
| 8MX546SG 8MX546SS 85MX546SG 85MX546SS 1066 x 610mm 8MX546MG 8MX556SG 8MX556SS 85MX556SS 1219 x 610mm 8MX556MG | Grid Shelves | Solid Shelves | Grid Shelves | Solid Shelves | |
| 8MX556SG 8MX556SS 85MX556SG 85MX556SS 1219 x 610mm 8MX556MG | 8MX536SG | 8MX536SS | 85MX536SG | 85MX536SS | 914 x 610mm |
| | 8MX546SG | 8MX546SS | 85MX546SG | 85MX546SS | 1066 x 610mm |
| 8MX566SG 8MX566SS 85MX566SG 85MX566SS 1524 x 610mm 8MX566MG | 8MX556SG | 8MX556SS | 85MX556SG | 85MX556SS | 1219 x 610mm |
| | 8MX566SG | 8MX566SS | 85MX566SG | 85MX566SS | 1524 x 610mm |



Polymer Racking - MetroMax I®

Features

- Corrosion proof polymer shelving & racking system
- Polymer shelves with stainless steel corners & a patented no tools shelf release mechanism, adjustable in 25mm increments
- All polymer posts
- Removable easy to clean polymer shelf mats (will fit most dishwashers) available as...
- Grid Shelves Open grid design allowing for the free circulation of air, improved visibility & light penetration
- Solid Shelves Solid design offering a tear free surface, ideal for the storage of sterilised DIN/ISO trays within theatres
- Built in Microban® antimicrobial product protection, to inhibit the growth of bacteria, mould, mildew & fungus that cause odours
- Available as stationary or mobile units
- Stationary Fitted with adjustable feet to allow levelling on uneven floors
- Mobile Fitted with 125mm swivel/braking castors
- Shelf Sizes
- Width 914, 1066, 1219 & 1524mm*
- Depth 610, 457mm*
- Shelf Capacity
- 914, 1066, 1219mm 363kg
- 1524mm 275kg
- Rack/Unit Capacity
- Stationary 907kg
- Mobile 339kg
- Label holders available for clear identification of items
- Fully configurable with dividers & ledges
- Lifetime anti corrosion warranty
- Delivered for self-assembly



Label Holders

Features

- Allows for clear identification of items
- Snaps onto shelf edge
- Accepts most commercial labels

Dimensions (w x h)

• 102 x 27mm

Standard Colours

Grey

| Item | Description | Weigh |
|---------|--------------|-------|
| 89989PX | Label Holder | 0.01k |

*other sizes available on request

914 x 610

1066 x 610

1219 x 610

1524 x 610

ves - 2035mm High Shelf Size (d x w)


85MX536MS

85MX546MS

85MX556MS

85MX566MS

Mobile Units





Shelf Divider

Features

- Corrosion proof polymer shelf divider (height 200mm)
- Use to create organised compartments or simply separate shelf contents
- Easily repositioned on grid & solid shelves
- Easy to clean

Dimensions (w x h)

• 610 x 200mm

Standard Colours

Grey

| Item | Description | Weight |
|----------|---------------|--------|
| 8MXD24-8 | Shelf Divider | 0.7kg |

Stackable Shelf Ledges

Features

- Corrosion proof stackable shelf ledge (height 51mm)
- Prevents larger items from protruding or falling
- Compatible with grid & solid shelves
- Easy to clean

Standard Colours

Chrome

| Item | Description | Weight |
|------------|---------------------|--------|
| Side Ledge | | |
| 8MXLS24-2S | Side Ledge - 610mm | 0.2kg |
| Rear Ledge | | |
| 8MXL36-2S | Rear Ledge - 914mm | 0.4kg |
| 8MXL42-2S | Rear Ledge - 1066mm | 0.6kg |
| 8MXL48-2S | Rear Ledge - 1219mm | 0.6kg |
| 8MXL60-2S | Rear Ledge - 1524mm | 0.8kg |







8Q566G2M

Polymer Racking - MetroMax Q[®]

Features

- Epoxy/polymer shelving & racking system
- Epoxy coated shelves & posts with a patented no tools shelf release mechanism, adjustable in 25mm increments
- Removable easy to clean polymer shelves (will fit most dishwashers)
- Open grid design allowing for the free circulation of air, improved visibility & light penetration
- Built in Microban® antimicrobial product protection, to inhibit the growth of bacteria, mould, mildew & fungus that cause odours
- Available as stationary or mobile units
- Stationary Fitted with adjustable feet to allow levelling on uneven floors
- Mobile Fitted with 125mm swivel/braking castors
- Shelf Sizes
- Width 914, 1066, 1219 & 1524mm*
- Depth 610, 457mm*
- Shelf Capacity
 - 914, 1066, 1219mm 363kg
 - 1524mm 275kg
- Rack/Unit Capacity
- Stationary 907kg
- Mobile 339kg
- Label holders available for clear identification of items
- Fully configurable with dividers & ledges
- Anti Corrosion warranty
- Frame & Posts Fifteen years
- Mats Lifetime
- · Delivered for self-assembly

| Stationary | / Units | | | Mobile Unit | ts *other siz | zes available | on request |
|------------|--------------------|--------|--------|-------------|--------------------|---------------|------------|
| 4 Shelves | Shelf Size (w x d) | Height | Weight | 4 Shelves | Shelf Size (w x d) | Height | Weight |
| 8Q536G2 | 914 x 610mm | 1590mm | 30.0kg | 8Q536G2M | 914 x 610mm | 1730mm | 34.6kg |
| 8Q546G2 | 1066 x 610mm | 1590mm | 33.6kg | 8Q546G2M | 1066 x 610mm | 1730mm | 38.2kg |
| 8Q556G2 | 1219 x 610mm | 1590mm | 35.6kg | 8Q556G2M | 1219 x 610mm | 1730mm | 40.2kg |
| 8Q566G2 | 1524 x 610mm | 1590mm | 44.4kg | 8Q566G2M | 1524 x 610mm | 1730mm | 49.0kg |
| 5 Shelves | Shelf Size (w x d) | Height | Weight | 5 Shelves | Shelf Size (w x d) | Height | Weight |
| 85Q537G2 | 914 x 610mm | 1895mm | 36.7kg | 85Q537G2M | 914 x 610mm | 2035mm | 41.3kg |
| 85Q547G2 | 1066 x 610mm | 1895mm | 41.2kg | 85Q547G2M | 1066 x 610mm | 2035mm | 45.8kg |
| 85Q557G2 | 1219 x 610mm | 1895mm | 43.7kg | 85Q557G2M | 1219 x 610mm | 2035mm | 48.3kg |
| 85Q567G2 | 1524 x 610mm | 1895mm | 54.7kg | 85Q567G2M | 1524 x 610mm | 2035mm | 59.3kg |



Label Holders

Features

- Allows for clear identification of items
- Snaps onto shelf edge
- Accepts most commercial labels

Dimensions (w x h)


• 102 x 27mm

Standard Colours

Grey

| Item | Description | Weight |
|--------|--------------|--------|
| 8Q04LH | Label Holder | 0.01kg |







Shelf Divider

Features

- Corrosion proof polymer shelf divider (height 200mm)
- Use to create organised compartments or simply separate shelf contents
- Easily repositioned on grid & solid shelves
- Easy to clean

Dimensions (w x h)

• 610 x 200mm

Standard Colours

Grey

| Item | Description | Weight |
|----------|---------------|--------|
| 8MQD24-8 | Shelf Divider | 0.7kg |

Stackable Shelf Ledges

Features

- Corrosion proof stackable shelf ledge (height 51mm)
- Prevents larger items from protruding or falling
- Compatible with grid & solid shelves
- Easy to clean

Standard Colours

Chrome

| Item | Description | Weight |
|------------|---------------------|--------|
| Side Ledge | | |
| 8MQLS24-2S | Side Ledge - 610mm | 0.2kg |
| Rear Ledge | | |
| 8MQL24-2S | Rear Ledge - 610mm | 0.2kg |
| 8MQL36-2S | Rear Ledge - 914mm | 0.4kg |
| 8MQL42-2S | Rear Ledge - 1066mm | 0.6kg |
| 8MQL48-2S | Rear Ledge - 1219mm | 0.6kg |
| 8MQL60-2S | Rear Ledge - 1524mm | 0.8kg |

OLYNPIC



Zinc Chrome Racking - Olympic[®]

Features

- Zinc Chrome shelving & racking system
- Open wire design shelves allowing for the free circulation of air, improved visibility & light penetration
- Chromate posts & shelves with a patented no tools tapered slit sleeve system that snaps together on each post, adjustable in 25mm increments
- Available as stationary or mobile units
- Stationary Fitted with adjustable feet to allow levelling on uneven floor
- Mobile Fitted with 125mm swivel/braking castors
- Shelf widths 914, 1066, 1219 & 1524mm*
- Shelf depths 610, 457mm*
- Shelf Capacities
- 914, 1066, 1219mm 364kg
- 1524mm 273kg
- Rack/Unit Capacities
- Stationary 907kg
- Mobile 339kg
- Label holders available for clear identification of items

4 Shelves Shelf Size (w x d) Height Weigh

8O966CM 1524 x 610mm 1730mm 49.0kg

5 Shelves Shelf Size (w x d) Height Weight

1730mm 34.6kg

1730mm 38.2kg

1730mm 40.2kg

2035mm 41.3kg

2035mm 45.8kg

2035mm 48.3kg

- Fully configurable with dividers & ledges
- 1 year performance warranty on shelves and posts
- Delivered for self-assembly

Mobile Units

1895mm 36.7kg 850537CM 914 x 610mm

8O536CM 914 x 610mm

8O546CM 1066 x 610mm

8O556CM 1219 x 610mm

85O547CM 1066 x 610mm

85O557CM 1219 x 610mm

1895mm 54.7kg 85O967CM 1524 x 610mm 2035mm 59.3kg

1590mm 30.0kg

1590mm 33.6ka

1590mm 35.6kg

1590mm 44.4kg

1895mm 41.2kg

1895mm 43.7kg



Label Holders

Features

- Allows for clear identification of items
- Snaps onto shelf edge
- Accepts most commercial labels

Dimensions (w x h)

• 75 x 32mm

Standard Colours

Grey

| Item | Description | Weight |
|--------|--------------|--------|
| 89990P | Label Holder | 0.1kg |

+44 (0)1258 484455 www.bristolmaid.com sales@bristolmaid.com

85O967C 1524 x 610mm



Shelf Divider

Features

- Keeps shelf contents orderly (height 203mm)
- Pressure fit to hold in place

Dimensions (w x h)

• 610 x 203mm

Standard Colours

Chrome

Shelf Ledges

Features

- Prevents items from protruding or falling (height 100mm)
- Pressure fit to hold in place
- Ledges can be stacked

Standard Colours

Chrome

| | | | Item | Shelf Length | Weight |
|--------|---------------|--------|----------|--------------|--------|
| | | | 8L24N-4C | 610mm | 0.90kg |
| | | | 8L36N-4C | 914mm | 1.25kg |
| | | | 8L42N-4C | 1066mm | 1.48kg |
| Item | Description | Weight | 8L48N-4C | 1219mm | 1.68kg |
| 8DD24C | Shelf Divider | 0.7kg | 8L60N-4C | 1524mm | 1.81kg |





High Density Storage - Top-Track®

Features

- Increases the storage capacity of a given area by up to 50%
- Maximises all available space
- Can eliminate the need to expand a given area through the use of optimised storage
- Creates moveable aisles
- Eliminates hard to reach or 'dead' spaces
- The guide track is positioned above the shelving system
- Aluminium overhead track system acts as a guide for mobile units
- Floors are easy to clean
- Carts & trolleys can easily be easily manoeuvred into & out of the active aisle
- The guide track compensates for uneven floor surfaces & keeps units in alignment
- Mobile units are designed to accommodate medium-duty applications
- The maximum weight capacity of a Top-Track® mobile unit is 410kg



Top Track Roller

Contact us to arrange your site-survey +44 (0)1258 484455



Sterile Container Storage



See Page 396









Sterile Services Equipment

- Preparation Tables Variable Height 370
 - Preparation Tables Fixed Height 371
 - Sterile Supplies Trolleys 372
 - Trolley Washer Compatible 374

Preparation Tables - Stainless Steel, Variable Height

Features

- Height adjustable preparation tables
- Uses a manual turn handle (folds away when not in use) • Static - 715-1015mm
- Mobile 812-1112mm
- Available as a static or mobile
- Static Fitted with adjustable feet for levelling on uneven floors
- Mobile Fitted with 80mm castors, 2 off braking
- 304 grade stainless steel frame & work surface Max capacity 70.0kg

Dimensions (w x d x h)

- Static 1500 x 750 x 715-1015mm • Mobile - 1500 x 750 x 812-1112mm

Standard Colours

Stainless Steel

Preparation Tables - Stainless Steel, Fixed Height

Features

- Fixed height preparation tables
- Worksurface height 950mm
- Available as a static or mobile
- Static Fitted with adjustable feet for levelling on uneven floors Mobile - Fitted with 100mm castors, 2 off braking
- 304 grade stainless steel frame, lower shelf & work surface
- Max capacity 70.0kg

Dimensions (w x d x h)

Overall Dimensions - 1500 x 750 x 950mm

Standard Colours

Stainless Steel



| Item | Description | Weight |
|-------|--|--------|
| PT050 | Preparation Table - Variable Height, Stainless Steel, Static | 49.0kg |
| PT110 | Preparation Table - Variable Height, Stainless Steel, Mobile | 53.0kg |

| Item | Description |
|-----------|--|
| PTS/115 | Preparation Table - Fixed Height, Stainles |
| PTS/115/M | Preparation Table - Fixed Height, Stainles |



Weight

ess Steel, Static 56.3kg ess Steel, Mobile 58.8kg



Sterile Supplies Trolleys - UN3291 Certified

Features

- Independently tested to UN3291 Certificate of Packaging Performance for the Transportation of Clinical Waste
- Suitable for the transportation of DIN & ISO instrument trays
- Mild steel frame
- Lightweight aluminium cladding, shelves & door/s
- Powder coated for durability
- Locking door/s…
 - Fitted with a flush fitting handle fitted with a CAM lock offering two-point locking (supplied with 2 keys)
 - When opened, doors fold back flat at either end of the trolley*
- Three adjustable shelves, (max 40kg per shelf)
- Shelf size...
- Small 534 x 614mm (w x d)
- Large 1128 x 614mm (w x d)
- Removable stainless steel drip tray
- Inset push handles at either end of trolley
- Low level bumper rail offers protection to both the trolley & fabric of the building
- Optional tow bar fitting...
- Tow bar fitted to the front (swivel castor)
- Antiluce fastener fitted to rear (fixed castor)
- Tow bar position approx. 360-380mm above floor
- Recommended max speed 7kph for towed trolleys
- 200mm heavy duty cushioned castors, 2 x fixed, 2 x swivel braking
- Maximum capacity..
 - Small 84k.0g per trolley
- Large 109.0kg per trolley

• Large - 1399 x 774 x 1257mm

Dimensions (w x d x h) • Small - 805 x 774 x 1257mm

Standard Colours

- Frame Grey
- Cladding / Shelves / Doors White



AGV Compatible Trolleys

Bristol Maid[™] is your ideal partner if you are looking for transport carts for Automated Guided Vehicle (AGV) systems. Through our experienced team we are able to offer solutions for Sterile Supplies, Laundry, Medical Records & Pharmacy distribution.

| Item | Description | Weight |
|-------------|--|---------|
| Small | | |
| ST/340/UN | Sterile Supplies Trolley - UN3291 Certified, Small | 74.5kg |
| ST/340/UN/T | Sterile Supplies Trolley - UN3291 Certified, Towing, Small | 78.5kg |
| Large | | |
| ST/690/UN | Sterile Supplies Trolley - UN3291 Certified, Large | 110.0kg |
| ST/690/UN/T | Sterile Supplies Trolley - UN3291 Certified, Towing, Large | 114.0kg |

Drip Tray



Sterile Supplies Trolleys - UN3291 Certified, Trolley Washer Disinfector Compatible



Features

- Independently tested to UN3291 Certificate of Packaging Performance for the Transportation of Clinical Waste
- Compatible with trolley washer disinfectors
- Suitable for the transportation of DIN & ISO instrument trays
- 304 grade stainless steel frame, cladding, shelves & doors
- Natural unpolished finish
- Locking door/s…
 - Fitted with a flush fitting handle fitted with a CAM lock offering two-point locking (supplied with 2 keys)
 - When opened, doors fold back flat at either end of the trolley*
 - Spring loaded plungers hold the doors open, essential when the trolley is being processed in a trolley washer disinfector
- One fixed, two adjustable shelves (max 50kg per shelf)
- Shelf size 1227 x 632mm (w x d)
- Stainless steel drip tray...
- Removable for cleaning
- Can be located at the rear of the trolley during cleaning, held in position using spring loaded plungers
- Inset push handles at either end of trolley
- High level corner bumpers offer protection to both the trolley & fabric of the building
- Optional tow bar fitting...
 - Tow bar fitted to the front (swivel castor)
 - Antiluce fastener fitted to rear (fixed castor)
 - Tow bar position approx. 360-380mm above floor
 - Recommended max speed 7kph for towed trolleys
- · 200mm heavy duty sealed bearing stainless steel castors,
- 2 x fixed, 2 x swivel braking
- Maximum capacity 136.0kg per trolley

Dimensions (w x d x h)

• 1407 x 762 x 1288mm

Standard Colours

Stainless Steel - Natural Unpolished





Drip Tray Locates at Rear of Trolley During Cleaning



High Level Corner Bumpers

AGV Compatible Trolleys Bristol Maid[™] is your ideal partner if you are

looking for transport carts for Automated Guided Vehicle (AGV) systems. Through our experienced team we are able to offer solutions for Sterile Supplies, Laundry, Medical Records & Pharmacy distribution.

| Item | Description | Weight |
|---------------|--|---------|
| STSS/690/UN | Sterile Supplies Trolley - Stainless Steel, Trolley Washer/Sterilisation Unit Compatible | 131.0kg |
| STSS/690/UN/T | Sterile Supplies Trolley - Stainless Steel, Trolley Washer/Sterilisation Unit Compatible, Towing | 137.0kg |











Storage Cupboards Cabinets & Lockers



- Storage Cabinets 380
- First Aid & Clinical Reagent Cabinets 381
 - Key Cabinets 382
- Storage / Personal Effects Lockers 384



SC/L/1870/7S5D

Storage Cupboards - Single Door, 1870mm High

Features

- Single door 1870mm high cupboards
- Single-point locking using a flush fitting CAM lock fitted with a handle (supplied with 2 keys)
- Available with either...
- Plastic Travs
- Shelves
- Shelves & Dividers
- 0.9mm (20 Gauge) mild steel welded construction
- Modular design allows for easy interchange with HTM 63 storage systems

Plastic Trays

- Interchangeable with our popular CT range (post 2016)
- Supplied with 7 x 100mm & 5 x 150mm trays
- Translucent high impact plastic trays feature an integral handle with plastic clip on label cover, (max load per tray 5kg)
- Each cupboard is supplied with perforated coloured label strips (white, blue, green, yellow & red)
- Trays can be subdivided in both length & width using simple slot in dividers (each tray supplied with 2 x long & 2 x short dividers)

Shelves

- Flat shelves adjustable at 50mm increments
- Supplied with 8 x shelves

Shelves & Dividers

- Flat shelves adjustable at 50mm increments
- Removable dividers 100mm high adjustable at 25mm increments
- Supplied with 8 x shelves & 16 x dividers

Dimensions (w x d x h)

Standard Colours

Cupboard • 500 x 536 x 1870mm Paint - White

870mm

- Trays
- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)

Options

SC/JK/1870 - Joining Kit

SC/DIV - Set of 16 Dividers

| /H Hinged | L/H Hinged | Description | Weigh |
|---------------|----------------|---|--------|
| C/R/1870/7S5D | SC/L/1870/7S5D | Single Door, 1870mm High, Plastic Trays | 83.6kg |
| C/R/1870/FS | SC/L/1870/FS | Single Door, 1870mm High, Flat Shelves | 64.0kg |
| C/R/1870/FSD | SC/L/1870/FSD | Single Door, 1870mm High, Flat Shelves & Dividers | 74.0kg |





SC/L/1870/FS

Joining Kits

Features

Allows two units to be joined together



Dividers

Features

- Spare / replacement dividers for SC cupboards supplied with flat shelves & dividers
- Removable 100mm high dividers locate into slots on the shelves (adjustable at 25mm increments)
- Compatible with 1870mm high units



| Item | Description |
|--------|--------------------|
| SC/DIV | Set of 16 Dividers |



Storage Cabinets - Single & Double Door

Features

- Adjustable & removable shelves
- Flush fitting handle fitted with a CAM lock offering single-point locking (supplied with 2 keys)
- 0.9mm (20 gauge) mild steel welded construction
- Mounting holes in back (fixings not included)
- Modular design allows for easy interchange with HTM 63 storage systems

Standard Colours

Paint - White



| Item | Shelves | Overall Dimension (w x d x h) | Weight | |
|-------------|---------|-------------------------------|--------|--|
| Single Door | | | | |
| SC020 | 2 | 600 x 300 x 600mm | 19.5kg | |
| SC025 | 3 | 600 x 300 x 900mm | 26.0kg | |
| Double Door | | | | |
| SC040 | 4 | 1000 x 300 x 600mm | 32.5kg | |
| SC045 | 6 | 1000 x 300 x 900mm | 45.0kg | |

Storage Cabinets - First Aid & Clinical Reagents

Features

- Removable shelves
- Flush fitting handle fitted with a CAM lock offering single-point locking (supplied with 2 keys)
- 0.9mm (20 gauge) mild steel welded construction
- Mounting holes in back (fixings not included)
- Supplied with 'First Aid' & 'Clinical Reagents' self-adhesive labels

Dimensions (w x d x h)

- Small 300 x 200 x 400mm
- Medium 400 x 200 x 400mm
- Large 400 x 200 x 600mm





Item Description

GC005 First Aid & Clinical Reagents Cabinet - Small GC010 First Aid & Clinical Reagents Cabinet - Medium GC015 First Aid & Clinical Reagents Cabinet - Large





How does it work?

600



Access User enters unique PIN code to gain access to the cabinet



Kev Removal Keys & keysets are individually locked An amber light indicates the return within the cabinet. Green LED's position for the key indicate those keys which the user can access. The user cannot remove keys highlighted by the Red LED's which continue to remain locked in place

Key Replacement

Storage Cabinets - Keys, Traka21

Features

- User, key & access rights administration
- Screen based audit trail & reporting capability or export via USB port
- 21 robust, long-life iFobs with security seals
- 21 locking receptors with LEDs
- Robust & compact aluminium & ABS housing with no exposed hinges or obvious access points
- Manual override & door release functions in case of emergency
- Audible alarms

Operation

- Stand-alone 'Plug & Play' solution with advanced RFID technology
- Capacity for up to 21 keys
- Touchscreen interface
- PIN access to designated keys or keysets which are individually locked in place
- · Keys are securely attached using special security seals
- Simple to configure through setup wizard
- No requirement for network connection or PC
- Concealed robust fixing points allow secure wall fixing
- Mains operated (optional battery back-up)

Standard Colours Dimensions (w x d x h)

· Aluminium & ABS - White

Traka21 iFob

Options

• 246 x 95 x 427mm

- SKR/25 Security Key Rings (x 25)
- SKR/100 Security Key Rings (x 100)
- SKB/BAT Battery Back-Up



Once attached to an iFob with a robust security seal, individual keys or keysets are assigned a specific receptor within the Traka cabinet & locked in place until released by an authorised user

| Item | Description | Weight |
|-------------|-----------------------|--------|
| SKC21/TRAKA | Key Cabinet - TRAKA21 | 4.0kg |

Storage Cabinets - Keys

Features

- Cabinets for the control & tracking of keys
- Available with either a CAM or electronic push button lock
- CAM Supplied with two keys
- Electronic Push Button Battery operated
- Adjustable hook bars accommodate different key lengths
- Removable control index for out of hours peace of mind
- Colour coded hook bars & tags ensure guick identification. (both numbered)
- · Provided with security seals & loops to secure keys that require evidence of tampering or removal

Standard Colours

Paint - Grev

Electronic Push Button Lock

- Uses a 10 button keypad and a # button. The # button is hidden beneath the handle when in the locked position
- Easy programming & code change
- Four operating codes with differing levels of authority Master code
- Sub-master code
- User code
- Technician code
- User codes are 4 digits
- Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code

CAM Lock

| ltem | Key Capacity | Overall Dimensions (w x d x |
|--------|--------------|-----------------------------|
| SKC20 | 20 Keys | 185 x 80 x 260mm |
| SKC30 | 30 Keys | 215 x 80 x 305mm |
| SKC35 | 35 Keys | 245 x 80 x 320mm |
| SKC48 | 48 Keys | 300 x 80 x 355mm |
| SKC64 | 64 Keys | 300 x 80 x 450mm |
| SKC100 | 100 Keys | 380 x 80 x 550mm |
| | | |





Electronic Push Button Lock

| x h) | Weight | Item | Key Capacity | Overall Dimensions (w x d x h) | Weight |
|------|--------|-----------|--------------|--------------------------------|--------|
| | 2.0kg | SKC20/CL | 20 Keys | 185 x 110 x 260mm | 2.0kg |
| | 3.0kg | SKC30/CL | 30 Keys | 215 x 110 x 305mm | 3.0kg |
| | 3.0kg | SKC35/CL | 35 Keys | 245 x 110 x 320mm | 3.0kg |
| | 4.0kg | SKC48/CL | 48 Keys | 300 x 110 x 355mm | 4.0kg |
| | 5.0kg | SKC64/CL | 64 Keys | 300 x 110 x 450mm | 5.0kg |
| | 8.0kg | SKC100/CL | 100 Keys | 380 x 110 x 550mm | 8.0kg |

Alternative colours available

BD

Blue Door

RD

Red Door

Replace suffix on the product code with any of the following to select your colour i.e. Single Bank with a Green Door - PL080/GD



GD Green Door



YD Yellow Door





Features

- Range of personal effects lockers available in single, double or triple banks in a number of compartment configurations
- Choice of locking options
- CAM lock (supplied with 2 keys)
- Hasp & staple (accepts a padlock/not supplied)
- Lock surround incorporates a 25mm round label holder
- Mounting holes provided in back (fixings not included)
- Optional stands & sloping tops

Dimensions (w x d x h)

- Single Bank 305 x 460 x 1780mm
- Double Bank 610 x 460 x 1780mm • Two - 255 x 440 x 815mm • Three - 255 x 440 x 525mm
- Triple Bank 915 x 460 x 1780mm

Standard Colours

- Body Grey Doors - Blue
- 6



Internal Compartments (w x d x h)


• One - 255 x 440 x 1680mm

• Four - 255 x 440 x 385mm • Five - 255 x 440 x 295mm

PL150/BD

Single Bank Compartments Yes No Coat Hook Yes Shelf Yes No No Hasp & Staple PL005/BD PL010/BD CAM Lock PL080/BD PL085/BD PL090/BD PL095/BD PL100/BD

| Double Bank | * * | 5 F | | | |
|---------------|----------|----------|----------|----------|----------|
| Compartments | 1 | 2 | 3 | 4 | 5 |
| Coat Hook | Yes | Yes | No | No | No |
| Shelf | Yes | No | No | No | No |
| Hasp & Staple | PL030/BD | PL035/BD | PL040/BD | PL045/BD | PL050/BD |
| CAM Lock | PL105/BD | PL110/BD | PL115/BD | PL120/BD | PL125/BD |
| | | | | | |






| • | | |
|-------------|--|--|
| | | |
| 5 | | |
| 2 | | |
| | | |
| 2 | | |
| | | |
| 5 | | |
| | | |
| ί. | | |
| 1 | | |
| ÷ | | |
| | | |
|) | | |
| 5 | | |
| έ. | | |
| | | |
| ק ק ק | | |
| 5 | | |
| 2 | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |

68

Hasp & Staple Lock



Stand

Features

- Compatible with single, double & triple bank 1780mm high personal effects lockers
- Keeps lockers clear of the floor assisting cleaning & aiding air-circulation
- Fitted with adjustable feet
- Increases the overall height of the locker by 150mm

Standard Colours

Stand - Grey



| ltem | Description |
|-------|---------------------|
| PL900 | Stand - Single Bank |
| PL905 | Stand - Double Bank |
| PL910 | Stand - Triple Bank |

Sloping Top

Features

- Compatible with single, double & triple bank 1780mm high personal effects lockers
- Prevents items being stored on the top of the locker & assists cleaning
- Increases the overall height of the locker by 150mm

Description

Sloping Top - Single Bank

Sloping Top - Triple Bank

Sloping Top - Double Bank

Standard Colours

Sloping Top - Grey

ltem

PL915

PL920

PL925





| 1 | |
|------|-------------|
| PL91 | 5 - Sloping |
| | |



Theatre Equipment

- Bowl Stands 388
- Bucket Stands 391
- Instrument Trolleys 392
- Mayo Tables & Operating Platforms 394
 - Suture Storage Trolleys 395
 - Sterile Container Storage 396

387



Bowl Stands - Fixed Height

Features

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Bowl retaining rings fitted with removable anti-static spiral buffers to support bowls
- Suitable for use with 360mm diameter stainless or polythene bowls (order separately)
- 75mm anti-static castors

Standard Colours

Stainless Steel

Options

- 2MECO079 360mm Diameter, Stainless Steel Bowl
- 2PLCO054 360mm Diameter, Polythene Bowl



2MECO079 Stainless Steel Bowl



| Item | Description | Overall Dimensions (w x d x h) | Weight |
|--------|--|--------------------------------|--------|
| SBS/2 | Bowl Stand - Fixed Height, Double Bowl, Side by Side | 735 x 400 x 855mm | 6.0kg |
| SBST/2 | Bowl Stand - Fixed Height, Double Bowl, Tiered | 400 x 400 x 855mm | 5.0kg |

Bowl Stands - Variable Height

Features

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Height adjustable upright using a handwheel (753-1100mm)
- Bowl retaining ring fitted with removable anti-static spiral buffers to support bowls
- Suitable for use with 360mm diameter stainless or polythene bowls (order separately)
- 50mm polypropylene/thermoplastic rubber twin-wheel castors (1 x anti-static)

Standard Colours

Stainless Steel

Options

- 2MECO079 360mm Diameter, Stainless Steel Bowl
- 2PLCO054 360mm Diameter, Polythene Bowl



2MECO079 Stainless Steel Bowl



2PLCO054 Polythene Bowl

| Item | Description | |
|--------|---|--|
| SBSA/1 | Bowl Stand - Variable Height, Single Bowl | |
| SBSA/2 | Bowl Stand - Variable Height, Double Bowl | |



Bowl Stands - Variable Height, Polyurethane, Easy Clean

Features

- · Five leg base with low centre of gravity to assist stability
- 25mm diameter scratch resistant stainless steel upright
- A seamless polyurethane over moulding encases the base providing...
- A flexible tear & water resistant finish
- Impact protection to the fabric of the building
- A shield for the twin wheel castors
- Smooth edges & contours, allowing for easy & effective cleaning (including steam cleaning)
- Coloured indicator caps provide department identification & segregation, reducing the risk of cross contamination
- Height adjustable upright using a handwheel
- Single Bowl 820 1125mm
- Double Bowl 820 1125mm
- Bowl retaining rings fitted with removable anti-static spiral buffers to support bowls
- Suitable for use with 360mm diameter stainless or polythene bowls (order separately)
- 50mm polypropylene/thermoplastic rubber twin wheel castors (1 x anti-static)

Dimensions (w x d x h)

Standard Colours

Bowl Retaining Rings - Stainless Steel

Base - Grev

- Single Bowl 650 x 650 x 820 1125mm Coloured Indicator Cap - Spearmint
- Double Bowl 746 x 650 x 820 1125mm

Options

- Coloured Indicator Cap Blue / Green / Red / Yellow
- 2MECO079 Stainless Steel Bowl
- 2PLCO054 Polythene Bowl



- Bowl retaining ring fitted with removable anti-static spiral buffers to support bowls
- Supplied with a 360mm diameter stainless steel bowl
- 75mm anti-static castors

Standard Colours Dimensions (w x d x h)

365 x 365 x 235mm

Options

Features

- 2MECO079 360mm Diameter, Stainless Steel Bowl
- 2PLCO054 360mm Diameter, Polythene Bowl



| tem | Description | Weight |
|------------|---|--------|
| BO/P/VS/SP | Bowl Stand - Variable Height, Polyurethane, Easy Clean, Single Bowl | 9.0kg |
| BO/P/VD/SP | Bowl Stand - Variable Height, Polyurethane, Easy Clean, Double Bowl | 9.8kg |

| Item | Description |
|--------|--|
| SBSL/1 | Bowl Stand - Kickabout, 360mm Diameter |





Bowl Stand - Kickabout, 360mm Diameter

Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning

Stainless Steel

7 litre stainless



Features

Options

- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Bucket retaining ring fitted with removable anti-static spiral buffers to support the bucket
- Model available fitted with a handle
- Supplied with a 15-litre graduated stainless steel bucket
- 75mm anti-static castors

Dimensions (w x d x h)

Standard Colours

- SBKSL/1 380 x 380 x 515mm
 Stainless Steel
- SBKSLH/1- 435 x 435 x 515mm

vith a 15-litre



| Item | Description | Weight |
|----------|--|--------|
| SBKSL/1 | Bucket Stand - 15 Litre Bucket | 5.0kg |
| SBKSLH/1 | Bucket Stand - 15 Litre Bucket, Handle | 5.5kg |

3.2kg



DTSF/450/2/SD

Fixed Shelves, Sides Down, Stainless Steel, 450mm Deep

Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength & ease of cleaning
- Fixed shelves with down turned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 75mm pressed steel non-marking castors (1 x anti-static, 3 x non-anti-static)

Options

- C001 80mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- C013 75mm, Pressed Steel Castors (4 x braking)
- Drawer Units Range of sizes / configurations (see page 107)
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers

Fixed Shelves, Sides Down, Stainless Steel, 600mm Deep

Features

- Meets BS4068:1977 Specification for hospital trolleys for instruments, dressings & for anaesthetists' use
- Seamlessly welded from 304 grade stainless steel for strength
 & ease of cleaning
- Fixed shelves with down turned sides
- Shelf heights 405 & 870mm
- Maximum capacity per shelf 40kg
- 100mm pressed steel anti-static non-marking castors

Options

- C002 100mm, Grey Polymer Castors (2 x braking / 2 x non-braking)
- C003 100mm, Pressed Steel Castors (2 x braking / 2 x non-braking)
- C015 100mm, Pressed Steel Castors (4 x braking)
- DP/C Infusion Pole & Bracket
- ZP005 Ring Buffers

| Item | Size | Shelf Size (w x d) | Overall Dimensions (w x d x h) | Weight |
|---------------|---------|--------------------|--------------------------------|--------|
| DTSF/450/2/SD | Small | 450 x 450mm | 470 x 470 x 870mm | 12.0kg |
| DTSF/600/2/SD | Medium | 600 x 450mm | 620 x 470 x 870mm | 14.5kg |
| DTSF/750/2/SD | Large | 750 x 450mm | 770 x 470 x 870mm | 18.0kg |
| DTSF/900/2/SD | X-Large | 900 x 450mm | 920 x 470 x 870mm | 21.0kg |

| Item | Size | Shelf Size (w x d) | Overa |
|---------------|--------|--------------------|---------|
| DTSF/66/2/SD | Small | 600 x 600mm | 620 x 6 |
| DTSF/96/2/SD | Medium | 900 x 600mm | 920 x 6 |
| DTSF/126/2/SD | Large | 1200 x 600mm | 1220 x |



Mayo Tables - Cantilever, Variable Height

Features

- Seamlessly welded frame from 304 grade stainless steel for strength & ease of cleaning
- Available with two castors & skid feet or four castors
- Height adjustment...
- SS/22/1 893-1480mm
- SS/22/2 858-1450mm
- Supplied with a removable stainless steel tray, 512 x 387 x 20mm (w x d x h)
- 50mm polypropylene/thermoplastic rubber twin-wheel anti-static castors

Dimensions (w x d x h)

Standard Colours

- SS/22/1 641 x 450 x 893-1480mm
- Stainless Steel



| Item | Description | Weight |
|---------|---|--------|
| SS/22/1 | Mayo Table - Cantilever, Variable Height, Two Castors & Skid Feet | 5.4kg |
| SS/22/2 | Mayo Table - Cantilever, Variable Height, Four Castors | 5.8kg |

Operating Platforms

Features

- Seamlessly welded frame from 304 grade stainless steel for strength & ease of cleaning
- Aluminium embossed non-slip tread
- Non-slip, anti-static rubber feet
- Max capacity 190kg

Dimensions (w x d x h)

- Small 400 x 300mm
- Large 600 x 300mm

Standard Colours

Stainless Steel



Aluminium embossed non-slip tread

| | T | |
|-------|--------|--|
| and a | | |
| | | |
| | | |
| | OP/432 | |
| | 017432 | |
| | | |

| ltem | Description | Weight |
|--------|--|--------|
| | | |
| OP/431 | Operating Platform - Small, 150mm High | 4.2kg |
| OP/432 | Operating Platform - Small, 200mm High | 4.4kg |
| Large | | |
| OP/631 | Operating Platform - Large, 150mm High | 6.0kg |
| OP/632 | Operating Platform - Large, 200mm High | 6.4kg |

Suture Storage Trolleys

Features

- · Angled shelves present contents for instant identification & easy access
- Allows logical organisation, eliminates confusion, allowing for easier & faster restocking
- Open wire design minimises dust accumulation, allowing for free circulation of air. improved visibility & light penetration
- No tools shelf adjustment in 25mm increments
- Each trolley supplied with 12 x shelf retainers to assist with segregation of contents
- Delivered flat pack for self-assembly
- 125mm pressed steel castors, 2 off braking

Options

- 89990P Label Holder, 75 x 32mm
- 8DCR17C Shelf Retainers

| Item | Description | Overall Dimensions (w x d x h) | Weight |
|-------------|---|--------------------------------|--------|
| 1 x Flat, 4 | x Angled Shelves | | |
| 8DC15EC | Trolley - Suture Storage, 1 x Flat / 4 x Angled Shelves, 610mm | 610 x 457 x 1524mm | 46.5kg |
| 8DC35EC | Trolley - Suture Storage, 1 x Flat / 4 x Angled Shelves, 914mm | 914 x 457 x 1524mm | 49.0kg |
| 8DC55EC | Trolley - Suture Storage, 1 x Flat / 4 x Angled Shelves, 1219mm | 1219 x 457 x 1524mm | 51.0kg |
| 5 x Angled | Shelves | | |
| 8DC16EC | Trolley - Suture Storage, 5 x Angled Shelves, 610mm | 610 x 457 x 1524mm | 50.0kg |
| 8DC36EC | Trolley - Suture Storage, 5 x Angled Shelves, 914mm | 914 x 457 x 1778mm | 51.5kg |
| 8DC56EC | Trolley - Suture Storage, 5 x Angled Shelves, 1219mm | 1219 x 457 x 1778mm | 52.5kg |





8DC56EC



Sterile Container Storage

Features

- Fully configurable freestanding modular racking suitable for the safe storage of sterile containers
- Compatible with both blue wrap & rigid containers, up to max 300mm x 600mm (width x depth)
- Containers stored on either...
- Telescopic shelves
- Wire arid shelves
- · Optional silicone pads offer additional cushioning for containers (suitable with telescopic shelves only)
- Anodised extruded aluminium uprights & cross members
- · Bolted construction using high quality moulded fixings provided rigidity & allows the system to be altered easily
- Available as either static or mobile units
- Static Adjustable feet allow levelling on uneven floors
- Mobile 50/80mm pressed steel castors
- Optional wall mounting brackets allow static units to be fixed securely to the wall. (when units are loaded over 100kg)
- Supplied for self-assembly

Dimensions (w x d x h)

Static

- Single Bay 394 x 642 x 1940mm
- Double Bay 737 x 642 x 1940mm

Mobile (50/80mm Castors)

- Single Bay 394 x 642 x 1996/2035mm
- Double Bay 737 x 642 x 1996/2035mm
- Standard Colours
- Uprights & Cross Members Anodised Aluminium
- Wire Grids Polished Stainless Steel
- Telescopic Shelves Anodised Aluminium
- Plastic Moulding Green / Grey

Options

- 2HH-7532293 Wall Mounting Bracket (x1)
- 2HH-7532277 10mm Plastic Cover Cap (x1)

Static

| Item | Description | Weight |
|------------------|---|--------|
| Single Bay | | |
| 4HH-SCS/S1/2W5T | Static, Single Bay, 2 x Wire Grid Shelves & 5 x Telescopic Shelves | 30.2kg |
| 4HH-SCS/S1/8W | Static, Single Bay, 8 x Wire Grid Shelves | 21.4kg |
| Double Bay | | |
| 4HH-SCS/S2/4W10T | Static, Double Bay, 4 x Wire Grid Shelves & 10 x Telescopic Shelves | 56.9kg |
| 4HH-SCS/S2/16W | Static, Double Bay, 16 x Wire Grid Shelves | 38.0kg |
| | | |

Mobile

| Item | Description | Weight |
|-------------------|---|--------|
| Single Bay | | |
| 4HH-SCS/M51/2W5T | Mobile, Single Bay, 2 x Wire Grid Shelves & 5 x Telescopic Shelves | 31.1kg |
| 4HH-SCS/M51/8W | Mobile, Single Bay, 8 x Wire Grid Shelves | 22.3kg |
| Double Bay | | |
| 4HH-SCS/M52/4W10T | Mobile, Double Bay, 4 x Wire Grid Shelves & 10 x Telescopic Shelves | 58.2kg |
| 4HH-SCS/M52/16W | Mobile, Double Bay, 16 x Wire Grid Shelves | 39.3kg |



2HH-7532255 - Wire Grid Shelf



2HH-7532335 - Telescopic Shelves



50mm Braking Pressed Steel Castor Adjustable Feet

Components

- 2HH-7532300-HH Base Unit (x 1)
- 2HH-7532310 Extension Unit (x 1)
- 2HH-7532335 Telescopic Shelf (x 1) • 2HH-7532255 - Wire Grid Shelf (x 1)
- 2HH-7532356 Runners (Pair)
- 2HH-7532320 Adjustable Foot Set (x 4)



Sterile Containers (not supplied)

• 2HH-7532321 - Adjustable Foot Set (x 2)

• 2HH-7532322 - 50mm Pressed Steel Castors (2 x Braking / 2 x Non-Braking) • 2HH-7532323 - 50mm Pressed Steel Castors (1 x Braking / 1 x Non-Braking) • 2HH-7532324 - 80mm Pressed Steel Castors (2 x Braking / 2 x Non-Braking) • 2HH-7532325 - 80mm Pressed Steel Castors (1 x Braking / 1 x Non-Braking)

Shown with BBraun Aesculap Sterile Containers (not supplied)

Room Design & Installation Service Available Contact us for more information



Workstations & Speciality Storage

Workstations 400

Storage Units 401

Catheter Storage Units 402





Other configurations

available. Contact our

customer services

team to discuss your

requirements

Workstations - Mobile, Polymer, Metro Starsys

Features

- Adaptable workstations
- Available as single, double & triple width units
- Moulded body, drawers & doors from a lightweight polymer that will not rust, chip, dent, peel or corrode
- Laminate worktop (height 914mm)
- 125mm swivel castors, 2 off braking

Standard Colours

- Body/Drawers/Doors Light Taupe
- Drawer/Door Handles Slate Blue



| Item | Description | Overall Dimensions (w x d x h) | Weight |
|-------------|---|--------------------------------|---------|
| 8SXRSENT1 | Single, Drawers, Overhead Unit | 600 x 726 x 2070mm | 121.1kg |
| 8SXRDENT2 | Double, Drawers & Door, Overhead Units | 1102 x 726 x 2070mm | 188.2kg |
| 8SXRD40MW10 | Double, Drawers, Overhead Units | 1102 x 726 x 2070mm | 189.0kg |
| 8SXRT38MW14 | Triple, Drawers & Doors, Overhead Units | 1606 x 726 x 2070mm | 219.8kg |

Storage Units - Mobile, Polymer, Metro Starsys

Features

Item

- Full height storage units
- Available as single, double & triple width units
- Moulded body, drawers & hinged/tambour doors from a lightweight polymer that will not rust, chip, dent, peel or corrode
- 125mm swivel castors, 2 off braking

Standard Colours

Options

- Body/Drawers/Doors Light Taupe Drawer/Door Handles - Slate Blue
- 8SXLABKIT Drawer/Door Label Holders (x 10) 8SXREHAN - Push Handle









| Overall Dimensions (w x d x h) | Weight |
|--------------------------------|---------|
| | |
| 578 x 632 x 1990mm | 187.7kg |
| 1080 x 709 x 1990mm | 159.2kg |
| 1669 x 709 x 1990mm | 187.8kg |
| | |
| 578 x 709 x 1990mm | 106.0kg |
| 1080 x 709 x 1990mm | 113.0kg |
| 1669 x 709 x 1990mm | 158.0kg |
| | |
| 578 x 709 x 2096mm | 108.9kg |
| 1080 x 709 x 2096mm | 152.9kg |









Features

- Full height catheter storage units available with either hinged or tambour door/s
- Catheters supported on retractable slides providing easy to access & safe storage
- Available as single & double width units
- Single 180 Catheters
- Double 360 Catheters
- Moulded body & hinged or tambour doors from a lightweight polymer that will not rust, chip, dent, peel or corrode

Standard Colours

- Body/Drawers/Doors Light Taupe
- Drawer/Door Handles Slate Blue

Options

- 8SXRCATHHK Catheter Hooks, (supports 12 catheters)
- SXRLABKIT Drawer/Door Label Holders. (x 10)
- 8SXREHAN Push Handle

Other configurations available. Contact our customer services team to discuss your requirements

125mm swivel castors, 2 off braking

Catheter Storage Units - Stationary, Polymer, Metro Starsys

Features

- Full height catheter storage units available with either hinged or tambour door/s
- Catheters supported on retractable slides providing easy to access & safe storage
- Available as single & double width units Single - 180 Catheters
- Double 360 Catheters
- Moulded body & hinged or tambour doors from a lightweight polymer that will not rust, chip, dent, peel or corrode

Standard Colours

- Body/Drawers/Doors Light Taupe
- Drawer/Door Handles Slate Blue

Options

- 8SXRCATHHK Catheter Hooks, (supports 12 catheters)
- 8SXRLABKIT Drawer/Door Label Holders, (x 10)
- 8SXR-FLR Floor Mounting Kit
- 8SXR72UUFLR Unit to Unit Filler Kit**
- 8SXR72BKFLR30 Back to Wall Kit (Half Depth Units)***
- 8SXR72BKFLR Back to Wall Kit (Full Depth Units)***
- **Covers the seam between two adjacent units
- **Fills the gap when a unit it is pulled away from the wall to be flush with the front of the worktop

| Item | Description | Overall Dimensions (w x d x h) | Weight |
|----------------|---|--------------------------------|---------|
| Hinged Doors | | | |
| 8SXRS72SCATH1 | Single, L/H Hinged Door, 180 x Catheters | 540 x 594 x 1842mm | 75.5kg |
| 8SXRD72SCATH1 | Double, R/H & L/H Hinged Door, 360 x Catheters | 1049 x 594 x 1842mm | 143.6kg |
| Tambour Doors | | | |
| 8SXRS80HCATHTD | Single, Half Depth, Tambour Door, 180 x Catheters | 540 x 457 x 2032mm | 60.3kg |
| 8SXRD80HCATHTD | Double, Half Depth, Tambour Door, 360 x Catheters | 360 1052 x 457 x 2032mm | 113.9kg |
| 8SXRS80SCATHTD | Single, Full Depth, Tambour Door, 180 x Catheters | 180 540 x 620 x 2032mm | 84.4kg |
| 8SXRD80SCATHTD | Double, Full Depth, Tambour Door, 360 x Catheters | 1052 x 620 x 2032mm | 140.2kg |

| Item | Description | Overall Dimensions (w x d x h) | Weight |
|----------------|--|--------------------------------|---------|
| Hinged Doors | | | |
| 8SXRS76MCATH1 | Single, L/H Hinged Door, 180 x Catheters | 578 x 632 x 1990mm | 100.1kg |
| 8SXRD76MCATH1 | Double, R/H & L/H Hinged Door, 360 x Catheters | 1080 x 632 x 1990mm | 136.3kg |
| Tambour Doors | | | |
| 8SXRS83MCATHTD | Single, Tambour Door, 180 x Catheters | 578 x 709 x 2096mm | 109.3kg |
| 8SXRD83MCATHTD | Double, Tambour Doors, 360 x Catheters | 1080 x 709 x 2096mm | 146.1kg |
| | | | |

8SXRD76MCATH1

5

ő

øð

N







| â | | Caratray Deaks | 40 40 50 50 | Cover | 05 | HTM71 | | Tiered | L Filing Trallave | 250 | Proposation Tables | | Dias & Desline Destars Armsheire | 202.204 | Courte | 154-173 |
|--|--------------------|-------------------------------------|------------------------------|---|----------------------|--|------------------|----------------------|-----------------------------------|------------------|-------------------------------------|--------------------|---|------------------|------------------------------------|---------------|
| A Abloy Protec Clig | 115, 140-141 | Caretray Racks Caretray Trolleys | 46-49, 56-59 40-45, 50-55 | Cover Fixed Height | 85 85 | Carts | 192-197 | | l Filing Trolleys fer Trolleys | 250 249-251 | Preparation Tables Fixed Height | 371 | Rise & Recline Porters Armchairs Saddle Stools | 323-324 306 | Couch Trolleys | 154-173 |
| Airways / Difficult Intubation Trolley | 287 | Chart | 72-73 | Variable Height | 84 | Distribution Trollevs | 197-198 | | Filing Cabinet | 251 | Variable Height | 370 | Soft Seating | 330-331, 334-335 | Airways / Difficult Intubation | 287 |
| Alcohol Gel Holders | 96 | Caretray Options | 60-69 | Controlled Drug Cabinets | 137-139 | Open Transfer/Transport Trolleys | 197-198 | | ne Bedside Cabinets | 22-29 | Pump Stands | 208-209 | Static Safe Tech Chairs | 344 | Aluminium Trays | 189 |
| Apron Dispensers | 94 | Carts | 00-03 | Couch | 107-100 | Racks | 200-201 | | ispensing Trolleys | 114 | R | 200-203 | Static Safe Tech Stools | 345 | Caretrav | 40-45, 50-55 |
| Arm / Leg Rest | 162 | Basix | 276-277 | Examination | 154-173 | Travs, Dividers & Accessories | 202-205 | | Infusion Stands | 208-214 | Rack | | Static Safe & Sterile Tech Chairs | 346 | Chart | 72 |
| R Rest | 102 | Flexline | 272-273, 286 | Fixed Height | 161 | Holder | 202-203 | Easy C | | 208-211 | Bedpan Drainage | 00 | Static Safe & Sterile Tech Stools | 347 | Clean Linen | 222, 224-227 |
| Baby | | HTM71 | 192-197 | Options | 181-183 | Alcohol Gel | 96 | Mild St | | 213 | Caretray | 46-49, 56-59 | Tech Chairs | 340 | Compact | 197 |
| Bath | 233 | Lifeline | 284-285 | Steps | 163 | Bag | 290-297 | Knock | | 214 | Catheter | 40-40, 00-00 Q8 | Tech Stools | 341 | Concealment | 84-85 |
| Cribs | 230-232 | Polymer | 272-275 | Therapy | 170-173 | Catheter Bag | 97 | | ess Steel | 212 | HTM71 | 200-201 | Twenty Four Hour Task & | | Cylinder | 88-90 |
| Mattress | 230-232 | Procedure | 272-277 | Treatment | 154-173 | Cylinder | 91, 266 | | Pump Stands | 208-209 | Urinal Drainage | 99 | Operators Chairs | 302-303 | Dispensing | 112-114 |
| Baby Cribs | 200 202 | Starsys | 274-275 | Variable Height | 154-173 | Document | 77 | | Storage Units | 401 | Racking | 00 | Visitor & Waiting Room Seating | 332-333 | Distribution | 129, 197, 198 |
| Variable Height | 230-231 | Catheter | 2 2. 0 | Couch Steps | 163 | Glove Box | 95 | | Workstations | 400 | HTM71 | 200-201 | Voyage Beam | 338-339 | Dressing | 102-106 |
| Fixed Height | 232 | Bag Holders | 97 | Curtain Screens | 80 | Mask | 94 | | r Seating | 330-331 | Polymer | 360-363 | Waiting Room Seating | 334-339 | General Purpose | 186-189 |
| Bag & Sackholders | 202 | Racks | 98 | Cylinder | | Sackholders | 290-297 | 0 | locality | 000 001 | Top Track High Density Storage | 366 | Zeta Modular Seating | 330-331 | Hamper | 223 |
| Free Standing | 295 | Storage Units | 402-403 | Holders | 91.266 | Sharps Box | 97 | Office, T | Task & Operators Chairs | 300-304 | Zinc Chrome | 364-365 | Security Clamp | 128 | Instrument | 102-106 |
| Hamper | 223 | Chairs | | Trolleys | 88-90 | Tab | 246 | Operatin | ng Platforms | 394 | Reclining Armchairs | 310-321 | Self-Administration Cabinets | 146-149 | Laminate Shelves | 189 |
| Removable Body | 290-294 | Bariatric | 15-16, 304-305 | Wall Brackets | 91 | Urine Bottle | 97 | | d & Chair Tables | | Reclining Bariatric Armchairs | 320-321 | Sepsis Trolleys | 352-355 | Medical Records | 236-251 |
| Plastic | 296-297 | Office, Task & Operators | 300-304 | Wall Plate | 91 | Wall Mounted Infusion | 215 | Econor | omy | 30 | Reclining Porters Armchair | 322-323 | Sharps Box Holder | 97 | Metal Dispensing | 114 |
| Soiled Linen | 218-223 | Low Back Patient Chairs | 308 | D | | I | | Lamina | ate | 32-33 | Refrigerators | 151 | Shelving | | Metal Distribution | 129 |
| Bariatric | | Mata Chairs | 332 | Document Holders | 77 | Infusion | | Polyme | ier | 31 | Residents Own Medication Cabinets | 146-149 | HTM71 | 200-201 | Metal Pharmacy | 116-125 |
| Commode | 11, 14 | Medical Grade Clean Room | 348 | Dispensers | | Clamps | 109, 215 | Overnigh | ht Chair | 238-329 | Rest | | Polymer | 360-363 | Mobile Workstation | 271 |
| Couches | 164-169 | Tech Chairs | 340 | Apron | 94 | Poles | 214 | P | | | Arm | 162 | Top Track High Density Storage | 366 | Open Transfer / Transport Trolleys | 199 |
| Folding Wheelchair | 8-9 | Medium & High Back | 309 | Glove | 95 | Stands | 208-214 | Paper Ta | ābs | 246 | Leg | 162 | Zinc Chrome | 364-365 | Pharmacy, Blister Pack (MDS) | 120-121 |
| Four Leg Chairs | 16, 305 | Patient Armchairs | | Syringe | 96 | Infusion Stands | | | Armchairs | 309 | Resuscitation Carts | 284-286 | Soft Seating | 330-331, 334-335 | Pharmacy, Biodose & Multimeds | 124 |
| High Back Chair | 15, 304 | Overnight Chair | 328 | Dispensing Trolleys | | Easy Clean | 208-211 | Patient (| | 308 | Resuscitation Trolleys | 280-283, 287 | Soiled Linen Bags | 223 | Pharmacy, Drawers | 118 |
| Phlebotomy Chair | 178, 257 | Phlebotomy | 178-180, 256-259 | High Capacity | 113 | Knockdown | 214 | | own Medication Cabinets | 146-149 | Rise & Recline Porters Armchair | 324-325 | Soiled Linen Trolleys | 218-222 | Pharmacy, Medinoxx | 125 |
| Portering Chairs | 10, 265 | Portering | 10, 262-267 | Security Clamp | 128 | Mild Steel | 213 | Perch St | | 307 | Room Dividers | 78-79 | Solid Screens | 81 | Pharmacy, Nomad | 122-123 |
| Rollator | 17 | Prima Chairs | 333 336-337 | Standard Capacity | 112, 114 | Stainless Steel | 212 | | al Effects Lockers | 384-385 | S | | Stand | 000.004 | Pharmacy, Patient Administration | 116-117 |
| Seating | 16, 304-305 | Prima Beam PU Saddle & Perch | 307 | Distribution Trolleys | 129 | Trolley Clamping | 109, 215 | | icy Trolleys | 101 | Sackholders | 290-297 | Bowl | 388-391 | Pharmacy, Trays | 119 |
| Shower Chair | 11-14 17 | PU Tech Chairs | 342 | Dressing Trolleys | 107 | Wall Mounted | 215 | | se® & MultiMeds | 124 | 20/50 Litre | 290-291 | Kickabout Bowl Kickabout Bucket | 391 391 | Phlebotomy | 254-255 |
| Training Suit | 17 | PU Tech Stools | 343 | Drawers Fixed Shelves | 107 104-106, 392-393 | Instrument Trolleys | 107 | | Pack (MDS) | 120-121 118 | 75/90 Litre | 292-293 | Mobile Infusion | 210-214 | Plaster | 270 188 |
| Treatment Chairs | 1/4-1/0 | Reclining Armchairs | 310-321 | | 104-106, 392-393 | Drawers Fixed Shelves | 107 | Drawe | | 126 | 75 Litre Front Opening | 294 295 | Mobile Pump | 208-209 | Polymer Shelves | 280-283, 287 |
| Walking Frame Beam Seating | 336-339 | Reclining Amerians | 320-321 | Options Removable Shelves | 102-103 | Options | 104-106, 392-393 | Medino | p Holders | 125 | Free Standing Fixed Body Plastic | 295 296-297 | Sterile Container Storage | 396-397 | Resuscitation Sepsis | 352-355 |
| Bedpan Drainage Rack | 330-339 | Reclining Porters Armchairs | 322-323 | Drug & Medicine | 102-103 | Removable Shelves | 102-103 | | n C5 Holder | 125 | Saddle Stools | 306-307 | Sterile Supplies Trolleys | 327-375 | Soiled Linen | 218-222 |
| Bedside Cabinets | 33 | Rise & Recline Porters Armchairs | 323-324 | Blister Pack Cabinets (MDS) | 136 | | 102-103 | | d Cassettes | 122-123 | ScanCell Carts | 192-197 | Storage | 521-515 | Stainless Steel Shelves | 186-187 |
| Melamine | 22-29 | Saddle Stools | 306 | Controlled Drug Cabinets | 137-139 | Key Cabinets | 382-383 | | nt Administration | 116-117 | Screens | 152-157 | Boxes | 150 | Sterile Supplies | 372-375 |
| Polymer | 20-21 | Soft Seating | 330-331, 334-335 | Cupboards | 142-145 | Kickabout Bowl Stand | 391 | | ity Clamp | 128 | Curtain | 80 | Cabinets | 380 | Suture Storage | 395 |
| Bins | 20 2 1 | Static Safe Tech Chairs | 344 | Deep Shelf Cabinets | 134-135 | Kickabout Bucket Stand | 391 | Trays | | 119 | Solid | 81 | Catheter | 402-403 | Wooden Dispensing | 112-113 |
| Freestanding | 295 | Static Safe Tech Stools | 345 | Shelves & Door Shelf Cabinets | 132-133 | L | | | cy & Vaccine Refrigerators | 151 | Seating | | Cupboards | 378-379 | U | |
| Removable Body | 290-294 | Static Safe & Sterile Tech Chairs | 346 | Refrigerators | 151 | Laminate Overbed & Chair Tables | 32-33 | | omy Trolleys | 254-255 | Bariatric | 15-16, 304-305 | Lockers | 384-385 | Urinal Drainage Rack | 99 |
| Plastic | 296-297 | Static Safe & Sterile Tech Stools | 347 | Drug & Medicine Cabinets | 132-136 | Laptop Holders | 126, 244 | Phleboto | omy Chair | | Office, Task & Operators | 300-304 | Mobile Units | 401 | Urine Bottle Holder | 97 |
| Body Tray | 85 | Tech Chairs | 340 | Drug & Medicine Cupboards | 142-145 | Leg Rest | 162 | Fixed I | Height | 180, 258-259 | Low Back Patient Chairs | 308 | Sterile Container | 396-397 | V | |
| Bowl Stands | 388-391 | Tech Stools | 341 | E | | Linen | | Variabl | ble Height | 178-179, 256-257 | Mata Chairs | 332 | Suture | 395 | Variable Height Examination | 154-183 |
| Bucket Stands | 391 | Twenty Four Hour Task & Opera- | 302-303 | Economy Overbed Table | 30 | Clean Linen Trolleys | 222, 224-227 | Plaster 1 | Trolleys | 270 | Medical Grade Clean Room | 348 | Suture Storage Trolleys | 395 | Couches & Chairs | |
| С | | tors Chairs | | Examination Chairs | 174-177 | Soiled Linen Bags | 223 | Plastic B | | 296-297 | Tech Chairs | 0-10 | Syringe Dispenser | 96 | Visitor & Waiting Room Seating | 332-333 |
| Cabinets | | Visitor & Waiting Room Seating | 332-333 | Examination Couches | 154-173 | Soiled Linen Trolleys | 218-222 | Plastic C | | 230-232 | Medium & High Back | 309 | Τ | | W | |
| Bedside, Melamine | 22-29 | Voyage Beam | 338-339 | F | | Locker | | | ure Carts | 272-277 | Patient Armchairs | | Tab Holders | 246 | Waiting Room Seating | 334-339 |
| Bedside, Polymer | 20-21 | Waiting Room Seating | 334-339 330-331 | Filing Pockets | 246 | Personal Effects | 384-385 | | Mask Holder | 94 | Overnight Chair | 328 | Table | | Wall Mounted Infusion Holders | 215 |
| Clinical Reagent | 381 | Zeta Modular Seating | 76 | Fixed Height Examination Couches | | Storage | 384-385 | | r Bedside Cabinets | 20-21 | Phlebotomy | 178-180, 256-259 | Laminate Overbed & Overchair | 32-33 | Waste Bins | 290-297 |
| Controlled Drug | 137-139 | Chart Boards Chart Trolleys | 70 72-73 | Fixed Height Phlebotomy Chairs | 180, 258-259 | Maak Halder | 0.4 | | r Overbed & Chair Tables | 31 | Portering Prima Chairs | 10, 262-267 333 | Mayo Belymer Overhed & Oversheir | 394 31 | Wooden Dispensing Trolleys | 112-113 |
| Drug & Medicine | 132-136 | Chart Workstations | 12-13 | First Aid Cabinets | 381 | Mask Holder | 94 | | r Shelving & Racking | 360-363 | Prima Chairs Prima Beam | 336-337 | Polymer Overbed & Overchair | 31 | X X Day Filing Cabinat | 251 |
| First Aid | 381 382-383 | Fixed Height | 74 | Free Standing Sackholders | 295 | Mayo Tables | 394 | | r General Purpose Trolleys | 188 | Prima Beam PU Saddle & Perch | 336-337 307 | Preparation Tablet Counter | 370-371 150 | X-Ray Filing Cabinet | 201 |
| Key Rationto Own | 382-383 146-149 | Variable Height | 74 | General Burnana Trallava | 186-189 | Medical Records Carefile Filing Pockets | 246 | Portering | | 262-267 | PU Tech Chairs | 342 | Technical Seating | 340-348 | Zinc Chrome Shelving & Racking | 364-365 |
| Patients Own Residents Own | 146-149 | Clamping Infusion Stand Trolley | 109, 215 | General Purpose Trolleys Glove Box Holders | 95 | Laptop Holders | 246 | Portering Bariatr | | 10, 265 | PU Tech Stools | 343 | Therapy Couch | 170-173 | Zine Gritome Sheiving & Racking | 304-303 |
| Self Administration | 146-149 | Clean Linen Trollevs | 222, 224-227 | Glove Box Holders Glove Box Dispenser | 95 95 | Motion C5 Holder | 244 245 | Front S | | 263-264 | Reclining Armchairs | 310-321 | Transfer & Transport Trollevs | 199 | | |
| Storage | 380 | Concealment Trolleys | , r <i></i> r | Н | 55 | Trolleys | 236-241 | Option | | 266-267 | Reclining Bariatric Armchairs | 320-321 | Treatment | 100 | | |
| X-Ray Filing | 251 | Body Tray | 85 | Hamper Bag Trolleys | 223 | Workstations | 242-243 | Rear S | | 262 | Reclining Porters Armchairs | 322-323 | Chair | 174-177 | | |
| st tay timig | 201 | , , | | Lamps. Dag Honoyo | 220 | | 2.2.2.10 | itear e | | | | | | | | |

E

| 0-9 | 4A2DR15 194 | 4B1VD10/E 195 | 5162 301 | 5RE/RRAE/ALWE/VVDEL/VPNB 313 | 5X-0299-000-010 17 | 8FLCOM-LT 272 | 8Q556G2 362 |
|------------------|---------------|-----------------------|------------------------------|---------------------------------|--------------------|--------------------|--------------------|
| 2FACO024 85 | 4A2DR18 194 | 4B1VD10/F 195 | 5163 301 | 5RE/RRAE/ASHI/VVDEL 316 | 5ZETA1 330 | 8FLCRP4 272 | 8Q556G2M 362 |
| | 4A2DSTM 198 | 4B1VD12/G 195 | 5201 307 | 5RE/RRAE/ASHI/VVDEL/VPNB 317 | 5ZETA1/B 330 | 8FLPED1 272 | 8Q566G2 362 |
| 2HH-100072 205 | | | | | | | |
| 2HH-240501 205 | 4A2DSTS 198 | 4B1VD12/H 195 | 5202 306 | 5RE/RRAE/BECK/VMULT 318 | 5ZETA1/CO 330 | 8FLTMENT2 272 | 8Q566G2M 362 |
| 2HH-240820 203 | 4A2OF09 192 | 4B1VD15/I 195 | 5205 307 | 5RE/RRAE/HADL610/VVDEL/VPNB 320 | 5ZETA2 330 | 8L24N-4C 365 | 8SXRD40MW10 400 |
| 2HH-240820 204 | 4A2OF10 192 | 4B1VD15/J 195 | 5260 302 | 5RE/RRAE/HADL710/VVDEL/VPNB 320 | 5ZETA2/2B 330 | 8L36N-4C 365 | 8SXRD72SCATH1 403 |
| | | | | | | | |
| 2HH-240821 203 | 4A2OF12 192 | 4B2DR09 194 | 5262 302 | 5RE/RRAE/SHAP/VVDEL 310 | 5ZETA2/B/LC 330 | 8L42N-4C 365 | 8SXRD76CM3 401 |
| 2HH-240821 204 | 4A2OF15 192 | 4B2DR10 194 | 5301 333 | 5RE/RRAE/SHAP/VVDEL/VPNB 311 | 5ZETA2/B/RC 330 | 8L48N-4C 365 | 8SXRD76CMHD1 401 |
| 2HH-241301 203 | 4A2OF18 192 | 4B2DR12 194 | 5301/A 333 | 5RE/RRAE/TWYF/VVDEL 314 | 5ZETA930C 330 | 8L60N-4C 365 | 8SXRD76MCATH1 402 |
| | 4A2OTTM 199 | 4B2DR15 194 | 5670 16 | 5RE/RRAE/TWYF/VVDEL/VPNB 315 | 5ZETAC45 330 | 8LECCRP2 284 | 8SXRD80HCATHTD 403 |
| 2HH-241301 204 | | | | | | | |
| 2HH-241302 203 | 4A2OTTS 199 | 4B2DR18 194 | 5670 305 | 5RE/RRPAH/DORC/VVDEL 324 | 5ZETAC45/B 330 | 8LECCRP3 284 | 8SXRD80SCATHTD 403 |
| 2HH-241302 204 | 4A2RS09 193 | 4B2OF09 192 | 5671 16 | 5RE/RRPAH/DORC/VVDEL/VPNB 325 | 5ZETAC90 330 | 8LECCRP4 284 | 8SXRD83MCATHTD 402 |
| 2HH-7011050 202 | 4A2RS10 193 | 4B2OF10 192 | 5671 305 | 5RLDF0119 151 | 5ZETAC90/B 330 | 8MBC1210DB 276 | 8SXRD83MGSTD 401 |
| | | | | | | | |
| 2HH-7011100 202 | 4A2RS12 193 | 4B2OF12 192 | 59003 303 | 5RLDF0119DIGILOCK 151 | 5ZETAH 330 | 8MBX1310ANES4 276 | 8SXRDENT2 400 |
| 2HH-7012100 202 | 4A2RS15 193 | 4B2OF15 192 | 5660XL 15 | 5RLDF0219 151 | 5ZETAH/B 330 | 8MBX3210TRAUMA 276 | 8SXRS72SCATH1 403 |
| 2HH-7012200 202 | 4A2RS18 193 | 4B2OF18 192 | 5660XL 304 | 5RLDF0219DIGILOCK 151 | 5ZETAQ 330 | 8MQD24-8 363 | 8SXRS76CM5 401 |
| 2HH-7014050 202 | 4A3DR09 194 | 4B2RS09 193 | 5AM3XRAY 251 | 5RLDF0519 151 | 85MX536MG 360 | 8MQL24-2S 363 | 8SXRS76MCATH1 402 |
| | | | | | | 0MQL24-20 303 | |
| 2HH-7014100 202 | 4A3DR10 194 | 4B2RS10 193 | 5FSH20WH 296 | 5RLDF0519DIGILOCK 151 | 85MX536MS 360 | 8MQL36-2S 363 | 8SXRS80HCATHTD 403 |
| 2HH-7014200 202 | 4A3DR12 194 | 4B2RS12 193 | 5HORI/1/A 335 | 5RLDG0519 151 | 85MX536SG 360 | 8MQL42-2S 363 | 8SXRS80SCATHTD 403 |
| 2HH-7810102 203 | 4A3DR15 194 | 4B2RS15 193 | 5HORI/2/A 335 | 5RLDG0519DIGILOCK 151 | 85MX536SS 360 | 8MQL48-2S 363 | 8SXRS83MCATHTD 402 |
| | 4A3DR18 194 | 4B2RS18 193 | 5HORI/3/A 335 | 5RSH42WH 297 | 85MX546MG 360 | 8MQL60-2S 363 | 8SXRS83MGSTD 401 |
| 2HH-7810103 203 | | | | | | | |
| 2HH-7810104 204 | 4A3OF09 192 | 4B2TD09/AA * 196 | 5ND/HAF/ALDE/VVDEL 309 | 5RSH70WH 297 | 85MX546MS 360 | 8MQLS24-2S 363 | 8SXRSANES 274 |
| 2HH-7810105 204 | 4A3OF10 192 | 4B2VD09/AA* 195 | 5ND/HAR/ALDE/VVDEL 309 | 5SL6 332 | 85MX546SG 360 | 8MX536MG 360 | 8SXRSARTH 274 |
| 2HH-7811204N 204 | 4A3OF12 192 | 4B2VD10/DD * 195 | 5ND/LC/ALDE/VVDEL 308 | 5SML6 332 | 85MX546SS 360 | 8MX536MS 360 | 8SXRSCAST 274 |
| | 4A30F15 192 | 4B3DR09 194 | 5ND/MAF/ALDE/VVDEL 309 | 5TC101 340 | 85MX556MG 360 | 8MX536SG 360 | 8SXRSCCU 274 |
| 2HH-7811205N 203 | | | | | | | |
| 2HH-7811205N 204 | 4A3OF18 192 | 4B3DR10 194 | 5ND/MAR/ALDE/VVDEL 309 | 5TC102 340 | 85MX556MS 360 | 8MX536SS 360 | 8SXRSDRS 274 |
| 2HH-7811206N 203 | 4A3OTTM 199 | 4B3DR12 194 | 5NOVA/1 334 | 5TC103 340 | 85MX556SG 360 | 8MX546MG 360 | 8SXRSENT1 400 |
| | 4A3OTTS 199 | 4B3DR15 194 | 5NOVA/2 334 | 5TC108 341 | 85MX556SS 360 | 8MX546MS 360 | 8SXRSGS2 401 |
| 2HH-7811304 204 | | | | | | | |
| 2HH-7811305 203 | 4A3RS09 193 | 4B3DR18 194 | 5PF332 336 | 5TC109 341 | 85MX566MG 360 | 8MX546SG 360 | 8SXRSIV 274 |
| 2HH-7811305 204 | 4A3RS10 193 | 4B3OF09 192 | 5PF333 336 | 5TC110 341 | 85MX566MS 360 | 8MX546SS 360 | 8SXRSTRAMA 274 |
| 2HH-7811306 203 | 4A3RS12 193 | 4B3OF10 192 | 5PF334 336 | 5TC201 344 | 85MX566SG 360 | 8MX556MG 360 | 8SXRT38MW14 400 |
| | | | 5PF335 336 | 5TC202 344 | 85MX566SS 360 | 8MX556MS 360 | |
| 2HH-7910201N 205 | 4A3RS15 193 | 4B3OF12 192 | | | | | 8SXRTGS2 401 |
| 2HH-7910204N 205 | 4A3RS18 193 | 4B3OF15 192 | 5PT336 337 | 5TC203 344 | 850537C 364 | 8MX556SG 360 | 8SXRTGS3 401 |
| 2HH-K100187 205 | 4B1DR09 194 | 4B30F18 192 | 5PT337 337 | 5TC208 345 | 850537CM 364 | 8MX556SS 360 | В |
| 2HH-K561039D 202 | 4B1DR10 194 | 4B3RS09 193 | 5PT338 337 | 5TC209 345 | 850547C 364 | 8MX566MG 360 | BA/3/B 232 |
| | | | | | | | |
| 2LOCK217 267 | 4B1DR12 194 | 4B3RS10 193 | 5PT339 337 | 5TC210 345 | 850547CM 364 | 8MX566MS 360 | BA/6 232 |
| 2LOCK218 267 | 4B1DR15 194 | 4B3RS12 193 | 5PU151 343 | 5TC304 346 | 850557C 364 | 8MX566SG 360 | BA/6/A 232 |
| 2LOKE011 267 | 4B1DR18 194 | 4B3RS15 193 | 5PU152 343 | 5TC305 346 | 850557CM 364 | 8MX566SS 360 | BA/VH 230 |
| 2LOKE134 21 | 4B1DSTM 198 | 4B3RS18 193 | 5PU153 343 | 5TC306 346 | 850967C 364 | 8MXD24-8 361 | BA/VH/NC 230 |
| | | | | | | | |
| 2MECO079 388 | 4B1DSTS 198 | 4B3VD09/AAA ** 195 | 5PU4 342 | 5TC308 347 | 850967CM 364 | 8MXL36-2S 361 | BA015 233 |
| 2PLCO026 281 | 4B10F09 192 | 4B3VD10/DDD ** 195 | 5PU5 342 | 5TC309 347 | 85Q537G2 362 | 8MXL42-2S 361 | BA030 231 |
| 2PLCO054 388 | 4B10F10 192 | 4D1DSTCM 197 | 5PU6 342 | 5TC310 347 | 85Q537G2M 362 | 8MXL48-2S 361 | BC1C/BE 22 |
| | 4B10F12 192 | 4D1DSTCS 197 | 5RE/OC/PORT/VVDEL 328 | 5TC311 348 | 85Q547G2 362 | 8MXL60-2S 361 | BC1C/GW 22 |
| 4A1DR09 194 | | | | | | | |
| 4A1DR10 194 | 4B10F15 192 | 4D1TD09/A 197 | 5RE/OC/PORT/VVLAVSAVL 328 | 5TC312 348 | 85Q547G2M 362 | 8MXLS24-2S 361 | BC1CPD/BE 22 |
| 4A1DR12 194 | 4B10F18 192 | 4D1TD09/B 197 | 5RE/RAE/ALWE/VVDEL 312 | 5TC313 348 | 85Q557G2 362 | 8MY1627-24BU 188 | BC1CPD/GW 22 |
| 4A1DR15 194 | 4B1RS09 193 | 4HH-HTM71/M51/400 201 | 5RE/RAE/ALWE/VVDEL/VPNB 313 | 5VOYAGE2 338 | 85Q557G2M 362 | 8MY1627-34BU 188 | BC1CPDWL/BE 27 |
| | | | | | | | |
| 4A1DR18 194 | 4B1RS10 193 | 4HH-HTM71/M51/600 201 | 5RE/RAE/ASHI/VVDEL 316 | 5VOYAGE3 338 | 85Q567G2 362 | 8MY2030-24BU 188 | BC1CPDWL/GW 27 |
| 4A1DSTM 198 | 4B1RS12 193 | 4HH-HTM71/M52/400 201 | 5RE/RAE/ASHI/VVDEL/VPNB 317 | 5VOYAGE4 338 | 85Q567G2M 362 | 8MY2030-34BU 188 | BC1CPDWR/GW 27 |
| 4A1DSTS 198 | 4B1RS15 193 | 4HH-HTM71/M52/600 201 | 5RE/RAE/BECK/VMULT 318 | 5X-0100-061-000/K 8 | 89989PX 361 | 8MY2636-25BU 188 | BC1CWL/BE 27 |
| 4A10F09 192 | 4B1RS18 193 | 4HH-HTM71/S1/400 200 | 5RE/RAE/SHAP/VVDEL 310 | 5X-0100-071-000/K 8 | 89990P 365 | 8MY2636-35BU 188 | BC1CWL/GW 27 |
| | | | | | | | |
| 4A10F10 192 | 4B1TD09/A 196 | 4HH-HTM71/S1/600 200 | 5RE/RAE/SHAP/VVDEL/VPNB 311 | 5X-0109-061-000/K 8 | 8BC2030-2DMB 188 | 80536C 364 | BC1CWR/BE 27 |
| 4A10F12 192 | 4B1TD09/B 196 | 4HH-HTM71/S2/400 200 | 5RE/RAE/TWYF/VVDEL 314 | 5X-0109-071-000/K 8 | 8BC2030-3DMB 188 | 80536CM 364 | BC1CWR/GW 27 |
| 4A10F15 192 | 4B1TD09/C 196 | 4HH-HTM71/S2/600 200 | 5RE/RAE/TWYF/VVDEL/VPNB 315 | 5X-0110-061-000/K 9 | 8BC2636-2DMB 188 | 80546C 364 | BC1CWRPD/BE 27 |
| 4A10F18 192 | 4B1TD10/D 196 | 4HH-SCS/M51/2W5T 396 | 5RE/RAM/ALWE/VVDEL 312 | 5X-0110-071-000/K 9 | 8BC2636-3DMB 188 | 80546CM 364 | BC1D/BE 22 |
| | | | | | | | |
| 4A10TTM 199 | 4B1TD10/E 196 | 4HH-SCS/M51/8W 396 | 5RE/RAM/ALWE/VVDEL/VPNB 313 | 5X-0142-061-000 14 | 8DC15EC 395 | 80556C 364 | BC1D/GW 22 |
| 4A10TTS 199 | 4B1TD10/F 196 | 4HH-SCS/M52/16W 396 | 5RE/RAM/ASHI/VVDEL 316 | 5X-0142-071-000 14 | 8DC16EC 395 | 80556CM 364 | BC1DPD/BE 22 |
| 4A1RS09 193 | 4B1TD10/G 196 | 4HH-SCS/M52/4W10T 396 | 5RE/RAM/ASHI/VVDEL/VPNB 317 | 5X-0150-061-000 11 | 8DC35EC 395 | 80966C 364 | BC1DPD/GW 22 |
| 4A1RS10 193 | 4B1TD12/H 196 | 4HH-SCS/S1/2W5T 396 | | 5X-0150-071-000 11 | 8DC36EC 395 | 80966CM 364 | BC1DPDWL/BE 27 |
| | | | 5RE/RAM/SHAP/VVDEL 310 | | | | |
| 4A1RS12 193 | 4B1TD12/I 196 | 4HH-SCS/S1/8W 396 | 5RE/RAM/SHAP/VVDEL/VPNB 311 | 5X-0151-061-000 11 | 8DC55EC 395 | 8Q04LH 363 | BC1DPDWL/GW 27 |
| 4A1RS15 193 | 4B1TD12/J 196 | 4HH-SCS/S2/16W 396 | 5RE/RAM/TWYF/VVDEL 314 | 5X-0151-071-000 11 | 8DC56EC 395 | 8Q04LHC 363 | BC1DPDWR/BE 27 |
| 4A1RS18 193 | 4B1VD09/A 195 | 4HH-SCS/S2/4W10T 396 | 5RE/RAM/TWYF/VVDEL/VPNB 315 | 5X-0155-066-000 12 | 8DD24C 365 | 8Q536G2 362 | BC1DPDWR/GW 27 |
| | 4B1VD09/B 195 | 5160 301 | | 5X-0156-066-000 13 | 8FL21P3369 286 | 8Q536G2M 362 | BC1DWL/BE 27 |
| 4A2DR09 194 | | | 5RE/RPAH/BATH/VVDEL 322 | | | | |
| 4A2DR10 194 | 4B1VD09/C 195 | 5160D 301 | 5RE/RPAH/BATH/VVDEL/VPNB 323 | 5X-0170-050-000 17 | 8FLCCU1 272 | 8Q546G2 362 | BC1DWL/GW 27 |
| 4A2DR12 194 | 4B1VD10/D 195 | 5161 301 | 5RE/RRAE/ALWE/VVDEL 312 | 5X-0180-076-000 17 | 8FLCLINE 272 | 8Q546G2M 362 | BC1DWR/BE 27 |
| | | | | | | | |

406

| 254 CRMS11/16/355D CTCHARM210/852D 73 284 CRMS11/16/35D S8 CTD1/082P 60 284 CRMS11/16/35D S8 CTD1/082P 60 284 CRMS11/16/054D S8 CTD1/102CA 60 284 CRMS21/61/1258D S9 CTD1/102CA 60 108 CRMS21/61/165 S9 CTD1/102SA 60 108 CRMS21/61/615 S9 CTD1/102SA 60 69 CRMS21/61/615 S9 CTD1/142PB 60 69 CRSM11/8/35D4 6 CTD2/08CA 60 108 CRSM11/8/35D4 6 CTD2/102B4 60 75 CRSM1/18/85D4 6 CTD2/102A60 60 75 CRSM1/18/80540 7 CTD2/142A60 60 76 CRSM2/18/105104 7 CTD2/142A60 60 76 CRSM2/18/105104 7 CTM1/108/152D 50 76 CRSM2/18/1682D47 CTM1/108/152D 50 | | | |
|--|-----|----------------------|--------------------------|
| 294 CRMS1/16/655D 58 CTD1/08CA 60 294 CRMS1/16/65D 58 CTD1/08CS 60 294 CRMS2/16/105105 59 CTD1/10CA 60 295 CRMS2/16/105105 59 CTD1/10CA 60 108 CRMS2/16/105105 59 CTD1/10CA 60 108 CRMS2/16/105105 59 CTD1/14CPB 60 69 CRSM1/18/256D 46 CTD2/08CA 60 69 CRSM1/18/256D 46 CTD2/14CPB 60 108 CRSM1/18/256D 46 CTD2/10CA 60 75 CRSM1/18/256D 46 CTD2/10CA 60 75 CRSM1/18/855D 46 CTD2/10CA 60 75 CRSM1/18/854D 46 CTD2/14CPB 60 76 CRSM1/18/854D 47 CTD2/14CPB 60 76 CRSM2/18/1051D 47 CTD2/14CPB 60 76 CRSM2/18/160 47 CTM1/08/152D 50 74 CRSM2/18/160 47 CTM1/08/152D 50 74 CRSM2/18/160 47 CTM1/08/152D 50 74 CRSM2/18/160 48 CTM1/10/152D 52 74 CRSM2/18/160 48 CTM1/10/152D 52 74 CRSS1/16/16/85D | 293 | CRMS1/16/3S6D 58 | CTCHARM2/10/6S4D 73 |
| 294 CRMS1/16/655D 58 CTD1/08CA 60 294 CRMS1/16/65D 58 CTD1/08CS 60 294 CRMS2/16/105105 59 CTD1/10CA 60 295 CRMS2/16/105105 59 CTD1/10CA 60 108 CRMS2/16/105105 59 CTD1/10CA 60 108 CRMS2/16/105105 59 CTD1/14CPB 60 69 CRSM1/18/256D 46 CTD2/08CA 60 69 CRSM1/18/256D 46 CTD2/14CPB 60 108 CRSM1/18/256D 46 CTD2/10CA 60 75 CRSM1/18/256D 46 CTD2/10CA 60 75 CRSM1/18/855D 46 CTD2/10CA 60 75 CRSM1/18/854D 46 CTD2/14CPB 60 76 CRSM1/18/854D 47 CTD2/14CPB 60 76 CRSM2/18/1051D 47 CTD2/14CPB 60 76 CRSM2/18/160 47 CTM1/08/152D 50 74 CRSM2/18/160 47 CTM1/08/152D 50 74 CRSM2/18/160 47 CTM1/08/152D 50 74 CRSM2/18/160 48 CTM1/10/152D 52 74 CRSM2/18/160 48 CTM1/10/152D 52 74 CRSS1/16/16/85D | 294 | CRMS1/16/4S5D 58 | CTCHARM2/10/8S2D 73 |
| 294 CRMS1/16/805 8 CTD1/08EPB 60 294 CRMS1/16/105100 59 CTD1/10CA 60 300 CRMS2/16/105100 59 CTD1/10EPB 60 108 CRMS2/16/105105 59 CTD1/10EPB 60 108 CRMS2/16/16512 59 CTD1/14EPB 60 69 CRMS2/16/65120 59 CTD1/14EPB 60 69 CRSM1/18/2560 46 CTD2/08EPB 60 108 CRSM1/18/2550 46 CTD2/08EPB 60 108 CRSM1/18/2550 46 CTD2/08EPB 60 108 CRSM1/18/2550 46 CTD2/10CA 60 75 CRSM1/18/2550 47 CTD2/14EPB 60 75 CRSM1/18/2580 47 CTD2/14EPB 60 76 CRSM2/18/15280 47 CTD2/14EPB 60 76 CRSM2/18/1520 47 CTM1/10/83D 50 74 CRSM2/18/1520 47 CTM1/10/83D 50 74 CRSS1/16/2560 48 CTM1/10/83D 50 74 CRSS1/16/650 48 CTM1/10/83D 52 74 CRSS1/16/650 48 CTM1/10/83D 52 74 CRSS1/16/650 48 CTM1/10/83D 52 74 CRSS1/16/65 | | | |
| 224 CRMS1/16/80 58 CTD1/10EX 60 30 CRMS2/16/105100 59 CTD1/10CA 60 108 CRMS2/16/15280 59 CTD1/10CA 60 108 CRMS2/16/15280 59 CTD1/14CA 60 108 CRMS2/16/65120 59 CTD1/14CA 60 69 CRSM1/18/2560 46 CTD2/08CA 60 69 CRSM1/18/2560 46 CTD2/08CA 60 75 CRSM1/18/2550 46 CTD2/14EPB 60 75 CRSM1/18/6540 46 CTD2/14CPB 60 75 CRSM1/18/0540 47 CTD2/14CPB 60 76 CRSM2/18/105104 7 CTD2/14EPB 60 76 CRSM2/18/161247 CTM1/08/1520 50 74 CTM1/08/3510 50 74 CRSM2/18/05104 CTM1/10/07320 52 138 CTM1/10/07320 52 138 CRSS1/16/052048 CTM1/10/07320 52 138 | | | |
| 30 CRMS2/16/10S10D 59 CTD1/10CA 60 CRMS2/16/16D 59 CTD1/10FPB 60 108 CRMS2/16/16D 59 CTD1/14CA 60 108 CRMS2/16/16D 59 CTD1/14CA 60 69 CRMS2/16/8510D 59 CTD1/14CFPB 60 69 CRSM11/8/256D 46 CTD2/08EPB 60 108 CRSM11/8/256D 46 CTD2/08EPB 60 108 CRSM11/8/255D 46 CTD2/102FPB 60 75 CRSM11/8/255D 46 CTD2/102FPB 60 75 CRSM11/8/255D 46 CTD2/102FPB 60 75 CRSM1/18/255D 47 CTD2/102FB 60 76 CRSM2/18/10510D 47 CTD2/14CA 60 76 CRSM2/18/10510D 47 CTD2/14CS 60 76 CRSM2/18/10510D 47 CTD1/14S31D 50 74 CRSM2/18/0512D 47 CTM1/08/132D 50 74 CRSS1/16/256D 48 CTM1/10/132D 52 74 CRSS1/16/256D 48 CTM1/10/132D 52 74 CRSS1/16/256D 48 CTM1/10/132D 52 75 CRSS1/16/256D 48 CTM1/10/132D 52 76 CRSS1/16/256D 48 | | | |
| 30 CRMS2/16/10S10D 59 CTD1/10CA 60 CRMS2/16/16D 59 CTD1/10FPB 60 108 CRMS2/16/16D 59 CTD1/14CA 60 108 CRMS2/16/16D 59 CTD1/14CA 60 69 CRMS2/16/8510D 59 CTD1/14CFPB 60 69 CRSM11/8/256D 46 CTD2/08EPB 60 108 CRSM11/8/256D 46 CTD2/08EPB 60 108 CRSM11/8/255D 46 CTD2/102FPB 60 75 CRSM11/8/255D 46 CTD2/102FPB 60 75 CRSM11/8/255D 46 CTD2/102FPB 60 75 CRSM1/18/255D 47 CTD2/102FB 60 76 CRSM2/18/10510D 47 CTD2/14CA 60 76 CRSM2/18/10510D 47 CTD2/14CS 60 76 CRSM2/18/10510D 47 CTD1/14S31D 50 74 CRSM2/18/0512D 47 CTM1/08/132D 50 74 CRSS1/16/256D 48 CTM1/10/132D 52 74 CRSS1/16/256D 48 CTM1/10/132D 52 74 CRSS1/16/256D 48 CTM1/10/132D 52 75 CRSS1/16/256D 48 CTM1/10/132D 52 76 CRSS1/16/256D 48 | 294 | CRMS1/16/8D 58 | CTD1/08SS 60 |
| CRMS2/16/12580 CTD1/10EPB C 108 CRMS2/16/145120 59 CTD1/10EPB 60 108 CRMS2/16/65120 59 CTD1/14EPB 60 69 CRSM1/18/2560 46 CTD2/08CA 60 69 CRSM1/18/2560 46 CTD2/08CA 60 75 CRSM1/18/2560 46 CTD2/0EPB 60 75 CRSM1/18/55D 46 CTD2/10EPB 60 75 CRSM1/18/85D 46 CTD2/14CA 60 75 CRSM1/18/85D 47 CTD2/14CPB 60 76 CRSM2/18/16D 47 CTD2/14CPB 60 76 CRSM2/18/16D 47 CTD2/14CPB 60 76 CRSM2/18/16D 47 CTM1/08/152D 50 76 CRSM2/18/16D 47 CTM1/08/152D 50 74 CRSS1/16/85D 48 CTM1/10/132D 52 738 CRSS1/16/85D 48 CTM1/10/132D 52 | | | |
| 108 CRMS2/16/461D 59 CTD1/14CS 60 108 CRMS2/16/6512D 59 CTD1/14CA 60 108 CRMS2/16/6512D 59 CTD1/14CA 60 69 CRSM1/18/256D 46 CTD2/08CP8 60 69 CRSM1/18/256D 46 CTD2/08CP8 60 108 CRSM1/18/255D 46 CTD2/08CP8 60 108 CRSM1/18/255D 46 CTD2/010CA 60 75 CRSM1/18/25D 47 CTD2/14CA 60 75 CRSM1/18/25D 47 CTD2/14CA 60 76 CRSM2/18/105100 47 CTD2/14CA 60 76 CRSM2/18/105100 47 CTD2/14CS 60 76 CRSM2/18/105100 47 CTD2/14SS 60 74 CRSM2/18/105100 47 CTM1/10/83D 50 74 CRSS1/16/65120 47 CTM1/10/83D 52 74 CRSS1/16/6510 48 CTM1/10/83D 52 75 CRSS1/16/65148 CTM1/10/83D 52 74 CRSS1/16/6510 48 CTM1/10/83D 52 75 CRSS1/16/6510 48 CTM1/10/83D 52 76 CRSS1/16/6510 48 CTM1/10/83D 54 77 CRSS2/16/16/510 | 30 | | |
| 108 CRMS2/16/4512D 59 CTD1/14.CR 60 108 CRMS2/16/6512D 59 CTD1/14EPB 60 69 CRSM1/18/25B 46 CTD2/08CPB 60 108 CRSM1/18/35D 46 CTD2/08CPB 60 108 CRSM1/18/35D 46 CTD2/08CPB 60 108 CRSM1/18/35D 46 CTD2/102PB 60 75 CRSM1/18/051D0 46 CTD2/14CA 60 76 CRSM2/18/1051D0 47 CTD2/14CA 60 76 CRSM2/18/1051D0 47 CTD2/14SPB 60 76 CRSM2/18/1651D0 47 CTM1/08/3D 50 74 CRSM2/18/1651D4 47 CTM1/08/3D 50 74 CRSM2/18/1651D48 CTM1/10/13D 52 78 CRSM1/18/35D48 CTM1/10/13D 52 74 CRSS1/16/35D48 CTM1/10/13D 52 75 CRSM2/18/16/31D 9 CTM1/10/10/13D 76 | | | |
| 108 CRMS2/16/6512D 59 CTD1/14EPB 60 69 CRMS1/16/8510D 59 CTD1/14SS 60 60 69 CRSM1/18/3S6D 46 CTD2/08EPS 60 108 CRSM1/18/3S5D 46 CTD2/01CPB 60 75 CRSM1/18/3S5D 46 CTD2/10CPB 60 75 CRSM1/18/3S5D 46 CTD2/10CPB 60 75 CRSM1/18/18/5D 46 CTD2/10CPB 60 76 CRSM2/18/10S10D 47 CTD2/14CPB 60 76 CRSM2/18/12D 47 CTM1/08/13D 50 74 CRSM2/18/13D 47 CTM1/08/13D 50 74 CRSM1/16/3SDD 48 CTM1/10/8/3D 52 74 CRSM1/16/3SDD 48 CTM1/10/8/3D 52 74 CRSM1/16/3SDD 48 CTM1/10/132D 52 74 CRSM1/16/3SDD 48 CTM1/10/132D 52 75 CRSM1/16/3SDD 48 CTM | 108 | CRMS2/16/16D 59 | CTD1/10SS 60 |
| 108 CRMS2/16/6512D 59 CTD1/14EPB 60 69 CRMS1/16/8510D 59 CTD1/14SS 60 60 69 CRSM1/18/3S6D 46 CTD2/08EPS 60 108 CRSM1/18/3S5D 46 CTD2/01CPB 60 75 CRSM1/18/3S5D 46 CTD2/10CPB 60 75 CRSM1/18/3S5D 46 CTD2/10CPB 60 75 CRSM1/18/18/5D 46 CTD2/10CPB 60 76 CRSM2/18/10S10D 47 CTD2/14CPB 60 76 CRSM2/18/12D 47 CTM1/08/13D 50 74 CRSM2/18/13D 47 CTM1/08/13D 50 74 CRSM1/16/3SDD 48 CTM1/10/8/3D 52 74 CRSM1/16/3SDD 48 CTM1/10/8/3D 52 74 CRSM1/16/3SDD 48 CTM1/10/132D 52 74 CRSM1/16/3SDD 48 CTM1/10/132D 52 75 CRSM1/16/3SDD 48 CTM | 108 | CRMS2/16/4S12D 59 | CTD1/14CA 60 |
| 69 CRMS2/16/8510D 59 CTD1/14/S5 60 69 CRSM1/18/2S6D 46 CTD2/08CA 60 69 CRSM1/18/2S6D 46 CTD2/08CPB 60 108 CRSM1/18/SSD 46 CTD2/10CA 60 75 CRSM1/18/SSD 46 CTD2/10CA 60 75 CRSM1/18/SSD 46 CTD2/10SCA 60 76 CRSM1/18/SSD 47 CTD2/10SCA 60 76 CRSM2/18/10S1DD 47 CTD2/14SCA 60 76 CRSM2/18/1651D 47 CTM1/10/8/3D 50 74 CRSM2/18/851D 47 CTM1/10/8/3D 52 74 CRSM1/16/SSD 48 CTM1/10/1/3D 52 74 CRSM1/16/SSD 48 CTM1/10/1/3D 52 75 CRSM1/16/SSD 48 CTM1/10/1/3D 52 74 CRSM2/18/16/SD 48 CTM1/10/1/3D 52 75 CRSM2/16/1/SD 48 CTM1/10/1/1/ | | | |
| 69 CRSM1/18/2S6D 46 CTD2/08CA 60 69 CRSM1/18/2S5D 46 CTD2/08EPB 60 108 CRSM1/18/S5D 46 CTD2/10EPB 80 75 CRSM1/18/S5D 46 CTD2/10EPB 80 75 CRSM1/18/S5D 46 CTD2/14CA 60 75 CRSM1/18/S5D 47 CTD2/14CA 60 76 CRSM2/18/15D 47 CTD2/14CA 80 76 CRSM2/18/15D 47 CTD2/14CA 80 76 CRSM2/18/15D 47 CTM1/08/152D 50 76 CRSM2/18/15D 47 CTM1/08/351D 50 74 CRSS1/16/2S6D 48 CTM1/10/8/351D 52 74 CRSS1/16/2S6D 48 CTM1/10/8/351D 52 78 CRSS1/16/2S6D 48 CTM1/10/8/351D 52 78 CRSS1/16/2S6D 48 CTM1/10/8/351D 52 74 CRSS1/16/2S6D 48 CTM1/10/8/351D 52 74 CRSS1/16/2S6D 48 CTM1/10/4S1D 52 75 CRSS1/16/2S6D 48 CTM1/10/4S1D 52 75 CRSS2/16/152D 49 CTM1/10/4S1D 52 76 CRSS2/16/152D 49 CTM1/10/4S1D 52 76 CRS2/16 | | | |
| 69 CRSM1/18/3SED 46 CTD2/08EP8 60 108 CRSM1/18/3SED 46 CTD2/10CA 60 75 CRSM1/18/SSD 46 CTD2/10CA 60 75 CRSM1/18/SSD 47 CTD2/10CA 60 75 CRSM1/18/SBD 47 CTD2/14EP8 60 76 CRSM2/18/10S10D 47 CTD2/14CA 60 76 CRSM2/18/1128D 47 CTD2/14CA 60 76 CRSM2/18/1128D 47 CTM1/10/8/13D 50 76 CRSM2/18/15D 47 CTM1/10/8/3D 50 74 CRSM2/18/15D 48 CTM1/10/8/3SD 50 74 CRSS1/16/2SD 48 CTM1/10/8/3D 52 738 CRSS1/16/2SD 48 CTM1/10/132D 52 74 CRSS1/16/6SD 48 CTM1/10/4SD 52 75 CRSS1/16/6SD 48 CTM1/10/4SD 52 76 CRSS1/16/6SD 48 CTM1/10/4SD 52 77 CRSS1/16/6SD 48 CTM1/10/4SD 52 78 CRSS1/16/6SD 48 CTM1/10/4SD 52 79 CRSS2/16/10SD 49 CTM1/10/4SD 52 79 CRSS2/16/16/15D 49 CTM1/11/4/2SD 54 79 CRSS2/16/16/16D | 69 | CRMS2/16/8S10D 59 | CTD1/14SS 60 |
| 69 CRSM1/18/4S5D 46 CTD2/08EP8 60 108 CRSM1/18/4S5D 46 CTD2/10CA 60 75 CRSM1/18/SSD 46 CTD2/10CA 60 75 CRSM1/18/SSD 47 CTD2/10CA 60 75 CRSM1/18/SAD 46 CTD2/10CA 60 76 CRSM2/18/10S10D 47 CTD2/14CP8 60 76 CRSM2/18/128D 47 CTD2/14CP8 60 76 CRSM2/18/128D 47 CTM1/10/8/13D 50 76 CRSM2/18/15D 47 CTM1/10/8/3D 50 74 CRSM2/18/15D 48 CTM1/10/8/3D 50 74 CRSS1/16/2SD 48 CTM1/10/8/3D 52 78 CRSS1/16/6SD 48 CTM1/10/8/3D 52 79 CRSS1/16/6SD 48 CTM1/10/8/3D 52 718 CRSS1/16/6SD 48 CTM1/10/8/3D 52 718 CRSS1/16/6SD 48 CTM1/10/8/3D 52 718 CRSS1/16/6SD 48 CTM1/10/8/3D 52 719 CRSS2/16/10SD 49 CTM1/10/8/3D 54 719 CRSS2/16/16/3D 49 CTM1/10/8/3D 54 719 CRSS2/16/16/3D 49 CTM1/11/4/3SD 54 719 CRSS | 69 | CRSM1/18/2S6D 46 | CTD2/08CA 60 |
| 108 CRSM1/18/4S5D 46 CTD2/0SS 60 108 CRSM1/18/5SD 46 CTD2/10EPB 60 75 CRSM1/18/6S4D 46 CTD2/10EPB 60 75 CRSM2/18/15D 47 CTD2/14CA 60 76 CRSM2/18/15D 47 CTD2/14CA 60 76 CRSM2/18/15D 47 CTD2/14CA 60 76 CRSM2/18/12D 47 CTM1/08/152D 50 76 CRSM2/18/812D 47 CTM1/08/3D 50 74 CRSM2/18/85D 48 CTM1/10/33D 52 74 CRSS1/16/256D 48 CTM1/10/33D 52 74 CRSS1/16/55D 48 CTM1/10/45D 52 75 CRSS1/16/55D 48 CTM1/10/45D 52 76 CRSS1/16/55D 48 CTM1/10/45D 52 77 CRSS1/16/55D 48 CTM1/10/45D 52 78 CRSS2/16/18D 49 CTM1/14/35D 54 79 CRSS2/16/18D 49 CTM1/14/35D 54 79 CRSS2/16/18D 49 CTM1/14/35D 54 79 CRSS2/16/18D 49 CTM1/14/45D 54 79 CRSS2/16/18D 49 CTM1/14/45D 54 718 CTM2/08/75D 51 <td< td=""><td></td><td></td><td></td></td<> | | | |
| 108 CRSM1/18/5SD 46 CTD2/10EPB 60 75 CRSM1/18/6S4D 46 CTD2/10EPB 60 75 CRSM218/10S10D 47 CTD2/10EPB 60 76 CRSM218/10S10D 47 CTD2/14CR 60 76 CRSM218/12SD 47 CTD2/14EPB 60 76 CRSM218/12SD 47 CTM1/08/13D 50 76 CRSM218/18/15D 47 CTM1/08/35 50 74 CRSM1/18/6SD 48 CTM1/10/13SD 52 74 CRSS1/16/2SED 48 CTM1/10/3SD 52 74 CRSS1/16/5SD 48 CTM1/10/13SD 52 78 CRSS1/16/5SD 48 CTM1/10/47D 52 79 CRSS1/16/6S4D 48 CTM1/10/47D 52 79 CRSS2/16/15SD 49 CTM1/10/47D 52 79 CRSS2/16/15SD 49 CTM1/10/47D 52 79 CRSS2/16/15D 49 CTM1/14/4S3D 54 79 CRSS2/16/15D 49 CTM1/14/4S3D 54 79 CRSS2/16/61D 49 CTM1/14/4S3D 54 719 CRSS2/16/61D 49 CTM1/14/4S3D 54 719 CRSS2/16/61D 49 CTM1/14/4S3D 54 718 CRS2/16/61D | | | |
| 75 CRSM1/18/6S4D 46 CTD2/10FPB 60 75 CRSM2/18/10510 47 CTD2/10FS 60 76 CRSM2/18/10510 47 CTD2/14CA 60 76 CRSM2/18/105 47 CTD2/14CA 60 76 CRSM2/18/150 47 CTM1/08/152D 50 76 CRSM2/18/150 47 CTM1/08/3D 50 74 CRSM1/18/2S6D 48 CTM1/10/13D 52 138 CRSS1/16/2S6D 48 CTM1/10/13D 52 138 CRSS1/16/4S5D 48 CTM1/10/13D 52 138 CRSS1/16/2S8D 48 CTM1/10/4S1D 52 139 CRSS2/16/12S8D 49 CTM1/14/3S4D 54 139 CRSS2/16/12S8D 49 CTM1/14/3S4D 54 138 CRSS2/16/12S8D 49 CTM1/14/3S4D 54 139 CRSS2/16/12S8D 49 CTM1/14/3S4D 54 138 CS005 163 CTM2 | 108 | | |
| 75 CRSM1/18/8D 46 CTD2/14CA 60 75 CRSM2/18/10S10D 47 CTD2/14CA 60 76 CRSM2/18/12SB 47 CTD2/14CA 60 76 CRSM2/18/12SB 47 CTD2/14CB 60 76 CRSM2/18/4S12D 47 CTM1/08/3D 50 74 CRSM2/18/4S12D 47 CTM1/08/3S1D 50 74 CRSM1/16/SED 48 CTM1/10/3S1D 50 74 CRSS1/16/SED 48 CTM1/10/3S2D 52 75 CRSS1/16/SED 48 CTM1/10/4S1D 52 78 CRSS1/16/SED 48 CTM1/10/4S1D 52 79 CRSS2/16/1S1D 49 CTM1/14/4S1D 54 79 CRSS2/16/1S1D 49 CTM1/14/4S3D 54 79 CRSS2/16/1S1D 49 CTM1/14/4S3D 54 79 CRSS2/16/1S1D 49 CTM1/14/4S3D 54 713 CRSS2/16/1S1D 49 CTM1/14 | 108 | CRSM1/18/5S5D 46 | CTD2/10CA 60 |
| 75 CRSM1/18/8D 46 CTD2/14CA 60 75 CRSM2/18/10S10D 47 CTD2/14CA 60 76 CRSM2/18/12SB 47 CTD2/14CA 60 76 CRSM2/18/12SB 47 CTD2/14CB 60 76 CRSM2/18/4S12D 47 CTM1/08/3D 50 74 CRSM2/18/4S12D 47 CTM1/08/3S1D 50 74 CRSM1/16/SED 48 CTM1/10/3S1D 50 74 CRSS1/16/SED 48 CTM1/10/3S2D 52 75 CRSS1/16/SED 48 CTM1/10/4S1D 52 78 CRSS1/16/SED 48 CTM1/10/4S1D 52 79 CRSS2/16/1S1D 49 CTM1/14/4S1D 54 79 CRSS2/16/1S1D 49 CTM1/14/4S3D 54 79 CRSS2/16/1S1D 49 CTM1/14/4S3D 54 79 CRSS2/16/1S1D 49 CTM1/14/4S3D 54 713 CRSS2/16/1S1D 49 CTM1/14 | 75 | CRSM1/18/6S4D 46 | |
| 75 CRSM2/18/10S10D 47 CTD2/14CA 60 76 CRSM2/18/12S8D 47 CTD2/14EPB 60 76 CRSM2/18/16D 47 CTD2/14SS 60 76 CRSM2/18/612D 47 CTM1/08/1S2D 50 74 CRSM2/18/6S12D 47 CTM1/08/1S2D 50 74 CRSS1/16/3S6D 48 CTM1/10/8/3S1D 50 74 CRSS1/16/3S6D 48 CTM1/10/13S2D 52 138 CRSS1/16/3S6D 48 CTM1/10/13S2D 52 138 CRSS1/16/13S6D 48 CTM1/10/15D 52 138 CRSS1/16/16S4D 48 CTM1/10/4S1D 52 138 CRSS1/16/16S1D 49 CTM1/14/14S5D 54 139 CRSS2/16/12S8D 49 CTM1/14/4S3D 54 139 CRSS2/16/16S12D 49 CTM1/14/4S3D 54 138 CRSS2/16/0812D 49 CTM1/14/4S3D 54 139 CS010 163 CTM2/08/6S2D 51 138 CS28/VBB 162 CTM2/08/085 51 139 CS010 163 | | | |
| 76 CRSM2/18/12S8D 47 CTD2/14EPB 60 76 CRSM2/18/16D 47 CTD2/14EPB 60 76 CRSM2/18/8512D 47 CTM1/08/152D 50 74 CRSM2/18/8512D 47 CTM1/08/152D 50 74 CRSS1/16/256D 48 CTM1/08/35D 52 138 CRSS1/16/256D 48 CTM1/10/32D 52 138 CRSS1/16/256D 48 CTM1/10/32D 52 138 CRSS1/16/55D 48 CTM1/10/451D 52 138 CRSS1/16/654D 48 CTM1/10/451D 54 139 CRSS2/16/10510D 49 CTM1/14/25D 54 139 CRSS2/16/161D 49 CTM1/14/453D 54 138 CRSS2/16/81DD 49 CTM1/14/453D 54 138 CRSS2/16/81DD 49 CTM1/14/453D 54 138 CS/28/VBB 162 CTM2/08/453D 51 138 CS/26/85163 CTM2/08/8515 | | | |
| 76 CRSM2/18/16D 47 CTD2/14SS 60 76 CRSM2/18/4S12D 47 CTM1/08/1S2D 50 74 CRSM2/18/6S12D 47 CTM1/08/3S1D 50 74 CRSS1/16/2S6D 48 CTM1/08/3S1D 50 74 CRSS1/16/3S6D 48 CTM1/10/13S2D 52 138 CRSS1/16/3S6D 48 CTM1/10/13S2D 52 138 CRSS1/16/3S6D 48 CTM1/10/14S1D 52 138 CRSS1/16/6S4D 48 CTM1/10/4S1D 52 139 CRSS2/16/10S1DD 49 CTM1/14/2S4D 54 139 CRSS2/16/8S12D 49 CTM1/14/S3D 54 138 CRSS2/16/8S12D 49 CTM1/14/S5D 54 138 CRSS2/16/8S12D 49 CTM1/14/S5D 54 138 CRSS2/16/8S12D 49 CTM1/14/S5D 54 138 CS/28/MB 162 CTM2/08/S3D 51 139 CS015 163 | 75 | CRSM2/18/10S10D 47 | CTD2/14CA 60 |
| 76 CRSM2/18/16D 47 CTD2/14SS 60 76 CRSM2/18/4S12D 47 CTM1/08/1S2D 50 74 CRSM2/18/6S12D 47 CTM1/08/3S1D 50 74 CRSS1/16/2S6D 48 CTM1/08/3S1D 50 74 CRSS1/16/3S6D 48 CTM1/10/13S2D 52 138 CRSS1/16/3S6D 48 CTM1/10/13S2D 52 138 CRSS1/16/3S6D 48 CTM1/10/13S2D 52 138 CRSS1/16/6S4D 48 CTM1/10/4S1D 52 139 CRSS2/16/10S1DD 49 CTM1/14/2S4D 54 139 CRSS2/16/12D 49 CTM1/14/4S3D 54 138 CRSS2/16/6S12D 49 CTM1/14/4S3D 54 138 CRSS2/16/6S12D 49 CTM1/14/4S3D 54 138 CRSS2/16/6S12D 49 CTM1/14/4S3D 54 138 CS/28/4B 162 CTM2/08/4S2D 51 139 CS015 163 | 76 | CRSM2/18/12S8D 47 | CTD2/14EPB 60 |
| 76 CRSM2/18/4S12D 47 CTM1/08/1S2D 50 74 CRSM2/18/4S12D 47 CTM1/08/3D 50 74 CRSM2/18/8S1DD 47 CTM1/08/3S1D 50 74 CRSS1/16/2S6D 48 CTM1/08/3S1D 50 138 CRSS1/16/2S6D 48 CTM1/10/1S3D 52 138 CRSS1/16/4S5D 48 CTM1/10/1S3D 52 138 CRSS1/16/6S4D 48 CTM1/10/4S5D 54 138 CRSS1/16/6S4D 48 CTM1/10/4S5D 54 139 CRSS2/16/10S10D 49 CTM1/14/2S4D 54 139 CRSS2/16/16D 49 CTM1/14/2S5D 54 139 CRSS2/16/16D 49 CTM1/14/2S5D 54 138 CRSS2/16/16D 49 CTM1/14/2S3D 54 139 CRSS2/16/16D 49 CTM1/14/0S1D 54 138 CRSS2/16/16D 19 CTM1/14/8S1D 54 138 CRSS2/16/0S12D 49 CTM1/14/0S1D 54 138 CRSS2/16/0S12D 49 CTM1/14/0S1D 54 139 CS015 163 CTM2/08/4S2D 51 139 CS015 163 CTM2/08/4S2D 51 139 CTA/AGH 65 CTM2/08/082D 53 144 | | | |
| 76 CRSM2/18/6S12D 47 CTM1/08/3D 50 74 CRSM2/18/8S10D 47 CTM1/08/3D 50 74 CRSS1/16/2S6D 48 CTM1/10/4S 50 138 CRSS1/16/2S6D 48 CTM1/10/4S1D 52 138 CRSS1/16/3S6D 48 CTM1/10/4S1D 52 138 CRSS1/16/5SD 48 CTM1/10/4D 52 138 CRSS1/16/6S4D 48 CTM1/10/4D 52 139 CRSS2/16/10S10D 49 CTM1/14/4S3D 54 139 CRSS2/16/12S8D 49 CTM1/14/4S3D 54 139 CRSS2/16/12S8D 49 CTM1/14/4S3D 54 139 CRSS2/16/6S12D 49 CTM1/14/4S3D 54 139 CRSS2/16/6S12D 49 CTM1/14/5S3D 54 138 CRSS2/16/6S12D 49 CTM1/14/5S3D 54 138 CS005 163 CTM2/08/3S4D 51 138 CS005 163 CTM2/08/6S2D 51 139 CS015 163 CTM2/08/6S2D 51 139 CT/ARH 65 CTM2/08/6S2D 51 139 CT/ARH 65 CTM2/08/6S2D 53 144 CT/CR/26 67 CTM2/08/6S2D 53 144 CT/CR/26 67 <td></td> <td></td> <td></td> | | | |
| 74 CRSM2/18/8510D 47 CTM1/08/351D 50 74 CRSS1/16/2S6D 48 CTM1/10/835D 50 138 CRSS1/16/3S6D 48 CTM1/10/13SD 52 138 CRSS1/16/4S5D 48 CTM1/10/13SD 52 138 CRSS1/16/6S4D 48 CTM1/10/4S1D 52 138 CRSS1/16/6S4D 48 CTM1/10/4S1D 52 139 CRSS2/16/10S1DD 49 CTM1/14/2S5D 54 139 CRSS2/16/16S12D 49 CTM1/14/2S5D 54 139 CRSS2/16/8S10D 49 CTM1/14/2S5D 54 138 CS025 163 CTM2/08/8S2D 51 139 CS015 163 CTM2/08/8S2D 51 139 CT/AGH 65 | 76 | | |
| 74 CRSM2/18/8510D 47 CTM1/08/351D 50 74 CRSS1/16/2S6D 48 CTM1/10/835D 50 138 CRSS1/16/3S6D 48 CTM1/10/13SD 52 138 CRSS1/16/4S5D 48 CTM1/10/13SD 52 138 CRSS1/16/6S4D 48 CTM1/10/4S1D 52 138 CRSS1/16/6S4D 48 CTM1/10/4S1D 52 139 CRSS2/16/10S1DD 49 CTM1/14/2S5D 54 139 CRSS2/16/16S12D 49 CTM1/14/2S5D 54 139 CRSS2/16/8S10D 49 CTM1/14/2S5D 54 138 CS025 163 CTM2/08/8S2D 51 139 CS015 163 CTM2/08/8S2D 51 139 CT/AGH 65 | 76 | CRSM2/18/6S12D 47 | CTM1/08/3D 50 |
| 74 CRSS1/16/2S6D 48 CTM1/08/4S 50 138 CRSS1/16/3S6D 48 CTM1/10/3S2D 52 138 CRSS1/16/3S6D 48 CTM1/10/3S2D 52 138 CRSS1/16/6S5D 48 CTM1/10/3S2D 52 138 CRSS1/16/6S5D 48 CTM1/10/4S1D 52 138 CRSS1/16/16D 48 CTM1/10/4S1D 52 139 CRSS2/16/10S10D 49 CTM1/14/3S4D 54 139 CRSS2/16/16S12D 49 CTM1/14/4S3D 54 138 CRSS2/16/6S12D 49 CTM1/14/4S3D 54 138 CRSS2/16/6S12D 49 CTM1/14/4S3D 54 138 CRSS2/16/6S12D 49 CTM1/14/8S1D 54 138 CRSS2/16/6S12D 49 CTM1/14/8S1D 54 138 CRS2/16/6S12D 49 CTM1/14/8S1D 54 138 CRS2/16/6S12D 49 CTM1/14/8S1D 54 139 CS015 G | | | |
| 138 CRSS1/16/3S6D 48 CTM1/10/IS3D 52 138 CRSS1/16/4S5D 48 CTM1/10/IS3D 52 138 CRSS1/16/5SD 48 CTM1/10/4D 52 138 CRSS1/16/6S4D 48 CTM1/10/4D 52 138 CRSS1/16/6S4D 48 CTM1/10/4S1D 52 139 CRSS2/16/10S10D 49 CTM1/14/2S4D 54 139 CRSS2/16/16D 49 CTM1/14/3S4D 54 139 CRSS2/16/6S12D 49 CTM1/14/4S3D 54 138 CS/28/VBB 162 CTM2/08/3S4D 51 139 CS010 163 CTM2/08/3S51 51 139 CT/AGH 65 CTM2/08/5D 53 144 CT/CPZ 7 CTM2 | | | |
| 138 CRSS1/16/4S5D 48 CTM1/10/4D 52 138 CRSS1/16/5S5D 48 CTM1/10/4D 52 138 CRSS1/16/6S4D 48 CTM1/10/4S1D 52 139 CRSS2/16/10S10D 49 CTM1/14/2S5D 54 139 CRSS2/16/12S8D 49 CTM1/14/2S5D 54 139 CRSS2/16/12D 49 CTM1/14/2S3D 54 139 CRSS2/16/8512D 49 CTM1/14/2S3D 54 138 CRSS2/16/8510D 49 CTM1/14/4S3D 54 138 CRSS2/16/8510D 49 CTM1/14/8S1D 54 138 CRSS2/16/8510D 49 CTM1/14/8S1D 54 138 CRSS2/16/8510D 49 CTM1/14/8S1D 54 138 CS0251663 CTM2/08/854D 51 139 CS015 163 CTM2/08/852D 51 144 CT/CR/2 67 CTM2/10/85D 53 144 CT/CR/2 67 CTM2/10/85D | | | |
| 138 CRSS1/16/6S4D 48 CTM1/10/4D 52 138 CRSS1/16/6S4D 48 CTM1/10/4S1D 52 138 CRSS1/16/8D 48 CTM1/10/4S1D 52 139 CRSS2/16/10S10D 49 CTM1/14/2S5D 54 139 CRSS2/16/12S4D 49 CTM1/14/4S3D 54 139 CRSS2/16/16S12D 49 CTM1/14/4S3D 54 138 CRSS2/16/6S1DD 49 CTM1/14/4S3D 54 138 CRS2/16/6S1DD 49 CTM1/14/6S1D 54 139 CS015 163 CTM2/08/3S4D 51 139 CT/AGH 65 CTM2/08/6S2D 53 144 CT/CR2 67 CTM2/10/08S4D 53 144 CT/CR1/6 CTM2/10/06S4D <td>138</td> <td>CRSS1/16/3S6D 48</td> <td>CTM1/10/1S3D 52</td> | 138 | CRSS1/16/3S6D 48 | CTM1/10/1S3D 52 |
| 138 CRSS1/16/6S4D 48 CTM1/10/4D 52 138 CRSS1/16/6S4D 48 CTM1/10/4S1D 52 138 CRSS1/16/8D 48 CTM1/10/4S1D 52 139 CRSS2/16/10S10D 49 CTM1/14/2S5D 54 139 CRSS2/16/12S4D 49 CTM1/14/4S3D 54 139 CRSS2/16/16S12D 49 CTM1/14/4S3D 54 138 CRSS2/16/6S1DD 49 CTM1/14/4S3D 54 138 CRS2/16/6S1DD 49 CTM1/14/6S1D 54 139 CS015 163 CTM2/08/3S4D 51 139 CT/AGH 65 CTM2/08/6S2D 53 144 CT/CR2 67 CTM2/10/08S4D 53 144 CT/CR1/6 CTM2/10/06S4D <td>138</td> <td>CRSS1/16/4S5D 48</td> <td>CTM1/10/3S2D 52</td> | 138 | CRSS1/16/4S5D 48 | CTM1/10/3S2D 52 |
| 138 CRSS1/16/6S4D 48 CTM1/10/4S1D 52 138 CRSS1/16/6D 48 CTM1/10/4S1D 52 139 CRSS2/16/10S1DD 49 CTM1/14/2S5D 54 139 CRSS2/16/12S8D 49 CTM1/14/2S5D 54 139 CRSS2/16/12S8D 49 CTM1/14/13S4D 54 139 CRSS2/16/6S12D 49 CTM1/14/4S3D 54 138 CRSS2/16/6S12D 49 CTM1/14/4S1D 54 138 CRSS2/16/6S12D 49 CTM1/14/8S1D 54 138 CS/28/VBB 162 CTM2/08/S4D 51 138 CS/28/VBB 162 CTM2/08/S5D 51 139 CS010 163 CTM2/08/S5D 51 139 CT/AGH 65 CTM2/08/S5D 53 144 CT/CR/26 CTM2/10/8SE 53 144 CT/CR/26 CTM2/10/8SE 53 144 CT/CYD 64 CTM2/10/8SE 53 <tr< td=""><td></td><td></td><td></td></tr<> | | | |
| 138 CRSS1/16/8D 48 CTM1/10/6S 52 139 CRSS2/16/10S10D 49 CTM1/14/2SBD 54 139 CRSS2/16/10S10D 49 CTM1/14/2SBD 54 139 CRSS2/16/16D 49 CTM1/14/2SBD 54 139 CRSS2/16/16D 49 CTM1/14/SSD 54 138 CRSS2/16/16S12D 49 CTM1/14/SSD 54 138 CRSS2/16/08S12D 49 CTM1/14/SSD 51 138 CS005 163 CTM2/08/3S4D 51 139 CS015 163 CTM2/08/08S2D 51 145 CT/IARH 65 CTM2/10/08S4D 53 144 CT/CR/2 67 CTM2/10/08S4D 53 144 CT/IDVI/D 69 CTM2/14/08S4D 55 <td></td> <td></td> <td></td> | | | |
| 139 CRSS2/16/10S10D 49 CTM1/14/2S5D 54 139 CRSS2/16/12S8D 49 CTM1/14/3S4D 54 139 CRSS2/16/12S8D 49 CTM1/14/4S3D 54 139 CRSS2/16/6S12D 49 CTM1/14/4S3D 54 138 CS028/16/8S1DD 49 CTM1/14/6S3D 51 139 CS010 163 CTM2/08/3S4D 51 139 CT/AGH 65 CTM2/08/6S2D 51 144 CT/CR/2 67 CTM2/10/8S4D 53 144 CT/CR/2 67 CTM2/10/05S4D 53 144 CT/CR/2 67 CTM2/10/05S4D 53 144 CT/DH/Ad CTM2/10/05S4D 53 144 | 138 | CRSS1/16/6S4D 48 | CTM1/10/4S1D 52 |
| 139 CRSS2/16/10S10D 49 CTM1/14/2S5D 54 139 CRSS2/16/12S8D 49 CTM1/14/3S4D 54 139 CRSS2/16/12S8D 49 CTM1/14/4S3D 54 139 CRSS2/16/6S12D 49 CTM1/14/4S3D 54 138 CS028/16/8S1DD 49 CTM1/14/6S3D 51 139 CS010 163 CTM2/08/3S4D 51 139 CT/AGH 65 CTM2/08/6S2D 51 144 CT/CR/2 67 CTM2/10/8S4D 53 144 CT/CR/2 67 CTM2/10/05S4D 53 144 CT/CR/2 67 CTM2/10/05S4D 53 144 CT/DH/Ad CTM2/10/05S4D 53 144 | 138 | CRSS1/16/8D 48 | CTM1/10/6S 52 |
| 139 CRSS2/16/12S8D 49 CTM1/14/3S4D 54 139 CRSS2/16/16D 49 CTM1/14/3S4D 54 139 CRSS2/16/16D 49 CTM1/14/3S3D 54 138 CRSS2/16/16S12D 49 CTM1/14/3S3D 54 138 CRSS2/16/0512D 49 CTM1/14/053D 54 138 CRSS2/16/0512D 49 CTM1/14/051D 54 138 CRS2/16/0512D 49 CTM1/14/051D 54 138 CRS2/16/051D 51 53 CTM2/08/052D 51 139 CS010 163 CTM2/08/052D 51 139 CT/AGH 65 CTM2/08/052D 51 145 CT/BIN 66 CTM2/10/08/05D 53 144 CT/CR/2 67 CTM2/10/08/05D 53 144 CT/CYD 64 CTM2/10/06/05D 53 144 CT/DI/X6 69 CTM2/10/06/05D 55 246 CT/LAB 69 CTM2/14/10 | | | |
| 139 CRSS2/16/16D 49 CTM1/14/4S3D 54 139 CRSS2/16/4512D 49 CTM1/14/4S3D 54 138 CRSS2/16/6512D 49 CTM1/14/4S3D 54 138 CRSS2/16/6512D 49 CTM1/14/4S5D 54 138 CS005 163 CTM2/08/3S4D 51 139 CS015 163 CTM2/08/4S2D 51 139 CT/AGH 65 CTM2/08/7S1D 51 145 CT/BIN 66 CTM2/10/08/2D 53 144 CT/CR/2 67 CTM2/10/08/2D 53 144 CT/DV/D 69 CTM2/10/08/2D 53 144 CT/DV/D 69 CTM2/10/08/2D 55 246 CT/DV/D 69 CTM2/14/10/08/2D 55 56 | | | |
| 138 CRSS2/16/4S12D 49 CTM1/14/SSD 54 138 CRSS2/16/6S1DD 49 CTM1/14/SSD 54 138 CRSS2/16/6S1DD 49 CTM1/14/SSD 54 138 CRSS2/16/6S1DD 49 CTM1/14/SSD 54 138 CSS2/16/6S1DD 49 CTM1/14/SSD 54 138 CSS2/16/6S1DD 49 CTM1/14/SSD 54 138 CS05 163 CTM2/08/SSD 51 139 CS010 163 CTM2/08/SD 51 139 CT/AGH 65 CTM2/08/SD 53 144 CT/CR12 67 CTM2/10/SSD 53 144 CT/CYD 64 CTM2/10/SSD 53 144 CT/OYD 69 CTM2/10/SSD 53 144 CT/OYD 69 CTM2/10/SSD 53 144 CT/DIV/D 69 CTM2/10/SSD 53 246 CT/DIV/S 69 CTM2/14/SSD 55 | | | |
| 138 CRSS2/16/4S12D 49 CTM1/14/SSD 54 138 CRSS2/16/6S1DD 49 CTM1/14/SSD 54 138 CRSS2/16/6S1DD 49 CTM1/14/SSD 54 138 CRSS2/16/6S1DD 49 CTM1/14/SSD 54 138 CSS2/16/6S1DD 49 CTM1/14/SSD 54 138 CSS2/16/6S1DD 49 CTM1/14/SSD 54 138 CS05 163 CTM2/08/SSD 51 139 CS010 163 CTM2/08/SD 51 139 CT/AGH 65 CTM2/08/SD 53 144 CT/CR12 67 CTM2/10/SSD 53 144 CT/CYD 64 CTM2/10/SSD 53 144 CT/OYD 69 CTM2/10/SSD 53 144 CT/OYD 69 CTM2/10/SSD 53 144 CT/DIV/D 69 CTM2/10/SSD 53 246 CT/DIV/S 69 CTM2/14/SSD 55 | 139 | CRSS2/16/16D 49 | CTM1/14/4S3D 54 |
| 138 CRSS2/16/6S12D 49 CTM1/14/6D 54 138 CRSS2/16/8S10D 49 CTM1/14/6D 54 138 CRSS2/16/8S10D 49 CTM1/14/6D 54 138 CS/28/VBB 162 CTM2/08/3S4D 51 138 CS005 163 CTM2/08/3S4D 51 139 CS015 163 CTM2/08/2S2D 51 139 CT/AGH 65 CTM2/08/7S1D 51 139 CT/AGH 65 CTM2/08/7S1D 51 144 CT/CR2 67 CTM2/10/4S5D 53 144 CT/CYD 64 CTM2/10/3S6D 53 144 CT/DI/XB 69 CTM2/10/8SD 53 144 CT/DI/XB 69 CTM2/10/8SD 53 144 CT/DI/XB 69 CTM2/14/108D 53 246 CT/ID/XB 69 CTM2/14/108D 55 56 CT/LHR/S 62 CTM2/14/108D 55 < | 139 | CRSS2/16/4S12D 49 | CTM1/14/5S3D 54 |
| 138 CRSS2/16/8510D 49 CTM1/14/851D 54 138 CS/28/VBB 162 CTM2/08/3S4D 51 138 CS005 163 CTM2/08/3S4D 51 139 CS015 163 CTM2/08/052D 51 139 CS015 163 CTM2/08/052D 51 139 CTACH 65 CTM2/08/751D 51 139 CT/ACH 65 CTM2/10/386D 53 144 CT/CR/2 67 CTM2/10/386D 53 144 CT/CR/5 67 CTM2/10/0584D 53 144 CT/DH/A4 67 CTM2/10/082D 53 144 CT/DH/A4 67 CTM2/10/082D 53 144 CT/DH/A 67 CTM2/10/082D 53 144 CT/DH/A 67 CTM2/10/082D 53 246 CT/DH/A 67 CTM2/14/082D 55 246 CT/LH/C5 63 CTM2/14/082D 55 56 CT/LH/C5 CTM2/14/082D 55 56 < | | | |
| 138 CS/28/VBB 162 CTM2/08/3S4D 51 138 CS005 163 CTM2/08/3S4D 51 139 CS010 163 CTM2/08/05 51 139 CS015 163 CTM2/08/05 51 139 CS015 163 CTM2/08/05 51 139 CT/AGH 65 CTM2/08/05 51 149 CT/AIR 65 CTM2/08/05 53 144 CT/CR/2 67 CTM2/10/03SeD 53 144 CT/CR/2 67 CTM2/10/05S4D 53 144 CT/CYD 64 CTM2/10/06S4D 53 144 CT/DI/VD 69 CTM2/10/06S2D 53 246 CT/DI/VD 69 CTM2/14/05S2D 55 246 CT/LH/25 63 CTM2/14/05S2D 55 78 CT/LAB 69 CTM2/14/05S2D 55 56 CT/LH/S 62 CTM2/14/05S2D 55 56 C | | | |
| 138 CS005 163 CTM2/08/4S3D 51 139 CS010 163 CTM2/08/4S2D 51 139 CS015 163 CTM2/08/6S2D 51 139 CT/AGH 65 CTM2/08/6S2D 51 139 CT/AGH 65 CTM2/08/6S2D 51 149 CT/AGH 65 CTM2/08/8S2D 51 144 CT/CR/2 67 CTM2/10/3S4D 53 144 CT/CYD 64 CTM2/10/4S4D 53 144 CT/CYD 64 CTM2/10/8S2D 53 144 CT/DI/VD 69 CTM2/10/8S2D 53 246 CT/IDV/D 69 CTM2/14/12D 55 246 CT/LAB 69 CTM2/14/12D 55 56 CT/LHR/S 62 CTM2/14/12D 55 56 CT/DI/MS 68 CTPHLEM1/10/032D 254 56 CT/DI/MS 68 CTPHLEM1/10/032D 255 <t< td=""><td>138</td><td>CRSS2/16/8S10D 49</td><td></td></t<> | 138 | CRSS2/16/8S10D 49 | |
| 138 CS005 163 CTM2/08/4S3D 51 139 CS010 163 CTM2/08/4S2D 51 139 CS015 163 CTM2/08/6S2D 51 139 CT/AGH 65 CTM2/08/6S2D 51 139 CT/AGH 65 CTM2/08/6S2D 51 149 CT/AGH 65 CTM2/08/8S2D 51 144 CT/CR/2 67 CTM2/10/3S4D 53 144 CT/CYD 64 CTM2/10/4S4D 53 144 CT/CYD 64 CTM2/10/8S2D 53 144 CT/DI/VD 69 CTM2/10/8S2D 53 246 CT/IDV/D 69 CTM2/14/12D 55 246 CT/LAB 69 CTM2/14/12D 55 56 CT/LHR/S 62 CTM2/14/12D 55 56 CT/DI/MS 68 CTPHLEM1/10/032D 254 56 CT/DI/MS 68 CTPHLEM1/10/032D 255 <t< td=""><td>138</td><td>CS/28/VBB 162</td><td>CTM2/08/3S4D 51</td></t<> | 138 | CS/28/VBB 162 | CTM2/08/3S4D 51 |
| 139 CS010 163 CTM2/08/6D 51 139 CS015 163 CTM2/08/6S2D 51 139 CT/AGH 65 CTM2/08/6S2D 51 139 CT/AGH 65 CTM2/08/7S1D 51 145 CT/BIN 66 CTM2/10/3S4D 53 144 CT/CR/2 67 CTM2/10/3S4D 53 144 CT/CR/5 67 CTM2/10/05S4D 53 144 CT/CR/5 67 CTM2/10/05S4D 53 144 CT/DH/A 67 CTM2/10/05S4D 53 144 CT/DH/DH/A4 67 CTM2/10/05S4D 53 246 CT/DH/M 69 CTM2/14/10SED 55 246 CT/LH 69 CTM2/14/10SED 55 56 CT/LH/C5 63 CTM2/14/4S8D 55 56 CT/LH/C5 63 CTM2/14/4S8D 55 56 CT/LH/S 62 CTPHLEM1/10/03S2D 254 | 138 | CS005 163 | |
| 139 CS015 163 CTM2/08/6S2D 51 139 CT/AGH 65 CTM2/08/6S2D 51 139 CT/AGH 65 CTM2/08/6S2D 51 145 CT/AIR 65 CTM2/08/6S2D 53 144 CT/CR/2 67 CTM2/10/3S6D 53 144 CT/CR/2 67 CTM2/10/4S5D 53 144 CT/CR/2 67 CTM2/10/4S5D 53 144 CT/CR/2 67 CTM2/10/4S5D 53 144 CT/DH/A 67 CTM2/10/8D 53 246 CT/DI/VD 69 CTM2/10/8D 55 246 CT/IVA 69 CTM2/14/108D 55 56 CT/LH/CS 63 CTM2/14/108D 55 56 CT/LH/RS 62 CTM2/14/148D 55 56 CT/LH/RS 62 CTM2/14/148D 55 56 CT/LH/RS 62 CTM2/14/14/1082D 55 56 | | | |
| 139 CT/AGH 65 CTM2/08/7S1D 51 139 CT/AR 65 CTM2/08/7S1D 51 145 CT/AIR 65 CTM2/08/7S1D 51 144 CT/CR/2 67 CTM2/10/3S0D 53 144 CT/CR/2 67 CTM2/10/4S5D 53 144 CT/CR/2 67 CTM2/10/5S4D 53 144 CT/CYD 64 CTM2/10/5S4D 53 144 CT/DI/D69 CTM2/10/8S2D 53 246 CT/DI/D69 CTM2/10/8S2D 53 246 CT/I/V 69 CTM2/14/10/5S2D 55 246 CT/LAB 69 CTM2/14/12D 55 56 CT/LHR/S 62 CTM2/14/8S4D 55 56 CT/DB/MS 68 CTPHLEM1/10/3S2D 254 56 CT/OB/MS 68 CTPHLEM1/10/4S1D 254 56 CT/OB/MS 88 CTPHLEM1/10/4S1D 254 57 CT/O | | | |
| 138 CT/AIR 65 CTM2/10/865 51 145 CT/BIN 66 CTM2/10/35D 53 144 CT/CR/2 67 CTM2/10/45D 53 144 CT/CR/5 67 CTM2/10/35D 53 144 CT/CR/5 67 CTM2/10/45D 53 144 CT/CR/5 67 CTM2/10/654D 53 144 CT/DH/A 67 CTM2/10/82D 53 246 CT/DIV/D 69 CTM2/10/82D 55 246 CT/LH/CS 67 CTM2/14/108D 55 246 CT/LH/CS 63 CTM2/14/108D 55 56 CT/LH/CS 63 CTM2/14/46SDD 55 56 CT/LH/CS 63 CTM2/14/46SDD 55 56 CT/LH/CS 63 CTM2/14/48SDD 55 56 CT/LH/CS 63 CTM2/14/48SDD 55 56 CT/LH/CS 63 CTPHLEM1/10/03S2D 254 | 139 | CS015 163 | CTM2/08/6S2D 51 |
| 138 CT/AIR 65 CTM2/10/865 51 145 CT/BIN 66 CTM2/10/35D 53 144 CT/CR/2 67 CTM2/10/45D 53 144 CT/CR/5 67 CTM2/10/35D 53 144 CT/CR/5 67 CTM2/10/45D 53 144 CT/CR/5 67 CTM2/10/654D 53 144 CT/DH/A 67 CTM2/10/82D 53 246 CT/DIV/D 69 CTM2/10/82D 55 246 CT/LH/CS 67 CTM2/14/108D 55 246 CT/LH/CS 63 CTM2/14/108D 55 56 CT/LH/CS 63 CTM2/14/46SDD 55 56 CT/LH/CS 63 CTM2/14/46SDD 55 56 CT/LH/CS 63 CTM2/14/48SDD 55 56 CT/LH/CS 63 CTM2/14/48SDD 55 56 CT/LH/CS 63 CTPHLEM1/10/03S2D 254 | 139 | CT/AGH 65 | CTM2/08/7S1D 51 |
| 145 CT/BIN 66 CTM2/10/386D 53 144 CT/CR/2 67 CTM2/10/385D 53 144 CT/CR/2 67 CTM2/10/485D 53 144 CT/CR/5 67 CTM2/10/584D 53 144 CT/CYD 64 CTM2/10/85D 53 144 CT/CYD 64 CTM2/10/85D 53 144 CT/DI/A69 CTM2/10/85D 53 246 CT/DI/VD 69 CTM2/10/85D 55 246 CT/ICN/D 69 CTM2/14/1085D 55 246 CT/IX 64 CTM2/14/1085D 55 246 CT/LAB 69 CTM2/14/1085D 55 78 CT/LAB 69 CTM2/14/1682D 55 56 CT/LH/CS 63 CTM2/14/1682D 55 56 CT/LH/RS 62 CTM2/14/1085D 55 56 CT/LH/RS 62 CTM2/14/1085D 55 56 CT/OB/MS 283 CTPHLEM1/10/352D 254 56 CT/OB/MS 283 CTPHLEM1/10/352D 255 57 CT/JPARM 255 CTPHLES1/10/085 255 57 CT/TRAY/D 69 CTRESUSM1/10/352D/SS 282 57 CT/TRAY/D 69 CTRESUSM1/10/352D/SS 282 | | | |
| 144 CT/CR/2 67 CTM2/10/4S5D 53 144 CT/CR/5 67 CTM2/10/4S5D 53 144 CT/CR/5 67 CTM2/10/4S5D 53 144 CT/CYD 64 CTM2/10/6S4D 53 144 CT/DYD 69 CTM2/10/8S2D 53 246 CT/DIV/D 69 CTM2/10/8S2D 55 246 CT/IV 69 CTM2/14/12D 55 246 CT/IV 69 CTM2/14/12D 55 246 CT/LAB 69 CTM2/14/16S2D 55 56 CT/LHR/S 62 CTM2/14/4S8D 55 56 CT/LHR/S 62 CTM2/14/8S6D 55 56 CT/OB/MS 283 CTPHLEM1/10/3S2D 254 56 CT/OB/MS 283 CTPHLEM1/10/4S1D 254 56 CT/OB/MS 283 CTPHLEM1/10/4S1D 255 57 CT/OB/MS 68 CTPHLES1/10/0S2D 255 | | | |
| 144 CT/CR/5 67 CTM2/10/5S4D 53 144 CT/CYD 64 CTM2/10/5S4D 53 144 CT/DYD 64 CTM2/10/5S4D 53 144 CT/DYD 64 CTM2/10/5S4D 53 246 CT/DIV/D 69 CTM2/10/8S2D 53 246 CT/DIV/D 69 CTM2/10/8S2D 55 246 CT/DIV/D 69 CTM2/14/10S2D 55 246 CT/LH/CS 63 CTM2/14/16S2D 55 78 CTLAB 69 CTM2/14/4S8D 55 56 CT/LH/CS 63 CTM2/14/4S8D 55 56 CT/DB/MS 68 CTPHLEM1/10/4S1D 254 56 CT/OB/MS 88 CTPHLEM1/10/0S2D 255 | | | |
| 144 CT/CYD 64 CTM2/10/684D 53 144 CT/DH/A4 67 CTM2/10/6D 53 246 CT/DIV/D 69 CTM2/10/8D 53 246 CT/DIV/D 69 CTM2/10/8D 55 246 CT/DIV/D 69 CTM2/14/10/8D 55 246 CT/IX 64 CTM2/14/10/8D 55 246 CT/LAB 69 CTM2/14/10/8D 55 78 CT/LAB 69 CTM2/14/168D 55 56 CT/LH/CS 63 CTM2/14/168D 55 56 CT/LH/RS 62 CTM2/14/10/8SD 25 56 CT/LH/RS 62 CTM2/14/10/8SD 25 56 CT/DB/MS 68 CTPHLEM1/10/8SD 25 56 CT/OB/MS 88 CTPHLEM1/10/8SD 25 57 CT/DARM 255 CTPHLES1/10/0S1D 255 57 CT/DARM 255 CTPHLES1/10/0S1D 255 | 144 | CT/CR/2 67 | CTM2/10/4S5D 53 |
| 144 CT/CYD 64 CTM2/10/684D 53 144 CT/DH/A4 67 CTM2/10/6D 53 246 CT/DIV/D 69 CTM2/10/8D 53 246 CT/DIV/D 69 CTM2/10/8D 55 246 CT/DIV/D 69 CTM2/14/10/8D 55 246 CT/IX 64 CTM2/14/10/8D 55 246 CT/LAB 69 CTM2/14/10/8D 55 78 CT/LAB 69 CTM2/14/168D 55 56 CT/LH/CS 63 CTM2/14/168D 55 56 CT/LH/RS 62 CTM2/14/10/8SD 25 56 CT/LH/RS 62 CTM2/14/10/8SD 25 56 CT/DB/MS 68 CTPHLEM1/10/8SD 25 56 CT/OB/MS 88 CTPHLEM1/10/8SD 25 57 CT/DARM 255 CTPHLES1/10/0S1D 255 57 CT/DARM 255 CTPHLES1/10/0S1D 255 | 144 | CT/CR/5 67 | CTM2/10/5S4D 53 |
| 144 CT/DH/A4 67 CTM2/10/8D 53 246 CT/DIV/D 69 CTM2/10/8D 55 246 CT/DIV/D 69 CTM2/14/10/8D 55 246 CT/I/V 69 CTM2/14/12D 55 246 CT/I/V 64 CTM2/14/12D 55 246 CT/I/V 64 CTM2/14/12D 55 246 CT/LAB 69 CTM2/14/14/82D 55 56 CT/LH/CS 63 CTM2/14/486D 55 56 CT/LH/R/S 62 CTM2/14/886D 55 56 CT/OB/MS 68 CTPHLEM1/10/3S2D 254 56 CT/OB/MS 88 CTPHLEM1/10/4S1D 254 56 CT/OB/MS 88 CTPHLEM1/10/4S1D 255 57 CT/DARM 255 CTPHLES1/10/052D 255 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/35 282 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/35 28 | | | |
| 246 CT/DIV/D 69 CTM2/14/1082D 53 246 CT/DIV/D 69 CTM2/14/1082D 55 246 CT/DIV/D 69 CTM2/14/102D 55 246 CT/I/E 67 CTM2/14/162D 55 246 CT/I/E 67 CTM2/14/162D 55 78 CT/LAB 69 CTM2/14/162D 55 56 CT/LH/C5 63 CTM2/14/48SDD 55 56 CT/LH/C5 62 CTM2/14/16S2D 254 56 CT/LH/C5 63 CTM2/14/16S2D 254 56 CT/LH/C5 62 CTM2/14/16S2D 254 56 CT/LH/C5 63 CTM2/14/16S2D 254 56 CT/DB/MS 68 CTPHLEM1/10/3S2D 254 56 CT/OB/MS 283 CTPHLEM1/10/3S2D 255 57 CT/SB/I5 68 CTPHLES1/10/04S/1D 255 57 CT/SB/I5 66 CTPHLES1/10/05S 282 <td></td> <td></td> <td></td> | | | |
| 246 CT/DIV/S 69 CTM2/14/10S6D 55 246 CT/FS 67 CTM2/14/10S6D 55 246 CT/FS 67 CTM2/14/10S6D 55 78 CT/LAB 69 CTM2/14/10S6D 55 56 CT/LH/C5 63 CTM2/14/10S2D 55 56 CT/LH/C5 62 CTM2/14/16S2D 254 56 CT/LH/R/S 62 CTM2/14/10/S2D 254 56 CT/LH/R/S 62 CTPHLEM1/10/S1D 254 56 CT/OB/MS 68 CTPHLEM1/10/S1D 254 56 CT/OB/MS 88 CTPHLEM1/10/S1D 255 57 CT/DB/MS 68 CTPHLES1/10/S2D 255 57 CT/JPARM 255 CTPHLES1/10/4S1D 255 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/S2 282 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/S2 282 57 CT/THARY/D CTRESUSM1/10/16S/S2 <td< td=""><td>144</td><td></td><td></td></td<> | 144 | | |
| 246 CT/DIV/S 69 CTM2/14/10S6D 55 246 CT/FS 67 CTM2/14/10S6D 55 246 CT/IV 64 CTM2/14/10S6D 55 78 CT/LAB 69 CTM2/14/16S2D 55 56 CT/LH/C5 63 CTM2/14/16S6D 55 56 CT/LH/NR/S 62 CTM2/14/16S6D 55 56 CT/LH/NR/S 62 CTM2/14/16S6D 55 56 CT/LH/NR/S 62 CTM2/14/10S1D 254 56 CT/LOB/MS 68 CTPHLEM1/10/3S2D 254 56 CT/OB/MS 283 CTPHLEM1/10/S1D 255 57 CT/DB/NS 68 CTPHLES1/10/052D 255 57 CT/DB/NS 66 CTPHLES1/10/052D 255 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/S2 282 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/S2 282 57 CT/THRAY/D 69 CTRESUSM1/10/3S2D/S2 282 57 CT/THRAY/D 69 CTRESUSM1/10/3S2D/S2 282 57 CT/THRAY/D 69 CTRESUSM1/10/3S2D/S2 282 57 CTCHARM1/10/3S2D 72 CTS1/08/1S2D 40 | 246 | CT/DIV/D 69 | CTM2/10/8S2D 53 |
| 246 CT/FS 67 CTM2/14/12D 55 246 CT/IV 64 CTM2/14/162D 55 246 CT/IV 64 CTM2/14/162D 55 78 CTLAB 69 CTM2/14/4S1D 55 56 CT/LHRS 63 CTM2/14/4S1D 55 56 CT/LHRS 62 CTM2/14/8S6D 55 56 CT/LHRS 62 CTPHLEM1/10/3S2D 254 56 CT/OB/MS 68 CTPHLEM1/10/4S1D 254 56 CT/OB/SS 68 CTPHLEM1/10/3S2D 255 57 CT/BARM 255 CTPHLEM1/10/3S2D 255 57 CT/SB/15 66 CTPHLES1/10/4S1D 255 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/58 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/58 57 CT/TRAY/D 69 CTRESUSM1/10/68/S 282 57 CTCHARM1/10/3S2D 72 CTS1/08/152D 40 | | CT/DIV/S 69 | CTM2/14/10S6D 55 |
| 246 CT/IV 64 CTM2/14/16S2D 55 78 CT/LAB 69 CTM2/14/16S2D 55 56 CT/LHR/G 69 CTM2/14/4S0D 55 56 CT/LHR/S 62 CTM2/14/4S6D 55 56 CT/LHR/S 62 CTM2/14/10/3S2D 254 56 CT/LHR/S 62 CTM2/14/10/3S2D 254 56 CT/OB/MS 68 CTPHLEM1/10/3S2D 254 56 CT/OB/MS 68 CTPHLEM1/10/S2D 255 57 CT/SBB/15 66 CTPHLES1/10/S2 255 57 CT/TRAY/D 69 CTRESUSM1/10/32D/SS 282 57 CT/THRAY/D 69 CTRESUSM1/10/SS 282 57 CTOHARM1/10/3S2D 72 CTS1/08/1S2 40 | | | |
| 78 CT/LAB 69 CTM2/14/4S10D 55 56 CT/LH/C5 63 CTM2/14/4S8D 55 56 CT/LH/C5 62 CTM2/14/4S8D 55 56 CT/LH/R/S 62 CTM2/14/4S8D 55 56 CT/LH/R/S 62 CTPHLEM1/10/3S2D 254 56 CT/OB/MS 68 CTPHLEM1/10/4S1D 254 56 CT/OB/MS 83 CTPHLEM1/10/4S1D 254 56 CT/OB/MS 88 CTPHLEM1/10/4S1D 255 57 CT/PARM 255 CTPHLES1/10/4S1D 255 57 CT/TRAY/D 69 CTRESUSM1/10/4S1D/2S2 282 57 CT/TRAY/D 69 CTRESUSM1/10/4S1D/SS 282 57 CT/TRAY/D 69 CTRESUSM1/10/4S1D/SS 282 57 CTCHARM1/10/3S2D 72 CTS1/08/1S2D 40 | | | |
| 56 CT/LH/C5 63 CTM2/14/858D 55 56 CT/LH/C5 63 CTM2/14/858D 55 56 CT/LH/C5 62 CTM2/14/850D 55 56 CT/LH/C5 62 CTPHLEM1/10/352D 254 56 CT/OB/MS 68 CTPHLEM1/10/4S1D 254 56 CT/OB/S5 68 CTPHLEM1/10/352D 255 57 CT/DB/S5 66 CTPHLES1/10/352D 255 57 CT/SB/15 66 CTPHLES1/10/4S1D 255 57 CT/TRAY/D 69 CTRESUSM1/10/32D/55 282 57 CT/TRAY/D 69 CTRESUSM1/10/4S1D/5S 282 57 CT/TRAY/D 69 CTRESUSM1/10/4S1D/5S 282 57 CT/TRAY/D 69 CTRESUSM1/10/4S1D/5S 282 57 CTCHARM.1/0/352D 72 CTS1/08/1S2D 40 | 246 | | |
| 56 CT/LH/C5 63 CTM2/14/858D 55 56 CT/LH/C5 63 CTM2/14/858D 55 56 CT/LH/C5 62 CTM2/14/850D 55 56 CT/LH/C5 62 CTPHLEM1/10/352D 254 56 CT/OB/MS 68 CTPHLEM1/10/4S1D 254 56 CT/OB/S5 68 CTPHLEM1/10/352D 255 57 CT/DB/S5 66 CTPHLES1/10/352D 255 57 CT/SB/15 66 CTPHLES1/10/4S1D 255 57 CT/TRAY/D 69 CTRESUSM1/10/32D/55 282 57 CT/TRAY/D 69 CTRESUSM1/10/4S1D/5S 282 57 CT/TRAY/D 69 CTRESUSM1/10/4S1D/5S 282 57 CT/TRAY/D 69 CTRESUSM1/10/4S1D/5S 282 57 CTCHARM.1/0/352D 72 CTS1/08/1S2D 40 | 78 | CT/LAB 69 | CTM2/14/4S10D 55 |
| 56 CT/LHNR/S 62 CTM2/14/8S6D 55 56 CT/LHR/S 62 CTPHLEM1/10/3S2D 254 56 CT/OB/MS 68 CTPHLEM1/10/4S1D 254 56 CT/OB/MS 88 CTPHLEM1/10/4S1D 254 56 CT/OB/MS 88 CTPHLEM1/10/4S1D 254 57 CT/SBB/15 66 CTPHLES1/10/4S1D 255 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/SS 282 57 CT/TRAY/S 69 CTRESUSM1/10/SS 20/SS 282 57 CT/TRAY/S 69 CTRESUSM1/10/SS 20/SS 282 57 CTCHARM1/10/3S2D 72 CTS1/08/1S2D 40 | | | |
| 56 CT/LHR/S 62 CTPHLEM1/10/3S2D 254 56 CT/OB/MS 68 CTPHLEM1/10/4S1D 254 56 CT/OB/MS 283 CTPHLEM1/10/4S1D 254 56 CT/OB/MS 283 CTPHLEM1/10/4S1D 255 57 CT/DB/SS 68 CTPHLES1/10/4S1D 255 57 CT/PARM 255 CTPHLES1/10/4S1D 255 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/SS 282 57 CT/TRAY/S 69 CTRESUSM1/10/4S1D/SS 282 57 CT/TRAY/S 69 CTRESUSM1/10/4S1D/SS 282 57 CTCHARM1/10/3S2D 72 CTS1/08/1S2D 40 | | | |
| 56 CT/OB/MS 68 CTPHLEM1/10/4S1D 254 56 CT/OB/MS 283 CTPHLEM1/10/6S 254 56 CT/OB/SS 68 CTPHLEM1/10/3S2D 255 57 CT/PARM 255 CTPHLES1/10/3S2D 255 57 CT/SB/JS 66 CTPHLES1/10/4S1D 255 57 CT/TRAY/D 69 CTRESUSM1/10/4S1D/SS 282 57 CTCHARM1/10/3S2D 72 CTS1/08/1S2D 40 | 56 | | |
| 56 CT/OB/MS 283 CTPHLEM1/10/6S 254 56 CT/OB/SS 68 CTPHLES1/10/032D 255 57 CT/PARM 255 CTPHLES1/10/4S1D 255 57 CT/SBB/15 66 CTPHLES1/10/4S1D 255 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/SS 282 57 CT/TRAY/S 69 CTRESUSM1/10/4S1D/SS 282 57 CT/TRAY/S 69 CTRESUSM1/10/4S1D/SS 282 57 CTCHARM1/10/3S2D 72 CTS1/08/1S2D 40 | 56 | | |
| 56 CT/OB/MS 283 CTPHLEM1/10/6S 254 56 CT/OB/SS 68 CTPHLES1/10/032D 255 57 CT/PARM 255 CTPHLES1/10/4S1D 255 57 CT/SBB/15 66 CTPHLES1/10/4S1D 255 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/SS 282 57 CT/TRAY/S 69 CTRESUSM1/10/4S1D/SS 282 57 CT/TRAY/S 69 CTRESUSM1/10/4S1D/SS 282 57 CTCHARM1/10/3S2D 72 CTS1/08/1S2D 40 | 56 | CT/OB/MS 68 | CTPHLEM1/10/4S1D 254 |
| 56 CT/OB/SS 68 CTPHLES1/10/3S2D 255 57 CT/PARM 255 CTPHLES1/10/4S1D 255 57 CT/SB/J5 CTPHLES1/10/4S1D 255 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/SS 282 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/SS 282 57 CT/WF 64 CTRESUSM1/10/4S1D/SS 282 57 CT/WF 64 CTRESUSM1/10/4S1D/SS 282 57 CTCHARM1/10/3S2D 72 CTS1/08/1S2D 40 | | | |
| 57 CT/PARM 255 CTPHLES1/10/4S1D 255 57 CT/SBB/15 66 CTPHLES1/10/6S 255 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/SS 282 57 CT/TRAY/S 69 CTRESUSM1/10/4S1D/SS 282 57 CT/TRAY/S 69 CTRESUSM1/10/4S1D/SS 282 57 CTWF 64 CTRESUSM1/10/6S/SS 282 57 CTCHARM1/10/3S2D 72 CTS1/08/1S2D 40 | | | |
| 57 CT/SBB/15 66 CTPHLES1/10/6S 255 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/SS 282 57 CT/TRAY/D 69 CTRESUSM1/10/4S1D/SS 282 57 CT/TRAY/D 69 CTRESUSM1/10/4S1D/SS 282 57 CTCHARM1/10/3S2D 72 CTS1/08/1S2D 40 | 56 | | |
| 57 CT/SBB/15 66 CTPHLES1/10/6S 255 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/SS 282 57 CT/TRAY/D 69 CTRESUSM1/10/4S1D/SS 282 57 CT/WF 64 CTRESUSM1/10/6S/SS 282 57 CTCHARM1/10/3S2D 72 CTS1/08/1S2D 40 | 57 | CT/PARM 255 | CTPHLES1/10/4S1D 255 |
| 57 CT/TRAY/D 69 CTRESUSM1/10/3S2D/SS 282 57 CT/TRAY/S 69 CTRESUSM1/10/3S2D/SS 282 57 CT/TRAY/S 69 CTRESUSM1/10/4S1D/SS 282 57 CT/WF 64 CTRESUSM1/10/4S2S 282 57 CTCHARM1/10/3S2D 72 CTS1/08/1S2D 40 | | | |
| 57 CT/TRAY/S 69 CTRESUSM1/10/4S1D/SS 282 57 CT/WF 64 CTRESUSM1/10/6S/SS 282 57 CTCHARM1/10/3S2D 72 CTS1/08/1S2D 40 | | | |
| 57 CT/WF 64 CTRESUSM1/10/6S/SS 282 57 CTCHARM1/10/3S2D 72 CTS1/08/1S2D 40 | 57 | | |
| 57 CT/WF 64 CTRESUSM1/10/6S/SS 282 57 CTCHARM1/10/3S2D 72 CTS1/08/1S2D 40 | 57 | CT/TRAY/S 69 | CTRESUSM1/10/4S1D/SS 282 |
| 57 CTCHARM1/10/3S2D 72 CTS1/08/1S2D 40 | | | CTRESUSM1/10/6S/SS 282 |
| | | | |
| 58 CTCHARM1/10/4S1D 72 CTS1/08/3D 40 | | | |
| | 58 | CTCHARM17/10/4S1D /2 | CTS1/08/3D 40 |

| BC1DWR/GW 27 | BR120/WHL 293 |
|--|---------------------------------------|
| BC2C/BE 23 | BR135/WHL 294 |
| BC2C/GW 23 | BR140/WHL 294 |
| BC2CWL/BE 28 | BR145/WHL 294 |
| BC2CWL/GW 28 | BR150/WHL 294 |
| BC2CWR/BE 28 | BR150/WHL 294 |
| BC2CWR/BC 28 | BT/12/A 30 |
| BC2CWR/GW 28 BC2D/GW 23 BC2D/GW 23 BC2DWL/GW 28 BC2DWL/GW 28 BC2DWR/BE 28 BC2DWR/BE 28 BC3C/BE 24 BC3C/GW 24 BC3CPD/GW 24 BC3CPD/GW 24 BC3CPDWL/BE 29 BC3CPDWL/BE 29 BC3CPDWR/GW 29 BC3CPDWR/GW 29 BC3CWR/GW 29 BC3CWR/GW 29 BC3CWR/GW 29 BC3CWR/GW 29 BC3CWR/GW 29 BC3DPD/GW 24 BC3DPD/GW 24 BC3DPD/GW 24 BC3DPD/GW 24 BC3DPD/GW 24 BC3DPD/GW 24 BC3DPD/GW 24 BC3DPD/GW 24 BC3DPD/GW 29 BC3DPD/GW 29 BC3DPWL/GW 29 BC3DPD/WR/GW 29 BC3DPDWL/GW 29 BC3DPDWR/GW 29 BC3DPWR/GW 29 BC3DPWR/GW 29 BC3DWR/GW 20 BC3DWR/GW 20 BC3DWR/G | C C C C C C C C C C C C C C C C C C C |

| CTS1/08/3S1D 40 | DH/A4/WF 77 | EBC900 21 | EC2430/VBB 169 | EPS015 208 | MR200 236 | OBL/GW1/SP/S 32 | PC539/FS 134 | PL050/BD 384 | Р |
|--------------------------------------|-------------------|----------------------------------|----------------------------------|------------------------|--------------------------|-----------------------|------------------------------|--------------------|------|
| CTS1/08/4S 40 | DP/2H 214 | EBC901 21 | EC2435/VBB 169 | G | MR210 236 | OBL/GW3/SP/S 32 | PC539/FS/LH 134 | PL055/BD 384 | P |
| CTS1/10/1S3D 42 | DP/4H 214 | EBC902 21 | EC2440/VBB 169 | G/16/TILT/VBB 162 | MR210CL 236 | OBP/SP/S 31 | PC539/FS/LH/LX 134 | PL060/BD 384 | F |
| CTS1/10/3S2D 42 | DP/C 109 | EBC910 21 | EC3105/VBB 172 | G/200/RS/VBB 262 | MR300 238 | OCL/BE0/GW/S 33 | PC539/FS/LX 134 | PL065/BD 384 | P |
| CTS1/10/4D 42 | DP/C 215 | EC005/VBB 161 | EC3110/VBB 172 | G/205/FS/VBB 264 | MR310 238 | OCL/BE1/GW/S 33 | PC539/LH 132 | PL070/BD 384 | |
| CTS1/10/4S1D 42 | DPKS/1 214 | EC010/VBB 161 | EC3115/VBB 172 | G/305/FS/VBB 263 | MR310CL 238 | OCL/BE3/GW/S 33 | PC539/LH/LX 132 | PL075/BD 384 | |
| CTS1/10/6S 42 | DPKS/1/4H 214 | | | | | | | | PTS |
| CTS1/14/2S5D 44 | DPWM/1 213 | EC015/VBB 160 | EC3120/VBB 173 | G/500/RS/CYL 266 | MR400 240 | OCL/GW0/SP/S 32 | PC539/LH/MDS 136 | PL080/BD 384 | PIS |
| CTS1/14/3S4D 44 | DPWM/1/4H 213 | EC020/VBB 160 | EC3125/VBB 173 | G/500/RS/VBB 10 | MR410 240 | OCL/GW1/SP/S 32 | PC539/LH/MDS/LX 136 | PL085/BD 384 | |
| CTS1/14/4S3D 44 | DPWM/1/W9 213 | EC1205/VBB 154 | EC3130/VBB 173 | G/500/RS/VBB 265 | MR410CL 240 | OCL/GW3/SP/S 32 | PC539/LX 132 | PL090/BD 384 | F |
| CTS1/14/5S3D 44 | DPWM/1/W9/4H 213 | EC1210/VBB 154 | EC3205/VBB 170 | G/CARD/LOCK 267 | MRF/L/LD/CAM 248 | OCP/SP/S 31 | PC539/MDS 136 | PL095/BD 384 | Р |
| CTS1/14/6D 44 | DPWS/1 212 | EC1215/VBB 154 | EC3210/VBB 170 | G/CARD/STATION 267 | MRF/L/LD/CL 248 | OP/431 394 | PC539/MDS/LX 136 | PL100/BD 384 | P |
| CTS1/14/8S1D 44 | DPWS/1/4H 212 | EC1220/VBB 154 | EC3215/VBB 170 | G/CLS 266 | MRF/L/OC 248 | OP/432 394 | PC636 132 | PL105/BD 384 | F |
| CTS2/08/3S4D 41 | DPWS/1/W9 212 | EC1225/VBB 154 | EC3220/VBB 170 | G/CLS/500 266 | MRF/S/LD/CAM 248 | OP/631 394 | PC636/FS 134 | PL110/BD 384 | |
| CTS2/08/4S3D 41 | DPWS/1/W9/4H 212 | EC1230/VBB 155 | EC3225/VBB 170 | G/DP 266 | MRF/S/LD/CL 248 | OP/632 394 | PC636/FS/LH 134 | PL115/BD 384 | |
| CTS2/08/6D 41 | DRW260 107 | EC1235/VBB 155 | EC3230/VBB 171 | G/RAIL-C 281 | MRF/S/OC 248 | Р | PC636/FS/LH/LX 134 | PL120/BD 384 | F |
| CTS2/08/6S2D 41 | DRW260/F 107 | EC1240/VBB 155 | EC3235/VBB 171 | G/SB 267 | MRW/D2DOC/CAM 243 | PC/1000/ST 144 | PC636/FS/LX 134 | PL125/BD 384 | |
| CTS2/08/7S1D 41 | DRW260/R 107 | EC1245/VBB 155 | EC3240/VBB 171 | GC005 381 | MRW/D2DOC/CL 243 | PC/300/ST 144 | PC636/LH 132 | PL130/BD 384 | R/E |
| CTS2/08/8S 41 | DRW350 107 | EC1250/VBB 155 | EC3245/VBB 171 | GC010 381 | MRW/D2DXP/CAM 243 | PC/500/ST 144 | PC636/LH/LX 132 | PL135/BD 384 | |
| CTS2/10/3S6D 43 | DRW350/F 107 | EC1305/VBB 156 | EC3250/VBB 171 | GC015 381 | MRW/D2DXP/CL 243 | PC/600/ST 144 | PC636/LX 132 | PL140/BD 384 | |
| CTS2/10/4S5D 43 | DRW350/R 107 | EC1310/VBB 156 | EC4305/VBB 174 | GH005 96 | MRW/D4D/CAM 243 | PC/800/ST 144 | PC639 132 | PL145/BD 384 | R/ |
| CTS2/10/5S4D 43 | DRW650 107 | EC1315/VBB 156 | EC4310/VBB 174 | GH015 96 | MRW/D4D/CL 243 | PC/RLBF 145 | PC639/FS 134 | PL150/BD 384 | 10 |
| CTS2/10/6S4D 43 | DRW650/F 107 | EC1320/VBB 156 | EC4315/VBB 174 | GMBH/15 65 | MRW/LH/C5 245 | PC/RLBS 145 | PC639/FS/LH 134 | PL900 385 | |
| CTS2/10/8D 43 CTS2/10/8S2D 43 | DRW650/R 107 | EC1325/VBB 156 | EC4320/VBB 174 | GMBH/WM 95 | MRW/LHNR/S 244 | PC1036 133 | PC639/FS/LH/LX 134 | PL905 385 | |
| CTS2/10/852D 43 CTS2/14/10S6D 45 | DTMR/450 103 | EC1320/VBB 150 | EC4325/VBB 174 | GP005 189 | MRW/LHNR/S 244 | PC1036/FS 135 | PC639/FS/LX 134 | PL903 385 | |
| CTS2/14/1056D 45 CTS2/14/12D 45 | DTMR/450/SR 103 | | | | | | | | |
| CTS2/14/12D 45 CTS2/14/16S2D 45 | DTMR/600 103 | EC1335/VBB 157 | EC4330/VBB 176 | GP010 189 | MRW/S1DOC/CAM 242 | PC1036/FS/IND 135 | PC639/LH 132 | PL915 385 | |
| CTS2/14/1032D 45 CTS2/14/4S10D 45 | DTMR/600/SR 103 | EC1340/VBB 157 | EC4335/VBB 176 | GP015 189 | MRW/S1DOC/CL 242 | PC1036/FS/IND/LX 135 | PC639/LH/LX 132 | PL920 385 | |
| CTS2/14/4310D 45 CTS2/14/6S8D 45 | | EC1345/VBB 157 | EC4340/VBB 176 | GP020 189 | MRW/S1DXP/CAM 242 | PC1036/FS/LX 135 | PC639/LH/MDS 136 | PL925 385 | |
| CTS2/14/038D 45 CTS2/14/8S6D 45 | DTSF/126/2/SD 106 | EC1350/VBB 157 | EC4345/VBB 176 | GP025 189 | MRW/S1DXP/CL 242 | PC1036/IND 133 | PC639/LH/MDS/LX 136 | PMH005 94 | |
| CYB/DE2 91 | DTSF/126/2/SD 393 | EC1405/VBB 158 | EC4405/VBB 175 | GP030 189 | MRW/S2D/CAM 242 | PC1036/IND/LX 133 | PC639/LX 132 | PT/LH/C5 127 | |
| CYB/F2 91 | DTSF/450 104 | EC1410/VBB 158 | EC4410/VBB 175 | GP035 189 | MRW/S2D/CL 242 | PC1036/LX 133 | PC639/MDS 136 | PT/LHNR/S 126 | |
| CYB/G2 91 | DTSF/450/2/SD 105 | EC1415/VBB 158 | EC4415/VBB 175 | GP040 189 | MS005 137 | PC1039 133 | PC639/MDS/LX 136 | PT/LHR/S 126 | |
| CYB/J2 91 | DTSF/450/2/SD 392 | EC1420/VBB 158 | EC4420/VBB 175 | 1 | MS010 137 | PC1039/FS 135 | PC836 133 | PT050 370 | |
| CYH/D1 91 | DTSF/600 104 | EC1425/VBB 158 | EC4425/VBB 175 | IS/P/2H/H/SP 210 | MS015 137 | PC1039/FS/IND 135 | PC836/FS 135 | PT110 370 | |
| CYH/E1 91 | DTSF/600/2/SD 105 | EC1430/VBB 159 | EC4430/VBB 177 | IS/P/2H/SP 210 | MS905 145 | PC1039/FS/IND/LX 135 | PC836/FS/IND 135 | PTD/EPB/A18C12 119 | |
| CYHP/D1 266 | DTSF/600/2/SD 392 | EC1435/VBB 159 | EC4435/VBB 177 | IS/P/4H/H/SP 210 | MS910 144 | PC1039/FS/LX 135 | PC836/FS/IND/LX 135 | PTD/EPB/BD69 124 | |
| CYHP/E1 266 | DTSF/66/2/SD 106 | EC1440/VBB 159 | EC4440/VBB 177 | IS/P/4H/SP 210 | MS911 144 | PC1039/IND 133 | PC836/FS/LX 135 | PTD/EPB/BP6 121 | |
| CYT/AF1 90 | DTSF/66/2/SD 393 | EC1445/VBB 159 | EC4445/VBB 177 | IV001 215 | MT/1/6F 128 | PC1039/IND/LX 133 | PC836/IND 133 | PTD/EPB/BP9 121 | |
| CYT/CD1 90 | DTSF/750 104 | EC1450/VBB 159 | EC5000/VBB 180 | IV002 215 | MT/1/6F 354 | PC1039/IND/MDS 136 | PC836/IND/LX 133 | PTD/EPB/DR/BP9 121 | |
| CYT/DE12 89 | DTSF/750/2/SD 105 | EC2205/VBB 164 | EC5000/VBB 259 | M | MT100 112 | PC1039/IND/MDS/LX 136 | PC836/LX 133 | PTD/EPB/LP24 119 | |
| CYT/E1 90 | DTSF/750/2/SD 392 | EC2210/VBB 164 | EC5305/VBB 178 | MC/DIV 145 | MT100CL 112 | PC1039/IND/MD3/LX 130 | PC839 133 | PTD/EPB/MX60 125 | |
| CYT/FG1 88 | DTSF/900 104 | EC2210/VBB 104 EC2215/VBB 164 | EC5305/VBB 256 | MC/JK/1870 145 | MT100CL 112 MT110 112 | PC1039/LX 133 | PC839/FS 135 | PTD/EPB/NC56 123 | |
| CYT/FG1/ST 88 | DTSF/900/2/SD 105 | | | | | | | | |
| CYT/FG2/ST 88 | DTSF/900/2/SD 392 | EC2220/VBB 164 | EC5310/VBB 178 | MC/JK/870 145 | MT110CL 112 | PC1039/MDS/LX 136 | PC839/FS/IND 135 | PTD/EPB/PA32 117 | |
| CYT/J1 89 | | EC2225/VBB 165 | EC5310/VBB 256 | MC/L/1870/7S5D 143 | MT120 112 | PC336/FS 134 | PC839/FS/IND/LX 135 | PTD/HSB/A18C12 119 | |
| CYT/J1/ST 89 | DTSF/96/2/SD 106 | EC2230/VBB 165 | EC5315/VBB 178 | MC/L/1870/FS 143 | MT120CL 112 | PC336/FS/LH 134 | PC839/FS/LX 135 | PTD/HSB/BD69 124 | |
| CYWP/1 91 | DTSF/96/2/SD 393 | EC2235/VBB 165 | EC5315/VBB 256 | MC/L/1870/FSD 143 | MT130 114 | PC336/FS/LH/LX 134 | PC839/IND 133 | PTD/HSB/BP6 121 | |
| D | DTSR/450 102 | EC2240/VBB 165 | EC5320/VBB 178 | MC/L/870/3S2D 142 | MT140 114 | PC336/FS/LX 134 | PC839/IND/LX 133 | PTD/HSB/BP9 121 | |
| D3 78 | DTSR/450/SR 102 | EC2305/VBB 166 | EC5320/VBB 256 | MC/L/870/FS 142 | MT150 114 | PC536 132 | PC839/LX 133 | PTD/HSB/DR/BP9 121 | |
| D5 78 | DTSR/600 102 | EC2310/VBB 166 | EC5325/VBB 178 | MC/R/1870/7S5D 143 | MT220 113 | PC536/FS 134 | PC905 145 | PTD/HSB/LP24 119 | SAG |
| DAF/1 94 | DTSR/600/SR 102 | EC2315/VBB 166 | EC5325/VBB 256 | MC/R/1870/FS 143 | MT220CL 113 | PC536/FS/LH 134 | PCB 145 | PTD/HSB/MX60 125 | |
| DAR/1 94 | DTSR/750 102 | EC2320/VBB 166 | EC5330/VBB 179 | MC/R/1870/FSD 143 | MT230 113 | PC536/FS/LH/LX 134 | PE/30 150 | PTD/HSB/NC56 123 | |
| DAR/3 94 | DTSR/750/SR 102 | EC2325/VBB 167 | EC5330/VBB 257 | MC/R/870/3S2D 142 | MT230CL 113 | PC536/FS/LX 134 | PL005/BD 384 | PTD/HSB/PA32 117 | |
| DD 78 | DTSR/900 102 | EC2330/VBB 167 | EC5335/VBB 179 | MC/R/870/FS 142 | MW005 271 | PC536/LH 132 | PL010/BD 384 | PTS/115 371 | SAG |
| DG/1 95 | DTSR/900/SR 102 | EC2335/VBB 167 | EC5335/VBB 257 | MC/ST 144 | MW010 287 | PC536/LH/LX 132 | PL015/BD 384 | PTS/115/M 371 | 0/10 |
| DG/1/30 95 | E | EC2340/VBB 167 | EC5340/VBB 179 | MR105 250 | MW015 287 | PC536/LH/MDS 136 | PL020/BD 384 | PTS/EPB/BD46 124 | |
| DG/3 95 | E/SOCKET/4 281 | EC2405/VBB 168 | EC5340/VBB 257 | MR103 250 | 0 | PC536/LH/MDS/LX 136 | PL025/BD 384 | PTS/EPB/BP4 120 | |
| DG/3/30 95 | EBC015 20 | EC2403/VBB 108 | EC5345/VBB 179 | MR110 250 | OBL/BE0/GW/S 33 | PC536/LX 130 | PL023/BD 384 | PTS/EPB/BP6 120 | SAG |
| DGL 145 | EBC020 20 | EC2410/VBB 168 | EC5345/VBB 179 EC5345/VBB 257 | MR120 251 MR130 249 | | | PL030/BD 384 PL035/BD 384 | PTS/EPB/DR/BP6 120 | SAU |
| DH/A3/50 77 | EBC025 20 | | | | OBL/BE1/GW/S 33 | PC536/MDS 136 | | | |
| DH/A3/WF 77 | EBC023 20 | EC2420/VBB 168 | EPS005 208 | MR140 249 | OBL/BE3/GW/S 33 | PC536/MDS/LX 136 | PL040/BD 384 | PTS/EPB/DR6 118 | |
| DH/A4/50 77 | EBC030 20 | EC2425/VBB 169 | EPS010 208 | MR140CL 249 | OBL/GW0/SP/S 32 | PC539 132 | PL045/BD 384 | PTS/EPB/LA16 119 | |
| | | | | | | | | | |

| WT/42 | 224 |
|------------------|-----|
| WT/42/FBD | 224 |
| WT/42/WT/COVER | 224 |
| WT/48 | 226 |
| WT/48/COVER | 226 |
| WT/48/COVER/TBAR | 226 |
| WT/48/FBD | 226 |
| WT/48/FBD/TBAR | 226 |
| WT/48/TBAR | 226 |
| WT/50 | 225 |
| WT/50/COVER | 225 |
| WT/50/COVER/TBAR | 225 |
| WT/50/FBD | 225 |
| WT/50/FBD/TBAR | 225 |
| WT/50/TBAR | 225 |
| WT/6/2/MS | 219 |
| WT/6/2/PL/MS | 219 |
| WT/6/2/PL/SS | 218 |
| WT/6/2/SS | 218 |
| WT/6/MS | 219 |
| WT/6/PL/MS | 219 |
| WT/6/PL/SS | 218 |
| WT/6/SS | 218 |
| WT/7/MS | 222 |
| WT/9/2/MS | 220 |
| WT/9/2/PL/MS | 220 |
| WT/9/MS | 220 |
| WT/9/PL/MS | 220 |
| WT/950/DS | 85 |
| WT/950/DS/LBT | 85 |
| WT/96/3/SS | 270 |
| WT/BODY/TRAY | 85 |
| WT/CBOARD | 281 |
| WT/CYL/D/E | 281 |
| WT/CYL/F | 281 |
| WT/DT | 281 |
| WT200 | 280 |
| WT210 | 280 |
| WT250 | 280 |
| V1230 | 200 |
| _ | 109 |
| 21 000 | |
| | |
| | |
| | |

| SAC/313/S/CL/VP 146 SKC48 383 SAC/313/S/VP 146 SKC44/32 383 SAC/413/B/CL 147 SKC64/333 382 SAC/413/B/CL/VP 147 SKC64/333 382 SAC/413/B/CL/VP 147 SK 108 383 SAC/413/B/CL/VP 146 SS/222 394 SAC/413/S/CL 146 SS/221 394 SAC/413/S/CL 146 SS/221 394 SAC/413/S/VP 146 SS/40/08 98 SAC/FS/2 148 SS/44/148 113 SAC/FS/2 148 SS/43/9 95 SAC/FS/2 148 SS/43/9 95 SAC/FS/2 148 SS/43/9 95 SAC/FS/2 148 SS/44/148 17340/UN 72 SB/221 150 ST/690/UN/T 72 SF/690/UN/T 74 SB/221 150 ST/690/UN/T 72 SF/690/UN/T 74 SB/221 150 ST/690/UN/T </th <th></th> <th></th> <th></th> <th></th> | | | | |
|--|-----------------|-----|------------------|-----|
| SAC/313/S/VP 146 SKC48/CL 383 SAC/413/B/CL/VP 147 SKC64/CL 383 SAC/413/B/CL/VP 147 SKC64/CL 383 SAC/413/B/CL/VP 147 SKC64/CL 383 SAC/413/B/CL/VP 146 SS/22/1 394 SAC/413/S/CL/VP 146 SS/22/2 394 SAC/413/S/CL/VP 146 SS/4/143 81 SAC/413/S/CL/VP 146 SS/4/143 81 SAC/A13/S/VP 146 SS/4/143 81 SAC/FS/2 148 SS/4/143 81 SAC/FS/2 148 ST/340/UN/T 372 SB/221 150 ST/690/UN/T 374 SB/21 150 ST/690/UN/T 374 SB/21 150 ST/690/UN/T 374 SB/21 150 ST/690/UN/T 374 SB/21 391 TBIN4/MS 68 SB/21 391 TBIN4/MS 68 SB/21 391 TBIN4/MS </td <td>SAC/313/S/CL/VP</td> <td>146</td> <td>SKC48</td> <td>383</td> | SAC/313/S/CL/VP | 146 | SKC48 | 383 |
| SAC/413/B 147 SKC64 383 SAC/413/B/CL/1 147 SKC64/CL 383 SAC/413/B/CL/V 147 SL 145 SAC/413/B/CL/V 146 SS/22/1 394 SAC/413/S/CL/V 146 SS/22/2 394 SAC/413/S/CL/V 146 SS/22/2 394 SAC/413/S/VP 146 SS/24/1 81 SAC/413/S/VP 146 SS/4/H3 81 SAC/FS/3 148 SS/4/H14 81 SAC/FS/2 145 SS 7340/UN/T SBS/2 38 SB/H15 66 <td< td=""><td></td><td></td><td></td><td></td></td<> | | | | |
| SAC/413/B/CL 147 SKC64/CL 383 SAC/413/B/CL/VP 147 SR SAC/413/B/CL/VP 147 SR SAC/413/B/VCL 146 SS/222 394 SAC/413/S/VP 146 SS/222 394 SAC/413/S/VP 146 SS/222 394 SAC/413/S/VP 146 SS/40/MS 98 SAC/FS/2 148 SS/44/H3 81 SAC/FS/2 148 SS/44/13 81 SAC/FS/2 148 SS/44/13 81 SAC/FS/2 148 SS/44/13 91 SAC/FS/2 148 ST/340/UN/ 372 SAC/FS/2 148 ST/340/UN/ 372 SAC/FS/2 150 ST/690/UN/ 372 SB/22 150 ST/690/UN/ 372 SB/21 150 STS/690/UN/ 374 SBH/14 STSS/690/UN/ 374 SBH/156 STSS/690/UN/ 374 SBH/1491 TBIN4/MS 68 SBS/2 388 TBIN4/MS 283 SBS/2 389 TBIN5/MS 68 SBS/2 389 TBIN5/MS 68 SBS/2 389 TBIN5/MS 68 SBS/2 389 TBIN5/MS 68 SC/J/K1870/TSD 378 UB00 SC/J/K1870/TSD 378 UB00 | | | | |
| SAC/413/B/CL/VP 147 SL 145 SAC/413/B/VP 147 SR 108 SAC/413/B/VP 147 SR 108 SAC/413/S/CL/VP 146 SS/22/1 394 SAC/413/S/CL/VP 146 SS/22/1 394 SAC/413/S/VC/VP 146 SS/44/143 81 SAC/FS/2 148 SS/44/143 81 SAC/FS/2 148 SS/44/143 81 SAC/FS/2 148 ST/340/UN/T 372 SAC/WP 148 ST/340/UN/T 372 SB/221 150 ST/690/UN/T 374 SB/21 50 ST/690/UN/T 374 SB/21 150 STSS/690/UN/T 374 SB/21 391 TBIN4/MS 68 SB/21 391 TBIN4/MS 68 SB/21 391 TBIN4/MS 68 SB/21 391 TBIN6/MS 68 SC/L/1870/FSD 378 U V | | | | |
| SAC/413/B/VP 147 SR 108 SAC/413/S/CL/VP 146 SS/221 1394 SAC/413/S/CL/VP 146 SS/222 234 SAC/413/S/CL/VP 146 SS/40/MS 98 SAC/413/S/CL/VP 146 SS/44/H3 81 SAC/FS/2 148 SS/44/H3 81 SAC/FS/3 148 SS/44/H3 81 SAC/FS/2 148 ST/340/UN 372 SAC/FS/3 148 ST/340/UN.T 372 SAC/FS/2 148 ST/340/UN.T 372 SB/221 150 ST/690/UN.T 372 SB/221 150 ST/690/UN.T 374 SB/21 150 ST/690/UN.T 374 SB/21 391 TBIN4/MS 68 SB/21 391 TBIN4/MS 68 SB/21 391 TBIN4/MS 68 SB/21 391 TBIN6/MS 68 SB/21 391 TBIN6/MS 68 | | | | |
| SAC/413/S 146 SS/22/2 394 SAC/413/S/CL 146 SS/22/2 394 SAC/413/S/CL 146 SS/22/2 394 SAC/413/S/CL 148 SS/40/MS 98 SAC/FS/2 148 SS/44/H3 81 SAC/FS/2 148 SS/44/14 81 SAC/FS/4 148 SS/44/3 95 SAC/LT 148 ST/340/UN/T 372 SAC/FS/4 148 ST/340/UN/T 372 SB/221 150 ST/690/UN/T 372 SB/221 150 STS/690/UN/T 372 SB/221 56 STSS/690/UN/T 374 SBH/15 66 STSS/690/UN/T 374 SBKSL/1 391 TBIN4/MS 68 SBS/2 386 TBIN5/MS 68 SBS/1/1 391 TBIN5/MS 68 SBS/2 383 TBIN6/MS 283 SBS/2 386 TBIN6/MS 283 | | | | |
| SAC/413/S/CL 146 SS/22/2 394 SAC/413/S/VC 146 SS/36/06 96 SAC/413/S/VC 148 SS/44/H3 81 SAC/FS/2 148 SS/44/H3 81 SAC/FS/2 148 SS/48/3 95 SAC/FS/1 148 SS/48/3 95 SAC/FS/1 148 ST/340/UN 372 SB/221 150 ST/690/UN 374 SB/221 150 ST/690/UN 374 SB/21 150 STSS/690/UN 374 SB/321 150 STSS/690/UN 374 SB/321 391 TBIN4/MS 68 SBKSL/1 391 TBIN4/MS 68 SBSA/2 388 TBIN5/SS 68 SBSA/2 389 TBIN5/SS 68 SC/L/1870/TSD 378 U U SC/L/1870/TSD 378 UB010 97 SC/L/1870/TSD 378 UB010 97 SC/L/1870/TSD 378 < | | | | |
| SAC/413/S/CL/VP 146 SS/36 96 SAC/413/S/VP 146 SS/40/MS 98 SAC/FS/2 148 SS/4/H4 81 SAC/FS/3 148 SS/4/H4 81 SAC/FS/3 148 SS/4/H4 81 SAC/FS/3 148 SS/4/H4 81 SAC/FS/3 148 ST/340/UN/T 372 SAC/WP 148 ST/340/UN/T 372 SB/221 150 ST/690/UN/T 372 SB/21 150 ST/5/690/UN/T 374 SB/21 150 ST/5/690/UN/T 374 SB/21 391 TBIN4/MS 68 SBKSL/1 391 TBIN4/MS 68 SBSA/2 388 TBIN5/MS 68 SBSA/1 391 TBIN6/MS 68 SC/L/1870/FSD 378 U U SC/L/1870/FSD 378 UB/9 97 SC/L/1870/FSD 378 UB/9 97 < | | | | |
| SAC/413/S/VP 146 SS/40/MS 98 SAC/FS/2 148 SS/44/H3 81 SAC/FS/2 148 SS/44/H4 15 SAC/FS/3 148 SS/44/H4 16 SAC/FS/4 148 SS/44/H4 17 SAC/FS/4 148 ST/340/UN/T 372 SAC/LT 148 ST/340/UN/T 372 SB/221 150 ST/690/UN/T 372 SB/221 150 STS/690/UN/T 374 SBH/15 66 STSS/690/UN/T 374 SBH/15 66 STSS/690/UN/T 374 SBKSL/1 391 TBIN4/MS 283 SBS/2 38 TBIN5/MS 68 SBSL/1 391 TBIN5/MS 68 SBSL/1 391 TBIN6/MS 283 SC/DIV 379 TBIN6/MS 283 SC/DIV 379 TBIN6/MS 68 SC/L/1870/FSD 378 UB00 97 | | | | |
| SAC/FS/2 148 SS/44/H3 81 SAC/FS/3 148 SS/44/H4 81 SAC/FS/4 148 SS/44/H4 81 SAC/FS/4 148 SS/44/H3 95 SAC/LT 148 ST/340/UN1 372 SB/221 150 ST/690/UN 374 SB/221 150 ST/690/UNT 374 SB/21 150 STSS/690/UNT 374 SBH/115 66 STSS/690/UNT 374 SBH/WM 97 T T SBKSL/1 391 TBIN4/MS 68 SBS/2 388 TBIN5/MS 68 SBSL 383 TBIN5/MS 68 SBSL/1 391 TBIN6/MS 283 SBSL/1 391 TBIN6/MS 68 SC/JK/1870/TSD 378 UB 97 SC/L/1870/TSD 378 UB010 97 SC/L/1870/TSD 378 UB010 97 SC/L/1870/FSD | | | | |
| SAC/FS/3 148 SS/48/3 95 SAC/FS/4 148 SS/48/3 95 SAC/LT 148 ST/34/UUNT 372 SAC/WP 148 ST/34/UUNT 372 SB/321 150 ST/690/UNT 372 SB/321 150 ST/690/UNT 374 SB/421 150 STSS/690/UNT 374 SB/41 301 TBIN4/MS 68 SBKSL/1 391 TBIN4/MS 68 SBKSL/1 391 TBIN5/MS 68 SBSA2 389 TBIN5/MS 68 SBSL/1 391 TBIN6/MS 68 SBSL/1 391 TBIN6/MS 68 SC/DIV 379 TBIN6/MS 68 68 SC/L/1870/FSD 378 U U SC/L/1870/FSD 378 UB00 97 SC/L/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/L/1870/ | | | | |
| SAC/FS/4 148 SS/48/3 95 SAC/LT 148 ST/340/UN 372 SAC/LT 148 ST/340/UN 372 SB/221 150 ST/690/UNT 372 SB/221 150 ST/690/UNT 372 SB/221 150 ST/S/690/UNT 374 SB/21 150 STS/690/UNT 374 SB/21 56 STS/690/UNT 374 SBH/15 66 STS/690/UNT 374 SBH/13 TBIN4/MS 283 SBS/2 SBS/2 38 TBIN5/MS 68 SBS/2 38 TBIN5/MS 68 SBS/2 38 TBIN5/MS 68 SBS/2 38 TBIN6/MS 283 SBS/2 38 TBIN5/MS 68 SBS/2 38 TBIN6/MS 283 SC/L/1870/FSD 378 UB'0 7 SC/L/1870/FSD 378 UB'0 97 SC/R/1870/FSD | | | | |
| SAC/LT 148 ST/340/UN 372 SAC/WP 148 ST/340/UNT 372 SB/221 150 ST/690/UN 372 SB/321 150 ST/690/UNT 372 SB/321 150 ST/690/UNT 374 SB/321 150 ST/690/UNT 374 SBH/15 66 STSS/690/UNT 374 SBH/14 91 TBIN4/MS 68 SBKSL/1 391 TBIN4/MS 68 SBS/2 389 TBIN4/MS 68 SBS/2 389 TBIN5/MS 68 SBS/2 389 TBIN5/MS 68 SBS/2 389 TBIN5/MS 68 SBS/1 391 TBIN6/MS 283 SC/L/1870/TSD 378 UB05 97 SC/L/1870/TSD 378 UB005 97 SC/L/1870/TSD 378 UB005 97 SC/L/1870/TSD 378 UB005 97 SC/L/1870/TS | | | | |
| SAC/WP 148 ST/340/UN/T 372 SB/221 150 ST/690/UN/T 372 SB/221 150 ST/690/UN/T 374 SB/21 150 ST/690/UN/T 374 SB/421 150 STSS/690/UN/T 374 SBH/TM 97 T 5555/690/UN/T 374 SBKSL/1 391 TBIN4/MS 68 585 58 58 58 58 58 58 58 58 58 58 58 56 58 56 58 56 56 56 58 56 56 56 56 56 56 56 56 56 56 56 56 56 56 56 57 56 56 57 | | | | |
| SB/221 150 ST/690/UN 372 SB/221 150 ST/690/UN 374 SB/21 150 ST/690/UN 374 SB/21 150 STSS/690/UN 374 SBH/15 66 STSS/690/UN 374 SBH/14 391 TBIN4/MS 68 SBKSL/1 391 TBIN4/MS 283 SBSA/2 388 TBIN5/MS 283 SBSA/2 389 TBIN5/MS 283 SBSL/1 391 TBIN6/MS 283 SBSL/1 391 TBIN6/MS 283 SC/L/1870/FSD 378 TBIN6/MS 283 SC/L/1870/FSD 378 UB09 97 SC/R/1870/FSD 378 UB00 97 | | | | |
| SB/321 150 ST/690/UN/T 372 SB/421 150 STSS/690/UN/T 374 SBH/15 G STSS/690/UN/T 374 SBH/15 G STSS/690/UN/T 374 SBH/15 G STSS/690/UN/T 374 SBK/L1 391 TBIN4/MS 68 SBS/2 388 TBIN4/MS 68 SBS/1 389 TBIN5/MS 68 SBS/2 389 TBIN5/MS 68 SBS/1 391 TBIN5/MS 68 SBS/12 388 TBIN6/MS 68 SC/DIV 379 TBIN6/MS 68 SC/L/1870/FSD 378 UB00 97 SC/L/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC | | | | |
| SB/421 150 STSS/690/UN 374 SBH/WM 97 T SBKSL/1 391 TBIN4/MS 68 SBKSL/1 391 TBIN4/MS 68 SBKSL/1 391 TBIN4/MS 68 SBKSL/1 391 TBIN4/MS 68 SBSA/2 389 TBIN5/MS 68 SBSA/2 389 TBIN5/MS 68 SBSL/1 391 TBIN5/MS 68 SBSL/1 391 TBIN5/MS 68 SC/JK1470 379 TBIN6/MS 283 SC/L/1870/FSD 378 UB010 97 SC/R/1870/FSD | | | | |
| SBH/15 66 STSS/690/UN/T 374 SBH/WM 97 T SBKSL/1 391 TBIN4/MS 283 SBKSL/1 391 TBIN4/MS 283 SBSA/2 388 TBIN4/MS 283 SBSA/1 390 TBIN5/MS 68 SBSA/2 389 TBIN5/MS 283 SBSL/1 391 TBIN6/MS 283 SBSL/2 388 TBIN6/MS 283 SC/JIV/379 TBIN6/MS 283 283 SC/JI/1870/FSD 378 UB/9 97 SC/L/1870/FSD 378 UB/9 97 SC/R/1870/FSD 378 UB/0 97 SC/R/1870/FSD < | | | | |
| SBH/WM 97 T SBKSL/1 391 TBIN4/MS 68 SBKSL/1 391 TBIN4/MS 68 SBKSL/1 391 TBIN4/MS 68 SBS/2 388 TBIN4/MS 68 SBS/2 389 TBIN5/MS 68 SBS/2 389 TBIN5/MS 68 SBS/1 391 TBIN5/MS 68 SBS/1/ 391 TBIN6/MS 68 SC/DIV 379 TBIN6/MS 68 SC/J/K1070 378 TH 25 SC/L/1870/FSD 378 UB005 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB005 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 | | | | |
| SBKSL/1 391 TBIN4/MS 68 SBKSL/1 391 TBIN4/MS 68 SBKSL/1 391 TBIN4/MS 68 SBS/2 388 TBIN4/MS 68 SBS/1 389 TBIN5/MS 68 SBS/1 389 TBIN5/MS 68 SBS/1 391 TBIN6/MS 68 SC/L/1870/TSD 379 TBIN6/MS 68 SC/L/1870/TSD 378 U V SC/L/1870/TSD 378 UB005 97 SC/R/1870/TSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FS 380 WK 84 SC0202 380 WK 8221 SC440 380 WK 8221 SEP/DID <t< td=""><td></td><td></td><td></td><td></td></t<> | | | | |
| SBKSLH/1 391 TBIN4/MS 283 SBS/2 388 TBIN4/SS 68 SBSA/1 389 TBIN5/MS 68 SBSA/2 389 TBIN5/MS 283 SBSA/2 389 TBIN5/MS 283 SBSL/1 391 TBIN6/MS 283 SBSL/2 388 TBIN6/MS 283 SC/DIV 379 TBIN6/MS 283 SC/L/1870/TSD 378 UB(9) 97 SC/L/1870/FSD 378 UB(0) 97 SC/R/1870/FSD 378 UB(0) 97 | | | | |
| SBS/2 388 TBIN4/SS 68 SBSA/1 389 TBIN5/MS 68 SBSA/2 389 TBIN5/MS 283 SBSL/1 391 TBIN5/MS 283 SBSL/1 391 TBIN5/MS 283 SBSL/1 391 TBIN6/MS 283 SC/DIV 379 TBIN6/MS 283 SC/L/1870/FSD 378 U 253 SC/L/1870/FSD 378 UB005 97 SC/R/1870/FSD 378 UB005 97 SC/R/1870/FSD 378 UB005 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB0010 97 < | | | | |
| SBSA/1 389 TBIN5/MS 68 SBSA/2 389 TBIN5/MS 68 SBSL/1 391 TBIN5/MS 68 SBSL/2 388 TBIN6/MS 68 SBSL/2 388 TBIN6/MS 68 SC/J/1870/TSD 379 TBIN6/SS 68 SC/L/1870/TSD 378 UB0 97 SC/L/1870/TSD 378 UB05 97 SC/L/1870/TSD 378 UB010 97 SC/R/1870/TSD 378 UB020 84 | | | | |
| SBSA/2 389 TBIN5/MS 283 SBSL/1 391 TBIN6/SS 68 SBST/2 388 TBIN6/MS 68 SC/DIV 379 TBIN6/MS 283 SC/L/1870 379 TBIN6/SS 68 SC/L/1870/TS5D 378 TBIN6/SS 68 SC/L/1870/FSD 378 UB00 97 SC/R/1870/FSD 378 UB010 82 | | | | |
| SBSL/1 391 TBIN5/SS 68 SBST/2 388 TBIN6/MS 68 SC/DIV 379 TBIN6/MS 283 SC/JK/1870 379 TBIN6/MS 283 SC/L/1870/FSD 378 TBIN6/MS 68 SC/L/1870/FSD 378 UB/9 7 SC/L/1870/FSD 378 UB/9 97 SC/R/1870/FSD 378 UB/09 97 SC/R/1870/FSD 378 UB/09 97 SC/R/1870/FSD 378 UB/00 97 SC/R/1870/FSD 378 WSA/4/C 80 SC045 380 WSA/4/C 80 <t< td=""><td></td><td></td><td></td><td></td></t<> | | | | |
| SBST/2 388 TBIN6/MS 68 SC/DIV 379 TBIN6/MS 283 SC/JK/1870 379 TBIN6/MS 68 SC/J/1870/FSD 378 TH 25 SC/J/1870/FSD 378 UB/9 97 SC/L/1870/FSD 378 UB005 97 SC/R/1870/FSD 378 UB005 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB005 84 SC025 380 WH/950/LB1 84 SC040 380 WSK/4/C 80 | | | | |
| SC/DIV 379 TBIN6/MS 283 SC/JK/1870 379 TBIN6/SS 68 SC/L/1870/TSSD 378 TH 25 SC/L/1870/TSSD 378 UB09 97 SC/L/1870/FSD 378 UB005 97 SC/L/1870/FSD 378 UB010 97 SC/R/1870/FSD 380 WK4/K 80 SC040 380 WK4/K 80 SEP/IDI 354 WT/10/MS12/P005 221 | | | | |
| SC/JK/1870 379 TBIN6/SS 68 SC/L/1870/FSD 378 UB/9 7 SC/L/1870/FSD 378 UB/9 97 SC/L/1870/FSD 378 UB/09 97 SC/L/1870/FSD 378 UB/00 97 SC/R/1870/FSD 378 UB/01 97 SC/201380 WU/1950/LBT 84 84 SC040 380 WSA/4/C 80 SEP/1DI 354 WT/10/MS 221 SEP/1L 355 WT/10/MS 223 SEP/ISB 352 WT/11/MS 223 SEP/1B 352 WT/14/MS 221 | SBST/2 | 388 | TBIN6/MS | 68 |
| SC/L/1870/7S5D 378 TH 25 SC/L/1870/FSD 378 UB/9 97 SC/R/1870/FSD 378 UB005 97 SC/R/1870/FSD 378 UB005 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 V V SC020 80 VH/950/LBT 84 SC040 380 W W SC45380 SC045 380 WSA/4/C 80 SEP/DID 354 WT/10/MS/2P005 221 SEP/ISBH 353 WT/10/SS/2P005 221 SEP/ISBH 353 WT/10/SS/2P005 221 SEP/SBH 352 WT/13/MS 222 SEP/ISBBM 352 WT/13/MS 221 SEPS/6/EPB 352 WT/14/MS 221< | SC/DIV | 379 | TBIN6/MS | 283 |
| SC/L/1870/FS 378 U SC/L/1870/FSD 378 UB09 97 SC/R/1870/FSD 378 UB005 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 380 W 84 SC040 380 W 84 SC045 380 WSK1/C 80 SEP/IDIS 54 WT/10/MS/2P005 221 SEP/IB 353 WT/10/SS/2P005 221 SEP/IB 352 WT/11/MS 223 SEP/IGHS | | | TBIN6/SS | 68 |
| SC/L/1870/FSD 378 UB/9 97 SC/R/1870/FSD 378 UB005 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 V V SC020 380 VH/950./BT 84 SC020 380 WK 84 SC040 380 WK 84 SC045 380 WSA/4/C 80 SEP/DID 354 WT/10/MS 221 SEP/LI 355 WT/10/MS 221 SEP/LB 353 WT/10/SS /2P005 221 SEP/SBH 353 WT/11/MS 223 SEP/IGE 352 WT/11/MS 223 SEPD/6/EPB 352 WT/14/MS 221 SEPS/6/HSB 352 WT/14/MS 221 SEPS/6/HSB 352 WT/14/SS 221 SKC100/CL< | SC/L/1870/7S5D | 378 | TH | 25 |
| SC/L/1870/FSD 378 UB/9 97 SC/R/1870/FSD 378 UB005 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 UB010 97 SC/R/1870/FSD 378 V V SC020 380 VH/950./BT 84 SC020 380 WK 84 SC040 380 WK 84 SC045 380 WSA/4/C 80 SEP/DID 354 WT/10/MS 221 SEP/LI 355 WT/10/MS 221 SEP/LB 353 WT/10/SS /2P005 221 SEP/SBH 353 WT/11/MS 223 SEP/IGE 352 WT/11/MS 223 SEPD/6/EPB 352 WT/14/MS 221 SEPS/6/HSB 352 WT/14/MS 221 SEPS/6/HSB 352 WT/14/SS 221 SKC100/CL< | | | | |
| SC/R/1870/7S5D 378 UB005 97 SC/R/1870/FS 378 UB010 97 SC/20 380 VH/950 84 SC020 380 WSA/4/C 80 SC403 380 WSA/4/C 80 SEP/DID 354 WT/10/MS/2P005 221 SEP/L1 355 WT/10/SS/2P005 221 SEP/ISB 353 WT/10/SS/2P005 221 SEP/ISB 352 WT/11/SS 223 SEP/IGHB 352 WT/14/MS 221 SEPS/6/EPB 352 WT/14/MS 221 SEPS/6/HSB 352 WT/14/MS 221 SKC100 383 WT/35/SS/2 186 SKC201/CL 383 WT/35/SS/3 186 | | | | |
| SC/R/1870/FS 378 UB010 97 SC/R/1870/FSD 378 V SC/R/1870/FSD 378 V SC/R/1870/FSD 378 VH/950 SC020 380 VH/950 SC025 380 VH/950/LBT SC040 380 W SC045 380 WSA/4/C SEP/DID 354 WT/10/MS 221 SEP/DID 354 WT/10/MS 221 SEP/LH 355 WT/10/SS/2P005 221 SEP/SH 353 WT/10/SS/2P005 221 SEP/SBBMP 353 WT/10/SS/2P005 221 SEP/ISBBMP 352 WT/13/MS 223 SEP/IGHEP 352 WT/14/SS 221 SEPS/6/EPB 352 WT/14/SS 221 SEPS/6/EPB 352 WT/14/SS 221 SKC100/CL 383 WT/33/SS/3 186 SKC2010/CL 383 WT/33/SS/3 186 SKC201/TRAKA 382 | | | | |
| SC/R/1870/FSD 378 V SC020 380 VH/950 84 SC020 380 VH/950 84 SC020 380 VH/950/LBT 84 SC040 380 W SC402 SC040 380 W SC47 SC045 380 WSA/4/C 80 SEP/DID 354 WT/10/MS 221 SEP/LI 355 WT/10/MS 221 SEP/LB 55 WT/10/SS/2P005 221 SEP/IG01 355 WT/10/SS/2P005 221 SEP/IG18 353 WT/10/SS 221 225 SEP/IG19 352 WT/11/MS 223 SEP5/6/EPB 352 WT/14/MS 221 SEP5/6/HSB 352 WT/14/SS 221 SKC100 383 WT/3/SS/212P005 186 SKC20/CL 383 WT/3/SS/31 186 SKC20/CL 383 WT/3/S/S3/3 186 SKC20/CL <t< td=""><td></td><td></td><td></td><td></td></t<> | | | | |
| SC020 380 VH/950 84 SC025 380 VH/950/LET 84 SC040 380 W 80 SC045 380 WSA/4/C 80 SEP/DID 354 WSK/4/C 80 SEP/DID 354 WSK/4/C 80 SEP/DID 354 WT/10/MS 221 221 SEP/LH 355 WT/10/S2/2P005 221 SEP/ISBH 353 WT/10/S2/2P005 221 SEP/ISBH 353 WT/11/S2 223 SEP/ISBH 352 WT/11/S2 223 SEPD/66/EPB 352 WT/14/MS 221 SEPS/6/HSB 352 WT/14/MS 221 SEPS/6/HSB 352 WT/14/MS 221 SKC100/CL 383 WT/35/S5/2 186 SKC200 383 WT/35/S5/3 186 SKC201/TRAKA 382 WT/35/S3/3 186 SKC30 383 WT/37/2 187 | | | | 5. |
| SC025 380 VH/950/LBT 84 SC040 380 W W SC040 380 WSA/4/C 80 SEP/DID 354 WSK/4/C 80 SEP/DID 354 WT/10/MS/2P005 221 SEP/L1 355 WT/10/SS/2P005 221 SEP/L1 355 WT/10/SS/2P005 221 SEP/SBH 353 WT/10/SS/2P005 221 SEP/SBH 353 WT/10/SS/2P005 221 SEP/SBH 353 WT/11/SS 223 SEP/SBH 352 WT/13/MS 222 SEPS/6/EPB 352 WT/14/SS 221 SEPS/6/EPB 352 WT/14/SS 221 SKC100 383 WT/33/SS/21 186 SKC200 383 WT/33/SS/21 186 SKC200 383 WT/33/SS/21 186 SKC200 383 WT/33/SS/21 186 SKC2010CL 383 WT/33/SS/31 186 | | | | 84 |
| SC040 380 W SC045 380 WSA/4/C 80 SEP/DID 354 WSK/4/C 80 SEP/DIS 354 WT/10/MS 221 SEP/LI 355 WT/10/MS/2P005 221 SEP/LH 355 WT/10/SS/2P005 221 SEP/SBH 353 WT/10/SS/2P005 221 SEP/SBH 353 WT/10/SS/2P005 221 SEP/ISBH 353 WT/10/SS/2P005 221 SEP/IG01 355 WT/11/MS 223 SEP/ID8BMP 353 WT/13/NS 222 SEP/IG6/EP8 352 WT/14/MS 221 SEPS/6/EPB 352 WT/14/MS 221 SEPS/6/HSB 352 WT/14/SS 221 SKC100 383 WT/3/SS/3/2P005 186 SKC20 383 WT/3/SS/31 186 SKC20/CL 383 WT/3/SS/31 186 SKC20/CL 383 WT/3/SS/31 186 | | | | |
| SC045 380 WSA/4/C 80 SEP/DID 354 WSK/4/C 80 SEP/DID 354 WT/10/MS 221 SEP/LI 355 WT/10/MS/ZP005 221 SEP/LB 355 WT/10/S/ZP005 221 SEP/ISBI 353 WT/10/SS/ZP005 221 SEP/ISBI 353 WT/10/SS/ZP005 221 SEP/ISBI 353 WT/11/SS 223 SEP/ISBI 352 WT/11/MS 223 SEPD/6/EPB 352 WT/14/MS 221 SEPS/6/EPB 352 WT/14/SS 221 SEPS/6/EPB 352 WT/14/SS 221 SEPS/6/EPB 352 WT/15/SS 221 SKC100/CL 383 WT/35/SS/2 186 SKC20/CL 383 WT/33/SS/3 186 SKC20/CL 383 WT/37/2 187 SKC30 383 WT/37/2 187 SKC30/CL 383 WT/37/3 187 | | | | |
| SEP/DID 354 WSK/4/C 80 SEP/DIS 354 WT/10/MS/22005 221 SEP/L 355 WT/10/MS/2P005 221 SEP/L 355 WT/10/SS/2P005 221 SEP/SBH 353 WT/10/SS/2P005 221 SEP/USBH 353 WT/10/SS/2P005 221 SEP/USBBM 353 WT/11/MS 223 SEPD/6/EPB 352 WT/13/SS 221 SEPD/6/HSB 352 WT/14/MS 221 SEPS/6/HSB 352 WT/14/MS 221 SEPS/6/HSB 352 WT/14/SS 221 SKC100 383 WT/33/SS/2 186 SKC200 383 WT/33/SS/2 186 SKC20/CL 383 WT/33/SS/3 186 SKC20/L 383 WT/33/SS/3/ZP005 186 SKC30/CL 383 WT/37/2/ZP010 187 SKC30/CL 383 WT/37/2/ZP010 186 SKC30/CL 383 WT/37/2/ | | | | |
| SEP/DIS 354 WT/10/MS 221 SEP/L 355 WT/10/MS/2/P005 221 SEP/L 355 WT/10/SS/2/P005 221 SEP/L5H 353 WT/10/SS/2/P005 221 SEP/S6H 353 WT/10/SS/2/P005 221 SEP/ISBBM 353 WT/10/SS/2/P005 223 SEP/ISBBM 353 WT/11/NS 223 SEP/IG6/EPB 352 WT/13/MS 222 SEPS/6/EPB 352 WT/14/SS 221 SEPS/6/EPB 352 WT/14/SS 221 SEPS/6/EPB 352 WT/14/SS 221 SKC100/CL 383 WT/33/SS/2 186 SKC200/CL 383 WT/33/SS/2 186 SKC201/CL 383 WT/33/SS/3 186 SKC201/TRAKA 382 WT/33/SS/3/ZP005 186 SKC30/CL 383 WT/37/2/ZP01 187 SKC30/CL 383 WT/37/2 187 SKC30/CL 383 | | | | |
| SEP/L 355 WT/10/MS/ZP005 221 SEP/LH 355 WT/10/SS 221 SEP/SBH 353 WT/10/SS/ZP005 221 SEP/SBH 353 WT/10/SS/ZP005 221 SEP/IB 355 WT/11/MS 223 SEP/USBBMP 353 WT/11/SS 223 SEPD/6/EPB 352 WT/14/MS 221 SEPS/6/EPB 352 WT/14/MS 221 SEPS/6/EPB 352 WT/14/MS 221 SEPS/6/EPB 352 WT/14/SS 221 SEPS/6/EPB 352 WT/14/SS 221 SKC100/CL 383 WT/33/SS/2 186 SKC20/CL 383 WT/33/SS/3 186 SKC20/L 383 WT/37/2 187 SKC30 383 WT/37/2 187 SKC30/CL 383 WT/37/2 187 SKC30/CL 383 WT/37/3 187 | | | | |
| SEP/LH 355 WT/10/SS 221 SEP/SBH 353 WT/10/SS/2P005 221 SEP/IG01 355 WT/11/SS 223 SEP/USBBMP 353 WT/11/SS 223 SEP/LSBBMP 352 WT/11/MS 221 SEPD/6/EPB 352 WT/14/MS 221 SEPS/6/HSB 352 WT/14/SS 221 SEPS/6/HSB 352 WT/15/SS 221 SKC100 383 WT/33/SS/2 186 SKC200 383 WT/33/SS/3 186 SKC20/L 383 WT/33/SS/3/2P005 186 SKC20/L 383 WT/33/SS/3/2P005 186 SKC30/S83 WT/37/2/2P005 186 SKC30/CL 383 WT/37/3/SS/3/2P005 186 SKC30/CL 383 WT/37/3/SS/3/2P005 186 SKC30/CL 383 WT/37/3/SS/3/2P005 186 SKC30/CL 383 WT/37/3/3/SS/3/2P005 186 SKC35 383 WT/3 | | | | |
| SEP/SBH 353 WT/10/SS/ZP005 221 SEP/IG01 355 WT/11/IMS 223 SEP/USBBMP 353 WT/11/IMS 223 SEPD/6/EPB 352 WT/13/MS 222 SEPD/6/HSB 352 WT/14/MS 221 SEPS/6/EPB 352 WT/14/MS 221 SEPS/6/EPB 352 WT/15/SS 221 SKC100 383 WT/15/SS 221 SKC200 383 WT/33/SS/21 186 SKC201/TRAKA 382 WT/33/SS/31 186 SKC201/TRAKA 382 WT/33/SS/31 186 SKC300 383 WT/37/2/ZP005 186 SKC301/TRAKA 382 WT/37/2/ZP015 186 SKC302 383 WT/37/2/ZP016 187 SKC30/CL 383 WT/37/2/ZP010 187 SKC30/CL 383 WT/37/3 187 | | | | |
| SEP/TIG01 355 WT/11/IMS 223 SEP/USBBMP 353 WT/11/S 223 SEPD/6/EPB 352 WT/13/MS 222 SEPD/6/EPB 352 WT/14/MS 221 SEPS/6/EPB 352 WT/14/MS 221 SEPS/6/EPB 352 WT/15/MS 221 SKC100 383 WT/15/SS 221 SKC100/CL 383 WT/33/SS/2/2P005 186 SKC20/CL 383 WT/33/SS/3/186 186 SKC20/L17RAKA 382 WT/33/SS/3/2P005 186 SKC30/283 WT/37/S1/2187 186 SKC30/21 187 SKC30/217RAKA 382 WT/37/21/2187 186 SKC30/21 187 SKC30/22 383 WT/37/21/2187 187 187 187 SKC30/23 383 WT/37/3 187 187 187 | | | | |
| SEP/USBBMP 353 WT/11/SS 223 SEPD/6/EPB 352 WT/13/MS 222 SEPD/6/EPB 352 WT/14/MS 221 SEPS/6/EPB 352 WT/14/SS 221 SEPS/6/HSB 352 WT/14/SS 221 SEPS/6/HSB 352 WT/14/SS 221 SKC100 383 WT/33/SS/2 186 SKC2002 383 WT/33/SS/3 186 SKC20/CL 383 WT/33/SS/3 186 SKC20/CL 383 WT/33/SS/3/ZP005 186 SKC30 383 WT/37/ZP101 187 SKC30/CL 383 WT/37/Z1010 187 SKC30/CL 383 WT/37/Z1010 187 | | | | |
| SEPD/6/EPB 352 WT/13/MS 222 SEPD/6/HSB 352 WT/14/MS 221 SEPS/6/EPB 352 WT/14/MS 221 SEPS/6/HSB 352 WT/15/MS 221 SKC100 383 WT/15/MS 221 SKC100/CL 383 WT/33/SS/2 186 SKC200 83 WT/33/SS/3 186 SKC201/TRAKA 382 WT/33/SS/3/ZP005 186 SKC300 383 WT/37/SS/3/ZP105 186 SKC30/CL 383 WT/37/SS/3/ZP105 186 SKC30/CL 383 WT/37/SS/3/ZP105 186 SKC30/CL 383 WT/37/S/ZP105 186 SKC30/CL 383 WT/37/2/ZP105 186 SKC30/CL 383 WT/37/3/ZP105 186 SKC30/CL 383 WT/37/3/ZP105 186 SKC35 383 WT/37/3/3 187 | | | | |
| SEPD/6/HSB 352 WT/14/IMS 221 SEPS/6/EPB 352 WT/14/ISS 221 SEPS/6/EPB 352 WT/15/ISS 221 SEPS/6/EB 352 WT/15/ISS 221 SKC100/CL 383 WT/33/SS/2 186 SKC200 83 WT/33/SS/3 186 SKC20/L 383 WT/33/SS/3/ZP005 186 SKC20/L 383 WT/33/SS/3/ZP005 186 SKC20/L 383 WT/37/Z/ZP005 186 SKC30/L 383 WT/37/Z/ZP010 187 SKC30/L 383 WT/37/Z/ZP10 187 SKC30/L 383 WT/37/3/ 187 | | | | |
| SEPS/6/EPB 352 WT/14/SS 221 SEPS/6/HSB 352 WT/15/SS 221 SKC100 383 WT/15/SS 221 SKC100/CL 383 WT/33/SS/2 186 SKC20/CL 383 WT/33/SS/2 186 SKC20/CL 383 WT/33/SS/3 186 SKC21/TRAKA 382 WT/33/SS/3/ZP005 186 SKC30 383 WT/37/SS/3/ZP105 186 SKC30 383 WT/37/S172 187 SKC30/CL 383 WT/37/2/ZP101 187 SKC30/CL 383 WT/37/2/ZP101 187 SKC30/S2 383 WT/37/2/ZP101 187 | | | | |
| SEPS/6/HSB 352 WT/15/MS 221 SKC100 383 WT/15/SS 221 SKC100/CL 383 WT/33/SS/2 186 SKC20 383 WT/33/SS/3 186 SKC20/CL 383 WT/33/SS/3 186 SKC20/CL 383 WT/33/SS/3/ZP005 186 SKC20/CL 383 WT/37/SS/3/ZP100 187 SKC30/CL 383 WT/37/2/ZP101 187 SKC30/CL 383 WT/37/2/ZP101 187 SKC30/SS 383 WT/37/3 187 | | | | |
| SKC100 383 WT/15/SS 221 SKC100/CL 383 WT/33/SS/2 186 SKC20 383 WT/33/SS/2 186 SKC20/CL 383 WT/33/SS/3 186 SKC20/L 383 WT/33/SS/3 186 SKC20/L 383 WT/33/SS/3/ZP005 186 SKC30 383 WT/37/ZP010 187 SKC30/CL 383 WT/37/2/ZP101 187 SKC33 383 WT/37/3 187 | | | | |
| SKC100/CL 383 WT/33/SS/2 186 SKC20 383 WT/33/SS/2/ZP005 186 SKC20/CL 383 WT/33/SS/3/186 186 SKC21/TRAKA 382 WT/33/SS/3/ZP005 186 SKC30 383 WT/37/Z12 187 SKC30/CL 383 WT/37/2/ZP010 187 SKC33/S33 WT/37/2/ZP010 187 SKC35 383 WT/37/3 | | | | |
| SKC20 383 WT/33/SS/2/ZP005 186 SKC20/CL 383 WT/33/SS/3/SS/3/ZP005 186 SKC21/TRAKA 382 WT/33/SS/3/ZP005 186 SKC30 383 WT/37/ZP101 187 SKC30/CL 383 WT/37/2/2101 187 SKC30/CL 383 WT/37/2/2101 187 SKC30/SC 383 WT/37/2/27010 187 | | | | |
| SKC20/CL 383 WT/33/SS/3 186 SKC21/TRAKA 382 WT/33/SS/3/2P005 186 SKC30 383 WT/37/2 187 SKC30/CL 383 WT/37/2/P010 187 SKC30/SC 383 WT/37/2/P010 187 SKC30/SC 383 WT/37/2/P010 187 | | | | |
| SKC21/TRAKA 382 WT/33/SS/3/ZP005 186 SKC30 383 WT/37/2 187 SKC30/CL 383 WT/37/2/ZP010 187 SKC33 383 WT/37/2/ZP101 187 SKC35 383 WT/37/3/3 187 | SKC20 | 383 | WT/33/SS/2/ZP005 | 186 |
| SKC30 383 WT/37/2 187 SKC30/CL 383 WT/37/2/ZP010 187 SKC35 383 WT/37/3 187 | SKC20/CL | 383 | WT/33/SS/3 | 186 |
| SKC30/CL 383 WT/37/2/ZP010 187 SKC35 383 WT/37/3 187 | SKC21/TRAKA | 382 | WT/33/SS/3/ZP005 | 186 |
| SKC30/CL 383 WT/37/2/ZP010 187 SKC35 383 WT/37/3 187 | SKC30 | 383 | WT/37/2 | 187 |
| | | | | |
| | SKC35 | 383 | WT/37/3 | 187 |
| | | | | |
| | | | | |

| PTS/EPB/MX40 PTS/EPB/NC30 PTS/EPB/PA24 PTS/HSB/BD4 PTS/HSB/BP6 PTS/HSB/DR/BP6 PTS/HSB/DR/BP6 PTS/HSB/LA16 PTS/HSB/LA16 PTS/HSB/LA10 PTS/HSB/MC30 PTS/HSB/NC30 PTS/HSB/PA24 | 125 122 116 124 120 120 120 118 119 125 122 116 | |
|---|--|--|
| R/BAG/BLUE R/BAG/GREEN R/BAG/GREY R/BAG/ORANGE R/BAG/ORANGE R/BAG/VELLOW R/C/TT R/C/145 RC/21/5 RC/21/5 RC/21/5 RC/51/5 RD/B5/SS RD/U5/SSS | 223 223 223 223 223 223 223 223 224 246 98 98 98 98 99 99 99 99 99 78 78 78 78 78 78 25 25 25 25 25 25 | |
| S SA040 SA040/CL SA045 SA045/CL SA045/CL SA045/CL SA0213/B/CL/VP SAC/213/B/CL/VP SAC/213/B/CL SAC/213/S/CL SAC/213/S/CL SAC/213/S/CL SAC/213/B/CL SAC/313/B/CL SAC/313/B/CL SAC/313/B/CL SAC/313/S/CL | 129 129 149 149 149 147 147 147 147 146 146 146 147 147 147 147 146 146 | |

Paint & Laminate

Vinvl

Paint

Painted surfaces should be cleaned using a soft damp cloth and polished completely dry. Do not use abrasive cloth. Persistent exposure to wet or humid conditions could result in rust marks or cause the paint to lift. Concentrated bleach, synthetic detergents, solvents, wax polishes or aerosol sprays should not be used.



Laminate & Melamine

Laminates and melamine should be cleaned using a soft damp cloth and polished completely dry. Do not use abrasive cloth. Concentrated bleach, synthetic detergents, solvents, wax polishes or aerosol sprays should not be used.



GW - Grey White LG - Light Grev WH - White BE - Beech

Ambla New Richmond in Colour - Spradling Valencia C5 - Memphis Faux Leather

Vinyl's are resistant to most mild acids, alkalis, drinks & household stains. Clean with a damp, soapy cloth and rinse well with clean water. Dry the surface thoroughly. A soft brush with soapy water may be used to remove stubborn stains. Ensure surfaces are then rinsed and dried. Do not use solvents, bleaches*, abrasives, synthetic detergents, wax polishes or aerosol sprays. Disinfectants may be used in dilution as specified by their manufacturers.

*Chlorine Releasing Agents: Levels of 1-10,000 parts per million should pose no problems as long as the material does not endure long soakings at the higher concentrations. Rinse well with clean water. These include Haz-Tabs & other disinfectants commonly used in hospital infection control programmes. Surfaces are NaDCC compatible.

Spradling Valencia C5

Ambla New Richmond in Colour



Memphis Faux Leather




VMVAN - Vanilla VMESP - Espresso VMGAR - Garden VMROY - Royal



Fabric & Plastic


Camira Advantage

Vacuum regularly. Wipe clean with a damp cloth or shampoo using proprietary upholstery shampoo.



Camira 24/7 Flax

Vacuum regularly. Professionally dry clean or wipe clean with a damp cloth.



VP Multi-Way Stretch

Wipe with mild detergent. Wash or dry clean using detergent up to 75°C. For disinfection use sodium hypochlorite at no more than 1,000 parts per million of active chlorine. Do not use biological or phenolic-based cleaners. Wash and dry completely before storage. Some surface wrinkling may occur.





Haz-Tabs or Chlor-Clean


SAVL - Savanna Lavender

Plastic Seat Shells

Plastic should be cleaned using a soft damp cloth and polished completely dry. Do not use abrasive cloth. Concentrated bleach, synthetic detergents, solvents, wax polishes or aerosol sprays should not be used.

Panaz Patternistas For Children

Vacuum regularly. Wipe clean with a damp cloth. Can be cleaned with water based proprietary cleaners and 5% sodium chlorite bleach solution,





Mobile Infusion Stands Polyurethane, Easy Clean

sales@bristolmaid.com +44 (0)1258 484455 www.bristolmaid.com

See Page 210

Technical Information

Dimensions & Weight

- Dimensions are stated in millimetres (width x depth x height)
- Weights are stated in kilograms
- Please ensure you are aware of the item size & weight prior to ordering

Castors

- All castors are fitted with non-marking tyres (with exception of some anti-static versions)
- Anti-static castors can be identified by either a yellow mark on the edge of the tyre, or where twin wheel castors are used a different colouration to the tyre (the result of the addition of a conductive material)
- Castors & wheel bearings do not require specific maintenance other than regular cleaning
- The use of oil will breakdown the grease applied to the bearings, reducing the life of the castor
- We recommend that fittings are checked regularly & re-tightened as required
- · Where anti-static castors are fitted, we recommend that the tyre is periodically cleaned to remove any dirt to ensure a good contact between the tyre & the floor

Locks & Keys

- · Locks should operate smoothly. If they do not please arrange urgent attention & where necessary, obtain a replacement
- Locks require only minimal maintenance and should only be lubricated using a dry lubricant such as graphite which can be purchased from most locksmiths (NEVER USE OIL)
- When ordering replacement locks or keys, please ensure that the product code & any letters/ numbers from the key & lock are available

Stainless Steel

- We use 304 & 430 grade stainless steel
- Grade 304 is non-magnetic & used on products which are not painted
- Grade 430 is magnetic & will be subsequently painted
- · Please note that the term 'stainless' is not entirely accurate, as certain chemicals can 'eat' into the surface. In particular, 'iodine' which can badly corrode/stain & as a result should not be used as a cleaning agent. Stainless steel products should always be dried after cleaning to reduce staining

Paint

- We only use lead free epoxy polyester powder paint
- Where possible colours are applied in special booths to allow for overspray to be recycled
- Alternative colours are available at additional cost

Wood, Melamine & Laminate

- We only use FSC accredited materials. This provides assurances that the material is produced from well-managed forests &/or recycled materials
- We use a combination of melamine faced MDF / chipboard & high pressure laminates bonded MDF of chipboard substrates
- Edges are protected with 2mm ABS. T-Barb edging or sealed with varnish
- Both melamine & laminate panels offer some resistance to water penetration under normal usage
- · We only use water based adhesives & varnishes

Upholstered Panels & Fabrics

- Vinvl upholstered panels are manufactured from 'Flame Retardant' materials. independently tested to BS 5852: 2006 Ignition Source 5, BS EN 1021-1: 2006, BS EN 1021-2: 2006 & BS 7176: 2007 Medium Hazard. Standard vinyl colours also contain an anti-microbial compound within its construction, which constantly migrates to the surface to replenish any material which is removed by cleaning or general usage
- · Fabric upholstered panels are manufactured from 'Flame Retardant' materials, independently tested to BS 5852: 2006 Ignition Source 5, BS EN 1021-1:2006, BS EN 1021-2:2006 & BS 7176: 2007 Low Hazard
- Ward screens are manufactured from 'Flame Retardant' materials, independently tested to BS 5438: Part 2

Electric Columns & Actuators

• Electric columns & actuators comply with EMC Directive 2014/30/EU. Low Voltage Directive 2014/35/EU & RoHS2 Directive 2011/65/EU

Preventative Maintenance, Spares & Service Support

Bristol Maid[™] products are designed & manufactured to cater for the demands of the healthcare environment, with the majority (under normal conditions) providing many years of trouble free service. Where annual maintenance or periodic inspections are required, we can offer comprehensive planned preventative maintenance through either a scheduled service contract or call out by our own trained Engineers. For your peace of mind & continued product reliability, please contact our customer service team for details on the following services.

- Preventative maintenance contracts & packages
- Supply of spare parts
- Site audit & equipment recommendations

Cleaning Guidance & Instructions

The maintenance of high standards of cleanliness on all surfaces and equipment is a crucial factor in the prevention and control of infection. Whilst all dust, dirt and moisture can harbour infection, the key risks are associated with contamination arising from contact with blood, body fluids, mucous membranes and damaged skin. In all cases surfaces, must be thoroughly cleaned and disinfected in order to destroy any pathogenic micro-organisms.

Bristol Maid products have been designed for use in demanding environments. Guidance & instructions can be found at - www.bristolmaid.com/cleaning

Technical Information

Medical Devices

Where applicable. Hospital Metalcraft product is marked with the CE symbol to demonstrate compliance with the essential requirements of any European Directives (& transposing regulations) that may apply & for which a Declaration of Conformity is held. Hospital Metalcraft is registered with the MHRA (Medicines & Healthcare Regulatory Authority) in connection with product classified as 'Class 1 Medical Devices' for purposes of the Medical Devices Directive (93/42/EEC) with the registration number CA001358.

WEEE Regulations

Hospital Metalcraft meets its obligations as a business to business producer of Electrical and Electronic Equipment (EEE) through:

- Membership of the Veolia Environmental Services WEEE, (Waste Electrical & Electronic Equipment), Producer Compliance Scheme (Registration No. WEE/KE0677RX)
- Marking EEE placed on the market by Hospital Metalcraft with the mandatory crossed out wheeled bin symbol & Bristol Maid™ logo
- Making information available, (for EEE placed on the market after April 2007), to assist those engaged in treatment, recovery & recycling activities
- · Arranging for the collection, treatment, recovery & environmentally sound disposal of obligated WEEE



Meet The Team

Here at Bristol Maid[™] our dedicated sales and customer service team are committed to finding a solution that works for you. Whether you require a full ward refurbishment or just simply some advice on our products, we are here to help. Refer to the map on the right for your regional sales representative or alternatively contact our friendly customer services team

Sales Director



an Farnfield 01258 451338 ifarnfield@bristolmaid.com

Technical Services Manager



416

Paul Open popen@bristolmaid.com


Sales Support



an Mills Sales Support Engineer (South) 07921 473241 mills@bristolmaid.com



Eric Hornby Sales Support Engineer (North) 07702 000953 ehornby@bristolmaid.com



Customer Services

Sophie Williams Customer Service Coordinator swilliams@bristolmaid.com

mchampion@bristolmaid.com

Sammy Foord

Mark Champion

Customer Service Team Leader

Assistant Customer Service Team Leader

sfoord@bristolmaid.com



Jessica Haines Customer Service Coordinator jhaines@bristolmaid.com

sales@bristolmaid.com +44 (0)1258 484455 www.bristolmaid.com



Kirsty Geoghegan Customer Service Team Leader kgeoghegan@bristolmaid.com



Sam Symes Customer Service Coordinator ssymes@bristolmaid.com



Stacey Finch sfinch@bristolmaid.com

Customer Service Coordinator

jbaines@bristolmaid.com

Export



Simon Haines Export Manager Tel: +44 (0)7917 207617 DDI: +44 (0)1258 484430 Email: shaines@bristolmaid.com

Bristol Maid's policy is to work together with selected partners who are able to support and advise in the selection of appropriate 'Healthcare Solutions', also providing a high level of customer care and portraying our traditional values. Our continued success and growth within international markets has created a demand for additional partners.

National Sales Manager



Steve Poller 07711 095899 spoller@bristolmaid.com

Grant Macpherson

External Sales Representatives



Scotland, N.Ireland & Northumberland 07768 057675 gmacpherson@bristolmaid.com



Chris Davis London 07725 543862 cdavis@bristolmaid.com



ohn Baines South & South West 07876 441477



Duncan Holder North 07790 916052 dholder@bristolmaid.com



Russell Sheard Midlands 07887 534755 rsheard@bristolmaid.com



Nick Moyster South East 07500 039301 nmoyster@bristolmaid.com



David Morton-Jones Export Coordinator +44 (0)1258 484468 dmortonjones@bristolmaid.com

Manufacturing Medical Furniture & Equipment Since 1953

With over 65 years of experience, Bristol Maid is a leading manufacturer and supplier of medical furniture & equipment. From design and manufacture to supply and installation, Bristol Maid has a reputation for providing high product quality, excellent customer service and value for money.

We take pride in working in close partnership with our customers applying innovative design and technology to deliver effective solutions for today while anticipating the needs of tomorrow.











