

Integrity | Professionalism | Cooperation | Innovation

**HEDY**  
Healthy Every Day



**HEDY**  
Healthy Every Day

## HealFusion Infusion Pump i7



HEDY Medical Device Co., Ltd.

Address: 10 th Floor, HEDY building, No 286, Science Avenue, Science City, Guangzhou, P.R. China

Tel: +86-20-83070896

Fax: +86-20-83070817

Email: overseas@hedymed.com

Website: www.hedymed.com





## Healfusion Infusion Pump


HealFusion Pump is a brand-new generation pumps developed by the elite team of HedyMed. It applies advanced modular design and it works with dock to provide multi -channel infusion. accurate date acquires automatically in real-time through powerful network transmission. It brings better care to patients and is very helpful for medical staff to get rid of burdensome record and operation. Sophisticated stepping motor and driving system make sure the accurate infusion result .Faulty tree analysis which optimizes every step assures the stability and safety of our pumps. Advanced touch screen technology makes all operations are convenient and intelligent. With HealFusion pumps launching in the market infusion treatment would become safer, more accurate and intelligent.

Infusion pump i7



## Accurate

- » All module uses high-precision stepper motor and system making the treatment more accurate.
- » i7 infusion pump has flow compensation function, which can avoid the error caused by pipeline collapse.



## Easy to Operate

- Infusion pump i7 applies 3.5 inch color touch screen and mobile style UI design which is easy to operation.
- Abundant infusion modes meet various clinic demands without complicated conversion and operation.
- The day/night mode could be switched automatically, and it is unnecessary to set one by one and noise pollution would be decreased effectively.



- Flexible stack meets demands of different bedside infusion. it is convenient for transportation of patients when it installs on the elaborate handle.



### >>> Modular Design

Infusion pump applies the same operation panel and every pump is an independent infusion module, all these are for the convenience of operation and users' training.

### >>> High Strength Shell

The shell is made from PBT plastic and it has well characteristics of resistance to corrosion. The corrosion risk of medicine leakage could be avoided.




- Stacked track  
Easy to stack pumps



- 3.5 inch large screen  
Observe screen information clearly



- Alarm light  
Easy to observe even in far distance



- Lock point  
Lock all pumps together firmly



- Anti-free flow button  
For safety



- Panel and touch screen  
Easy to operation






Safe



Dural CPU design assures the safety of system.

Dynamic Pressure Surveillance (DPS )

11 levels of alarm range is adjustable, it not only helps medical staff to set reasonable pressure alarm range, but also helps to find out pipeline off, leakage of injection and occlusion in time.



Intelligent

It can work with our patent designed dock cabin for one-stop power supply and pipeline management. It makes the ward more comfortable and brings you to an intelligent infusion management generation.



Integrity | Professionalism | Cooperation | Innovation

**HEDY**  
Healthy Every Day



## HealFusion Syringe Pump s7



HEDY Medical Device Co., Ltd.

Address: 8th Floor, HEDY Innovative Industry Park, No. 63, Punan Road, Huangpu District, Guangzhou, P.R. China

Tel: +86-20-83070896

Fax: +86-20-83070817

Email: [overseas@hedymed.com](mailto:overseas@hedymed.com)

Website: [www.hedymed.com](http://www.hedymed.com)





## Healfusion Syringe Pump


HealFusion Pump is a brand-new generation pumps developed by the elite team of HedyMed. It applies advanced modular design and it works with dock to provide multi-channel infusion. It acquires accurate date automatically in real-time through powerful network transmission. It brings better care to patients and is very helpful for medical staff to get rid of burdensome record and operation. Sophisticated stepping motor and driving system make sure the accurate infusion result. Faulty tree analysis which optimizes every step assures the stability and safety of our pumps. Advanced touch screen technology makes all operations are convenient and intelligent. With HealFusion pumps launching in the market infusion treatment would become safer, more accurate and intelligent.

Syringe pump s7



## Accurate

- »» All module uses high-precision stepper motor and driving system making the treatment more accurate.
- »» The minimum resolution of this high-precision stepper motor is 36.6nm/step, and the velocity of driving system is 0.01ml/h, which meet the need of high accuracy injection for toxic/anesthetic drug or neonate.




S7 syringe pump has fast start functions, which can eliminate error caused by starting delay.

## Easy to Operate



- Syringe pump s7 applies 3.5 inch color touch screen and mobile style UI design which is easy to operation.
- Abundant infusion modes meet various clinic demands without complicated conversion and operation.
- The day/night mode could be switched automatically, and it is unnecessary to set one by one and noise pollution would be decreased effectively.



- Flexible stack meets demands of different bedside syringe. it is convenient for transportation of patients when it installs on the elaborate handle.



### >>> Modular Design

Syringe pump applies the same operation panel and every pump is an independent infusion module, all these are for the convenience of operation and users' training.

### >>> High Strength Shell

The shell is made from PBT plastic and it has well characteristics of resistance to corrosion. The corrosion risk of medicine leakage could be avoided.



- 3.5 inch large screen  
Observe screen information clearly



- Alarm light  
Easy to observe even in far distance



- Piston clamp  
Convenient to load pumps



- Stacked track  
Easy to stack pumps

- Lock point  
Lock all pumps together firmly

- Panel and touch screen  
Easy to operation

- Injector size  
5/10/20/30/50/60ml




Safe



Dural CPU design assures the safety of system.

Dynamic Pressure Surveillance (DPS )

11 levels of alarm range is adjustable, it not only helps medical staff to set reasonable pressure alarm range, but also helps to find out pipeline off, leakage of injection and occlusion in time.



Intelligent

It can work with our patent designed dock cabin for one-stop power supply and pipeline management. It makes the ward more comfortable and brings you to an intelligent infusion management generation.

