

**Product Manual** 

#### **VICTOR MEDICAL INSTRUMENTS CO., LTD.**

Add: No.21 Huafeng Road, Dingyan Sub-district, Wujin District, Changzhou, Jiangsu, China Zip code:213018 Tel: +86 519 88411338 Fax: +86 519 88411358







## **Company Profile**

Victor Medical Instruments Co., Ltd is a professional firm dedicated to developing, manufacturing and marketing for auto surgical staplers and Minimally Invasive Surgery related products. By providing surgeons with comprehensive surgical solutions, and offering our customers with advanced surgical devices, professional medical education, as well as superior customer service, Victor Medical has been constantly practicing the creed of being "professional, focused and concentrated."

By strictly comply to the Quality System operational requirements of ISO13485, GMP and Brazil ANVISA, the Victor Medical products have been successively certified by CFDA, The US FDA, CE, The Korean FDA, India CDSCO, Saudi Arabia SFDA, Iran and relevant administrations from other nations. Our products have so far been marketed in 16 Countries around the globe and hundreds of Medical Institutions in China. The company has created a strong R&D team. As a result, we obtained 35 Chinese patents, including 14 invention patents, and 9 patents that are being evaluated. Victor Medical has been successively awarded as Jiangsu Province High-tech Enterprise, Small & Medium Sized Technology Enterprise in Jiangsu Province, Private Technological Enterprise in Jiangsu Province, Engineering Technology Research Center in Changzhou City. It also receives awards of National Creative Funds, Jiangsu Province Renowned Trademark and Technological Creation.

By upholding the principles of striving for perfection and continuous creations, Victor Medical is aspired to become the leader in the field of Chinese surgical devices.

# VICTOR MEDICAL PROFILE

•Leaded and drafted the national industry standard for Endoscopic Linear Cutter and Reload 2019 ● To be recognized as Industrial Design and Research Center in Changzhou City • Established experimental base for post-doctoral innovation in Changzhou City Established a postgraduate workstation in cooperation with Southeast University Explored new way of industrial-academic-research cooperation and promoted patent transfer by 2018 cooperating with clinician • Got approval of the registration of India, Saudi Arabia and Iran markets Merger and reorganization with Canopus, and had a sucessful IPO 2017 at China New Third Board Market with Trade No 834478 Obtained 2 national patents Annual sales outpaced 20 million US Dollars • The 3<sup>rd</sup> generation staplers for open surgery launched into market 2016 Obtained US FDA approvalsObtained 6 national patents ● To be awarded national innovation fund by new generation Endoscopic Linear Cutter 2015 ● To be recognized as Engineering Technology Research Center in Changzhou City Obtained 3 national patents • New generation Endoscopic Linear Cutter launched into market 2014 • Got the approval of Brazil ANVISA and KFDA GMP Inspection Obtained 3 national patents ● To be recognized as Small & Medium Sized Technology Enterprise in Jiangsu Province ●To be awarded national innovation fund
●Obtained 1 national patent ● To be recognized as Jiangsu Province High-tech Enterprise • Research & Develop the 1st generation Endoscopic Linear Cutter Obtained 2 national patents ● To be recognized as Science and Technology Enterprise in Jiangsu Province Obtained 9 national patents • The 1st generation stapler for open surgery launched into market Obtained 6 national patents

● Foundation of Victor Medical ● Acquired CE and ISO 13485 certification

2013

2012

2011

2010



# Single Use Endoscopic Linear Cutter












## The staple height reaches up to 5.0mm for green reload

The staple height reaches up to 5.0mm for green reload
Ensure suturing extra—thick tissue
To effectively prevent the staples opening by extra—thick tissue





#### Different staple heights from inner to outer rows

Adapt to different tissue thickness Excellent hemostatic effect To provide optimal blood supply



## Grasping surface design

Providing the least tissue slippage during firing



## Reinforced and one-piece I-beam design

Effectively control of the space between reload and anvil, consistently force control of the parallel clamping, reducing tissue damage in second closure, ensure optimal B-shape formation



## Extra-thin reinforced and one-piece anvil

Stable performance during firing. Optimal staple formation. To be easier passing the narrow tissue gap with the extra—thin anvil





#### oriented design normal design

## Drop-shaped oriented anvil design

Ensure consistent B-shape formation, excellent hemostasis





Victor Medical

## 6 level angle adjustment from left to right

The articulating angle could be adjusted to 15 degrees, 30 degrees and 45 degrees. 6 level angle adjustment could easily reach into the anastomosis position





## One-handed Operation

Easy and convenient to close and reopen the jaw by one-handed operation system



## Single Use Endoscopic Linear Cutter



#### Articulating Reload


Figure	Item Code	Suture Length	Staple Height	Compressed Tissue Thickness of Application	Color	Color	Packaging
	PECFR30Y	30mm	2.3mm/2.3mm/2.3mm	0.70-1.20mm	Grey		6pcs/box
	PECFR30W	30mm	2.0mm/2.5mm/3.0mm	0.88-1.50mm	White	$\bigcirc$	6pcs/box
	PECFR30B	30mm	3.0mm/3.5mm/4.0mm	1.50-2.25mm	Blue		6pcs/box
	PECFR45Y	45mm	2.3mm/2.3mm/2.3mm	0.70-1.20mm	Grey		6pcs/box
	PECFR45W	45mm	2.0mm/2.5mm/3.0mm	0.88-1.50mm	White	$\bigcirc$	6pcs/box
	PECFR45B	45mm	3.0mm/3.5mm/4.0mm	1.50-2.25mm	Blue		6pcs/box
	PECFR45G	45mm	4.0mm/4.5mm/5.0mm	2.25-3.00mm	Green		6pcs/box
- Frag 401	PECFR60Y	60mm	2.3mm/2.3mm/2.3mm	0.70-1.20mm	Grey		6pcs/box
	PECFR60W	60mm	2.0mm/2.5mm/3.0mm	0.88-1.50mm	White	$\bigcirc$	6pcs/box
	PECFR60B	60mm	3.0mm/3.5mm/4.0mm	1.50-2.25mm	Blue		6pcs/box
	PECFR60A	60mm	4.2mm/4.2mm/4.2mm	1.75-2.50mm	Gold		6pcs/box
	PECFR60G	60mm	4.0mm/4.5mm/5.0mm	2.25-3.00mm	Green		6pcs/box

#### Straight Reload

Figure	Item Code	Suture Length	Staple Height	Compressed Tissue Thickness of Application	Color	Color	Packaging
	PECR30Y	30mm	2.3mm/2.3mm/2.3mm	0.70-1.20mm	Grey		6pcs/box
	PECR30W	30mm	2.0mm/2.5mm/3.0mm	0.88-1.50mm	White	$\bigcirc$	6pcs/box
	PECR30B	30mm	3.0mm/3.5mm/4.0mm	1.50-2.25mm	Blue		6pcs/box
_	PECR45Y	45mm	2.3mm/2.3mm/2.3mm	0.70-1.20mm	Grey		6pcs/box
	PECR45W	45mm	2.0mm/2.5mm/3.0mm	0.88-1.50mm	White	$\bigcirc$	6pcs/box
	PECR45B	45mm	3.0mm/3.5mm/4.0mm	1.50-2.25mm	Blue		6pcs/box
	PECR45G	45mm	4.0mm/4.5mm/5.0mm	2.25-3.00mm	Green		6pcs/box
	PECR60Y	60mm	2.3mm/2.3mm/2.3mm	0.70-1.20mm	Grey		6pcs/box
	PECR60W	60mm	2.0mm/2.5mm/3.0mm	0.88-1.50mm	White	0	6pcs/box
	PECR60B	60mm	3.0mm/3.5mm/4.0mm	1.50-2.25mm	Blue		6pcs/box
	PECR60A	60mm	4.2mm/4.2mm/4.2mm	1.75-2.50mm	Gold		6pcs/box
_	PECR60G	60mm	4.0mm/4.5mm/5.0mm	2.25-3.00mm	Green		6pcs/box

#### Single Use Endoscopic Linear Cutter Parameter Chart





#### Solution for Bariatric Procedure **O**



The green reload, which is applied to extra-thick tissue of pylorus



Single Use Endoscopic Linear Cutter—Extra Long



Single Use Bladeless Trocar—Extra Long

### Solution for Thoracic Procedure •



The green reload, which is applied to extra thick and hard tissue of lung



Single Use Endoscopic Linear Cutter—Short




The grey reload, which provides excellent hemostatic effect



## Single Use Endoscopic Linear Cutter







## Enhanced compression space control system

Ensure space between anvil and reload keeping the same from proximal to distal side. During firing, to provide consistent B-sharp formation



## Real one-handed operation

Easy and convenient to open and close the jaw by one-handed operation



## Auxiliary grasping force of reload surface

Avoid the tissue slipping out during firing



## Safety block mechanism

Ensure the knife can be returned to the original position anytime in case of emergency



## Multiple reload options meet for different requirement of tissue thickness

Provide more options to suit for different tissues



## Adaptive angle adjustment

Adjust the reload angle freely

#### Single Use Endoscopic Linear Cutter Parameter Chart

Description	Item Code	Shaft Length	Item Code	Shaft Length	Item Code	Shaft Length
Short	SEC30S/SEFC30S	250mm	SEC45S/SEFC45S	265mm	SEC60S/SEFC60S	280mm
Standard	SEC30/SEFC30	310mm	SEC45/SEFC45	325mm	SEC60/SEFC60	340mm
Extra long	SEC30L/SEFC30L	410mm	SEC45L/SEFC45L	425mm	SEC60L/SEFC60L	440mm

### Single Use Endoscopic Linear Cutter Reload Parameter Chart

Item Code	Suture Length	Staple Height	Staple Formation Height
ECR30Y		2.0mm/2.0mm/2.0mm	0.75mm
ECR30W		2.5mm/2.5mm/2.5mm	1.0mm
ECR30B(H)	30mm	3.5mm/3.5mm/3.5mm	1.5mm
ECR30G(H)		4.2mm/4.2mm/4.2mm	2.3mm
ECR30A(H)		3.8mm/3.8mm/3.8mm	1.8mm
ECR45Y		2.0mm/2.0mm/2.0mm	0.75mm
ECR45W		2.5mm/2.5mm/2.5mm	1.0mm
ECR45B(H)	47mm	3.5mm/3.5mm/3.5mm	1.5mm
ECR45G(H)		4.2mm/4.2mm/4.2mm	2.3mm
ECR45A(H)		3.8mm/3.8mm/3.8mm	1.8mm
ECR60Y		2.0mm/2.0mm/2.0mm	0.75mm
ECR60W		2.5mm/2.5mm/2.5mm	1.0mm
ECR60B(H)	59mm	3.5mm/3.5mm/3.5mm	1.5mm
ECR60G(H)		4.2mm/4.2mm/4.2mm	2.3mm
ECR60A(H)		3.8mm/3.8mm/3.8mm	1.8mm

# Single Use Bladeless Trocar & Veress Needle

#### Wide product range

Standard trocar, pediatric trocar, extra-long trocar, thyroid trocar

#### One-way thread design

Smooth insertion and provide sufficient abdominal wall fixation during penetration

## High quality medical grade PC material

Solid and flexible obturator to ensure safe insertion

#### Abundant kit models

Multiple combinations, multiple options Applicable for different departments and surgical procedure







#### Patented bladeless olecranon design

Avoid risk of injury to internal organs, fascia suture is not required



#### Bivalved elastic sealing

Excellent leakage prevention
Wear-resisitance to support frequent
instrument exchanges



#### Universal seal design

Perfect wrapping instrument rod to provide better prevention of leakage. Three-layers structure and special materials to ensure gas proofness performance. Accommodate instruments from 5mm to 12mm, no need converter



#### 12mm optical trocar

Visual access to ensure the insertion control

## Extra-long trocar



Applicable for bariatric surgery and obese patients



Extra-long trocar models 5mm、10mm、12mm、15mm



Simplified valve body upper cover design Provide sufficient operation space

#### Standard Trocar & Extra-long Trocar & Veress Needle Parameter Chart

Item Code	Description	Specification
VT05	5mm Small head	5*92mm
VTB05	5mm	5*96mm
VTB10	10mm	10*96mm
VTB12	12mm	12*96mm
VTB15	15mm	15*110mm
VT05L	5mm extra-long	5*137mm
VTB10L	10mm extra-long	10*147mm
VTB12L	12mm extra-long	12*147mm
VTB15L	15mm extra-long	15*147mm
VTN12	12mm veress needle	120mm
VTN15	15mm veress needle	150mm

#### Trocar Kit Specification Chart

Name								
VTK02         3         2 pcs         1 pc (Small head)         2 pcs (Extra-long)         1 pc (Extra-long)           VTBK10-1         3         1 pc         2 pcs         1 pc           VTBK10-2         4         2 pcs         2 pcs         1 pc           VTBK10-2-1         4         2 pcs (Small head)         2 pcs         1 pc           VTBK10-3         2         1 pc         1 pc           VTBK10-4         3         2 pcs         1 pc         1 pc           VTBK10-4-1         3         2 pcs (Small head)         1 pc         1 pc           VTBK10-5         2         1 pc         1 pc         1 pc           VTBK10-6-1         2         1 pc         1 pc         1 pc           VTBK12-1         2         1 pc         1 pc         1 pc           VTBK12-2         4         2 pcs         1 pc         1 pc           VTBK12-3-1         4         2 pcs (Small head)         2 pcs         1 pc           VTBK12-4-1         3         1 pc (Small head)         2 pcs         1 pc           VTBK12-5-1         3         2 pcs (Small head)         1 pc         1 pc	Item Code	Qty for Kit	3mm	5mm	10mm	12mm	15mm	VTN12
VTK04         3         2 pcs (Extra-long)         1 pc (Extra-long)           VTBK10-1         3         1 pc         2 pcs         1 pc           VTBK10-2         4         2 pcs         2 pcs         1 pc           VTBK10-2-1         4         2 pcs (Small head)         2 pcs         1 pc           VTBK10-3         2         1 pc         1 pc           VTBK10-4         3         2 pcs         1 pc         1 pc           VTBK10-4-1         3         2 pcs (Small head)         1 pc         1 pc           VTBK10-5         2         1 pc         1 pc         1 pc           VTBK10-6-1         2         1 pc         1 pc         1 pc           VTBK12-1         2         2 pcs         1 pc           VTBK12-2         4         2 pcs         1 pc         1 pc           VTBK12-3-1         4         2 pcs (Small head)         2 pcs         1 pc           VTBK12-4-1         3         1 pc (Small head)         2 pcs         1 pc           VTBK12-5         3         2 pcs (Small head)         1 pc         1 pc           VTBK12-5-1         3         2 pcs (Small head)         1 pc         1 pc	VTK01	4		2 pcs (Extra-long)	1 pc (Extra-long)	1 pc (Extra-long)		
VTBK10-1         3         1 pc         2 pcs         1 pc           VTBK10-2         4         2 pcs         2 pcs         1 pc           VTBK10-2-1         4         2 pcs (Small head)         2 pcs         1 pc           VTBK10-3         2         1 pc         1 pc           VTBK10-4         3         2 pcs         1 pc         1 pc           VTBK10-4-1         3         2 pcs (Small head)         1 pc         1 pc           VTBK10-5         2         1 pc         1 pc         1 pc           VTBK10-6-1         2         1 pc         1 pc         1 pc           VTBK12-1         2         2 pcs         1 pc         1 pc           VTBK12-2         4         2 pcs         1 pc         1 pc           VTBK12-3-1         4         2 pcs         2 pcs         1 pc           VTBK12-4-1         3         1 pc (Small head)         2 pcs         1 pc           VTBK12-5         3         2 pcs         1 pc         1 pc           VTBK12-5-1         3         2 pcs (Small head)         1 pc         1 pc	VTK02	3	2 pcs	1 pc (Small head)				
VTBK10-2         4         2 pcs         2 pcs         1 pc           VTBK10-2-1         4         2 pcs (Small head)         2 pcs         1 pc           VTBK10-3         2         1 pc         1 pc           VTBK10-4         3         2 pcs         1 pc         1 pc           VTBK10-4-1         3         2 pcs (Small head)         1 pc         1 pc           VTBK10-5         2         1 pc         1 pc         1 pc           VTBK10-6-1         2         1 pc         1 pc         1 pc           VTBK12-1         2         2 pcs         1 pc         1 pc           VTBK12-2         4         2 pcs         1 pc         1 pc           VTBK12-3         4         2 pcs         2 pcs         1 pc           VTBK12-3-1         4         2 pcs (Small head)         2 pcs         1 pc           VTBK12-4-1         3         1 pc (Small head)         2 pcs         1 pc           VTBK12-5         3         2 pcs (Small head)         1 pc         1 pc           VTBK12-5-1         3         2 pcs (Small head)         1 pc         1 pc	VTK04	3		2 pcs (Extra-long)	1 pc (Extra-long)			
VTBK10-2-1       4       2 pcs (Small head)       2 pcs       1 pc         VTBK10-3       2       1 pc       1 pc         VTBK10-4       3       2 pcs       1 pc       1 pc         VTBK10-4-1       3       2 pcs (Small head)       1 pc       1 pc         VTBK10-5       2       1 pc       1 pc       1 pc         VTBK10-6-1       2       1 pc       1 pc       1 pc         VTBK12-1       2       2 pcs       1 pc       1 pc         VTBK12-2       4       2 pcs       1 pc       1 pc         VTBK12-3       4       2 pcs       2 pcs       1 pc         VTBK12-3-1       4       2 pcs (Small head)       2 pcs       1 pc         VTBK12-4-1       3       1 pc (Small head)       2 pcs       1 pc         VTBK12-5       3       2 pcs (Small head)       1 pc       1 pc         VTBK12-5-1       3       2 pcs (Small head)       1 pc       1 pc	VTBK10-1	3		1 pc	2 pcs			1 pc
VTBK10-3       2       1 pc       1 pc         VTBK10-4       3       2 pcs       1 pc       1 pc         VTBK10-4-1       3       2 pcs (Small head)       1 pc       1 pc         VTBK10-5       2       1 pc       1 pc       1 pc         VTBK10-6-1       2       1 pc       1 pc       1 pc         VTBK12-1       2       2 pcs       1 pc         VTBK12-2       4       2 pcs       1 pc         VTBK12-3       4       2 pcs       2 pcs       1 pc         VTBK12-3-1       4       2 pcs (Small head)       2 pcs       1 pc         VTBK12-4-1       3       1 pc (Small head)       2 pcs       1 pc         VTBK12-5       3       2 pcs       1 pc       1 pc         VTBK12-5-1       3       2 pcs (Small head)       1 pc       1 pc	VTBK10-2	4		2 pcs	2 pcs			1 pc
VTBK10-4       3       2 pcs       1 pc       1 pc         VTBK10-4-1       3       2 pcs (Small head)       1 pc       1 pc         VTBK10-5       2       1 pc       1 pc       1 pc         VTBK10-6-1       2       1 pc       1 pc       1 pc         VTBK12-1       2       2 pcs       1 pc       1 pc         VTBK12-2       4       2 pcs       1 pc       1 pc         VTBK12-3       4       2 pcs       2 pcs       1 pc         VTBK12-3-1       4       2 pcs (Small head)       2 pcs       1 pc         VTBK12-4-1       3       1 pc (Small head)       2 pcs       1 pc         VTBK12-5       3       2 pcs       1 pc       1 pc         VTBK12-5-1       3       2 pcs (Small head)       1 pc       1 pc	VTBK10-2-1	4		2 pcs (Small head)	2 pcs			1 pc
VTBK10-4-1       3       2 pcs (Small head)       1 pc       1 pc         VTBK10-5       2       1 pc       1 pc       1 pc         VTBK10-6-1       2       1 pc       1 pc       1 pc         VTBK12-1       2       2 pcs       1 pc         VTBK12-2       4       2 pcs       1 pc         VTBK12-3       4       2 pcs       2 pcs       1 pc         VTBK12-3-1       4       2 pcs (Small head)       2 pcs       1 pc         VTBK12-4-1       3       1 pc (Small head)       2 pcs       1 pc         VTBK12-5       3       2 pcs       1 pc       1 pc         VTBK12-5-1       3       2 pcs (Small head)       1 pc       1 pc	VTBK10-3	2		1 pc	1 pc			
VTBK10-5       2       1 pc       1 pc       1 pc         VTBK10-6-1       2       1 pc       1 pc       1 pc         VTBK12-1       2       2 pcs       1 pc       1 pc         VTBK12-2       4       2 pcs       1 pc       1 pc         VTBK12-3       4       2 pcs       2 pcs       1 pc         VTBK12-3-1       4       2 pcs (Small head)       2 pcs       1 pc         VTBK12-4-1       3       1 pc (Small head)       2 pcs       1 pc         VTBK12-5       3       2 pcs       1 pc       1 pc         VTBK12-5-1       3       2 pcs (Small head)       1 pc       1 pc	VTBK10-4	3		2 pcs	1 pc			1 pc
VTBK10-6-1       2       1 pc       1 pc       1 pc         VTBK12-1       2       2 pcs       1 pc       1 pc         VTBK12-2       4       2 pcs       1 pc       1 pc         VTBK12-3       4       2 pcs       2 pcs       1 pc         VTBK12-3-1       4       2 pcs (Small head)       2 pcs       1 pc         VTBK12-4-1       3       1 pc (Small head)       2 pcs       1 pc         VTBK12-5       3       2 pcs       1 pc       1 pc         VTBK12-5-1       3       2 pcs (Small head)       1 pc       1 pc	VTBK10-4-1	3		2 pcs (Small head)	1 pc			1 pc
VTBK12-1       2       2 pcs       1 pc         VTBK12-2       4       2 pcs       1 pc         VTBK12-3       4       2 pcs       2 pcs       1 pc         VTBK12-3-1       4       2 pcs (Small head)       2 pcs       1 pc         VTBK12-4-1       3       1 pc (Small head)       2 pcs       1 pc         VTBK12-5       3       2 pcs       1 pc       1 pc         VTBK12-5-1       3       2 pcs (Small head)       1 pc       1 pc	VTBK10-5	2			1 pc			1 pc
VTBK12-2       4       2 pcs       1 pc       1 pc         VTBK12-3       4       2 pcs       2 pcs       1 pc         VTBK12-3-1       4       2 pcs (Small head)       2 pcs       1 pc         VTBK12-4-1       3       1 pc (Small head)       2 pcs       1 pc         VTBK12-5       3       2 pcs       1 pc       1 pc         VTBK12-5-1       3       2 pcs (Small head)       1 pc       1 pc	VTBK10-6-1	2		1 pc	1 pc			1 pc
VTBK12-3       4       2 pcs       2 pcs       1 pc         VTBK12-3-1       4       2 pcs (Small head)       2 pcs       1 pc         VTBK12-4-1       3       1 pc (Small head)       2 pcs       1 pc         VTBK12-5       3       2 pcs       1 pc       1 pc         VTBK12-5-1       3       2 pcs (Small head)       1 pc       1 pc	VTBK12-1	2				2 pcs		1 pc
VTBK12-3-1       4       2 pcs (Small head)       2 pcs       1 pc         VTBK12-4-1       3       1 pc (Small head)       2 pcs       1 pc         VTBK12-5       3       2 pcs       1 pc       1 pc         VTBK12-5-1       3       2 pcs (Small head)       1 pc       1 pc	VTBK12-2	4		2 pcs	1 pc	1 pc		
VTBK12-4-1         3         1 pc (Small head)         2 pcs         1 pc           VTBK12-5         3         2 pcs         1 pc         1 pc           VTBK12-5-1         3         2 pcs (Small head)         1 pc         1 pc	VTBK12-3	4		2 pcs		2 pcs		1 pc
VTBK12-5         3         2 pcs         1 pc         1 pc           VTBK12-5-1         3         2 pcs (Small head)         1 pc         1 pc	VTBK12-3-1	4		2 pcs (Small head)		2 pcs		1 pc
VTBK12-5-1 3 2 pcs (Small head) 1 pc 1 pc	VTBK12-4-1	3		1 pc (Small head)		2 pcs		1 pc
	VTBK12-5	3		2 pcs		1 pc		1 pc
VTBK15-1 5 2 pcs 2 pcs 1 pc 1 pc	VTBK12-5-1	3		2 pcs (Small head)		1 pc		1 pc
	VTBK15-1	5		2 pcs	2 pcs		1 pc	1 pc
VTBK15-1-1 5 2 pcs (Small head) 2 pcs 1 pc 1 pc	VTBK15-1-1	5		2 pcs (Small head)	2 pcs		1 pc	1 pc

Note: For all kits, if 2pcs of same code included in one box, one of them is sleeve.







Patented small head design

Provide sufficient operation space





#### Cupped oriented design

Easy for instruments to insert and prevent from being bent and deformed





#### Thread design


Strengthen abdominal wall fixation, no need suture fixation





#### Stainless steel obturator

Provide reliable insertion




#### Minimally Invasive Surgery

The diameter of crochet is 1.2mm, can reduce one incision and one sleeve to lower operating expense.

#### Save Operation Time

Unique double crochet design, no need to completely withdraw from the body to insert again, reducing trauma and effectively saving time.

Item Code	Descripiton	Specification
VT03	3mm	3*52mm
VT05	Small head 5mm	5*92mm
VTK02	VTK02 Trocar Kit	Trocar 3mm*1 、small head trocar 5mm*1、trocar sleeve 3mm*1







#### Independent visual stripper

Ergonomic handle, easier to peel
Visible peeling under lens guidance to ensure safe operation





#### Widened thread

Meet the special requirement of thyroid surgery Provide better fixation and reduce the damage





#### Small head design

More suitable for approaching from areola and underarm Easy to operate

Item Code	Descripiton	Specification
VTK01	VTK01 Trocar Kit	Trocar 5mm*1, trocar 10mm*1, trocar sleeve 5mm*1, obturator 12mm*1
VTK04	VTK04 Trocar Kit	Trocar extra-long 5mm*1, trocar extra-long 10mm*1,trocar extra-long sleeve 5mm*1





One-piece anvil provide strong support

Ensure the "B" shape formation



5.8mm extra-thin anvil

Easy to be placed into purse-string



Variable speed tail

Protect stomas while removing stapler, uniformly squeeze tissue when closing



Single edged circular knife

ensure sharper cutting



Indicated window design
Tissue compression scale

Fire until the indicators within the green zone, ensure the safe suture



Patented clamping slot design

Ensure better connection and smooth placement, reduces the damage when placing into the tissue



### Single Use Circular Stapler Parameter Chart

Figure	Item Code	Outer Dia	Inner Dia	Staple Qty	Staple Height	Staple Formation Height
	(E)SCS21	21mm	13mm	18	4.5mm	1.9mm
	SCSA21	21mm	13mm	18	4.5mm	1.9mm
	(E)SCS23	23mm	15mm	20	4.5mm	1.9mm
	SCSA23	23mm	15mm	20	4.5mm	1.9mm
	(E)SCS25	25mm	17mm	22	4.8mm	2.0mm
	SCSA25	25mm	17mm	22	4.8mm	2.0mm
	(E)SCS27	27mm	19mm	24	4.8mm	2.0mm
	SCSA27	27mm	19mm	24	4.8mm	2.0mm
	(E)SCS29	29mm	21mm	26	4.8mm	2.0mm
	SCSA29	29mm	21mm	26	4.8mm	2.0mm
	(E)SCS31	31mm	23mm	28	5.0mm	2.2mm
	SCSA31	31mm	23mm	28	5.0mm	2.2mm
	(E)SCS33	33mm	25mm	32	5.2mm	2.2mm
	SCSA33	33mm	25mm	32	5.2mm	2.2mm

## Single Use Linear Cutter





#### Auto bumped knife protector

Avoid misfiring from empty cartridge and protect the surgeons and nurse in operation



#### Floating reload

Cannot be fired if the instrument don't close properly



Adjustable bidirectional firing knob for easy control of firing

Applicable for different operating habits



#### Reinforced instrument body

Ensure good staple formation after multiple firings

Figure	Item Code	Suture Length	Staple Qty	Staple Height	Formation Height
	SLC60B	63mm	64	3.8mm	1.5mm
	SLC80B	83mm	84	3.8mm	1.5mm
In .	SLC100B	103mm	104	3.8mm	1.5mm
Twee	SLC60G	63mm	64	4.8mm	2.0mm
	SLC80G	83mm	84	4.8mm	2.0mm
	SLC100G	103mm	104	4.8mm	2.0mm
	LCR60B	63mm	64	3.8mm	1.5mm
	LCR80B	83mm	84	3.8mm	1.5mm
	LCR100B	103mm	104	3.8mm	1.5mm
	LCR60G	63mm	64	4.8mm	2.0mm
The second second	LCR80G	83mm	84	4.8mm	2.0mm
	LCR100G	103mm	104	4.8mm	2.0mm







Patented 3 rows of staple line design

Excellent hemostatic effect, better blood supply

Accurate fixe



#### Alignment pin design

Accurate fixation, ensure good staple formation



#### Segmental firing

More safe in second close

Item Code	Suture Length	Staple Line	Staple Qty	Staple Height	Formation Height	Color	Color
PSLS30B2/G2	30mm	2	13	3.8/4.5mm	1.5/2.0mm	Blue/Green	••
PSLS45B2/G2	45mm	2	17	3.8/4.5mm	1.5/2.0mm	Blue/Green	••
PSLS60B2/G2	60mm	2	23	3.8/4.5mm	1.5/2.0mm	Blue/Green	••
PSLS100B2/G2	100mm	2	41	3.8/4.5mm	1.5/2.0mm	Blue/Green	••
PSLS30W3/B3/G3	30mm	3	20	2.5/3.8/4.5mm	1.0/1.5/2.0mm	White/Blue/Green	$\bigcirc \bullet \bullet$
PSLS45B3/G3	45mm	3	26	3.8/4.5mm	1.5/2.0mm	Blue/Green	••
PSLS60B3/G3	60mm	3	35	3.8/4.5mm	1.5/2.0mm	Blue/Green	••

## Single Use Hemorrhoidal Stapler & Accessories





24.1ml ultra-large capacity One-piece central rod design
Accommodate more tissue Ensure excellent concentricity



Ensure excellent concentricity and verticality. Ensure consistent of the staple formation



53mm central rod design
Providing better
operation view



Single edged circular knife
Ensure sharper cutting



Obvious feedback when firing
High quality gasket provide
obvious audible and tactile
feedback



Variable speed tail

Protect stomas while removing stapler, uniformly squeeze tissue when closing stapler



More staples
Providing excellent
hemostatic effect



Tissue compression scale

Fire until the indicators within the green zone, ensure the safe suture

Indicate window design

Item Code	Outer Dia	Inner Dia	Staple Qty	Staple Height	Formation Height	Color	Color
SHS32	32mm	24mm	32	4.0mm	2.0mm	Green	
SHS34	34mm	26mm	34	4.0mm	2.0mm	White	





Circular Anal Dilator



Suture Threader

## Single Use Anal Hemorrhoid Ligator Device & Accessories



#### Characteristics

- Disposable device, prevent cross infection.
- Mild pain, barely bleeding.
- No scar after operation, no stricture of anus or any other complications, convenient for the follow-up treatment.
- Just short stay or even no hospitalization, faster recovery, cost-effective, little impact on patients' life.
- Simple operation, 5–10 minutes, saving the cost.

Parameter Item Code	Outer Dia of Air Suction Pipe	Tolerance	Inner Dia of Loop	Outer Dia of Loop	Tolerance	Inner Dia of Anoscope	Outer Dia of Barrel	Tolerance
VMTZQ-12x1								
VMTZQ-12x2						20	40 5	. 4 5
VMTZQ-12x3			4.5	0	105	20mm	19.5mm	± 1.5mm
VMTZQ-12x4	40	I 4 Frans						
VMTQ-12x1	12mm ± 1.5mm	1.5mm	6mm	± 0.5mm				
VMTQ-12x2						1	/	/
VMTQ-12x3								
VMTQ-12x4								

● 1-2 pcs of ligation loop subassembly will be included in each ligator device package.



## Single Use Circumcision Stapler



#### Indication

Clinical operation for foreskin cutting and suturing.

#### Product Characteristics

- Humanized design of balanus cover, protect the prepuce bridle effective, Prevent excessive excision of the prepuce bridle.
- Accurate control of the scope of excision, regular cutting edge, beautiful, small scar, fast recovery.
- Orange safety block makes it safer, prevent mistaken firing, more convenient operation.
- Different types of rings are available according to individual differences.
- Disposable sterilizing product. Prevent various infections, no complications after operation.

#### Single Use Circumcision Stapler Parameter Chart

Item Code	Outer Dia of Blade	Staple Height	Staple Qty
VMBQA-11	11mm		9
VMBQA-13	13mm		10
VMBQA-15	15mm		11
VMBQA-17	16.5mm		12
VMBQA-19	18mm	2.5mm	13
VMBQA-22	22mm		16
VMBQA-25	25mm		18
VMBQA-29	29mm		20
VMBQA-34	34mm		24

## Single Use Endoscopic Retriever











6 sizes provide more options for different sizes specimen

Medical grade nontoxic TPU material Retriever bag could open automatically and smoothly easy to operation

Excellent performance of elasticity and intensity to ensure reliable and convenient operation

Figure	Item Code	Bag Height	Bag Diameter	Bag Capacity
	VER09	150mm	90mm	150ml
	VER12	200mm	120mm	500ml
	VERB60	175mm	60mm	100ml
<del></del>	VERB80	175mm	80mm	300ml
	VERB100 175mm 100m	100mm	550ml	
	VERB130	175mm	130mm	750ml


## Single Use Suction & Irrigation Cannula



Suction & Irrigation Cannula				
Item Code				
VESI05				
Description				
Suction & Irrigation Cannula				
Details				
ф 5mm × 320mm				

- Visualized suction and irrigating button.
- Ergonomic design handle, easy to operate.

## Single Use Skin Stapler






Rapid and convenient operation to benefit for wound healing

Single Use device, prevent cross infection

Regular/wide staple options, applicable for different wounds

## Single Use Wound Protector and Retractor

#### Open incision without damage





Squeeze cutting edge to prevent leakage •

#### Protect incision

- Prevent tumour implantation metastasis in incision
  - Prevent infection of incision •
- Prevent incision from burn due to misoperation •

Item Code	Sizes	Channel Diameter	Applicative Incision Size
VIPA	60/70-60/150	6cm	2-4cm
	80/90-80/150	8cm	2.5-6cm
	120/130-120/250	12cm	5-9cm
	150/160-150/250	15cm	7-12cm
	180/190-180/250	18cm	9-14cm
	220/230-220/250	22cm	11-17cm
	270/280-270/250	27cm	20-22cm
	320/330-320/250	32cm	22-27cm
VIPB	40/30-15/20	1.5cm	2cm
	50/40-25/25	2.5cm	3cm
	70/70-35/25	3.5cm	4cm
	120/100-60/40	6.0cm	6cm



### **CERTIFICATES**













